


Lifetime Limited Warranty

If you have a question about your MOBB Healthcare device or this warranty, please contact an authorized MOBB Healthcare dealer.

MOBB Healthcare warrants to the original consumer purchaser for as long as they own their MOBB product that it will be free from defects in material and workmanship under normal use. This warranty is only applicable to products used in residential applications within the United States and Canada and is not transferable to subsequent owners. Should defects in material or workmanship exist during the applicable warranty period in any MOBB Healthcare Product, MOBB Healthcare will elect at its discretion to either repair or replace the defective material. This warranty covers only parts and materials of the products supplied by MOBB Healthcare. Not covered under this warranty are wear and tear items or abuse of the product in a manner that was not intended for the product to be used. This warranty shall become void if the product is in any way modified, improperly adjusted or abused.

MHKWHF



990000000042



MOBB Healthcare Ltd.

116A Viceroy Rd. Unit 5-6
Concord, Ontario, Canada
L4K 2M1
www.mobbhhc.com



Knee Walker

300 lbs weight capacity



CAUTION

TO REDUCE THE RISK OF ACCIDENT:

- Before using the Knee Walker, you should be trained by a healthcare professional or dealer.
- Engage the brake lock and practice bending, reaching and transfers on and off the Knee Walker. DO NOT use the Knee Walker by yourself until you have been properly trained, you have practiced with someone and you feel confident you can maneuver by yourself.
- The use of the Knee Walker is individual to the person. You should develop your own methods for use based on your personal level of function and ability.
- NEVER attempt a maneuver that has not been practiced.
- Be aware of your surroundings when operating the Knee Walker. Look for hazards and avoid them.



operating the knee walker


The Knee Walker is a mobility device that provides strength and comfort to anyone who may have an injury below the knee. It allows an even distribution of body weight by supporting half of the weight on the device and the other half on the non-injured leg. The device has a cushioned platform to support both the injured leg and ankle to make certain that no weight will be placed on the lower leg.

The device is designed for easy steering and stopping. Brakes are also featured on this device for controlled movement and safety.

- Placement of injured leg on pad should be centered side to side and positioned forward to cover the full length of pad.
- With injured leg on pad, stand as straight as possible (adjust the height of the Knee Walker as necessary).
- With injured foot pointing down, the propulsion leg should be maintained as close as possible to the pad.
- While on your toes, start out with small steps with propulsion leg. Concentrate on keeping that leg next to the pad. With practice, you should be able to move forward without deviating to the opposite side.

assembly instructions

1. Straighten the tiller to the upright position. Slide locking mechanism to the side to allow tiller to reach full upright position. Release locking mechanism to lock tiller in position.



2. Slide locking mechanism into groove on tiller, depress locking handle to secure tiller lock.



3. Insert pad post into receptacle on frame. Secure pad in place by inserting hexagonal screw into hole and tighten with star shaped rotary knob.



4. Attach basket by sliding basket receptacles down over hooks on the front of the tiller.



using the hand brake

The Knee Walker's brake is similar to a bicycle and also has a locking feature.


Brake Operation as shown in figure below:

- **D1:** Normal use state, brake disengaged.
- **D2:** Engage the brake as needed in the use process of the Walker
- **D3:** Press down the brake lever to keep braking state



adjusting the hand brake

- For the brake adjustment turn the adjuster counterclockwise to tighten the brake, or clockwise to loosen the brake. Keep the cable adjuster nut against the brake.



adjusting handle height

To adjust the handle height, loosen the adjustment screw on the handle by turning it counter clockwise until it can be removed. Adjust the hand grip to the desired height and reinstall the adjustment screw by turning it clockwise until it is tightened.

adjusting pad height

To adjust the pad height loosen rotary knob under pad (shown in step 3 of assembly), and then remove hexagonal screw. Adjust pad to the desired height then reinstall hexagonal screw and tighten with star shaped rotary knob.

to fold tiller

To fold the tiller down, release the tiller locking handle (as shown in step 2 of assembly) and slide assembly to the side until the tiller is disengaged. Fold the tiller down until it rests on the pad and release the locking mechanism.

additional cautions

Always consult with your health care practitioner or authorized dealer on the proper fitment and usage of your MOBB Knee Walker. The instructions here are a guide and should be used in conjunction with professional assistance.

- The MOBB Knee Walker is designed to replace the use of crutches for long term lower limb injury recovery. It should be used at slow and controlled speeds. If not, tipping can occur with too much speed while turning. While becoming accustomed to it's operation in the learning phase, apply the brake cautiously when starting, stopping and turning.
- Always check the fasteners on the tiller, wheels, brakes, and pad for tightness before use.
- Do not sit, stand or lean on steering mechanism.
- Steering mechanism is not weight bearing.
- Operate at your own risk.