

Gravity 7 Long Term Care Pressure Redistribution Mattress

HCPCS: E0184-1,E0184-1





FIGURE A

- High quality therapeutic foam mattress provides pressure redistribution and shear/friction reduction
- Drive's most popular therapeutic foam mattress maximizes value and comfort with 3 full layers of latex free foam
- Full width top layer uses unique die-cutting to optimize zoned pressure redistribution
- Bottom layer provides full length horizontally scored articulation cuts extending the durability and life of the mattress
- High density foam perimeters provide solid edge support and facilitate safe patient transfers
- Available with or without raised side rails
- Visco foam heel section offloads pressure from vulnerable heels
- Removable multi-stretch cover is fluid resistant, low shear and vapor permeable helping protect the skin from friction and moisture
- Meets flammability standard CFR 1632 and CFR 1633

PRODUCTS	Description
15770	GRAVITY 7 MATTRESS, 80",1/EA
15777	GRAVITY 7 MATTRESS W/ELEV PERIM,80",1/EA
15784	GRAVITY 7 MATTRESS 84",1/EA
15785	GRAVITY 7 MATTRESS W/ELEV PERIM 84",1/EA
15876	GRAVITY 7 MATTRESS 76", 1/EA
15877	GRAVITY 7 MATTRESS W/ELEV PERIM 76",1/EA

SPECIFICATIONS	Description
Warranty	7 Year Limited Foam 3 Year Cover
Weight Capacity	350 lb
Side Rails	Yes
Cover Material	Stretch, Fluid-Resistant Nylon Top Waterproof Vinyl Bottom
Mattress Layers of Foam	3
CA Prop 65	[WARNING: This product can expose you to chemicals that are known to the State of California to cause cancer and birth defects or other reproductive harm. For more information go to www.P65Warnings.ca.gov .]

ACCESSORIES	Description
15770C	Mattress Cover 15770C,1/ea, EA
15777C	Mattress Cover 15777C,1/ea, EA
15784C	Mattress Cover 15784C,1/ea, EA
15785C	Mattress Cover 15785C,1/ea, EA
15877C	Mattress Cover 15877,1/ea, EA
15876C	Mattress Cover 15876,1/ea, EA

ACCESSORIES**Description**

14427QF	Quick Fill Pump w/adapters, 1 EA/12 CV, EA
---------	--