



Pride Mobility Products Co.

5096 South Service Rd, Beamsville, ON, L3J 1V4

Phone: 888-570-1113 | Fax: 866-514-1303 | Email: productorders@pridemobility.com

Q6 Edge® 3 Stretto Order Form with Contoured Seating

Account Number: _____ Date: _____

Provider Name: _____

Contact: _____

Phone: _____ Fax: _____

Email: _____

PO Number: _____

Marked for: _____

Ship to Address: _____

City: _____ Prov: _____ PC: _____



INTRODUCTION

This form is interactive when viewed with [Adobe Acrobat Reader](#) and may not function correctly if opened with applications other than Acrobat. Complete the form by placing checks in the desired boxes and provide information in the interactive fields. Buttons shown at the bottom of the form may be utilized to print or submit the order form through a desktop email application. To email via a web-based application, please 'Save As' and attach the PDF to your email.


Thank you for choosing Quantum!

PATIENT INFORMATION

Weight: _____ Height: _____ A. Shoulder Width: _____ D. Max Sitting Height: _____ G. Thigh Depth: _____

B. Chest Width: _____ E. Shoulder Height: _____ H. Lower Leg Length: _____

C. Hip Width: _____ F. Axilla Height: _____ I. Elbow Height: _____



1. BASE MODEL

REQUIRED. Select a group 3 power base. 300 lb. weight capacity with standard 6mph programmable motors and off board charger. Dimensions of base without legrests: 20.5"Wx35.5"L

Quantum Q6 Edge 3 Stretto\$8,450

Part: Q6 EDGE 3 STRETTO

2. DRIVE MOTORS

Optional. Requires Q-Logic 3 electronics. The standard 6 mph motor can be upgraded to Accu-Trac motors that delivers enhanced traction and tracking, especially with switched drive controls.

Accu-Trac Motor Technology\$995
Part: XRF000268. HCPCS: K0108

3. BASE COLOR

REQUIRED. Select one shroud color. Colors shown here may differ from actual product depending on monitor/printer calibration used.

3.1a. Standard Shrouds

 Red Corvette Part: DGN1900199-PRD1	 Deep Purple Part: DGN1900199-PPR5	 Midnight Blue Part: DGN1900199-PBU10	 Green Machine Part: DGN1900199-PGN11
 Back in Black Part: DGN1900199-PBK1	 Raspberry Beret Part: DGN1900199-PPK4	 Electric Blue Part: DGN1900199-PBU9	 Military Green Part: DGN1900199-PGN8
 Night Train Part: DGN1900199-PBK3	 Pretty in Pink Part: DGN1900199-PPK3	 Blue Jean Part: DGN1900199-PBU11	
 Grey Street Part: DGN1900199-PSL3	 Evolution Orange Part: DGN1900199-POR7	 Lemon Crush Part: DGN1900199-PYL4	

3.1b. Hydrographic Shrouds

Military Digital Camo.....\$325
Part: DGN1801101-WCM7



Faded Glory.....\$325
Part: DGN1801101-WFL2

4. POWER ADJUSTABLE SEAT HEIGHT

Optional. The seat to floor height will be increased. Lift/elevate will be operated through the joystick. Power elevating seat includes a 1/4 speed inhibit. ^Please verify the appropriate billing code to be used for this item with all third-party payors.

12" Power Lift for Captain Seat with iLevel..... \$2,995

Part: PTOINDV3857. HCPCS: E2298/E2300^

Available with swing away legrests or manually adjustable center mount foot platform.

5. BASE MOUNTED LIGHTS

Optional. Select one option.

Switch Fender Lights..... No Charge

Part: DGN1810041

Fender lights operate by pressing included single switch button. Fender lighting is not available with Single Actuator Control through Toggle Box. Fender lights come standard with track mount hardware.



Fender Lights

6. QUANTUM BACKUP CAMERA

Optional. The backup camera may interfere with armrest mounted accessories such as media device holder.

Quantum Backup Camera\$475
Part: MEC1910365

1. Optional Camera Mounting Bracket

Backup Camera Captain's Headrest Bracket\$35
Part: MEC1901231

Backup Camera Track Mount Bracket\$35
Part: MEC2104227

2. Optional Display Mounting Bracket

Backup Camera Q-Logic Joystick Mount Bracket\$145
Part: MEC2106123



Backup Camera

7. DRIVE WHEELS

REQUIRED. Select a type of drive wheel. Accent colors are available for 14" tires.

7.1. Drive Wheels

3"x14" Pneumatic Tires..... No Charge
Part: WHL1713925

3"x14" Flat Free Tires..... No Charge
Part: WHL1704058

4"x14" Knobby, Flat Free Black Tires..... \$120
Part: XRF000368

2"x12.5" Pneumatic Tires..... No Charge
Part: WHL1810666

2"x12.5" Flat Free Tires..... No Charge
Part: WHL1810665

4"x14" Knobby, Pneumatic Black Tires..... \$120
Part: XRF000367

Drive wheels section continued on following page.

7. DRIVE WHEELS CONTINUED**7.2. Drive Wheel Accent Color for 14" Tires**

Not available with 12.5" tires.

Back in Black..... No Charge
Part: DGN164245-PBK1 x2Red Corvette..... No Charge
Part: DGN164245-PRD1 x2Night Train..... No Charge
Part: DGN164245-PBK3 x2Grey Street..... No Charge
Part: DGN164245-PSL3 x2Deep Purple..... No Charge
Part: DGN164245-PPR5 x2Raspberry Beret..... No Charge
Part: DGN164245-PPK4 x2Pretty in Pink..... No Charge
Part: DGN164245-PPK3 x2Evolution Orange.... No Charge
Part: DGN164245-POR7 x2Midnight Blue No Charge
Part: DGN164245-PBU10 x2Electric Blue No Charge
Part: DGN164245-PBU9 x2Blue Jean No Charge
Part: DGN164245-PBU11 x2Lemon Crush No Charge
Part: DGN164245-PYL4 x2Green Machine No Charge
Part: DGN164245-PGN11 x2Military Green No Charge
Part: DGN164245-PGN8 x2Military Digital Camo... \$50 ea.
Part: DGN164245-WCM7 x2Faded Glory \$50 ea.
Part: DGN164245-WFL2 x2**8. BATTERIES****Optional.** The powerchair requires 2 batteries and will be installed if ordered.NF-22, Gel Cell, Interceptor 2 Battery \$300 ea.
Part: BAT2004376. HCPCS: E2361
Uses 8 AMP chargerGel Cell, 40AH, Interceptor 2 Battery..... \$200 ea.
Part: BAT2100316. HCPCS: K0108
Uses 5 AMP charger**9. PATIENT'S WEIGHT****REQUIRED.** The patient's weight is required to quote and order the power chair.Up to 150 lb User Weight
Part: PTOINDV3977151-250 lb User Weight
Part: PTOINDV3985250-300 lb User Weight
Part: PTOINDV3984**10. EXTENDED SEAT TOWERS****Optional.**Extended Seat Towers No Charge
Part: MEC2105559**11. ELECTRONICS****REQUIRED.** Select one type of drive control from the following sections. The No Joystick option is used to omit the joystick and utilize another control, such as a specialty control found on the [Alternative Controls and Electronics Order Form](#).**11.1a. Non-Expandable Controls**

Select a joystick below or omit selection if another controller type is desired. Non-expandable does not allow supplementary drive controllers.

QLNES 6 Key Joystick Standard
Part: CTL2205379
Controls up to 2 actuators.QLNES 10 Key Joystick No Charge
Part: CTL2205380
Controls up to 2 actuators through joystick and seat mounted, LED, full light package.

QLNES Joystick

11.1b. Expandable Controls

Select expandable controls or omit selection if another controller type is desired. Expandable controls allow multiple or specialized controllers.

2. Q-Logic 3 Joystick

Q-Logic 3 EX Joystick includes XLR USB charger.

Q-Logic 3 EX Joystick \$900
Part: CTL167724Q-Logic 3e Joystick \$620
Part: CTL1910600
Can operate full light package. Not available with alternative or specialty controls.Q-Logic 3 EX Joystick for Full Lights \$900
Part: CTL167725No Joystick, Setup for Q-Logic 3
Part: XRF000010
Requires selection of Enhanced Display

Q-Logic 3 EX Joystick

11.2. Q-Logic Magnetic Charger

Available with QLNES. Included with Q-Logic 3e Joystick.

Magnetic to XLR Charger Adapter with Wall Holder \$100
Part: ELE2103635

Electronics section continued on following page.

11. ELECTRONICS CONTINUED**11.3. Q-Logic 3 Joystick Armor**

REQUIRED with Q-Logic 3 EX Joystick selections above. Not available with Q-Logic 3e Joystick. Reference colors shown on page 2. *Night train joystick armor will be the black, soft touch.

Black, Soft TouchStandard
Part: DGN124307

Also matches Night Train shrouds.

Red Corvette No Charge
Part: DGN124336-PRD1

Back in Black No Charge
Part: DGN124336-PBK1

Grey Street No Charge
Part: DGN124336-PSL3

Deep Purple No Charge
Part: DGN124336-PPR5

Raspberry Beret No Charge
Part: DGN124336-PPK4

Pretty in Pink No Charge
Part: DGN124336-PPK3

Evolution Orange No Charge
Part: DGN124336-POR7

Midnight Blue No Charge
Part: DGN124336-PBU10

Electric Blue No Charge
Part: DGN124336-PBU9

Blue Jean No Charge
Part: DGN124336-PBU11

Lemon Crush No Charge
Part: DGN124336-PYL4

Green Machine No Charge
Part: DGN124336-PGN11

Military Green No Charge
Part: DGN124336-PGN8

Military Digital Camo No Charge
Part: DGN124336-WCM7

Faded Glory No Charge
Part: DGN124336-WFL2

12. JOYSTICK HANDLES

Optional. The optional handle will be mounted to the joystick and the standard handle will be shipped with the unit.

Mushroom Handle (B) \$120
Part: ACCASMB1860.

Chin Cup (C) \$85
Part: ACCASMB186190

T-Bar (D) \$160
Part: ACCASMB1858.

Medium Ball Handle (E) \$70
Part: ACCASMB1859.

Soft Ball Handle (F) \$120
Part: ACCASMB1856.

Extended I-Handle (G) \$140
Part: ACCASMB1857.

Bodypoint 3" Goal Post Handle \$180
Part: PLSKNCP1057.

Bodypoint 4" Goal Post Handle \$180
Part: PLSKNCP1058.



Joystick Handles

13. JOYSTICK MOUNTING BRACKETS

Optional. Select a mount to be used with a joystick selected in previous section.

13.1a. Fixed, In-line Joystick Bracket

Fixed, In-line, Bracket for Q-Logic 3 Standard
Part: FRMASMB10617

13.1b. Swing-Away (Retractable) Joystick Bracket

Allows the joystick to retract to the side of the armrest while maintaining the same joystick orientation.

Swing-away, Q-Logic 3 \$265
Part: ACC1908434. HCPCS: E1028

13.1c. Basic Swing Away Joystick Mount

Basic Swing Away, Q-Logic, Left Side \$265
Part: ACC2201041. HCPCS: E1028

Basic Swing Away, Q-Logic, Right Side \$265
Part: ACC2201042. HCPCS: E1028

13.1d. Specialty Joystick Mounts

Select a bracket and side to mount it.

Retract4 Mount \$450
Part: XRF000003. HCPCS: E1028


Full 180° of movement that translates into multi-position use, straight, right offset, left offset, or fold away. Joystick orientation varies as bracket is adjusted.

Multi-Axis Mount \$100
Part: XRF000004
Allows height adjustment of joystick.


13.2. Joystick Mount Location

Left Side Mount
Part: PTOINDV1225

Right Side Mount
Part: PTOINDV1226



Swing Away Joystick



Retract4



Multi-Axis Joystick Mount

14. CONTOURED SEATING**Optional.** Select only one seating type.**14.1. Occupied Transit, WC-19****Optional.**

WC-19 Occupied Transit Option for Contoured Seat \$400

Part: KIT122609

To determine if your chair is WC-19 compliant (occupied), please refer to this website: <https://wc-transportation-safety.umtri.umich.edu/>**14.2a. Q-Captain Seat**

16"W x 16"D No Charge

Part: SET2005145

18"W x 16"D No Charge

Part: SET2005149

20"W x 16"D No Charge

Part: SET2005152

16"W x 18"D No Charge

Part: SET2005146

18"W x 18"D No Charge

Part: SET2005150

20"W x 18"D No Charge

Part: SET2005153

16"W x 20"D No Charge

Part: SET2005147

18"W x 20"D No Charge

Part: SET2005151

20"W x 20"D No Charge

Part: SET2005154

14.2b. Q-Captain Seat with Solid Seat Pan

16"W x 16"D No Charge

Part: SET2005394

18"W x 16"D No Charge

Part: SET2005397

20"W x 16"D No Charge

Part: SET2005400

16"W x 18"D No Charge

Part: SET2005395

18"W x 18"D No Charge

Part: SET2005398

20"W x 18"D No Charge

Part: SET2005401

16"W x 20"D No Charge

Part: SET2005396

18"W x 20"D No Charge

Part: SET2005399

20"W x 20"D No Charge

Part: SET2005402

14.3. Seat to Floor Height**REQUIRED.**

Lowest Seat to Floor

Part: PTOINDV3558

Middle Seat to Floor

Part: PTOINDV3559


Highest Seat to Floor

Part: PTOINDV3560

	Cushion	SSP
High	22.5"	18.5"
Mid	21.5"	17.5"
Low	20.5"	16.5"
Q-Captain Seat to Floor Height		

14.4. Stealth Cushion

Optional. Specific cushion/component sizes available, part numbers, weight capacities, MSRP, HCPCS code, and special order cushions will be found on the respective complete [Stealth Cushion Order Forms](#). Part numbers in this section are for reference only. This section allows a quick selection of a single cushion that matches the seat pan width and depth to alleviate listing all part numbers and codes on this power base form. The complete cushion part number and MSRP is a product of the cushion model, top foam, liner, cover, width, and depth as shown below. Reference the subsequent chart for available top foam, liner, and cover options for each cushion type. See specific cushion order form for all sizes, top foam, liners, covers, parts numbers, and MSRP.

**1. Cushion Type**

TRU-Comfort 2 SPP
Reference part: XRF000116.
Part series: TC2-###XXYY.

Essence SPP
Reference part: XRF000118.
Part series: ESS-###XXYY.

Essence Plus SPP
Part series: ESSP-###XXYY.

Solution SPP
Reference part: XRF000113.
Part series: SOL-###XXYY.

Spectrum Foam SPP
Reference part: XRF000112.
Part series: SPCF-###XXYY.

Spectrum Gel SPP
Reference part: XRF000117.
Part series: SPCG-###XXYY.

Premiere P
Reference part: XRF000122.
Part series: PREM-###XXYY.

Glacial SP
Reference part: XRF000114.
Part series: GLA-###XXYY.

Zen SP
Reference part: XRF000115.
Part series: ZEN-###XXYY.

Simplicity G
Reference part: XRF000111.
Part series: SIM-###XXYY.

Contoured seating section continued on following pages.

14. CONTOURED SEATING CONTINUED**2. Essence Cushion Wedge Kit**

Only available with Essence cushion selection above.

Small Wedge Positioning Kit\$45
 Part: POP1909006-S
 14"W and smaller

Small Wedge Positioning Kit\$45
 Part: POP1712340-S
 15"W and larger

Medium Wedge Positioning Kit.....\$85
 Part: POP1909006-M
 14"W and smaller

Medium Wedge Positioning Kit.....\$85
 Part: POP1712340-M
 15"W and larger

Large Wedge Positioning Kit.....\$90
 Part: POP1909006-L
 14"W and smaller

Large Wedge Positioning Kit.....\$90
 Part: POP1712340-L
 15"W and larger

3. Cushion Top Foam

Visco Top Foam
 Reference part: XRF000380

FRG Top Foam
 Reference part: XRF000381

No Top Foam
 Reference part: XRF000382

4. Cushion Liner

Stealth Urethane Liner
 Reference part: XRF000383

Stealth Polyurethane Liner
 Reference part: XRF000384

No Liner
 Reference part: XRF000385

5. Cushion Cover

CoolCore Cover
 Reference part: XRF000386

Mesh Cover
 Reference part: XRF000387

Stealth Tek Cover
 Reference part: XRF000388

	None	Visco/V.Elastic	FRG	None	Urethane	Polyurethane	CoolCore	Mesh	Stealth Tek
	Top Foam			Liner			Cover		
TRU-Comfort 2	-	✓	-	✓	-	✓	✓	-	-
Essence	-	✓	-	✓	-	✓	✓	✓	-
Solution	-	✓	-	-	✓	✓	✓	✓	✓
Spectrum Foam	-	✓	✓	-	✓	✓	✓	✓	✓
Spectrum Gel	-	✓	✓	-	✓	✓	✓	✓	✓
Premiere	-	✓	✓	✓	✓	✓	✓	-	-
Glacial	✓	-	-	✓	✓	✓	✓	✓	✓
Zen	-	✓	✓	✓	✓	✓	✓	✓	✓
Simplicity	✓	-	-	✓	-	-	-	-	✓

6. Optional Base Insert

Option selection available with Glacial, Premiere, Simplicity, Solution, Spectrum Foam, Spectrum Gel, and Zen cushions.

ABS Flat Rigidizer
 Reference part: XRF000389
 MSRP: \$48-\$69

7. Additional Cushion Liner

Up to 5 may be selected below. See respective cushion order form for MSRP and part numbers.

Additional Stealth Urethane Liner
 Reference Part: XRF000390
 QTY _____

Additional Stealth Polyurethane Liner
 Reference Part: XRF000391
 QTY _____

8. Additional Cushion Cover

Up to 5 may be selected below. See respective cushion order form for MSRP and part numbers.

Additional CoolCore Cover
 Reference Part: XRF000392
 QTY _____

Additional Stealth Tek Cover
 Reference Part: XRF000394
 QTY _____

Additional Mesh Cover
 Reference Part: XRF000393
 QTY _____

14.5. Armrests**Optional.** Configured in pairs. Left (L) and right (R) boxes may be checked to configure different size arms within one section arm type.**1a. Q-Captain, Flip-Up, Height Adjustable Armrests**

Includes full length, flat arm pads.

Full Length, Left Armrest No Charge
 Part: MEC2201304. HCPCS: K0020

Full Length, Right Armrest..... No Charge
 Part: MEC2201303. HCPCS: K0020

1b. Drop-In, Height Adjustable, Armrests (T-Arm)

Includes flat arm pads.

☐ ☐ Full Length, Drop-In Armrests..... No Charge
 L R Left part: XRF000298. Right part: XRF000299. HCPCS: E0973

☐ ☐ Desk Length, Drop-In Armrests No Charge
 L R Left part: XRF000300. Right part: XRF000301. HCPCS: E0973

1c. Single Post, Height Adjustable, Removable Armrests

Height adjustment from 11"-17" in half-inch increments.

☐ ☐ Full Length, Single Post Armrests\$145 ea.
 L R Left part: ARM2000906. Right part: ARM2000907. HCPCS: E0973

☐ ☐ Desk Length, Single Post Armrests.....\$145 ea.
 L R Left part: ARM2000908. Right part: ARM2000909. HCPCS: E0973

1d. Bariatric, Height Adjustable, Removable Armrests

Select armrests and arm pads. Armrests are sold in same-size pairs.

☐ ☐ Full Length, Bariatric Armrests, 8"-12" H\$145 ea.
 L R Left part: FRMASMB6728. Right part: FRMASMB6729. HCPCS: E0973

☐ ☐ Full Length, Bariatric Armrests, 10"-14"H\$145 ea.
 L R Left part: FRMASMB6730. Right part: FRMASMB6731. HCPCS: E0973

1e. Heavy Duty, Height Adjustable, Removable Armrests

Select armrests and arm pads. Armrests are sold in same-size pairs.

☐ ☐ Full Length, Heavy Duty Armrests, 8"-12" H\$130 ea.
 L R Left part: XRF000253. Right part: XRF000254. HCPCS: E0973

☐ ☐ Full Length, Heavy Duty Armrests, 10"-14"H.....\$130 ea.
 L R Left Part: XRF000257. Right part: XRF000258. HCPCS: E0973

Contoured seating section continued on following page.

14. CONTOURED SEATING CONTINUED**1f. Pediatric, Removable, Height Adjustable Armrests**

Height adjustment from 7.75" to 8.25". Includes desk length, straight armpads.

L R Pediatric, Removable Armrests \$145 ea.
 Left Part: ARM142876. Right part: ARM142877. HCPCS: E0973

14.7. Armrest Height**REQUIRED with Bariatric, Heavy Duty, and Pediatric Height Adjustable armrests.****1. Left Armrest Height**

7.75" H PTOINDV3722	8.25" H PTOINDV3724	10" H PTOINDV3537	12" H PTOINDV3539	14" H PTOINDV3541	16" H PTOINDV3543
8" H PTOINDV3535	9" H PTOINDV3536	11" H PTOINDV3538	13" H PTOINDV3540	15" H PTOINDV3542	

2. Right Armrest Height

7.75" H PTOINDV3721	8.25" H PTOINDV3723	10" H PTOINDV3452	12" H PTOINDV3454	14" H PTOINDV3456	16" H PTOINDV3458
8" H PTOINDV3450	9" H PTOINDV3451	11" H PTOINDV3453	13" H PTOINDV3455	15" H PTOINDV3457	

14.8. Armpads**REQUIRED except with Flip-Up, Drop-In, and Pediatric armrests.****1a. Straight Armpads for Single Post Armrests**

Straight, Full Length Left and Right No Charge Part: ARM134275	Straight, Desk Length Left, Full Length Right No Charge Part: ARM134280
Straight, Full Length Left, Desk Length Right No Charge Part: ARM134281	Straight, Desk Length Left and Right No Charge Part: ARM134254

1b. Waterfall Armpads for Single Post Armrests

Waterfall, Full Length Left and Right No Charge Part: ARM134285	Waterfall, Desk Length Left, Full Length Right No Charge Part: ARM134288
Waterfall, Full Length Left, Desk Length Right No Charge Part: ARM134289	Waterfall, Desk Length Left and Right No Charge Part: ARM134282

1c. Falcon Gel Armpads for Single Post Armrests

Falcon Gel, Full Length Left and Right \$300 Part: ARM136438	Falcon Gel, Desk Length Left, Full Length Right \$250 Part: ARM136439
Falcon Gel, Full Length Left, Desk Length Right \$250 Part: ARM136442	Falcon Gel, Desk Length Left and Right \$200 Part: ARM136437

1d. Armpads for Bariatric or Heavy Duty Armrests

Bariatric Armpads \$80 Part: UPSASMB1027 Only for use with Bariatric Armrests.	Waterfall Armpads No Charge Part: UPSASMB1024 Only for use with Heavy Duty Armrests
Straight Armpads No Charge Part: UPSASMB1021 Only for use with Heavy Duty Armrests	

14.9. Quad Release Lever**Optional. Available with Single Post and Pediatric armrests.**

Quad Release Lever, Pair No Charge Part: FRMASMB10517
--

15. LAP BELTS**REQUIRED. Select one lap belt below.****15.1a. Standard Lap Belts**

38" Lap Belt Part: ACCBELT1006	50" Lap Belt Part: ACCBELT1000	60" Lap Belt Part: ACCBELT1002	70" Lap Belt Part: ACCBELT1003	80" Lap Belt Part: ACCBELT1004	90" Lap Belt Part: ACCBELT1001
-----------------------------------	-----------------------------------	-----------------------------------	-----------------------------------	-----------------------------------	-----------------------------------

15.1b. Optional Lap Belts

Stealth Retractable Lap Belt (IFRB5350) \$175 Part: ST-IFRB5350	Stealth Positioning, 60" Single Pull, Push Button, \$134 Extra Large Padded (11.5"), 2" Belt (IFPB5315) Part: ST-IFPB5315-Q
Stealth Positioning, 60" Single Pull, Push Button, Medium \$134 (7.5") Padded, 1.5" Belt (IFPB5213) Part: ST-IFPB5213-Q	Stealth Positioning, 60" Single Pull, Push Button, \$134 Non-padded, 2" Belt (IFPB5310) Part: ST-IFPB5310-Q
Stealth Positioning, 60" Single Pull, Push Button, Large \$134 Padded (9.5"), 2" Belt (IFPB5314) Part: ST-IFPB5314-Q	

16. LEGRESTS

Optional. All efforts will be made to accommodate legrest length selection, but minor changes may need to be made to meet shipping guidelines. Some adjustments may be needed upon receipt of the unit. Reference the section's compatibility matrix for legrest interference.

16.1a. Basic Center Mount Foot Platform

Select legrest length in following section.

1. Center Mount Foot Platform Type

Basic Center Mount \$200
Part: XRF000323

2. Basic Center Mount Foot Platform Extension and Angle Adjustable Footplates

Short Extension, 5"x9" Footplates No Charge Medium Extension, 5"x9" Footplates No Charge
Part: MEC1906675 Part: MEC1906680

3. Basic Center Mount Calf Pads

Calf pads are standard with a power articulating foot platform.

Small Calf Pads, Basic Center Mount \$210
Part: RIG2104369

4. Basic Center Mount Additional Accessories

Basic Center Mount Footplate Connector \$55
Part: RIG2104240

16.1b. Center Mount Foot Platform

Select legrest length in following section.

1. Center Mount Foot Platform Assembly

Power Articulating Foot Platform \$3,145 Quick Adjustable Center Mount Foot Platform \$265
Part: PTOINDV3754. HCPCS: E1012+K0040 Part: XRF000311. HCPCS: K0040
Only available with power lift option.

2. Center Mount Lower Extension and Angle Adjustable Footplates

Pediatric extensions are not compatible with Power Articulating Foot Platform.

Pediatric Extension, 5"x6" Footplates Part: RIG152577	Medium Extension, 5"x8" Footplates Part: RIG2404016	Long Extension, 5"x8" Footplates Part: RIG2404017
Pediatric Extension, 5"x8" Footplates Part: RIG2404015	Medium Extension, 5"x11" Footplates Part: RIG2404018	Long Extension, 5"x11" Footplates Part: RIG2404019
Medium Extension, 5"x6" Footplates Part: RIG150206	Long Extension, 5"x6" Footplates Part: RIG151505	

3a. Center Mount Calf Pads

Calf pads are standard and required with a power articulating foot platform.

Small Calf Pads, AFP Only, 6"Wx6"H (pair) Standard Small Calf Pads, Quick Adjustable Center Mount Only, \$210
Part: RIGASMB7120018 6"Wx6"H (pair)
Part: XRF000156

3b. Center Mount Calf Panel

Small Calf Panel, 18"Wx8"H \$265
Part: FRMASMB7110012

4. Center Mount Modular Footplate Extension Components

Optional. Assembly required. May require seating configuration adjustments to avoid front caster interference.

2" Bell-Shaped Length Extension \$55 Part: RIG152519 Only available with the 5x6 footplates.	Outside Edge \$55 Part: RIG2404921 Available with 5"x8" footplates.
Width Extension for 8" \$55 Part: RIG2404919	Front and Rear Cap \$55 Part: RIG2404923 Available with 5"x8" footplates.
Width Extension for 11" \$55 Part: RIG2404920	

5. Center Mount Additional Accessories

Small Foam Calf Pad Cover (pair) \$150 Part: ST-CSCS0775-Q	Footplate Riser \$150 Part: XRF000348 Select legrest length between 2" to 9" at 1" increments.
Small Foam Panel Cover (pair) \$150 Part: ST-CSCP0920-Q	Footplate Riser Calf Panel No Charge Part: FRMASMB7110011 Only available and required with Footplate Riser selection and extension lengths of 6" to 9".
Footplate Connector \$55 Part: RIGASMB7120023	

6. Stealth Gel Footplate Covers

5"x8" Tapered Gel Footplate Covers \$275 Part: ST-FPG5008T-Q. HCPCS: K0108	5"x8" Tapered Gel Footplate Covers for Width Extension \$275 Part: ST-AFPG5008E-Q. HCPCS: K0108
5"x11" Tapered Gel Footplate Covers \$275 Part: ST-AFPG5011T-Q. HCPCS: K0108	5"x11" Tapered Gel Footplate Cover for Width Extension \$275 Part: ST-AFPG5011E-Q. HCPCS: K0108

Legrests seating section continued on following page.

16. LEGRESTS CONTINUED**16.1c. Swing-Away Legrests**

Select a legrest upper, a legrest extension, and type of footplate. Pair (P), left (L), and right (R) configurations are selected by checking the appropriate boxes.

1. 60°-90° Swing-Away Legrest Upper

Adjustable in 5° increments.

P L R

60°-90° Swing-Away Legrest, Static Seat \$142.50 ea.
Pair part: XRF000082. Left: XRF000083. Right: XRF000084
Not available with 11"-13" extensions.

P L R

60°-90° Swing-Away Legrest, Short, Static Seat \$142.50 ea.
Pair part: XRF000085. Left: XRF000086. Right: XRF000087
For use only with 11"-13" lower extensions.**2. Lower Extension**

P L R

11"-13" Length
Pair part: XRF000088.
Left: XRF000089.
Right: XRF000090.

P L R

13"-16.5" Length
Pair part: XRF000094.
Left: XRF000095.
Right: XRF000096.

P L R

15.5"-19.5" Length
Pair part: XRF000097.
Left: XRF000098.
Right: XRF000099.**3b. Angle Adjustable Footplate(s)**¹The 6"x11" footplates may cause interference with seat widths less than 18". ²The 8"x11" footplates may cause interference with seat widths less than 20".

P L R

5"Wx8"D Footplates with Heel Loops \$87.50 ea.
Pair part: FRMASMB4735. Left: RIG137791. Right: RIG137792
HCPCS: K0040+E0951

P L R

6"Wx11"D¹ Footplates with Heel Loops \$87.50 ea.
Pair part: FRMASMB4741. Left: RIG137796. Right: RIG137797
HCPCS: K0040+E0951

P L R

6"Wx8"D¹ Footplates with Heel Loops \$87.50 ea.
Pair part: FRMASMB4738. Left: RIG137793. Right: RIG137795
HCPCS: K0040+E0951

P L R

8"Wx11"D² Footplates with Heel Loops \$87.50 ea.
Pair part: FRMASMB4744. Left: RIG137798. Right: RIG137799
HCPCS: K0040+E0951**16.1d. Elevating Legrests**

Select a legrest upper, a legrest extension, and type of footplate. Pair (P), left (L), and right (R) configurations are selected by checking the appropriate boxes. Select legrest length in following section.

1. Elevating Legrests Upper

P L R

Elevating Legrest Assembly \$275 ea.
Pair part: FRMASMB6483. Left: RIG138012. Right: RIG138013. HCPCS: E0990
Select lower extension and footplates**2a. Lower Extension**

P L R

11"-14" Length
Pair: FRMASMB4755. Left: RIG138026. Right: RIG138030
For reinforced only

P L R

14"-17" Length
Pair: FRMASMB4756. Left: RIG138033. Right: RIG138035
For reinforced only

P L R

8.5"-12" HD Drop-Ins/12.5"-16.5" ELR, Ext Length
Pair part: XRF000101. Left: XRF000102. Right: XRF000103
For angle adjustable only

P L R

13"-18" HD Drop-Ins/15.5"-19.5" ELR, Ext Length
Pair part: XRF000104. Left: XRF000105. Right: XRF000106
For angle adjustable only**3a. Reinforced Cast Footplate(s)**

P L R

7.25"Wx6"D Reinforced Footplate..... Standard
Pair part: FRMASMB4748. Left: RIG138055. Right: RIG138056

P L R

8.5"Wx6"D Reinforced Footplate Standard
Pair part: FRMASMB4747. Left: RIG138060. Right: RIG138062
May cause interference on seat widths less than 18."**3b. Angle Adjustable Footplate(s)**¹The 6"x11" footplates may cause interference with seat widths less than 18". ²The 8"x11" footplates may cause interference with seat widths less than 20".

P L R

5"Wx8"D Footplates with Heel Loops \$87.50 ea.
Pair part: FRMASMB4735. Left: RIG137791. Right: RIG137792
HCPCS: K0040+E0951

P L R

6"Wx11"D¹ Footplates with Heel Loops \$87.50 ea.
Pair part: FRMASMB4741. Left: RIG137796. Right: RIG137797
HCPCS: K0040+E0951

P L R

6"Wx8"D¹ Footplates with Heel Loops \$87.50 ea.
Pair part: FRMASMB4738. Left: RIG137793. Right: RIG137795
HCPCS: K0040+E0951

P L R

8"Wx11"D² Footplates with Heel Loops \$87.50 ea.
Pair part: FRMASMB4744. Left: RIG137798. Right: RIG137799
HCPCS: K0040+E0951**16.2. Legrest Accessories**Heel Loops (pair)..... \$35
Part: FRMASMB7873. HCPCS: E0951
Standard on angle adjustable and multi axis footplates when noted.Wheel Bumper \$15 ea.
Part: XRF000100

Legrests seating section continued on following page.

16. LEGRESTS CONTINUED**16.3. Set Legrest Extension Length****REQUIRED with center mounts and matching left and right rigging configurations. Part number is PTOINDV plus number following # below.****1a. Centermount Extension Length**

2" L (#3989)	4" L (#3991)	6" L (#3993)	8" L (#3680)	9" L (#3682)	10" L (#3684)	11" L (#3686)	12" L (#3688)	13" L (#3690)	14" L (#3692)	15" L (#3694)
3" L (#3990)	5" L (#3992)	7" L (#3994)	8.5" L (#3681)	9.5" L (#3683)	10.5" L (#3685)	11.5" L (#3687)	12.5" L (#3689)	13.5" L (#3691)	14.5" L (#3693)	

1b. Front Rigging Extension Length

Non-matching or individual legrest configurations lengths are not set by the manufacturer. See previous section for center mount extension lengths.

7" L (#3470)	8" L (#3472)	9" L (#3474)	10" L (#3476)	11" L (#3478)	12" L (#3480)	13" L (#3482)	14" L (#3484)	15" L (#3486)
7.5" L (#3471)	8.5" L (#3473)	9.5" L (#3475)	10.5" L (#3477)	11.5" L (#3479)	12.5" L (#3481)	13.5" L (#3483)	14.5" L (#3485)	

17. POWER ELECTRONICS**REQUIRED with actuator seating.**

Single Actuator Through Controller and Interface (AM3) . No Charge

Part: XRF000335. HCPCS: E2310

Single Actuator Through Controller and Interface (AAM)..... \$1,675

Part: XRF000068. HCPCS: E2310

18. SEAT MOUNTED LED LIGHTS**Optional.** Requires Q-Logic 3 Joystick for Full Light Package or Q-Logic 3e Joystick.

Seat Mounted, LED, Full Light Package\$770

Part: XRF000130

Not available with Base Mounted, LED, Full Light Package.

19. ACCESSORIES**Optional.** Multiple accessories can be mounted along the various tracks but some interference may occur. Accessories may need to be installed by the provider.

Cup Holder, Q-Captain\$20

Part: ACC2006355

For use with Flip-Up armrests.

Cup Holder, Bariatric Arms or Gel Armpads\$20

Part: ACCASMB2632

Cup Holder, Tube Style Armrests.....\$20

Part: XRF000221

For HD Removable, Single Post, or Pediatric Arms.

Oxygen Holder.....\$110

Part: ACC2006353. HCPCS: E2208

Not available with Power Lift option

Oxygen Holder.....\$110

Part: ACC2400456. HCPCS: E2208

Select with Power Lift option.

Cane/Crutch Holder\$195

Part: ACC2006352. HCPCS: E2207

Not available with Power Lift option.

Cane/Crutch Holder\$195

Part: ACC2400455. HCPCS: E2207

Select with Power Lift option.

Rear Basket\$150

Part: ACC2006351

Not available with Power Lift option.

Rear Basket\$150

Part: ACC2400454

Select with Power Lift option.

Saddle Bag\$15

Part: DWR1020D010

Large, Black, Weather Cover\$240

Part: ACCCOVR1021

Quantum Backpack\$200

Part: ACC1905853

Grip Lap Board (20005).....\$100

Part: ACC152992

Grip Drink Holder (20085)\$100

Part: ACC152993

Reflector Kit\$35

Part: ACC124670

Clothing Guard (Pair).....\$275

Part: ACC125007

Available with Flip-Up armrests only. Not available with Elevating Legrests.

Transfer Handles\$152.50 ea.

Part: MEC2205680. HCPCS: E0705

Available with Flip-Up armrests only. Not available with Elevating Legrests.

Privacy Bag\$80

Part: ACC2005723

Personal Item Pouch\$175

Part: ACC2101017

Remote Stop Switch\$1,495

Part: XRF000379

Aftermarket Accessory Mounting Bracket

Select a quantity up to 5.

Aftermarket Accessory Mounting Bracket, Slotted\$100 ea.

Part: ACC2407046

QTY

Aftermarket Accessory Mounting Bracket, Non-slotted \$100 ea.

Part: ACC2405496

QTY

22. LAP TRAY**22.1. Simple Solutions Lap Tray****Mounting bracket is required with lap tray selection. See chart below for compatibility.****1. Lap Tray**Simple Solutions Small Tray with Cut Out \$360
Part: ST-SSLXS010001171-Q. HCPCS: E0950Simple Solutions Medium Tray with Cut Out..... \$360
Part: ST-SSLXM010001171-Q. HCPCS: E0950Simple Solutions Large Tray with Cut Out..... \$360
Part: ST-SSLXL010001171-Q. HCPCS: E0950**2. Mounting Hardware**Swing Away Mounting Hardware.....No Charge
Part: ST-SSLXHD10000060-Q

Simple Solutions Lap Tray Armrest Compatibility	Width											
	10	11	12	13	14	16	17	18	19	20	21	22
Flip-Back	S	S	S	S	S	S	SM	ML	ML	ML	L	L
Flip-Back, Cane Mounted	S	S	S	S	S	SM	SM	ML	ML	L	L	
Drop-In, T-Arm	S	S	S	S	S	SM	SM	ML	ML	L	L	
Two-Post, Flip-Back	S	S	S	S	S	SM	SM	ML	ML	L	L	
Cantilever	S	S	S	S	S	M	M	ML	ML	L	L	
Q-Captain, Flip-Up						SM		ML		L		
Single Post/ Heavy Duty/ Bariatric	S	S	S	S	S	SM	SM	S	ML	ML	L	L

S:ST-SSLXS010001171-Q M:ST-SSLXM010001171-Q L:ST-SSLXL010001171-Q

:Compatible

:Armrest configuration may not hold tray position

*Some armrests listed may not be available on this order form.