



Pride Mobility Products Co.

5096 South Service Rd, Beamsville, ON, L3J 1V4

Phone: 888-570-1113 | Fax: 866-514-1303 | Email: productorders@pridemobility.com

Q6 Edge® 3 Stretto Order Form TRU Balance® 3 Positioning

Account Number: _____ Date: _____

Provider Name: _____

Contact: _____

Phone: _____ Fax: _____

Email: _____

PO Number: _____

Marked for: _____

Ship to Address: _____

City: _____ Prov: _____ PC: _____



INTRODUCTION

This form is interactive when viewed with [Adobe Acrobat Reader](#) and may not function correctly if opened with applications other than Acrobat. Complete the form by placing checks in the desired boxes and provide information in the interactive fields. Buttons shown at the bottom of the form may be utilized to print or submit the order form through a desktop email application. To email via a web-based application, please 'Save As' and attach the PDF to your email.


Thank you for choosing Quantum!

PATIENT INFORMATION

Weight: _____ Height: _____ A. Shoulder Width: _____ D. Max Sitting Height: _____ G. Thigh Depth: _____

B. Chest Width: _____ E. Shoulder Height: _____ H. Lower Leg Length: _____

C. Hip Width: _____ F. Axilla Height: _____ I. Elbow Height: _____



1. BASE MODEL

REQUIRED. Select a group 3 power base. 300 lb. weight capacity with standard 6mph programmable motors and off board charger. Dimensions of base without legrests: 20.5"Wx35.5"L

Quantum Q6 Edge 3 Stretto\$8,450

Part: Q6 EDGE 3 STRETTO

2. DRIVE MOTORS

Optional. Requires Q-Logic 3 electronics. The standard 6 mph motor can be upgraded to Accu-Trac motors that delivers enhanced traction and tracking, especially with switched drive controls.

Accu-Trac Motor Technology\$995
Part: XRF000268, HCPCS: K0108

3. BASE COLOR

REQUIRED. Select one shroud color. Colors shown here may differ from actual product depending on monitor/printer calibration used.

3.1a. Standard Shrouds

 Red Corvette Part: DGN1900199-PRD1	 Deep Purple Part: DGN1900199-PPR5	 Midnight Blue Part: DGN1900199-PBU10	 Green Machine Part: DGN1900199-PGN11
 Back in Black Part: DGN1900199-PBK1	 Raspberry Beret Part: DGN1900199-PPK4	 Electric Blue Part: DGN1900199-PBU9	 Military Green Part: DGN1900199-PGN8
 Night Train Part: DGN1900199-PBK3	 Pretty in Pink Part: DGN1900199-PPK3	 Blue Jean Part: DGN1900199-PBU11	
 Grey Street Part: DGN1900199-PSL3	 Evolution Orange Part: DGN1900199-POR7	 Lemon Crush Part: DGN1900199-PYL4	

3.1b. Hydrographic Shrouds

 Military Digital Camo\$325 Part: DGN1801101-WCM7	 Faded Glory\$325 Part: DGN1801101-WFL2
--	--

4. BASE MOUNTED LIGHTS

Optional. Select one option.

Switch Fender Lights..... No Charge
Part: DGN1810041
Fender lights operate by pressing included single switch button. Fender lighting is not available with Single Actuator Control through Toggle Box.

**5. QUANTUM BACKUP CAMERA**

Optional. The backup camera may interfere with armrest mounted accessories such as media device holder.

Quantum Backup Camera\$475
Part: MEC1910365

1. Optional Camera Mounting Bracket

Backup Camera Headrest Bracket..... \$35 Part: MEC1910373 Works with .5" x .5" vertical supports only.	Backup Camera Track Mount Bracket \$35 Part: MEC2104227
--	--

2. Optional Display Mounting Bracket

Backup Camera Q-Logic Joystick Mount Bracket\$145
Part: MEC2106123

**6. DRIVE WHEELS**

REQUIRED. Select a type of drive wheel. Accent colors are available for 14" tires.

6.1. Drive Wheels

3"x14" Pneumatic Tires..... No Charge Part: WHL1713925	2"x12.5" Pneumatic Tires No Charge Part: WHL1810666
3"x14" Flat Free Tires..... No Charge Part: WHL1704058	2"x12.5" Flat Free Tires No Charge Part: WHL1810665
4"x14" Knobby, Flat Free Black Tires..... \$120 Part: XRF000368	4"x14" Knobby, Pneumatic Black Tires..... \$120 Part: XRF000367

6.2. Drive Wheel Accent Color for 14" Tires

Not available with 12.5" tires.

Back in Black..... No Charge Part: DGN164245-PBK1 x2	Deep Purple..... No Charge Part: DGN164245-PPR5 x2	Midnight Blue No Charge Part: DGN164245-PBU10 x2	Green Machine No Charge Part: DGN164245-PGN11 x2
Red Corvette..... No Charge Part: DGN164245-PRD1 x2	Raspberry Beret..... No Charge Part: DGN164245-PPK4 x2	Electric Blue No Charge Part: DGN164245-PBU9 x2	Military Green No Charge Part: DGN164245-PGN8 x2
Night Train..... No Charge Part: DGN164245-PBK3 x2	Pretty in Pink No Charge Part: DGN164245-PPK3 x2	Blue Jean No Charge Part: DGN164245-PBU11 x2	Military Digital Camo \$50 ea. Part: DGN164245-WCM7 x2
Grey Street No Charge Part: DGN164245-PSL3 x2	Evolution Orange..... No Charge Part: DGN164245-POR7 x2	Lemon Crush..... No Charge Part: DGN164245-PYL4 x2	Faded Glory.....\$50 ea. Part: DGN164245-WFL2 x2

7. BATTERIES

Optional. The powerchair requires 2 batteries and will be installed if ordered.

NF-22, Gel Cell, Interceptor 2 Battery \$300 ea.

Part: BAT2004376. HCPCS: E2361

Uses 8 AMP charger

Gel Cell, 40AH, Interceptor 2 Battery \$200 ea.

Part: BAT2100316. HCPCS: K0108

Uses 5 AMP charger

8. PATIENT'S WEIGHT

REQUIRED. The patient's weight is required to quote and order the power chair.

Up to 150 lb User Weight

Part: PTOINDV3977

151-250 lb User Weight

Part: PTOINDV3985

250-300 lb User Weight

Part: PTOINDV3984

9. EXTENDED SEAT TOWERS

Optional. Only available with TRU Balance Static Seat or Power Recline Only.

Extended Seat Towers No Charge

Part: MEC2107220

10. ELECTRONICS

REQUIRED. Select a type of drive control from the following sections. The No Joystick option is used to omit the joystick and utilize another control, such as a specialty control at the end of this form or found on the [Alternative Controls and Electronics Order Form](#).

10.1a. Non-Expandable Controls

Select a joystick below or omit selection if another controller type is desired. Non-expandable does not allow supplementary drive controllers.

QLNES 6 Key Joystick Standard

Part: CTL2205379

Controls up to 2 actuators.

QLNES 10 Key Joystick No Charge

Part: CTL2205380

Controls up to 2 actuators through joystick and seat mounted, LED, full light package.

10.1b. Expandable Controls

Select expandable controls or omit selection if another controller type is desired. Expandable controls allow multiple or specialized controllers.



2. Q-Logic 3 Joystick

Q-Logic 3 EX Joystick includes XLR USB charger.

Q-Logic 3 EX Joystick \$900

Part: CTL167724

Q-Logic 3 EX Joystick for Full Lights \$900

Part: CTL167725

Q-Logic 3e Joystick \$620

Part: CTL1910600

Can operate full light package. Not available with alternative or specialty controls.

No Joystick, Setup for Q-Logic 3

Part: XRF000010



10.2. Q-Logic Magnetic Charger

Available with QLNES. Included with Q-Logic 3e Joystick.

Magnetic to XLR Charger Adapter with Wall Holder \$100

Part: ELE2103635

10.3. Q-Logic 3 Joystick Armor

REQUIRED with Q-Logic 3 EX Joystick selections above. Not available with Q-Logic 3e Joystick. Reference colors shown on page 2. *Night train joystick armor will be the black, soft touch.

Black, Soft Touch Standard

Part: DGN124307

Also matches Night Train shrouds.

Red Corvette No Charge

Part: DGN124336-PRD1

Back in Black No Charge

Part: DGN124336-PBK1

Grey Street No Charge

Part: DGN124336-PSL3

Deep Purple No Charge

Part: DGN124336-PPR5

Raspberry Beret No Charge

Part: DGN124336-PPK3

Pretty in Pink No Charge

Part: DGN124336-PPK3

Evolution Orange No Charge

Part: DGN124336-POR7

Midnight Blue No Charge

Part: DGN124336-PBU10

Electric Blue No Charge

Part: DGN124336-PBU9

Blue Jean No Charge

Part: DGN124336-PBU11

Lemon Crush No Charge

Part: DGN124336-PYL4

Green Machine No Charge

Part: DGN124336-PGN11

Military Green No Charge

Part: DGN124336-PGN8

Military Digital Camo No Charge

Part: DGN124336-WCM7

Faded Glory No Charge

Part: DGN124336-WFL2

11. JOYSTICK HANDLES**Optional.** The optional handle will be mounted to the joystick and the standard handle will be shipped with the unit.Mushroom Handle (B) \$120
Part: ACCASMB1860.Chin Cup (C) \$85
Part: ACCASMB186190T-Bar (D) \$160
Part: ACCASMB1858.Medium Ball Handle (E) \$70
Part: ACCASMB1859.

Soft Ball Handle (F) \$120

Part: ACCASMB1856.

Extended I-Handle (G) \$140
Part: ACCASMB1857.Bodypoint 3" Goal Post Handle \$180
Part: PLSKNCP1057.Bodypoint 4" Goal Post Handle \$180
Part: PLSKNCP1058.**12. JOYSTICK MOUNTING BRACKETS****Optional.** Select a mount to be used with a joystick selected in previous section.**12.1a. Fixed, In-line Joystick Bracket**Fixed, In-line, Bracket for Q-Logic Standard
Part: FRMASMB10617**12.1b. Swing-Away (Retractable) Joystick Bracket**

Allows the joystick to retract to the side of the armrest while maintaining the same joystick orientation.

Swing-away, Q-Logic \$265
Part: ACC1908434. HCPCS: E1032**12.1c. Basic Swing Away Joystick Mount**Basic Swing Away, Q-Logic, Left Side \$265
Part: ACC2201041. HCPCS: E1032Basic Swing Away, Q-Logic, Right Side \$265
Part: ACC2201042. HCPCS: E1032**12.1d. Offset Swing Away Joystick Mount**Offset Swing Away, Q-Logic, Left Side \$265
Part: ACC142742. HCPCS: E1028Offset Swing Away, Q-Logic, Right Side \$265
Part: ACC142743. HCPCS: E1028**12.1e. Specialty Joystick Mounts**



Select a bracket and side to mount it.

Retract4 Mount \$450
Part: XRF000003. HCPCS: E1032

Full 180° of movement that translates into multi-position use, straight, right offset, left offset, or fold away. Joystick orientation varies as bracket is adjusted.

Multi-Axis Mount \$100
Part: XRF000004

Allows height adjustment of joystick.

**12.2. Joystick Mount Location**☐ Left Side Mount
Part: PTOINDV1225☐ Right Side Mount
Part: PTOINDV1226**13. TRU BALANCE® 3 SEATING - QUANTUM INTELLIGENT SYSTEMS™****REQUIRED.** Maximum seat size is 20"W x 20"D.**13.1. TRU Balance 3 Configuration**TRU Balance 3 Static Seat \$1,045
Part: PTOINDV3442
20° of limited recline adjustment (-5/ +15) via gas strut within back cover.TRU Balance 3 Power Tilt Only \$6,245
Part: PTOINDV3443. HCPCS: E1002
50° Center of Gravity (CG) power tilt seating system.TRU Balance 3 Power Recline Only \$7,875
Part: PTOINDV3444. HCPCS: E1004
175° of power recline.TRU Balance 3 Power Tilt & Recline \$13,125
Part: PTOINDV3446. HCPCS: E1007
50° CG power tilt with 175° of power recline. Full drive lockout at 20°.**13.2. 12" Power Lift****Optional.** ^Please verify the appropriate billing code to be used for this item with all third-party payors.TRU Balance 3 12" Power Lift with iLevel® \$2,995
Part: PTOINDV3773. HCPCS: E2298/E2300^
12" of scissor lift with iLevel System.**13.3. Occupied Transit Kit****Optional.** The power base is standard with Unoccupied Transit Loops. 60" Sure-Lok Lap Belt is required with Occupied Transit Option.WC-19 Occupied Transit Option \$400
Part: KIT152047To determine if your chair is WC-19 compliant (occupied), please refer to this website: <https://wc-transportation-safety.umtri.umich.edu/>**13.4. Seat Width****REQUIRED with seat selection.** Maximum seat width is 20". Part number is SETWDTH plus number following # below. 12" Width is not available with 16" and 17" Depth. 13" and 15" widths are only available with non-matching back width.

12" W (#1002)	13" W (#1003)	14" W (#1004)	15" W (#1005)	16" W (#1006)	17" W (#1007)	18" W (#1008)	19" W (#1009)	20" W (#1010)	21" W (#1011)	22" W (#1012)
------------------	------------------	------------------	------------------	------------------	------------------	------------------	------------------	------------------	------------------	------------------

TRU Balance 3 Seating section continued on following pages.

13. TRU BALANCE 3 SEATING CONTINUED**13.5. Seat Width Package**

REQUIRED with seat selection.

12"-16" Width Package
Part: XRF00000614"-16" Width Package
Part: XRF00000716"-20" Width Package
Part: XRF00000820"-24" Width Package
Part: XRF000158
Maximum width is 22".**13.6. Seat Depth**

REQUIRED with seat selection. Maximum seat depth is 20". Part number is SETDPTH plus number following # below.

12" D
(#1002)13" D
(#1003)14" D
(#1004)15" D
(#1005)16" D
(#1006)17" D
(#1007)18" D
(#1008)19" D
(#1009)20" D
(#1010)21" D
(#1011)22" D
(#1012)**13.7. Seat Depth Package**

REQUIRED with seat selection.

12"-16" Depth Package
Part: XRF00001114"-16" Depth Package
Part: XRF00001216"-20" Depth Package
Part: XRF00001320"-24" Depth Package
Part: XRF000159
Maximum depth is 22".**13.8. Seat to Floor Height**


REQUIRED.

Lowest Seat to Floor
Part: PTOINDV3558Middle Seat to Floor
Part: PTOINDV3559Highest Seat to Floor
Part: PTOINDV3560

Low	Med	High
Static Seat/ Recline Only		
16"	17"	18"
Power Tilt/Tilt & Recline		
17.5"	18.5"	19.5"
Any Power Lift		
18"	19"	20"
TRU Balance 3 Seat to Floor Heights		

13.9. Stealth Seat Cushion

Optional. Specific cushion/component sizes available, part numbers, weight capacities, MSRP, HCPCS code, and special order cushions will be found on the respective complete order forms. Part numbers in this section are for reference only. This section allows a quick selection of a single cushion that matches the seat pan width and depth to alleviate listing all part numbers and codes on this power base form. The complete cushion part number and MSRP is a product of the cushion model, top foam, liner, cover, width, and depth as shown below. Reference the subsequent chart for available top foam, liner, and cover options for each cushion type. See specific cushion order form for all sizes, top foam, liners, covers, parts numbers, and MSRP.

**1. Cushion Type**

TRU-Comfort 2 SPP
Reference part: XRF000116.
Part series: TC2-###XXYY.

Essence SPP
Reference part: XRF000118.
Part series: ESS-###XXYY.

Essence Plus SPP
Part series: ESSP-###XXYY.

Solution SPP
Reference part: XRF000113.
Part series: SOL-###XXYY.

Spectrum Foam SPP
Reference part: XRF000112.
Part series: SPCF-###XXYY.

Spectrum Gel SPP
Reference part: XRF000117.
Part series: SPCG-###XXYY.

Premiere P
Reference part: XRF000122.
Part series: PREM-###XXYY.

Glacial SP
Reference part: XRF000114.
Part series: GLA-###XXYY.

Zen SP
Reference part: XRF000115.
Part series: ZEN-###XXYY.

Simplicity G
Reference part: XRF000111.
Part series: SIM-###XXYY.

2. Essence Cushion Wedge Kit

Only available with Essence cushion selection above.

Small Wedge Positioning Kit.....\$45
Part: POP1909006-S
14"W and smaller

Medium Wedge Positioning Kit.....\$85
Part: POP1909006-M
14"W and smaller

Large Wedge Positioning Kit.....\$90
Part: POP1909006-L
14"W and smaller

Small Wedge Positioning Kit.....\$45
Part: POP1712340-S
15"W and larger

Medium Wedge Positioning Kit.....\$85
Part: POP1712340-M
15"W and larger

Large Wedge Positioning Kit.....\$90
Part: POP1712340-L
15"W and larger

TRU Balance 3 Seating section continued on following pages.

13. TRU BALANCE 3 SEATING CONTINUED**3. Cushion Top Foam**

Visco Top Foam
Reference part: XRF000380

FRG Top Foam
Reference part: XRF000381

No Top Foam
Reference part: XRF000382

4. Cushion Liner

Stealth Urethane Liner
Reference part: XRF000383

Stealth Polyurethane Liner
Reference part: XRF000384

No Liner
Reference part: XRF000385

5. Cushion Cover

CoolCore Cover
Reference part: XRF000386

Mesh Cover
Reference part: XRF000387

Stealth Tek Cover
Reference part: XRF000388

	None	Visco/V.Elastic	FRG	None	Urethane	Polyurethane	CoolCore	Mesh	Stealth Tek
	Top Foam			Liner					
TRU-Comfort 2	—	✓	—	✓	—	✓	✓	—	—
Essence	—	✓	—	✓	—	✓	✓	✓	—
Solution	—	✓	—	—	✓	✓	✓	✓	✓
Spectrum Foam	—	✓	✓	—	✓	✓	✓	✓	✓
Spectrum Gel	—	✓	✓	—	✓	✓	✓	✓	✓
Premiere	—	✓	✓	—	✓	✓	✓	—	—
Glacial	✓	—	—	✓	✓	✓	✓	✓	✓
Zen	—	✓	✓	✓	✓	✓	✓	✓	✓
Simplicity	✓	—	—	✓	—	—	—	—	✓

6. Optional Base Insert

Option selection available with Glacial, Premiere, Simplicity, Solution, Spectrum Foam, Spectrum Gel, and Zen cushions.

ABS Flat Rigidizer
Reference part: XRF000389

7. Additional Cushion Liner

Up to 5 may be selected below. See respective cushion order form for MSRP and part numbers.

Additional Stealth Polyurethane Liner
Reference Part: XRF000391

8. Additional Cushion Cover

Up to 5 may be selected below. See respective cushion order form for MSRP and part numbers.

Additional CoolCore Cover
Reference Part: XRF000392

Additional Stealth Tek Cover
Reference Part: XRF000394

Additional Mesh Cover
Reference Part: XRF000393

13.10. Back Type Selection

REQUIRED. Match a choice in this section with the options in the following back cushion configuration section.

Redesigned, Static Type Back..... \$735
Part: XRF000014
12" and 24" back width will use the TRU Balance 3 Back assembly. Select a back cushion configuration in the subsequent section.

Back Cane Option
Part: XRF000015
Available only with TRU Balance 3 Static Seat and Power Tilt (1 actuator) selections. Select cane type, back cane height, and back cushion configuration below. Widths less than 15" must select Omit Back Pan/Cushion.

Redesigned, Flip Forward Back..... \$735
Part: XRF000016
12" and 24" back width will use the TRU Balance 3 Back assembly. Available with TRU Balance 3 Static Seat and Power Tilt (1 actuator). Only available with 2-post or pediatric armrests. Not available with recline configurations. Select a back cushion type in the subsequent section.

13.11. Back Cane Selection (Back Cane Option Only)

REQUIRED when Back Cane Option is selected above. Select type and height of back canes.

1. Back Cane Type

Angle Adjustable Back Cane..... No Charge
Part: MEC140266

Flip Forward Back Cane No Charge
Part: XRF000273

2" Endomorph (1" per side) Back Cane..... No Charge
Part: XRF000269

4" Endomorph (2" per side) Back Cane..... No Charge
Part: XRF000271

2. Back Cane Height

16"H Cane
Part: FRMASMB15834

17"H Cane
Part: FRMASMB15835

20"H Cane
Part: FRMASMB15838

21"H Cane
Part: FRMASMB15839

24"H Cane
Part: FRMASMB15842

25"H Cane
Part: FRMASMB15843

TRU Balance 3 Seating section continued on following pages.

13. TRU BALANCE 3 SEATING CONTINUED**13.12. Back Cushion Configuration**

REQUIRED. Sport Back and TRU-Comfort Back selections up to 22" back width will be the Sport 4 or TRU-Comfort 4 back cushions. Sport back and TRU-Comfort back selections of 24" back width will be the Sport or TRU-Comfort 2 back cushions. Select a configuration type and applicable dimensions. See the Sport, Q-Back, and TRU-Comfort back availability chart for stock order back size ranges. Other size backs may be available with special order and extended lead time.

1a. Sport Back/Omit/After-Market Interface

Sport Back No Charge
Part: XRF000017

Complete following width package and size.

Omit Back Pan/Cushion No Charge
Part: PTOINDV1287

Select a following width package.

After-market Back Interface No Charge

Part: ST-SSAPT4-Q. Select a set width in following section, range of 12"-21". Overall seat depth may need to be adjusted to accommodate certain backs. May interfere with gel/waterfall arm pads and arm through. Does not include brackets to mount after-market backs.

1b. TRU-Comfort Back

TRU-Comfort Back (Static/Flip Forward), 12"-15" Width \$835
Part: SET123622. HCPCS: E2620

Complete following width package and size.

TRU-Comfort Back (Static/Flip Forward), 16" to 20" Width ... \$1,035
Part: TC4V10W21W. HCPCS: E2620

Complete following width package and size.

TRU-Comfort Back, Back Cane Option \$835
Part: TC4V10W21WBC. HCPCS: E2620

Complete following width package and size.

TRU-Comfort Back (Static/Flip Forward), 22" Width \$1,235
Part: TC4V22W32W. HCPCS: E2621

Complete following width package and size.

Overall Height	Tilt or Static Width												Recline Width											
	12"	14"	15"	16"	17"	18"	19"	20"	21"	22"	24"	12"	14"	15"	16"	17"	18"	19"	20"	21"	22"	24"		
	14"	T	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	14"	
15"	T	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	15"	
16"	S	S	T	S	T	S	T	S	T	S	T	S	-	-	-	-	-	-	-	-	-	-	16"	
17"	S	S	T	S	T	S	T	S	T	S	T	S	-	-	-	-	-	-	-	-	-	-	17"	
18"	-	S	T	S	T	S	T	S	T	S	T	S	-	-	-	-	-	-	-	-	-	-	18"	
19"	-	S	T	S	T	S	T	S	T	S	T	S	S	S	T	S	T	S	T	S	T	S	19"	
20"	-	S	T	S	T	S	T	S	T	S	T	S	S	S	T	S	T	S	T	S	T	S	20"	
21"	-	S	T	S	T	S	T	S	T	S	T	S	-	S	T	S	T	S	T	S	T	S	21"	
22"	-	S	T	S	T	S	T	S	T	S	T	S	-	S	T	S	T	S	T	S	T	S	22"	
23"	-	S	T	S	T	S	T	S	T	S	T	S	-	S	T	S	T	S	T	S	T	S	23"	
24"	-	S	T	S	T	S	T	S	T	S	T	S	-	S	T	S	T	S	T	S	T	S	24"	
25"	-	S	T	S	T	S	T	S	T	S	T	S	-	S	T	S	T	S	T	S	T	S	25"	
26"	-	S	T	S	T	S	T	S	T	S	T	S	-	S	T	S	T	S	T	S	T	S	26"	
27"	-	S	T	S	T	S	T	S	T	S	T	S	-	S	T	S	T	S	T	S	T	S	27"	
28"	-	S	T	S	T	S	T	S	T	S	T	S	-	S	T	S	T	S	T	S	T	S	28"	
29"	-	S	T	S	T	S	T	S	T	S	T	S	-	S	T	S	T	S	T	S	T	S	29"	
30"	-	-	-	-	-	-	-	-	-	-	-	-	-	S	T	S	T	S	T	S	T	S	30"	
31"	-	-	-	-	-	-	-	-	-	-	-	-	-	S	T	S	T	S	T	S	T	S	31"	
32"	-	-	-	-	-	-	-	-	-	-	-	-	-	S	T	S	T	S	T	S	T	S	32"	

S

 Sport Back

Q

 Q-Back

T

 TRU-Comfort Back

Sport, Q-Back, and TRU-Comfort Back Availability Chart

S Sport Back Q Q-Back T TRU-Comfort Back

Sport, Q-Back, and TRU-Comfort Back Availability Chart

TRU Balance 3 Seating section continued on following pages.

13. TRU-BALANCE 3 SEATING CONTINUED

1c. ADI2.0 Backrest

Specify back type, width, height, and contour in the part number configurator below. Previously selected back cane height must be equal to the greater of overall height range. **Available SOF for Provider installation only.**

For more information please see [ADI2.0 Backrest order form](#)

ADI2.0 Backrest 10"W - 20"W \$1,495
 For use with back canes only.
 Indicate part number below.
 ADP: SE0002152

ADI2.0 Backrest 22"W - 24"W \$1,695
 For use with back canes only.
 Indicate part number below.
 ADP: SE0002153

Specify Part Number for Back Selection

OPTIONS INDEX	ADI2	Back Type: AL AFT	Width: 10"-20" (1" increments), 22", 24"	Height: N = Nino L = Low M = Medium S = Standard T = Tall X = Extra Tall	Contour: A = Active D = Deep E = Extra Deep	2 = FRG Foam	P = Quick-Release
ENTER SELECTION	ADI2					2	P

BACK CONTOUR*:



BACK SERIES TYPES*:

(AL) Aluminum Series

(*ALL BACKS SHOWN ARE 18" WIDE)

(AFT) Aluminum Flat Top Series



BACK HEIGHTS AND WIDTH:

(Available for AL Back Shell Only)

Nino: 8" inch (10"W - 14"W)

Low: 10 inch (10"W - 24"W)

Medium: 13 inch (12"W - 24"W)

(Available for AL and AFT Shell Types)

Standard: 16 inch (13"W - 24"W)

Tall: 18 inch (16"W - 24"W)

Extra Tall: 20 inch (16"W - 24"W)

13. TRU BALANCE 3 SEATING CONTINUED**1f. Q-Back 2.0, AL Series, Active Contour****Not available with Back Cane Option.** Select a back type and the specific overall back height in subsequent section.

Active Contour, 14"Wx19-20"H Tilt/22"-23" H Recline..... \$985 Part: QB2AL14AS. HCPCS: E2613	Active Contour, 18"Wx21-22"H Tilt/24"-25" H Recline..... \$1,145 Part: QB2AL18AT. HCPCS: E2613
Active Contour, 14"Wx21-22"H Tilt/24"-25" H Recline..... \$985 Part: QB2AL14AT. HCPCS: E2613	Active Contour, 18"Wx23-24"H Tilt/26"-27" H Recline..... \$1,145 Part: QB2AL18AX. HCPCS: E2613
Active Contour, 16"Wx19-20"H Tilt/22"-23" H Recline..... \$1,145 Part: QB2AL16AS. HCPCS: E2613	Active Contour, 18"Wx25-26"H Tilt/28"-29" H Recline..... \$1,145 Part: QB2AL18AXX. HCPCS: E2613
Active Contour, 16"Wx21-22"H Tilt/24"-25" H Recline..... \$1,145 Part: QB2AL16AT. HCPCS: E2613	Active Contour, 20"Wx19-20"H Tilt/22"-23" H Recline..... \$1,145 Part: QB2AL20AS. HCPCS: E2613
Active Contour, 16"Wx23-24"H Tilt/26"-27" H Recline..... \$1,145 Part: QB2AL16AX. HCPCS: E2613	Active Contour, 20"Wx21-22"H Tilt/24"-25" H Recline..... \$1,145 Part: QB2AL20AT. HCPCS: E2613
Active Contour, 16"Wx25-26"H Tilt/28"-29" H Recline..... \$1,145 Part: QB2AL16AXX. HCPCS: E2613	Active Contour, 20"Wx23-24"H Tilt/26"-27" H Recline..... \$1,145 Part: QB2AL20AX. HCPCS: E2613
Active Contour, 18"Wx19-20"H Tilt/22"-23" H Recline..... \$1,145 Part: QB2AL18AS. HCPCS: E2613	Active Contour, 20"Wx25-26"H Tilt/28"-29" H Recline..... \$1,145 Part: QB2AL20AXX. HCPCS: E2613

1g. Q-Back 2.0, AL Series, Deep Contour**Not available with Back Cane Option.** Select a back type and the specific overall back height in subsequent section.

Deep Contour, 14"Wx19-20"H Tilt/22"-23" H Recline \$1,145 Part: QB2AL14DS. HCPCS: E2620	Deep Contour, 18"Wx21-22"H Tilt/24"-25" H Recline \$1,325 Part: QB2AL18DT. HCPCS: E2620
Deep Contour, 14"Wx21-22"H Tilt/24"-25" H Recline \$1,145 Part: QB2AL14DT. HCPCS: E2620	Deep Contour, 18"Wx23-24"H Tilt/26"-27" H Recline \$1,325 Part: QB2AL18DX. HCPCS: E2620
Deep Contour, 16"Wx19-20"H Tilt/22"-23" H Recline \$1,325 Part: QB2AL16DS. HCPCS: E2620	Deep Contour, 18"Wx25-26"H Tilt/28"-29" H Recline \$1,325 Part: QB2AL18DXX. HCPCS: E2620
Deep Contour, 16"Wx21-22"H Tilt/24"-25" H Recline \$1,325 Part: QB2AL16DT. HCPCS: E2620	Deep Contour, 20"Wx19-20"H Tilt/22"-23" H Recline \$1,325 Part: QB2AL20DS. HCPCS: E2620
Deep Contour, 16"Wx23-24"H Tilt/26"-27" H Recline \$1,325 Part: QB2AL16DX. HCPCS: E2620	Deep Contour, 20"Wx21-22"H Tilt/24"-25" H Recline \$1,325 Part: QB2AL20DT. HCPCS: E2620
Deep Contour, 16"Wx25-26"H Tilt/28"-29" H Recline \$1,325 Part: QB2AL16DXX. HCPCS: E2620	Deep Contour, 20"Wx23-24"H Tilt/26"-27" H Recline \$1,325 Part: QB2AL20DX. HCPCS: E2620
Deep Contour, 18"Wx19-20"H Tilt/22"-23" H Recline \$1,325 Part: QB2AL18DS. HCPCS: E2620	Deep Contour, 20"Wx25-26"H Tilt/28"-29" H Recline \$1,325 Part: QB2AL20DXX. HCPCS: E2620

TRU Balance 3 Seating section continued on following pages.

13. TRU BALANCE 3 SEATING CONTINUED**1h. Q-Back 2.0, AL Series, Extra Deep Contour****Not available with Back Cane Option.** Select a back type and the specific overall back height in subsequent section.Extra Deep Contour, 14"Wx19-20"H Tilt/22"-23" H Recline \$1,225
Part: QB2AL14ES. HCPCS: E2620Extra Deep Contour, 14"Wx21-22"H Tilt/24"-25" H Recline \$1,225
Part: QB2AL14ET. HCPCS: E2620Extra Deep Contour, 16"Wx19-20"H Tilt/22"-23" H Recline..... \$1,410
Part: QB2AL16ES. HCPCS: E2620Extra Deep Contour, 16"Wx21-22"H Tilt/24"-25" H Recline..... \$1,410
Part: QB2AL16ET. HCPCS: E2620Extra Deep Contour, 16"Wx23-24"H Tilt/26"-27" H Recline \$1,410
Part: QB2AL16EX. HCPCS: E2620Extra Deep Contour, 16"Wx25-26"H Tilt/28"-29" H Recline \$1,410
Part: QB2AL16EXX. HCPCS: E2620Extra Deep Contour, 18"Wx19-20"H Tilt/22"-23" H Recline \$1,410
Part: QB2AL18ES. HCPCS: E2620Extra Deep Contour, 18"Wx21-22"H Tilt/24"-25" H Recline \$1,410
Part: QB2AL18ET. HCPCS: E2620Extra Deep Contour, 18"Wx23-24"H Tilt/26"-27" H Recline \$1,410
Part: QB2AL18EX. HCPCS: E2620Extra Deep Contour, 18"Wx25-26"H Tilt/28"-29" H Recline \$1,410
Part: QB2AL18EXX. HCPCS: E2620Extra Deep Contour, 20"Wx19-20"H Tilt/22"-23" H Recline \$1,410
Part: QB2AL20ES. HCPCS: E2620Extra Deep Contour, 20"Wx21-22"H Tilt/24"-25" H Recline \$1,410
Part: QB2AL20ET. HCPCS: E2620Extra Deep Contour, 20"Wx23-24"H Tilt/26"-27" H Recline \$1,410
Part: QB2AL20EX. HCPCS: E2620Extra Deep Contour, 20"Wx25-26"H Tilt/28"-29" H Recline \$1,410
Part: QB2AL20EXX. HCPCS: E2620**1i. Q-Back 2.0, AFT Series, Active Contour****Not available with Back Cane Option.** Select a back type and the specific overall back height in subsequent section.AFT, Active, 16"Wx19-20"H Tilt/22"-23" H Recline \$1,145
Part: QB2AFT16AS. HCPCS: E2613AFT, Active, 18"Wx19-20"H Tilt/22"-23" H Recline \$1,145
Part: QB2AFT18AS. HCPCS: E2613**1j. Q-Back 2.0, AFT Series, Deep Contour****Not available with Back Cane Option.** Select a back type and the specific overall back height in subsequent section.AFT, Deep Contour, 16"Wx19-20"H Tilt/22"-23" H Recline..... \$1,325
Part: QB2AFT16DS. HCPCS: E2620AFT, Deep Contour, 16"Wx21-22"H Tilt/24"-25" H Recline..... \$1,325
Part: QB2AFT16DT. HCPCS: E2620AFT, Deep Contour, 18"Wx19-20"H Tilt/22"-23" H Recline \$1,325
Part: QB2AFT18DS. HCPCS: E2620AFT, Deep Contour, 18"Wx21-22"H Tilt/24"-25" H Recline..... \$1,325
Part: QB2AFT18DT. HCPCS: E2620AFT, Deep Contour, 20"Wx21-22"H Tilt/24"-25" H Recline \$1,325
Part: QB2AFT20DT. HCPCS: E2620**1k. Q-Back 2.0, AFT Series, Extra Deep Contour****Not available with Back Cane Option.** Select a back type and the specific overall back height in subsequent section.AFT, Extra Deep Contour, 16"Wx25-26"H Tilt/28"-29" H Recline \$1,410
Part: QB2AFT16EXX. HCPCS: E2620**13.13. TRU Balance 3 Back Shroud****REQUIRED.** The back shroud is a feature of the TRU Balance 3 seating system and is available in matching shroud colors. The back shroud is not available when Omit Back Pan/Cushion is selected. Reference colors shown on page 2.Red Corvette
Part: XRF000019Grey Street
Part: XRF000023Pretty in Pink
Part: XRF000027Electric Blue
Part: XRF000030Green Machine
Part: XRF000033Back in Black
Part: XRF000021Deep Purple
Part: XRF000025Evolution Orange
Part: XRF000028Blue Jean
Part: XRF000031Military Green
Part: XRF000034Night Train
Part: XRF000022Raspberry Beret
Part: XRF000026Midnight Blue
Part: XRF000029Lemon Crush
Part: XRF000032**13.14. Lateral Positioning Wedge****Optional.**Lateral Positioning Wedge, pair \$15 ea.
Part: SET123465

Available with Sport and TRU-Comfort backs. Not available on power tilt backs 16"-17"H with limited adjustability on 18" H. Not available on power recline back 19"-20"H with limited adjustability on 21"H.

13.15. Back Width Package**REQUIRED** except when Omit Back Pan/Cushion is selected. Available with +/-" of seat width.Match Seat Width Selections
Part: PTOINDV346814"W-16"W Package
Part: XRF00003718"W-22"W Package
Part: XRF00003910"W-14"W Package
Part: XRF00003616"W-20"W Package
Part: XRF00003820"W-24"W Package
Part: XRF000040**13.16. Back Width****REQUIRED** except when Omit Back Pan/Cushion is selected. Available with +/-" of seat width.12" Back Width
Part: PTOINDV356216" Back Width
Part: PTOINDV356419" Back Width
Part: PTOINDV356722" Back Width
Part: PTOINDV357014" Back Width
Part: PTOINDV356317" Back Width
Part: PTOINDV356520" Back Width
Part: PTOINDV356824" Back Width
Part: PTOINDV357115" Back Width
Part: PTOINDV385618" Back Width
Part: PTOINDV356621" Back Width
Part: PTOINDV3569

TRU Balance 3 Seating section continued on following pages.

13. TRU BALANCE 3 SEATING CONTINUED**13.17. Back Cushion Configuration - Overall Back Height**

Static Seat and Power Tilt available height range: 14"-29"H. *Odd heights include 4" gap and even heights include 3" gap between seat pan and bottom of back cushion.*
Power Recline available height range: 19"-32"H. *Odd heights include 6" gap and even heights include 7" gap. See the Sport and TRU-Comfort back availability chart for stock order back size ranges.*

14" Height Part: SETHGHT1027	19" Height Part: SETHGHT1032	24" Height Part: SETHGHT1037	29" Height Part: SETHGHT1042
15" Height Part: SETHGHT1028	20" Height Part: SETHGHT1033	25" Height Part: SETHGHT1038	30" Height Part: SETHGHT1043
16" Height Part: SETHGHT1029	21" Height Part: SETHGHT1034	26" Height Part: SETHGHT1039	31" Height Part: SETHGHT1102
17" Height Part: SETHGHT1030	22" Height Part: SETHGHT1035	27" Height Part: SETHGHT1040	32" Height Part: SETHGHT1103
18" Height Part: SETHGHT1031	23" Height Part: SETHGHT1036	28" Height Part: SETHGHT1041	

13.18. Power Tilt Privacy Flap

Optional. A privacy flap is standard with power recline systems. This selection is to add a flap for power tilt only systems.

Privacy Flap \$40
 Part: XRF000351

14. TRU BALANCE 3 ARMRESTS

REQUIRED. Select the desired type of armrest below. Please read notes to ensure seating compatibility. Be sure to also select an armrest size and armpad type.

14.1a. Flip-Back, 8"-16" H Armrests

Choose a left and right armrest and set height in the following section. Not available with waterfall armpads. Not available with Back Cane Option.

☐ ☐ Flip-Back, 8"-16"H Armrests\$145 ea. Left part: ARM1806662. Right part: ARM1806663 Not available with TRU Balance 3 Static Seat.	☐ ☐ Flip-Back, 8"-16"H Armrests for Static Seat\$145 ea. Left part: XRF000320. Right part: XRF000321
---	--

Armrest Support

Optional. Available up to 2" difference between seat and back width via endomorph or mesomorph configuration. Not available less than 16" width seat or back.

☐ ☐ Flip-Back Armrest Downpost \$115 ea. Left part: XRF000265. Right part: XRF000266
--

14.1b. Drop-In, Height Adjustable, Removable, T-Arm Armrests

Only available with TRU Balance 3 Static Seat. Includes flat armpads.

☐ ☐ Full Length, Drop-In Armrests No Charge Left part: ARM1803253. Right part: ARM156742	☐ ☐ Desk Length, Drop-In Armrests No Charge Left part: ARM1803254. Right part: ARM156743
--	--

14.1c. Two-Post, Height Adjustable, Removable Armrests

Select the armrest and set the left and right heights in the following section. These armrests may interfere with offset backs/ back canes.

1a. Two-Post Armrest for Power Seating

☐ ☐ Full Length, Two-Post Armrest\$145 ea. Left part: XRF000047. Right part: XRF000048	☐ ☐ Desk Length, Two-Post Armrest\$145 ea. Left part: XRF000049. Right part: XRF000050
--	--

14.1d. Height Adjustable, Removable, Heavy Duty Armrests

Select a type of HD removable armrests and the set height in the following section. Only one height package is selectable as a normal configuration. Payer coding, coverage, and bundling guidelines may apply.

8"-12"H, Single Post, Height Adjustable, Removable,\$130 ea. Heavy Duty (Sold in pairs) Part: XRF000055. HCPCS w/ Static Seat: E0973	10"-14"H, Single Post, Height Adjustable, Removable,\$130 ea. Heavy Duty (Sold in pairs) Part: XRF000056. HCPCS w/ Static Seat: E0973
--	---

14.2. Armrest Height

REQUIRED. Set left and right armrest height. Note the available height range in the option description.

1. Left Armrest Height

8" H PTOINDV3535	10" H PTOINDV3537	12" H PTOINDV3539	14" H PTOINDV3541	16" H PTOINDV3543
9" H PTOINDV3536	11" H PTOINDV3538	13" H PTOINDV3540	15" H PTOINDV3542	

TRU Balance 3 armrests section continued on following pages.

14. TRU BALANCE 3 ARMRESTS CONTINUED**2. Right Armrest Height**

8" H PTOINDV3450	10" H PTOINDV3452	12" H PTOINDV3454	14" H PTOINDV3456	16" H PTOINDV3458
9" H PTOINDV3451	11" H PTOINDV3453	13" H PTOINDV3455	15" H PTOINDV3457	

14.3. Armrest Size**REQUIRED.** Select a left (L) and right (R) armrest size. This also affects armpad selection.

☐ L ☐ R 10" Length L part: XRF000057. R part: XRF000123	☐ L ☐ R 14" Length L part: XRF000058. R part: XRF000124	☐ L ☐ R 17" Length L part: XRF000059. R part: XRF000125 Only available with Flip-Back, 8"-16" H Armrests and straight arm pads.
--	--	---

15. TRU BALANCE 3 ARMPADS**REQUIRED.** Select to omit an armpad selection or select a left and right armpad type. Please note that waterfall and gel armpads may cause interference with the seat back.

No Armpad Selection..... No Charge Part: XRF000060 Not available with Flip-back, 8"-16" H armrests.	Interface for Aftermarket Pad No Charge Part: XRF000308 Only available with Flip-back, 8"-16" armrests.
---	---

15.1a. Foam Armpads

☐ L ☐ R Straight, Armpad Standard Left part: XRF000061. Right part: XRF000126	☐ L ☐ R Waterfall, Armpad No Charge Left part: XRF000062. Right part: XRF000127 Not available with Flip-back, 8"-16" armrests.
--	--

15.1b. Gel Armpads

Sold only in pairs. Armpads are 3" wide. Manufacturer part number for full length is 314G and desk length is 310G.

a. Falcon Gel Armpads for Flip-Back Armrests

Falcon Gel, Full Length Left and Right \$300 Part: ARM1806807	Falcon Gel, Desk Length Left, Full Length Right..... \$250 Part: ARM1806810
Falcon Gel, Full Length Left, Desk Length Right..... \$250 Part: ARM1806813	Falcon Gel, Desk Length Left and Right \$250 Part: ARM1806806

b. Falcon Gel Armpads for Tube Style Armrests

Falcon Gel, Full Length Left and Right \$300 Part: ARM136438	Falcon Gel, Desk Length Left, Full Length Right..... \$250 Part: ARM136439
Falcon Gel, Full Length Left, Desk Length Right..... \$250 Part: ARM136442	Falcon Gel, Desk Length Left and Right \$200 Part: ARM136437

16. TRU BALANCE 3 ARMREST BRACKETS**16.1. Joystick Bracket****REQUIRED.**

Standard Fixed Receiver Standard Part: XRF000063	Stealth Height Adjustable Receiver..... \$194 Part: XRF000064. HCPCS: K0108
---	--

17. TRU BALANCE 3 HEADRESTS**Optional.** Select a headrest pad and bracket.**17.1. Headrests, Comfort Plus with Cool Core**

Only available with 14" to 22" back widths of Static Type or Flip Forward back, without occupied transit.

1. Headrest Pads

Stealth 8" Comfort Plus, Cool Core (CP08CCFPM) \$279 Part: ST-CP08CCFPM-Q. HCPCS: E0955	Stealth 14" Comfort Plus, Cool Core (CP14CCFPM)..... \$279 Part: ST-CP14CCFPM-Q. HCPCS: E0955
Stealth 10" Comfort Plus, Cool Core (CP10CCFPM)..... \$279 Part: ST-CP10CCFPM-Q. HCPCS: E0955	Stealth 6", Comfort Plus, Cool Core (CP06CCFPM)..... \$279 Part: ST-CP06CCFPM-Q. HCPCS: E0955

2. Headrest Brackets

Comfort Plus Headrest Mounting Bracket, UniLink..... \$485 Part: ST-UHB2D33200-Q. HCPCS: E1033	Multi-Axis Removable Headrest Bracket..... \$271 Part: ST-HMO475-17M-Q. HCPCS: E1033
---	---

17.2. Q Cam Lock**Optional.**

Q Cam Lock..... No Charge Part: XRF000376
--

Q Cam Lock Availability (Overall height)	TRU-Comfort/Sport		Q-Back 2	
	Standard	\$50	Standard	\$50
Tilt/Static/Back Canes	14"-21"H	22"-29"H	19"-24"H	25"-26"H
Power Recline Back	19"-21"H	22"-32"H	22"-25"H	26"-29"H

Account Number: _____ Date: _____ Quote Name/Marked For: _____

18. LAP BELTS

REQUIRED. Select a lap belt below.

18.1a. Standard Lap Belts

- ☐ **38" Lap Belt** Part: ACCBELT1006 ☐ **50" Lap Belt** Part: ACCBELT1000 ☐ **60" Lap Belt** Part: ACCBELT1002 ☐ **70" Lap Belt** Part: ACCBELT1003 ☐ **80" Lap Belt** Part: ACCBELT1004 ☐ **90" Lap Belt** Part: ACCBELT1001

18.1b. Optional Lap Belts

Stealth Retractable Lap Belt (IFRB5350).....\$175
Part: ST-IFRB5350

Stealth Positioning, 60" Single Pull, Push Button, Medium..... \$134
(7.5") Padded, 1.5" Belt (IFPB5213)
Part: ST-IFPB5213-Q

Stealth Positioning, 60" Single Pull, Push Button, Large \$134
Padded (9.5"), 2" Belt (IFPB5314)
Part: ST-IFPB5314-Q

Stealth Positioning, 60" Single Pull, Push Button, Extra Large. \$134
Padded (11.5"), 2" Belt (IFPB5315)
Part: ST-IFPB5315-Q

60" Sure-Lok Lap Belt (Non-retractable)..... \$200
Required with ISO 7176-19 Compliant Transit Option.
Part: ACCBELT1016

19. TRU BALANCE 3 LEGREST CONTROLS

Optional. Specific power legrest options are selected later in this form. Select the control type here, if power legrests are desired, to assist the selection of the power positioning electronics.

- ☐ **Independent Legrests (adds 2 actuator control requirement)** Part: XRF000108 ☐ **Articulating Foot Platform (adds 1 actuator control requirement)** Part: XRF000110
- ☐ **Combined Legrests (adds 2 actuator control requirement)** Part: XRF000109

20. LEGRESTS

Optional. All efforts will be made to accommodate legrest length selection, but minor changes may need to be made to meet shipping guidelines. Some adjustments may be needed upon receipt of the unit. Reference the section's compatibility matrix for legrest interference.

20.1a. Basic Center Mount Foot Platform

Select legrest length in following section.

1. Center Mount Foot Platform Type

Basic Center Mount.....\$200
Part: XRF000323

Short	Medium
8"-10"	10"-14"
Basic Center Mount Extension Range	

2. Basic Center Mount Foot Platform Extension and Angle Adjustable Footplates

Short Extension, 5"x9" FootplatesNo Charge
Part: MEC1906675

Medium Extension, 5"x9" FootplatesNo Charge
Part: MEC1906680

3. Basic Center Mount Calf Pads

Calf pads are standard with a power articulating foot platform.

Small Calf Pads, Basic Center Mount.....\$210
Part: RIG2104369

4. Basic Center Mount Additional Accessories

Basic Center Mount Footplate Connector \$55
Part: RIG2104240

20.1b. Center Mount Foot Platform

If power articulating foot platform is selected, be sure power legrest control has been selected in earlier section. Select legrest length in following section.

1. Center Mount Foot Platform Type

Power Articulating Foot Platform \$3,145
Part: PTOINDV3754. HCPCS: E1012+K0040

Quick Adjustable Center Mount Foot Platform \$265
Part: XRF000311. HCPCS: K0040

2. Extension and Angle Adjustable Footplates

Pediatric extensions are not compatible with Power Articulating Foot Platform.

Pediatric Extension, 5"x6" Footplates
Part: RIG152577

Medium Extension, 5"x11" Footplates
Part: RIG2404018

Pediatric Extension, 5"x8" Footplates
Part: RIG2404015

Long Extension, 5"x6" Footplates
Part: RIG151505

Medium Extension, 5"x6" Footplates
Part: RIG150206

Long Extension, 5"x8" Footplates
Part: RIG2404017

Medium Extension, 5"x8" Footplates
Part: RIG2404016

Long Extension, 5"x11" Footplates
Part: RIG2404019

Medium	Long
11"-17"	14"-20"
Power AFP Extension Range	

Pediatric	Medium	Long
8.5"-12.5"	11"-16.5"	14"-20"
Quick Adj Center Mount Extension Range		

3a. Center Mount Calf Pads

Calf pads are standard and required with a power articulating foot platform.

Small Calf Pads, AFP Only, 6"Wx6"H (pair) Standard
Part: RIGASMB7120018

Small Calf Pads, Quick Adjustable Center Mount Only, \$210
6"Wx6"H (pair)
Part: XRF000156

3b. Center Mount Calf Panel

Small Calf Panel, 18"Wx8"H.....\$210
Part: FRMASMB7110012

Legrests section continued on following pages.

20. LEGRESTS CONTINUED**4. Center Mount Modular Footplate Extension Components**

Optional. Assembly required. May require seating configuration adjustments to avoid front caster interference.

- 2" Bell-Shaped Length Extension \$55
 Part: RIG152519
 Only available with the 5x6 footplates.
- Width Extension for 8" \$55
 Part: RIG2404919
- Width Extension for 11" \$55
 Part: RIG2404920

- Outside Edge \$55
 Part: RIG2404921
 Available with 5"x8" or 5"x11" footplates.
- Front and Rear Cap \$55
 Part: RIG2404923
 Available with 5"x8" or 5"x11" footplates.

5. Center Mount Additional Accessories

- Small Foam Calf Pad Cover (pair) \$150
 Part: ST-CSCS0775-Q
- Small Foam Panel Cover (pair) \$150
 Part: ST-CSCP0920-Q
- Footplate Connector \$55
 Part: RIGASMB7120023

- Footplate Riser \$150
 Part: XRF000348
 Select legrest length between 2" to 9" at 1" increments.
- Footplate Riser Calf Panel No Charge
 Part: FRMASMB7110011
 Only available and required with Footplate Riser selection and extension lengths of 6" to 9".

6. Stealth Gel Footplate Covers

- 5"x8" Tapered Gel Footplate Covers \$275
 Part: ST-FPG5008T-Q. HCPCS: K0108
- 5"x11" Tapered Gel Footplate Covers \$275
 Part: ST-AFPG5011T-Q. HCPCS: K0108

- 5"x8" Tapered Gel Footplate Covers for Width Extension \$275
 Part: ST-AFPG5008E-Q. HCPCS: K0108
- 5"x11" Tapered Gel Footplate Cover for Width Extension \$275
 Part: ST-AFPG5011E-Q. HCPCS: K0108

20.1c. Swing-Away Legrests

Select a legrest upper, a legrest extension, and type of footplate. Pair (P), left (L), and right (R) configurations are selected by checking the appropriate boxes.

1. 60°-90° Swing-Away Legrest Upper

Adjustable in 5° increments.

- ☐ P ☐ L ☐ R 60°-90° Swing-Away Legrest \$285
 Pair part: RIG120062. Left: FRMASMB17558.
 Right: FRMASMB17559
 Not available with 11"-11" or 11"-16.5" extensions.
- ☐ P ☐ L ☐ R 60°-90° Swing-Away Legrest, Short \$285
 Pair part: RIG136481. Left: RIG132607.
 Right: RIG132609.
 For use only with 11"-13" or 11"-16.5" lower extensions.

60°-90° Swing-Away Legrest Compatibility Matrix	Power Recline or Static Seat			Power Tilt			Power Lift		
	Low	Med.	High	Low	Med.	High	Low	Med.	High
Extruded, 11"-13"	!	✓	✓	✓	✓	✓	✓	✓	✓
Extruded, 13"-16.5"	!	✓	✓	✓	✓	✓	✓	✓	✓
Extruded, 15.5"-19.5"	!	!	!	!	!	!	!	!	✓
Angle Adjustable, 13"-16.5"	!	✓	✓	✓	✓	✓	✓	✓	✓
Angle Adjustable, 15.5"-19.5"	!	!	!	!	!	✓	!	✓	✓
Multi-Axis, 11"-16.5"	!	✓	✓	✓	✓	✓	✓	✓	✓
Multi-Axis, 15.5"-19.5"	!	!	!	!	!	!	!	!	✓
11"-13" extension causes no interference at all seat-to-floor heights.									

- ☒ This rigging configuration causes no interference issues at this seat-to-floor height.
- ☐ This rigging configuration may have less than 1.75" clearance when fully extended and may interfere with the ground or reduce obstacle climbing ability.

2. Lower Extension

- ☐ P ☐ L ☐ R 11"-13" Length
 Pair part: XRF000088.
 Left: XRF000089.
 Right: XRF000090.
- ☐ P ☐ L ☐ R 13"-16.5" Length
 Pair part: XRF000094.
 Left: XRF000095.
 Right: XRF000096.
- ☐ P ☐ L ☐ R 15.5"-19.5" Length
 Pair part: XRF000097.
 Left: XRF000098.
 Right: XRF000099.

3. Angle Adjustable Footplate(s)¹The 6"x11" footplates may cause interference with seat widths less than 18". ²The 8"x11" footplates may cause interference with seat widths less than 20".

- ☐ P ☐ L ☐ R 5"Wx8"D Footplates with Heel Loops \$87.50 ea.
 Pair part: FRMASMB4735. Left: RIG137791. Right: RIG137792
 HCPCS: E0951
- ☐ P ☐ L ☐ R 6"Wx11"D¹ Footplates with Heel Loops \$87.50 ea.
 Pair part: FRMASMB4741. Left: RIG137796. Right: RIG137797
 HCPCS: E0951
- ☐ P ☐ L ☐ R 6"Wx8"D¹ Footplates with Heel Loops \$87.50 ea.
 Pair part: FRMASMB4738. Left: RIG137793. Right: RIG137795
 HCPCS: E0951
- ☐ P ☐ L ☐ R 8"Wx11"D² Footplates with Heel Loops \$87.50 ea.
 Pair part: FRMASMB4744. Left: RIG137798. Right: RIG137799
 HCPCS: E0951

Legrests section continued on following pages.

20. LEGRESTS CONTINUED**20.1d. Power Elevating Legrests**

Be sure power legrest control has been selected in previous section. Select an upper extension type and footplate. Pair (P), left (L), and right (R) configurations are selected by checking the appropriate boxes. Select legrest length in following section.

1. Power Elevating Legrest Upper Assembly

Includes calf pads.

Power Elevating Legrests, Pair \$2,865

Part: XRF000341. HCPCS: E1010+K0040+E0951

Legrest Length Matrix	Seat to Floor Option		
	Low	Middle	High
Power ELR: Static Seat, Recline Only			
8.5"-14.5" Lower Length	8.5"- 12.5"	8.5"- 13.5"	8.5"- 14.5"

Legrest Length Matrix	Seat to Floor Option		
	Low	Middle	High
Power ELR: Tilt Only, Tilt and Recline, Any Power Lift			
8.5"-14.5" Lower Length	8.5"- 14"	8.5"- 14.5"	8.5"- 14.5"

2. Lower Extension

8.5" to 14.5" Lower Length

Part: XRF000344

3. Angle Adjustable Footplate with Heel Loop and Wheel Bumper

*The 6"W footplates may cause interference with seat widths less than 16". *The 8"W footplates may cause interference with seat widths less than 19".

4"Wx6"D, Pair

Part: FRMASMB7448

5"Wx8"D, Pair

Part: FRMASMB7450

6"Wx11"D¹, Pair

Part: FRMASMB7452

5"Wx6"D, Pair

Part: FRMASMB7449

6"Wx8"D¹, Pair

Part: FRMASMB7451

8"Wx11"D², Pair

Part: FRMASMB7453

20.1e. Elevating Legrests

Select a legrest upper, a legrest extension, and type of footplate. Pair (P), left (L), and right (R) configurations are selected by checking the appropriate boxes. Select legrest length in following section.

1. Elevating Legrests Upper

☐ ☐ ☐ Elevating Legrest Assembly \$275 ea.

Pair part: FRMASMB6483. Left: RIG138012

Right: RIG138013. HCPCS: E0990

Select lower extension and footplates

2a. Lower Extension

☐ ☐ ☐ 11"-14" Length
Pair: FRMASMB4755
Left: RIG138026
Right: RIG138030
For reinforced only

☐ ☐ ☐ 14"-17" Length
Pair: FRMASMB4756
Left: RIG138033
Right: RIG138035
For reinforced only

☐ ☐ ☐ 8.5"-12" HD Drop-Ins/
12.5"-16.5" ELR,
Ext Length
Pair part: XRF000101.
Left: XRF000102.
Right: XRF000103.
For angle adjustable only

☐ ☐ ☐ 13"-18" HD Drop-Ins/
15.5"-19.5" ELR,
Ext Length
Pair part: XRF000104.
Left: XRF000105.
Right: XRF000106
For angle adjustable only

Elevating Legrest Compatibility Matrix	Power Recline or Static Seat			Power Tilt			Power Lift		
	Low	Med.	High	Low	Med.	High	Low	Med.	High
Reinforced, 11"-14"	✓	✓	✓	✓	✓	✓	✓	✓	✓
Reinforced, 14"-17"	!	!	✓	!	✓	✓	✓	✓	✓
Angle Adjustable, 12"-16.5"	✓	✓	✓	✓	✓	✓	✓	✓	✓
Angle Adjustable, 15.5"-19.5"	!	!	!	!	!	✓	!	✓	✓
 This rigging configuration causes no interference issues at this seat-to-floor height.									
 This rigging configuration may have less than 1.75" clearance when fully extended and may interfere with the ground or reduce obstacle climbing ability.									

3a. Reinforced Cast Footplate(s)

☐ ☐ ☐ 7.25"Wx6"D Reinforced Footplate..... Standard
Pair part: FRMASMB4748. Left: RIG138055. Right: RIG138056

☐ ☐ ☐ 8.5"Wx6"D Reinforced Footplate..... Standard
Pair part: FRMASMB4747. Left: RIG138060. Right: RIG138062
May cause interference on seat widths less than 18."

3b. Angle Adjustable Footplate(s)

*The 6"x11" footplates may cause interference with seat widths less than 18". *The 8"x11" footplates may cause interference with seat widths less than 20".

☐ ☐ ☐ 5"Wx8"D Footplates with Heel Loops \$87.50 ea.
Pair part: FRMASMB4735. Left: RIG137791. Right: RIG137792
HCPCS: K0040+E0951

☐ ☐ ☐ 6"Wx11"D¹ Footplates with Heel Loops \$87.50 ea.
Pair part: FRMASMB4741. Left: RIG137796. Right: RIG137797
HCPCS: K0040+E0951

☐ ☐ ☐ 6"Wx8"D¹ Footplates with Heel Loops \$87.50 ea.
Pair part: FRMASMB4738. Left: RIG137793. Right: RIG137795
HCPCS: K0040+E0951

☐ ☐ ☐ 8"Wx11"D² Footplates with Heel Loops \$87.50 ea.
Pair part: FRMASMB4744. Left: RIG137798. Right: RIG137799
HCPCS: K0040+E0951

20.2. Legrest Accessories

Heel Loops (pair)..... \$35

Part: FRMASMB7873. HCPCS: E0951

Standard on angle adjustable and multi axis footplates when noted.

Wheel Bumper..... \$15 ea.

Part: XRF000100

Legrests section continued on following pages.

20. LEGRESTS CONTINUED**20.3. Set Legrest Extension Length****REQUIRED** with center mounts and matching left and right rigging configurations. Part number is PTOINDV plus number following # below.**1a. Centermount Extension Length**

2" L (#3989)	4" L (#3991)	6" L (#3993)	8" L (#3680)	9" L (#3682)	10" L (#3684)	11" L (#3686)	12" L (#3688)	13" L (#3690)	14" L (#3692)	15" L (#3694)
3" L (#3990)	5" L (#3992)	7" L (#3994)	8.5" L (#3681)	9.5" L (#3683)	10.5" L (#3685)	11.5" L (#3687)	12.5" L (#3689)	13.5" L (#3691)	14.5" L (#3693)	

1b. Front Rigging Extension Length

Non-matching or individual legrest configurations lengths are not set by the manufacturer. See previous section for center mount extension lengths.

7" L (#3470)	8" L (#3472)	9" L (#3474)	10" L (#3476)	11" L (#3478)	12" L (#3480)	13" L (#3482)	14" L (#3484)	15" L (#3486)
7.5" L (#3471)	8.5" L (#3473)	9.5" L (#3475)	10.5" L (#3477)	11.5" L (#3479)	12.5" L (#3481)	13.5" L (#3483)	14.5" L (#3485)	

21. TRU BALANCE 3 POWER POSITIONING ELECTRONICS**REQUIRED** with power seating or legrest options. Select one option from section 21.1 or 21.3.**21.1a. Tilt Only through Toggle**

Available with QLNES only.

Single Actuator Control, Power Tilt Only, Through Toggle No Charge

Part: ELE1905907

Select toggle side. May require extended lead time when selected with Fender Lights.

Left Side Toggle

Part: PTOINDV1839

Right Side Toggle

Part: PTOINDV1838

21.1b. Actuator Control through Joystick**1. Actuator Control Selection**

Single Actuator Through Q-Logic 3 Controller and Interface (AM3) No Charge

2 Actuators Through Q-Logic 3 Controller and Interface (AM3) \$1,250

Not funded by ADP

3 Actuators Through Q-Logic 3 Controller and Interface (AM3) \$1,500

Not funded by ADP

Up to 5 Actuators Through Q-Logic 3 Controller and Interface (AAM) \$1,675

ADP: SECX15005

2. Q-Logic 3 Optional Switch Control

Only available with Q-Logic 3 AAM, multiple actuator configurations.

4-Function Push Button \$225

Part: CTL123218

See next section for switch mounting.

Only available with AAM.

4-Way Toggle \$575

Part: ELEASMB5925

Only available with AAM.

4-Way Button Toggle \$575

Part: ELEASMB5926

Only available with AAM.

3. User Controls Side Selection**Required with 4-Way Toggle and 4-Way Button Toggle.**

Mount User Controls on Left

Part: PTOINDV3996

Mount User Controls on Right

Part: PTOINDV3997

21.1c. Actuator Control with iAccess Switch

iAccess Flat Button \$575

Part: CTL144009

See next section for switch mounting.

iAccess Toggle Button \$575

Part: CTL144010

See next section for switch mounting.

21.2. 4-Function Push Button/iAccess Mount**REQUIRED** with 4-Function Push Button or iAccess options. Under arm mount is not available with TRU Balance 3 HD armrests.

Joystick Mount No Charge

Part: MEC1907500

Attendant Mount \$35

Part: MEC1906564

Under Arm Mount, Flip-Back Armrest No Charge

Part: MEC1906531

Under Arm Mount, Flip-Back Armrest, plus iAccess Guard \$35

Part: MEC2001392

Under Arm Mount, Tube Style Armrest No Charge

Part: MEC1906530

Under Arm Mount, Tube Style Armrest, plus iAccess Guard \$35

Part: MEC2001391

1. Under Arm Mount Side Selection

Left Side Under Arm Mount

Part: PTOINDV3678

Right Side Under Arm Mount

Part: PTOINDV3679



22. TRU BALANCE 3 ATTENDANT CONTROL**Optional.** *Attendant control will be at no charge when ordered with Tilt and Recline. Mounting bracket is included.

Attendant Control \$950
 Part: XRF000078. HCPCS: E2331

23. SEAT MOUNTED LED LIGHTS**Optional.** Requires Q-Logic 3 Joystick for Full Light Package or Q-Logic 3e Joystick.

Seat Mounted, LED, Full Light Package \$770
 Part: XRF000130
 Not available with Base Mounted, LED, Full Light Package.

**24. PROVIDER INSTALLED ACCESSORIES****Optional.** Multiple accessories can be mounted along the various tracks but some interference may occur.

Accessory Bar \$90
 Part: XRF000128
 Horizontal rail attached to rear of seat. Standard with the Trilogy vent tray and oxygen holder.

Backpack Holder (Pair) (A) \$42
 Part: ACC147774

Cane/Crutch Holder \$195
 Part: ACC152515. HCPCS: E2207

Cell Phone/ Media Holder, X-Grip (E) \$95
 Part: ACC155835
 Not available with Heavy Duty armrests.

Clothing Guard (Pair) \$275
 Part: ACC125007

Compact Vent Tray \$1,395
 Part: ACC133889. HCPCS: E1029

Cup Holder \$20
 Part: ACC125003
 For use with track-style armrests.

Cup Holder (Tube-style Armrest) \$20
 Part: ACCASMB2634
 For use with tube-style armrests.

Flag Kit \$125
 Part: ACC1700443

Flip-Back Armrest, Side Track Accessory Bracket .. \$25
 Part: MEC1809178

Glove Box (C) \$85
 Part: ACC125011

Grip Drink Holder (20085) \$100
 Part: ACC152993

Grip Lap Board (20005) \$100
 Part: ACC152992

Hydration System (B) \$328
 Part: ACC125010

Large, Black, Weather Cover \$240
 Part: ACCCOVR1021

Oxygen Holder \$110
 Part: ACC125006. HCPCS: E2208
 Includes accessory bar, a credit will be applied PTOIN-DV3769.

Personal Item Hook \$20
 Part: ACC125008

Personal Item Pouch \$175
 Part: ACC2101017

Privacy Bag \$80
 Part: ACC2005723

Push Handles (pair) \$140
 Part: ACC2005027
 Not available with oxygen holder or vent tray.

Q-Logic EX/NE Programmer \$625
 Part: ELEASMB6976

Quantum Backpack \$200
 Part: ACC1905853

Rear View Mirror \$35
 Part: ACC155833


Reflector Kit \$35
 Part: ACC124670

Stealth Large Tablet Holder \$565
 Part: ACC141818

Transfer Handle (D) \$152.50 ea.
 Part: MEC2205680 x2. HCPCS: E0705

Depth Adjustable Vent Tray \$1,000
 Part: ACC2302500. HCPCS: E1029
 Includes accessory bar, a credit will be applied PTOIN-DV3769.

Remote Stop Switch \$1,495
 Part: XRF000379

**Aftermarket Accessory Mounting Bracket**

Select a quantity up to 5.

Aftermarket Accessory Mounting Bracket, Slotted \$100 ea.
 Part: ACC2407046

QTY

Aftermarket Accessory Mounting Bracket, Non-slotted \$100 ea.
 Part: ACC2405496

QTY

26. POSITIONING COMPONENTS**Optional.** More options are available on the [Positioning Components Order Form](#). Please provide quantity next to check box where applicable.**26.1. Stealth Thoracic Laterals**

Swing-away laterals may interfere with TRU Balance 3 armrests and 2-post, flip up armrests. These laterals are not available with the back cane option.

1. Bracket

Swing-away, UniLink, Redesign Back, pair (ULS4042000-Q)... \$1032
 Part: ST-ULS4042000-Q. HCPCS: E1034
 Only available with Redesigned back assemblies.

Swing-away, UniLink, 125° w/Track, pair (ULPK404108B00) ... \$966
 Part: ST-ULPK404108B00-Q. HCPCS: E1034
 Only available with Sport or TRU Comfort backs.

Swing-away, UniLink, Back Cane, pair (ULPK404104B00) \$966
 Part: ST-ULPK404104B00-Q. HCPCS: E1034
 Only available with back cane and omit back configurations.

Swing-away, UniLink, Q-Back 2, Left (ULPK404103L00) \$483
 Part: ST-ULPK404103L00-Q. HCPCS: E1034
 Only available with Active Contour Q-Backs.

Swing-away, UniLink, Q-Back 2, Right (ULPK404103R00) \$483
 Part: ST-ULPK404103R00-Q. HCPCS: E1034
 Only available with Active Contour Q-Backs.

2. Pads

Curved Pads, 3"x5", pair (TWBL3X5)..... \$164
 Part: INDPART2663. HCPCS: E0956

Curved Pads, 4"x6", pair (TWBL4X6)..... \$
 Part: INDPART2664. HCPCS: E0956

3. Installation

Install Laterals
 Part: PTOINDV3499

Do Not Install
 Part: PTOINDV3500

26.2. Stealth UniLink Pelvic/Thigh Guide**1. UniLink Pelvic/Thigh Hardware**

^ Please verify the appropriate billing code depending upon the placement of the support pad.

UniLink, 5" Height Adjustable, \$405 ea.
 Removable Mounting Hardware
 QTY Part: ST-UTTP5-Q. HCPCS: E1028/E1034^

UniLink, 7" Height Adjustable, \$423 ea.
 Removable Mounting Hardware
 QTY Part: ST-UTTP7-Q. HCPCS: E1028/E1034^

UniLink Multi-Angle Thigh Guide, Left \$585 ea.
 Part: ST-UTT8A11602-L-Q. HCPCS: E1028/E1034^

UniLink Multi-Angle Thigh Guide, Right..... \$585 ea.
 Part: ST-UTT8A11602-R-Q. HCPCS: E1028/E1034^

UniLink Flip Down Hip/Thigh Support, Left..... \$483 ea.
 Part: ST-UTS25100-L-Q. HCPCS: E1028/E1034^

UniLink Flip Down Hip/Thigh Support, Right \$483 ea.
 Part: ST-UTS25100-R-Q. HCPCS: E1028/E1034^

Flip Down with LEFT Full Surface Contact \$459 ea.
 Feature (TWBADD-L & LPHW407)
 Part: POS140428. HCPCS: E1028/E1034^

Flip Down with RIGHT Full Surface Contact..... \$459 ea.
 Feature (TWBADD-R & LPHW407)
 Part: POS140429. HCPCS: E1028/E1034^

2. UniLink Pelvic/Thigh Pads

^ Please verify the appropriate billing code to be used with all third-party payors.

UniLink, 4" x 4" Pad \$164 ea.
 QTY Part: ST-UTP544056100-Q. HCPCS: E0953/E0956^

UniLink, 4" x 8" Pad..... \$164 ea.
 QTY Part: ST-UTP548056100-Q. HCPCS: E0953/E0956^

UniLink, 4" x 12" Pad \$164 ea.
 QTY Part: ST-UTP54D056100-Q. HCPCS: E0953/E0956^

3. UniLink Pelvic/Thigh Additional Links

UniLink, 1.5" Link \$121 ea.
 QTY Part: ST-UIL001150A000-Q

UniLink, 2" Link..... \$121 ea.
 QTY Part: ST-UIL001200A000-Q

4. UniLink Pelvic/Thigh Guide Installation

Install on Left
 Part: PTOINDV3585

Install on Right
 Part: PTOINDV3586

Do Not Install
 Part: PTOINDV3719

26.3. Stealth Positioning Chest Harness

Upper Classic, Medium (IFPS1103) \$221
 Part: ST-IFPS1103-Q. HCPCS: E0960

Upper Classic, Large (IFPS1104) \$221
 Part: ST-IFPS1104-Q. HCPCS: E0960

Upper Contour, Medium (IFPS2103)..... \$221
 Part: ST-IFPS2103-Q. HCPCS: E0960

Upper Contour, Large (IFPS2104) \$221
 Part: ST-IFPS2104-Q. HCPCS: E0960

**26.4. Stealth Positioning Chest Straps**

Chest Strap, Medium w/ Belt Anchor, \$140
 (IFCS4113)
 Part: ST-IFCS4113-Q. HCPCS: E0960

Chest Strap, Large w/ Belt Anchor, \$140
 (IFCS4114)
 Part: ST-IFCS4114-Q. HCPCS: E0960



Positioning components section continued on following page.

26. POSITIONING COMPONENTS CONTINUED**26.5. Stealth Swing-Away Headrest Mounts**

These swing-away headrest mounts are only available with the Stealth Comfort Series headrests.

1. Facial Lateral Hardware

UniLink Facial Lateral, Swing Away, Right.....\$496
 Part: ST-SAR2A-R-Q. HCPCS: E1028

UniLink Facial Lateral, Swing Away, Left\$496
 Part: ST-SAR2A-L-Q. HCPCS: E1028

UniLink Facial Lateral, Dual Swing Away \$765
 2 Upper & 1 Lower Right
 Part: ST-SAR4A-DBU-RL-Q. HCPCS: E1028

UniLink Facial Lateral, Dual Swing Away \$765
 2 Upper & 1 Lower Left
 Part: ST-SAR4A-DBU-LL-Q. HCPCS: E1028

UniLink Facial Lateral, Dual Swing Away \$674
 2 Upper
 Part: ST-SAR4A-DBU-Q. HCPCS: E1028

**2. Facial Lateral Pads and Accessories**


Egg Switch Mount (ESM) (A) \$70 ea.
 Part: INDPART2815

Figure 8 Bracket with 13" Rod.....\$366 ea.
 and Egg Switch Mount (B)
 Part: ST-F8675-13-Q

Facial Pad, 4".5x2", Cool Core (900M-Q) \$96 ea.
 Part: ST-900M-Q. HCPCS: K0108

Facial Pad, 6"x2", Cool Core (920M-Q)..... \$96 ea.
 Part: ST-920M-Q. HCPCS: K0108

Swing-Away Joystick Mount, 18" Rod (C).....\$100 ea.
 Part: ST-SUS18-Q
 Allows mounting of facial pads and attaches to the headrest

**26.6. Stealth Elbow Support****1. Elbow Support Bracket**

Left, Swing-Away, Elbow Support\$279
 Bracket Tube Style Arms (ARMS-ES-L)
 Part: POS124064. HCPCS: E1028

Left, UniLink, Fixed, Elbow Support Bracket..... \$309
 for Flip-Back, 8"-16"H Armrests (UETB4021FL-Q)
 Part: ST-UETB4021FL-Q

Right, Swing-Away, Elbow Support.....\$279
 Bracket Tube Style Arms (ARMS-ES-R)
 Part: POS124066. HCPCS: E1028

Right, UniLink, Fixed, Elbow Support Bracket..... \$309
 for Flip-Back, 8"-16"H Armrests (UETB4021FR-Q)
 Part: ST-UETB4021FR-Q

2. Elbow Support Pad

3"Wx3"H Elbow Pad (ARMS-ES-3X3) \$85 ea.
 Part: POS124065. HCPCS: K0108

4"Wx4"H Elbow Pad (ARMS-ES-4X4)..... \$85 ea.
 Part: POS124067. HCPCS: K0108

**26.7. Stealth Foot Boxes**

Stealth Small Foot Boxes, 5"Wx8"Dx7"H..... \$725
 (SPLE107) (pair)
 Part: INDPART2786. HCPCS: E0954

Stealth Medium Foot Boxes, 6"Wx10"Dx9"H \$725
 (SPLE108) (pair)
 Part: INDPART2787. HCPCS: E0954

Stealth Large Foot Boxes, 7"Wx12"Dx11"H..... \$725
 (SPLE109) (pair)
 Part: INDPART2788. HCPCS: E0954

26.8. Stealth Positioning Shoe Holder

Small Shoe Holder, 3/4" Webbing (IFSH4211)..... \$93 ea.
 Part: ST-IFSH4211-Q. HCPCS: E0951+E0952

Medium Shoe Holder, 3/4" Webbing (IFSH4212) \$93 ea.
 Part: ST-IFSH4212-Q. HCPCS: E0951+E0952

Large Shoe Holder, 3/4" Webbing (IFSH4213) \$93 ea.
 Part: ST-IFSH4213-Q. HCPCS: E0951+E0952

27. LAP TRAY**27.1. Quantum Lap Tray**

Mounting bracket is required with lap tray selection. These trays are compatible with Flip-Back armrests.

1a. Black Lap Tray

Black Lap Tray, 11" W \$360
 Part: MEC1905328-11B. HCPCS: E0950
 For use with 12"-13" back width.

Black Lap Tray, 13" W \$360
 Part: MEC1905328-13B. HCPCS: E0950
 For use with 14"-15" back width.

Black Lap Tray, 15" W \$360
 Part: MEC1905328-15B. HCPCS: E0950
 For use with 16"-17" back width.

Black Lap Tray, 17" W \$360
 Part: MEC1905328-17B. HCPCS: E0950
 For use with 18"-19" back width.

Black Lap Tray, 19" W \$360
 Part: MEC1905328-19B. HCPCS: E0950
 For use with 20"-22" back width.

1b. Clear Lap Tray

Clear Lap Tray, 11" W \$360
 Part: MEC1905328-11C. HCPCS: E0950
 For use with 12"-13" back width.

Clear Lap Tray, 13" W \$360
 Part: MEC1905328-13C. HCPCS: E0950
 For use with 14"-15" back width.

Clear Lap Tray, 15" W \$360
 Part: MEC1905328-15C. HCPCS: E0950
 For use with 16"-17" back width.

Clear Lap Tray, 17" W \$360
 Part: MEC1905328-17C. HCPCS: E0950
 For use with 18"-19" back width.

Clear Lap Tray, 19" W \$360
 Part: MEC1905328-19C. HCPCS: E0950
 For use with 20"-22" back width.

2. Mounting Hardware

Mounting Bracket, Chest Bar/Lap Tray No Charge
 Part: MEC2104886

27.2. Simple Solutions Lap Tray

Mounting bracket is required with lap tray selection. See chart below for compatibility.

1. Lap Tray

Simple Solutions Small Tray with Cut Out \$750
 Part: ST-SSLXS010001171-Q. HCPCS: E0950

Simple Solutions Medium Tray with Cut Out \$750
 Part: ST-SSLXM010001171-Q. HCPCS: E0950

Simple Solutions Large Tray with Cut Out \$750
 Part: ST-SSLXL010001171-Q. HCPCS: E0950

2. Mounting Hardware

Swing Away Mounting Hardware No Charge
 Part: ST-SSLXHD10000060-Q

Simple Solutions Lap Tray Armrest Compatibility	Width											
	10	11	12	13	14	16	17	18	19	20	21	22
Flip-Back	S	S	S	S	S	S	SM	ML	ML	ML	L	L
Flip-Back, Cane Mounted	S	S	S	S	S	SM	SM	ML	ML	L	L	
Drop-In, T-Arm	S	S	S	S	S	SM	SM	ML	ML	L	L	
Two-Post, Flip-Back	S	S	S	S	S	SM	SM	ML	ML	L	L	
Cantilever	S	S	S	S	S	M	M	ML	ML	L	L	
Q-Captain, Flip-Up						SM		ML		L		
Single Post/ Heavy Duty/ Bariatric	S	S	S	S	S	SM	SM	ML	ML	ML	L	L

S:ST-SSLXS010001171-Q M:ST-SSLXM010001171-Q L:ST-SSLXL010001171-Q

S:Compatible

M:Armrest configuration may not hold tray position

*Some armrests listed may not be available on this order form.

28. ALTERNATIVE CONTROLS AND ELECTRONICS**Optional.** Select a type of specialty controls below. More options are available on the [Alternative Controls and Electronics order form](#).**28.1. Enhanced Display****Q-Logic 3 Enhanced Display Kit**..... \$2,045

Part: CTL167726

Select a color option below. Includes enhanced display, flexible gooseneck and mounting bracket. Features: 3.5" LCD display; up to 5 profiles to control drive, seat, and auxiliary functions; Bluetooth connectivity for mobile devices and mouse emulation; integrated single switch scanner; IR controls.

Q-Logic 3 Enhanced Display Shroud

Required with Q-Logic 3 Enhanced Display selection above. Reference colors shown on page 2.

Black, Soft Touch Standard

Part: DGN137728

Deep Purple

Part: DGN137729-PPR5

Electric Blue

Part: DGN137729-PBU9

Military Digital Camo \$50

Part: DGN137729-WCM7

Red Corvette

Part: DGN137729-PRD1

Raspberry Beret

Part: DGN137729-PPK4

Blue Jean

Part: DGN137729-PBU11

Faded Glory \$50

Part: DGN137729-WFL2

Back in Black

Part: DGN137729-PBK1

Pretty in Pink

Part: DGN137729-PPK3

Lemon Crush

Part: DGN137729-PYL4

Carbon Fiber \$50

Part: DGN137729-WCB1

Night Train

Part: DGN137729-PBK3

Evolution Orange

Part: DGN137729-POR7

Green Machine

Part: DGN137729-PGN11

Grey Street

Part: DGN137729-PSL3

Midnight Blue

Part: DGN137729-PBU10

Military Green

Part: DGN137729-PGN8

28.2. SCIM Module**Optional.**☐ **SCIM Module for Q-Logic 3** \$850

Part: CTL2001963

Allows use of joystick display for specialty devices. Available with Q-Logic 3 EX Joystick selection above. Up to 5 profiles to control drive, seat, and auxiliary functions, excluding sip-n-puff and single switch functions. A switch is required to operate the input device attached to the SCIM module.

28.3. Q-Logic 3 SAJ**Q-Logic 3 EX SAJ** \$850

Part: CTL143980

Must be used in conjunction with Enhanced Display or Q-Logic 3 Joystick.

28.4. Head Array**i-Drive 4.0 Bluetooth Tri-array Series**..... \$7,406**3 Switch with Egg Switch**

Part: ST-IDHBT200-U-Q. HCPCS: E2330

Adult, Curved, i-Drive Tri-Array with Mounting Hardware \$7,406

Part: ST-IDHBT400-U-Q. HCPCS: E2330

i-Drive 4.0 Bluetooth Pro Series..... \$7,486**3 Switch with Egg Switch**

Part: ST-IDHBT300-U-Q. HCPCS: E2330

Junior, Curved, i-Drive Tri-Array with Mounting Hardware \$7,406

Part: ST-IDHBT450-U-Q. HCPCS: E2330

28.5. Sip and Puff Drive**1. Sip and Puff Module****Sip and Puff Module for Q-Logic 3**..... \$1,695

Part: ST-SP625M-U-Q. HCPCS: E2325

Includes module and mounting hardware.

2. Sip and Puff Straw and Mount**The straw and mount do not include electrical components to operate the Sip&Puff with a head array. This component requires the selection of a Q-Logic Sip&Puff module or additional module from Stealth Products for operation.****Stealth Sip and Puff for Comfort Series, Left (SP625M-CL)** \$496

Part: ST-SP625M-CL-Q. HCPCS: E2326

Includes straw and mounting hardware. Use with Stealth Comfort Plus headrests.

Stealth Sip and Puff for Comfort Series, Right (SP625M-CR) \$496

Part: ST-SP625M-CR-Q. HCPCS: E2326

Includes straw and mounting hardware. Use with Stealth Comfort Plus headrests.

Part: ST-SP625-U-Q. HCPCS: E2326

Includes straw and mounting hardware. Use with Stealth i-Drive Head Arrays.

3. Sip & Puff Accessories**Sip & Puff Straws, 12 pack (SP625-S12)**..... \$47

Part: ST-SP625-S12-Q



Enhanced Display



SCIM Module



EX SAJ



Adult Pro-Series



Sip & Puff Module



Stealth Sip & Puff

Alternative controls and electronics section continued on following page.

28. ALTERNATIVE CONTROLS AND ELECTRONICS CONTINUED**28.6. Egg Switches**

Red Egg Switch (ESRED)..... \$116
Part: SWTMCRO1077

Blue Egg Switch (ESBLU)..... \$116
Part: SWTMCRO1079

Black Egg Switch (ESBLK)..... \$116
Part: SWTMCRO1081

Green Egg Switch (ESGRN) \$116
Part: SWTMCRO1078

Yellow Egg Switch (ESYEL) \$116
Part: SWTMCRO1080

**28.7. Switch Adapter**

2 to 1 Switch Adapter..... \$135
Part: DWR1361H037
Connects 2 single switches to a stereo jack to act as a dual switch.

28.8. Environmental Controls

Q-Logic 3 ECU Module Kit \$1,780
Part: CTL168057

28.9. USB Mobile Device Charger

USB Mobile Device Charger \$185
Part: ELE168069

Can charge 1 device. A multiplier harness is required and must be added when a Static Seat or Power Tilt Only seating system is chosen. The selection of more than 3 modules requires an additional multiplier.

