



Pride Mobility Products Co.

5096 South Service Rd, Beamsville, ON, L3J 1V4

Phone: 888-570-1113 | Fax: 866-514-1303 | Email: productorders@pridemobility.com

Q6 Edge® 3 Stretto Order Form TRU Balance® 4 Positioning

Account Number: _____ Date: _____

Provider Name: _____

Contact: _____

Phone: _____ Fax: _____

Email: _____

PO Number: _____

Marked for: _____

Ship to Address: _____

City: _____ Prov: _____ PC: _____



INTRODUCTION

This form is interactive when viewed with [Adobe Acrobat Reader](#). Complete the form by placing checks in the desired boxes and provide information in the interactive fields. Buttons shown at the bottom of the form may be utilized to print or submit the order form through a desktop email application. To email via a web-based application, please 'Save As' and attach the PDF to your email.

This order form contains a large variety of options to fulfill various patient needs. Descriptions and section notes, such as optional or required, are included to help you complete your order. Other part numbers can be found on the quote acknowledgment after this order form is submitted. Please contact Quantum® Sales at 888-570-1113 if assistance is needed. Send the completed order form by fax (866-514-1303) or email (productorders@pridemobility.com). Incomplete forms may delay the quote or order. Customer service will contact you if the order is incomplete or if there are compatibility issues. If special order requests are needed, be sure that the Patient Information section is completed

All prices are MSRP. Prices, specifications, part numbers, and availability are subject to change without notice. Prices and part numbers as shown are only available when configured as a complete power chair. Please contact technical service for accurate parts ordering.

Thank you for choosing Quantum!


PATIENT INFORMATION

Weight: _____ Height: _____

A. Shoulder Width: _____ D. Max Sitting Height: _____ G. Thigh Depth: _____

B. Chest Width: _____ E. Shoulder Height: _____ H. Lower Leg Length: _____

C. Hip Width: _____ F. Axilla Height: _____ I. Elbow Height: _____



1. BASE MODEL

REQUIRED. 300 lb. weight capacity with standard 6mph programmable motors and off board charger. XLR USB charger is included. Dimensions of base without legrests: 20.5"Wx35.5"L

Quantum Q6 Edge 3 Stretto 3SP-SS\$8,450

Part: Q6 EDGE 3 STRETTO 3

2. DRIVE MOTORS

Optional. Requires Q-Logic 3 electronics. The standard 6 mph motor can be upgraded to Accu-Trac motors that delivers enhanced traction and tracking, especially with switched drive controls.

Accu-Trac Motor Technology\$995

Part: XRF000268. HCPCS: K0108

Account Number: _____ Date: _____ Marked for: _____

3. BASE COLOR

REQUIRED. Select one shroud color. Colors shown here may differ from actual product depending on monitor/printer calibration used.

3.1a. Standard Shrouds

 Red Corvette Part: DGN1900199-PRD1	 Deep Purple Part: DGN1900199-PPR5	 Midnight Blue Part: DGN1900199-PBU10	 Green Machine Part: DGN1900199-PGN11
 Back in Black Part: DGN1900199-PBK1	 Raspberry Beret Part: DGN1900199-PPK4	 Electric Blue Part: DGN1900199-PBU9	 Military Green Part: DGN1900199-PGN8
 Night Train Part: DGN1900199-PBK3	 Pretty in Pink Part: DGN1900199-PPK3	 Blue Jean Part: DGN1900199-PBU11	
 Grey Street Part: DGN1900199-PSL3	 Evolution Orange Part: DGN1900199-POR7	 Lemon Crush Part: DGN1900199-PYL4	

3.1b. Hydrographic Shrouds



Military Digital Camo \$325
Part: DGN1801101-WCM7



Faded Glory..... \$325
Part: DGN1801101-WFL2

4. BASE MOUNTED LIGHTS

Optional. Select one option.

Switch Fender Lights..... No Charge
Part: DGN1810041

Fender lights operate by pressing included single switch button. Fender lighting is not available with Single Actuator Control through Toggle Box.

5. QUANTUM BACKUP CAMERA

Optional. The backup camera may interfere with armrest mounted accessories such as media device holder.

Quantum Backup Camera \$475
Part: MEC1910365

1. Optional Mounting Bracket

Backup Camera Headrest Bracket..... \$30
Part: MEC1910373

Works with .5" x .5" vertical supports only.

Backup Camera Track Mount Bracket \$35
Part: MEC2104227

2. Optional Display Mounting Bracket

Backup Camera Q-Logic Joystick Mount Bracket \$145
Part: MEC2106123



6. DRIVE WHEELS

REQUIRED. Select a type of drive wheel. Accent colors are available for 14" tires.

6.1. Drive Wheels

3"x14" Pneumatic Tires..... No Charge
Part: WHL1713925

2"x12.5" Pneumatic Tires No Charge
Part: WHL1810666

3"x14" Flat Free Tires..... No Charge
Part: WHL1704058

2"x12.5" Flat Free Tires No Charge
Part: WHL1810665

6.2. Drive Wheel Accent Color for 14" Tires

Not available with 12.5" tires.

Back in Black..... No Charge
Part: DGN164245-PBK1 x2

Pretty in Pink No Charge
Part: DGN164245-PPK3 x2

Green Machine No Charge
Part: DGN164245-PGN11 x2

Red Corvette..... No Charge
Part: DGN164245-PRD1 x2

Evolution Orange..... No Charge
Part: DGN164245-POR7 x2

Military Green No Charge
Part: DGN164245-PGN8 x2

Night Train..... No Charge
Part: DGN164245-PBK3 x2

Midnight Blue No Charge
Part: DGN164245-PBU10 x2

Military Digital Camo..... \$50 ea.
Part: DGN164245-WCM7 x2

Grey Street..... No Charge
Part: DGN164245-PSL3 x2

Electric Blue No Charge
Part: DGN164245-PBU9 x2

Faded Glory \$50 ea.
Part: DGN164245-WFL2 x2

Deep Purple..... No Charge
Part: DGN164245-PPR5 x2

Blue Jean No Charge
Part: DGN164245-PBU11 x2

Raspberry Beret..... No Charge
Part: DGN164245-PPK4 x2

Lemon Crush..... No Charge
Part: DGN164245-PYL4 x2

7. BATTERIES

Optional. The powerchair requires 2 batteries and will be installed if ordered.

NF-22, Gel Cell, Interceptor 2 Battery \$300 ea.
Part: BAT2004376 x2. HCPCS: E2361
Uses 8 AMP charger

8. PATIENT'S WEIGHT

REQUIRED. The patient's weight is required to quote and order the power chair.

Up to 150 lb User Weight
Part: PTOINDV3977

151 lb-250 lb User Weight
Part: PTOINDV3985

251 lb-300 lb User Weight
Part: PTOINDV3984

9. ELECTRONICS

REQUIRED. Select a type of drive control from the following sections. The No Joystick option is used to omit the joystick and utilize another control, such as a specialty control at the end of this form or found on the [Alternative Controls and Electronics Order Form](#).

9.1. Expandable Controls - Quantum Intelligent Systems™

Select expandable controls or omit selection if another controller type is desired. Expandable controls allow multiple or specialized controllers.

2. Q-Logic 3 Joystick

Q-Logic 3 EX Joystick\$900
Part: CTL167724

No Joystick, Setup for Q-Logic 3
Part: XRF000010

Q-Logic 3e Joystick\$620
Part: CTL1910600

Can operate full light package. Not available with alternative or specialty controls.



Q-Logic 3 EX Joystick

9.2. Q-Logic 3 Joystick Armor

REQUIRED with Q-Logic 3 Joystick selections above. Not available with Q-Logic 3e Joystick. Reference colors shown on page 2.

Red Corvette
Part: DGN124336-PRD1

Deep Purple
Part: DGN124336-PPR5

Midnight Blue
Part: DGN124336-PBU10

Green Machine
Part: DGN124336-PGN11

Carbon Fiber
Part: DGN124336-WCB1

Back in Black
Part: DGN124336-PBK1

Raspberry Beret
Part: DGN124336-PPK4

Electric Blue
Part: DGN124336-PBU9

Military Green
Part: DGN124336-PGN8

☐ Night Train*
See Black, Soft Touch

Pretty in Pink
Part: DGN124336-PPK3

Blue Jean
Part: DGN124336-PBU11

Military Digital Camo
Part: DGN124336-WCM7

Grey Street
Part: DGN124336-PSL3

Evolution Orange
Part: DGN124336-POR7

Lemon Crush
Part: DGN124336-PYL4

Faded Glory
Part: DGN124336-WFL2

10. JOYSTICK HANDLES

Optional. The optional handle will be mounted to the joystick and the standard handle will be shipped with the unit.

Mushroom Handle (B)\$120
Part: ACCASMB1860.

Extended I-Handle (G)\$140
Part: ACCASMB1857.

Chin Cup (C)\$85
Part: ACCASMB186190

Bodypoint 3" Goal Post Handle\$180
Part: PLSKNCP1057.

T-Bar (D)\$160
Part: ACCASMB1858.

Bodypoint 4" Goal Post Handle\$180
Part: PLSKNCP1058.

Medium Ball Handle (E)\$70
Part: ACCASMB1859.

Soft Ball Handle (F)\$120
Part: ACCASMB1856.



Joystick Handles

11. JOYSTICK MOUNTING BRACKETS

Optional. Select a mount to be used with a joystick selected in previous section.

11.1a. Fixed, In-line Joystick Bracket

Fixed, In-line, Bracket for Q-Logic 3Standard
Part: FRMASMB10617

11.1b. Swing-Away (Retractable) Joystick Bracket

Allows the joystick to retract to the side of the armrest while maintaining the same joystick orientation.

Swing-away, Q-Logic 3\$265
Part: ACC1908434. HCPCS: E1028

11.1c. Basic Swing Away Joystick Mount


Basic Swing Away, Q-Logic, Left Side\$265
Part: ACC2201041. HCPCS: E1028

Basic Swing Away, Q-Logic, Right Side\$265
Part: ACC2201042. HCPCS: E1028


11.1d. Offset Swing Away Joystick Mount

Offset Swing Away, Q-Logic, Left Side\$265
Part: ACC142742. HCPCS: E1028

Offset Swing Away, Q-Logic, Right Side\$265
Part: ACC142743. HCPCS: E1028



Swing Away Joystick



Retract4

Joystick mounting brackets section continued on following page.

11. JOYSTICK MOUNTING BRACKETS CONTINUED**11.1d. Specialty Joystick Mounts**

Select a bracket and side to mount it.

Retract4 Mount\$450
Part: XRF000003. HCPCS: E1028

Full 180° of movement that translates into multi-position use, straight, right offset, left offset, or fold away. Joystick orientation varies as bracket is adjusted.

Multi-Axis Mount\$100

Part: XRF000004

Allows height adjustment of joystick.

Stealth Swing-Away/Flip Down Bracket\$300
Part: XRF000005. HCPCS: E1028

Allows the joystick to pivot to side of armrest, or flip-down if bracket is adjusted.

**11.2. Joystick Mount Location**

Left Side Mount

Part: PTOINDV1225

Right Side Mount

Part: PTOINDV1226

12. TRU BALANCE® 4 SEATING

Maximum seat size is 22"W x 22"D.

12.1. TRU Balance 4 Configuration**REQUIRED.**

TRU-Balance 4 Power Tilt Only\$8,495

Part: PTOINDV4014. HCPCS: E1002

50° Center of Gravity (CG) power tilt seating system.

TRU-Balance 4 Power Tilt and Recline\$14,995

Part: PTOINDV3954. HCPCS: E1007

50° of posterior tilt and 165° of recline.

TRU-Balance 4 Power iLevel Lift, Tilt and Recline\$17,990

Part: PTOINDV3954. MSRP: \$14,995 (T&R) + \$2,995 (iLevel)

12" Power Lift, 50° of posterior tilt and 165° of recline. Power Articulating Foot Platform is mandatory with this selection.

12.3. Anterior Tilt

Chest support is required with anterior tilt 20 and 30 degree functions. See chest support section. Omitting a chest support selection requires the provider to facilitate aftermarket chest supports that best fits consumer needs.

1. Anterior Tilt Angle

Only available with Power Lift.

10° Anterior Tilt Angle No Charge

Part: MEC2009356

20° Anterior Tilt Angle\$2,695

Part: MEC2009399. HCPCS: K0108

Requires chest support and knee block.

30° Anterior Tilt Angle\$3,195

Part: MEC2009400. HCPCS: K0108

Requires chest support and knee block.

0° Anterior Tilt No Charge

Part: PTOINDV3980

12.4. Knee Block Assembly**REQUIRED with 20° or 30° anterior tilt. Optional with 10° and 0° anterior tilt.**

Anterior Knee Block Assembly (ST-KBATB4ATC-1-Q)\$700

Part: ST-KBATB4ATC-1-Q. HCPCS: K0108

Anterior Knee Block Assembly, Curved (ST-KBATB4ATC-U-Q) ..\$700

Part: ST-KBATB4ATC-U-Q. HCPCS: K0108

12.5. Occupied Transit Kit**Optional.** The power base is standard with Unoccupied Transit Loops. 60" Sure-Lok Lap Belt is required with Occupied Transit Option.

WC-19 Occupied Transit Option\$400

Part: KIT152047

To determine if your chair is WC-19 compliant (occupied), please refer to this website: <http://wc-transportation-safety.umtri.umich.edu/crash-tested-product-lists/wheelchairs>**12.6. Seat Width****REQUIRED with seat selection. 12"-15" widths not available with Power Lift. Maximum seat width is 22". Part number is SETWDTH plus number following # below.**

12" W (#1002)	13" W (#1003)	14" W (#1004)	15" W (#1005)	16" W (#1006)	17" W (#1007)	18" W (#1008)	19" W (#1009)	20" W (#1010)	21" W (#1011)	22" W (#1012)
------------------	------------------	------------------	------------------	------------------	------------------	------------------	------------------	------------------	------------------	------------------

12.7. Seat Width Package**REQUIRED with seat selection.**

12"-16" Width Package

Part: XRF000006

Available with Power Tilt Only.

16"-20" Width Package

Part: XRF000008

20"-24" Width Package

Part: XRF000158

Maximum width is 22".





12.8. Seat Depth**REQUIRED with seat selection. 12"-13" depth not available with Power Lift. Maximum seat depth is 22". Part number is SETDPTH plus number following # below.**

12" D (#1002)	13" D (#1003)	14" D (#1004)	15" D (#1005)	16" D (#1006)	17" D (#1007)	18" D (#1008)	19" D (#1009)	20" D (#1010)	21" D (#1011)	22" D (#1012)
------------------	------------------	------------------	------------------	------------------	------------------	------------------	------------------	------------------	------------------	------------------




TRU Balance 4 Seating section continued on following pages.

12. TRU BALANCE 4 SEATING CONTINUED**12.9. Seat Depth Package****REQUIRED** with seat selection.**1a. Power Tilt Depth Package**12"-16" Depth Package
Part: XRF00001116"-20" Depth Package
Part: XRF00001320"-24" Depth Package
Part: XRF000159
Maximum depth is 22".**1b. Power Lift Depth Package**16"-22" Depth Package
Part: XRF000309**12.10. Seat to Floor Height****REQUIRED.**Low Seat to Floor
Part: PTOINDV3976Middle Seat to Floor
Part: PTOINDV3559Highest Seat to Floor
Part: PTOINDV3560

Low	Med	High
Q6 Edge 3 Stretto		
18.0"	18.5"	19.0"
TRU Balance 4 Seat to Floor Heights		

12.11. Stealth Seat Cushion - Quick Selection**Optional.** This section allows a quick selection of a single cushion that matches the seat pan width and depth to alleviate listing all part numbers and codes on this power base form. See specific order form for non-matching sizes, alternate covers and liners, or additional covers and liners.**Stealth TRU-Comfort 2 SPP Cushion**Skin protection and positioning with Coolcore, 4-way stretch, breathable cover.
Part series: ST-TC2SPPXXYY. HCPCS: E2607/E2608TRU-Comfort 2 SPP, Matching Seat Size
Reference part: XRF000116.**Stealth Solution SPP Cushion**Skin protection and positioning with Coolcore, 4-way stretch, breathable cover.
Includes urethane liner.
Part series: ST-SOLSPPXXYY. HCPCS: E2607/E2608Solution SPP, Matching Seat Size
Reference part: XRF000113.**Stealth Spectrum Foam SPP Cushion**Skin protection and positioning with Coolcore, 4-way stretch, breathable cover.
Includes urethane liner.
Part series: ST-SPCFSPXXYY. HCPCS: E2607/E2608Spectrum Foam SPP, Matching Seat Size
Reference part: XRF000112.**Stealth Premiere P Cushion**Positioning with Coolcore, 4-way stretch, breathable cover.
Part series: ST-PREMPXXYY. HCPCS: E2605/E2606Premiere P, Matching Seat Size
Reference part: XRF000122.**Stealth Zen SP Cushion**Skin protection with Coolcore, 4-way stretch, breathable cover.
Part series: ST-ZENSPXXYY. HCPCS: E2603/E2604Zen SP, Matching Seat Size
Reference part: XRF000115.**Stealth Spectrum Gel SPP Cushion**Skin protection and positioning with Coolcore, 4-way stretch, breathable cover.
Includes urethane liner.

Part series: ST-SPCGSPPXXYY. HCPCS: E2607/E2608

Spectrum Gel SPP, Matching Seat Size
Reference part: XRF000117. MSRP: \$507/\$585**Stealth Essence SPP Cushion**Skin protection and positioning with Coolcore, 4-way stretch, breathable cover.
Part series: ST-ESSSPPXXYY. HCPCS: E2607/E2608Essence SPP, No Wedge, Matching Seat Size
Reference part: XRF000118.Essence SPP, Small Wedge, Matching Seat Size
Reference part: XRF000119.Essence SPP, Medium Wedge, Matching Seat Size
Reference part: XRF000120.Essence SPP, Large Wedge, Matching Seat Size
Reference part: XRF000121.Essence Plus SPP, Large Wedge, Matching Seat Size
Reference part: XRF000121.**Stealth Glacial SP Cushion**Skin protection with Coolcore, 4-way stretch, breathable cover.
Part series: ST-GLSPXXYY. HCPCS: E2603/E2604Glacial SP, Matching Seat Size
Reference part: XRF000114.**Stealth Simplicity G Cushion**General use cushion with Stealth Tek 2, 2-way stretch, breathable cover.
Part series: ST-SIMGXXYY. HCPCS: E2601/E2602Simplicity G, Matching Seat Size
Reference part: XRF000111.

TRU Balance 4 Seating section continued on following pages.

12. TRU BALANCE 4 SEATING CONTINUED**12.12. Back Type Selection****REQUIRED.** Match a choice in this section with the options in the following back cushion configuration section.Redesigned, Static Type Back \$735
Part: XRF000014

Select a back cushion configuration in the subsequent section.

Back Cane Option No Charge
Part: XRF000015Available only with TRU-Balance 3 Power Tilt (1 actuator) selections.
Select cane type, back cane height, and back cushion configuration below.Redesigned, Flip Forward Back \$735
Part: XRF000016

Not available with power recline systems. Only available with 2-post, flip back or pediatric armrests. Not available with 12"-14" width. Select a back cushion type in the subsequent section.

12.13. Back Cane Selection**REQUIRED when Back Cane Option is selected above.** Select type and height of back canes.**1. Back Cane Type**Angle Adjustable Back Cane No Charge
Part: MEC140266Flip Forward Back Cane No Charge
Part: XRF0002732" Endomorph (1" per side) Back Cane No Charge
Part: XRF0002694" Endomorph (2" per side) Back Cane No Charge
Part: XRF000271**2. Back Cane Height**16"H Cane
Part: FRMASMB1583420"H Cane
Part: FRMASMB1583824"H Cane
Part: FRMASMB1584217"H Cane
Part: FRMASMB1583521"H Cane
Part: FRMASMB1583925"H Cane
Part: FRMASMB15843**12.14. Back Cushion Configuration****REQUIRED.** Select a configuration type and applicable dimensions. See the Sport, Q-Back, and TRU-Comfort back availability chart for stock order back size ranges. Other size backs may be available with special order and extended lead time.**1a. Sport Back/Omit/After-Market Interface**Sport Back No Charge
Part: XRF000017

Complete following width package and size.

Omit Back Pan/Cushion No Charge
Part: PTOINDV1287

Select a following width package.

After-market Back Interface No Charge
Part: ST-SSAPT4-Q. Select a set width in following section, range of 12"-21". Overall seat depth may need to be adjusted to accommodate certain backs. May interfere with gel/waterfall arm pads and arm through. Does not include brackets to mount after-market backs.**1b. TRU-Comfort Back**TRU-Comfort Back (Static/Flip Forward), Up to 15" Width \$835
Part: SET123622. HCPCS: E2620

Complete following width package and size.

TRU-Comfort Back (Static/Flip Forward), 16" to 20" Width ... \$1,035
Part: TC4V10W21W. HCPCS: E2620

Complete following width package and size.

TRU-Comfort Back (Static/Flip Forward), 22" Width \$1,235
Part: TC4V22W32W. HCPCS: E2621

Complete following width package and size.

TRU-Comfort Back, Back Cane Option \$835
Part: TC4V10W21WBC. HCPCS: E2620

Complete following width package and size.

1c. Q-Back 2.0, AL Series, Active Contour

Select a back type and the specific overall back height in subsequent section.

Active Contour, 14"Wx19-20"H Tilt/22"-23" H Recline \$9,85
Part: QB2AL14AS. HCPCS: E2613Active Contour, 14"Wx21-22"H Tilt/24"-25" H Recline \$985
Part: QB2AL14AT. HCPCS: E2613Active Contour, 16"Wx19-20"H Tilt/22"-23" H Recline \$1,145
Part: QB2AL16AS. HCPCS: E2613Active Contour, 16"Wx21-22"H Tilt/24"-25" H Recline \$1,145
Part: QB2AL16AT. HCPCS: E2613Active Contour, 16"Wx23-24"H Tilt/26"-27" H Recline \$1,145
Part: QB2AL16AX. HCPCS: E2613Active Contour, 16"Wx25-26"H Tilt/28"-29" H Recline \$1,145
Part: QB2AL16AXX. HCPCS: E2613Active Contour, 18"Wx19-20"H Tilt/22"-23" H Recline \$1,145
Part: QB2AL18AS. HCPCS: E2613Active Contour, 18"Wx21-22"H Tilt/24"-25" H Recline \$1,145
Part: QB2AL18AT. HCPCS: E2613Active Contour, 18"Wx23-24"H Tilt/26"-27" H Recline \$1,145
Part: QB2AL18AX. HCPCS: E2613Active Contour, 18"Wx25-26"H Tilt/28"-29" H Recline \$1,145
Part: QB2AL18AXX. HCPCS: E2613Active Contour, 20"Wx19-20"H Tilt/22"-23" H Recline \$1,145
Part: QB2AL20AS. HCPCS: E2613Active Contour, 20"Wx21-22"H Tilt/24"-25" H Recline \$1,145
Part: QB2AL20AT. HCPCS: E2613Active Contour, 20"Wx23-24"H Tilt/26"-27" H Recline \$1,145
Part: QB2AL20AX. HCPCS: E2613Active Contour, 20"Wx25-26"H Tilt/28"-29" H Recline \$1,145
Part: QB2AL20AXX. HCPCS: E2613

TRU Balance 4 Seating section continued on following pages.

12. TRU BALANCE 4 SEATING CONTINUED**1d. Q-Back 2.0, AL Series, Deep Contour**

Select a back type and the specific overall back height in subsequent section.

Deep Contour, 14"Wx19-20"H Tilt/22"-23" H Recline	\$1,145
Part: QB2AL14DS. HCPCS: E2620	
Deep Contour, 14"Wx21-22"H Tilt/24"-25" H Recline	\$1,145
Part: QB2AL14DT. HCPCS: E2620	
Deep Contour, 16"Wx19-20"H Tilt/22"-23" H Recline	\$1,325
Part: QB2AL16DS. HCPCS: E2620	
Deep Contour, 16"Wx21-22"H Tilt/24"-25" H Recline	\$1,325
Part: QB2AL16DT. HCPCS: E2620	
Deep Contour, 16"Wx23-24"H Tilt/26"-27" H Recline	\$1,325
Part: QB2AL16DX. HCPCS: E2620	
Deep Contour, 16"Wx25-26"H Tilt/28"-29" H Recline	\$1,325
Part: QB2AL16DXX. HCPCS: E2620	
Deep Contour, 18"Wx19-20"H Tilt/22"-23" H Recline	\$1,325
Part: QB2AL18DS. HCPCS: E2620	

Deep Contour, 18"Wx21-22"H Tilt/24"-25" H Recline	\$1,325
Part: QB2AL18DT. HCPCS: E2620	
Deep Contour, 18"Wx23-24"H Tilt/26"-27" H Recline	\$1,325
Part: QB2AL18DX. HCPCS: E2620	
Deep Contour, 18"Wx25-26"H Tilt/28"-29" H Recline	\$1,325
Part: QB2AL18DXX. HCPCS: E2620	
Deep Contour, 20"Wx19-20"H Tilt/22"-23" H Recline	\$1,325
Part: QB2AL20DS. HCPCS: E2620	
Deep Contour, 20"Wx21-22"H Tilt/24"-25" H Recline	\$1,325
Part: QB2AL20DT. HCPCS: E2620	
Deep Contour, 20"Wx23-24"H Tilt/26"-27" H Recline	\$1,325
Part: QB2AL20DX. HCPCS: E2620	
Deep Contour, 20"Wx25-26"H Tilt/28"-29" H Recline	\$1,325
Part: QB2AL20DXX. HCPCS: E2620	

1e. Q-Back 2.0, AL Series, Extra Deep Contour

Select a back type and the specific overall back height in subsequent section.

Extra Deep Contour, 14"Wx19-20"H Tilt/22"-23" H Recline	\$1,225
Part: QB2AL14ES. HCPCS: E2620	
Extra Deep Contour, 14"Wx21-22"H Tilt/24"-25" H Recline	\$1,225
Part: QB2AL14ET. HCPCS: E2620	
Extra Deep Contour, 16"Wx19-20"H Tilt/22"-23" H Recline	\$1,410
Part: QB2AL16ES. HCPCS: E2620	
Extra Deep Contour, 16"Wx21-22"H Tilt/24"-25" H Recline	\$1,410
Part: QB2AL16ET. HCPCS: E2620	
Extra Deep Contour, 16"Wx23-24"H Tilt/26"-27" H Recline	\$1,410
Part: QB2AL16EX. HCPCS: E2620	
Extra Deep Contour, 16"Wx25-26"H Tilt/28"-29" H Recline	\$1,410
Part: QB2AL16EXX. HCPCS: E2620	
Extra Deep Contour, 18"Wx19-20"H Tilt/22"-23" H Recline	\$1,410
Part: QB2AL18ES. HCPCS: E2620	

Extra Deep Contour, 18"Wx21-22"H Tilt/24"-25" H Recline	\$1,410
Part: QB2AL18ET. HCPCS: E2620	
Extra Deep Contour, 18"Wx23-24"H Tilt/26"-27" H Recline	\$1,410
Part: QB2AL18EX. HCPCS: E2620	
Extra Deep Contour, 18"Wx25-26"H Tilt/28"-29" H Recline	\$1,410
Part: QB2AL18EXX. HCPCS: E2620	
Extra Deep Contour, 20"Wx19-20"H Tilt/22"-23" H Recline	\$1,410
Part: QB2AL20ES. HCPCS: E2620	
Extra Deep Contour, 20"Wx21-22"H Tilt/24"-25" H Recline	\$1,410
Part: QB2AL20ET. HCPCS: E2620	
Extra Deep Contour, 20"Wx23-24"H Tilt/26"-27" H Recline	\$1,410
Part: QB2AL20EX. HCPCS: E2620	
Extra Deep Contour, 20"Wx25-26"H Tilt/28"-29" H Recline	\$1,410
Part: QB2AL20EXX. HCPCS: E2620	

12.15. TRU Balance 4 Back Shroud**REQUIRED.** The back shroud is a feature of the TRU Balance 4 seating system and is available in matching shroud colors. The back shroud is not available when Omit Back Pan/Cushion is selected. Reference colors shown on page 2.

Red Corvette Part: XRF000019	Grey Street Part: XRF000023	Pretty in Pink Part: XRF000027	Electric Blue Part: XRF000030	Green Machine Part: XRF000033
Back in Black Part: XRF000021	Deep Purple Part: XRF000025	Evolution Orange Part: XRF000028	Blue Jean Part: XRF000031	Military Green Part: XRF000034
Night Train Part: XRF000022	Raspberry Beret Part: XRF000026	Midnight Blue Part: XRF000029	Lemon Crush Part: XRF000032	

12.16. Lateral Positioning Wedge**Optional.** Available with Sport and TRU-Comfort backs. Not available on power recline back 19"-20"H with limited adjustability on 21"H.

Lateral Positioning Wedge	\$15 ea.
Part: SET123465 x2	

TRU Balance 4 Seating section continued on following pages.

12. TRU BALANCE 4 SEATING CONTINUED**12.17. Back Width Package****REQUIRED** except when Omit Back Pan/Cushion is selected. Available with +/-4" of seat width.Match Seat Width Selections
Part: PTOINDV346814"W-16"W Package
Part: XRF00003716"W-20"W Package
Part: XRF00003818"W-22"W Package
Part: XRF000039**12.18. Back Width****REQUIRED** except when Omit Back Pan/Cushion is selected. Available with +/-4" of seat width.12" Back Width
Part: PTOINDV356216" Back Width
Part: PTOINDV356418" Back Width
Part: PTOINDV356620" Back Width
Part: PTOINDV356822" Back Width
Part: PTOINDV357014" Back Width
Part: PTOINDV356317" Back Width
Part: PTOINDV356519" Back Width
Part: PTOINDV356721" Back Width
Part: PTOINDV3569

	Tilt Only Width											Recline Width										
	12"	14"	15"	16"	17"	18"	19"	20"	21"	22"		12"	14"	15"	16"	17"	18"	19"	20"	21"	22"	
14"	S	T	S	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	14"
15"	S	T	S	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	15"
16"	S	T	S	T	S	T	S	T	S	T	S	T	S	T	S	T	S	T	S	T	S	16"
17"	S	T	S	T	S	T	S	T	S	T	S	T	S	T	S	T	S	T	S	T	S	17"
18"	-	S	T	S	T	S	T	S	T	S	T	S	T	S	T	S	T	S	T	S	T	18"
19"	-	S	T	S	T	S	T	S	T	S	T	S	T	S	T	S	T	S	T	S	T	19"
20"	-	S	T	S	T	S	T	S	T	S	T	S	T	S	T	S	T	S	T	S	T	20"
21"	-	S	T	S	T	S	T	S	T	S	T	S	T	S	T	S	T	S	T	S	T	21"
22"	-	S	T	S	T	S	T	S	T	S	T	S	T	S	T	S	T	S	T	S	T	22"
23"	-	S	T	S	T	S	T	S	T	S	T	S	T	S	T	S	T	S	T	S	T	23"
24"	-	S	T	S	T	S	T	S	T	S	T	S	T	S	T	S	T	S	T	S	T	24"
25"	-	S	T	S	T	S	T	S	T	S	T	S	T	S	T	S	T	S	T	S	T	25"
26"	-	S	T	S	T	S	T	S	T	S	T	S	T	S	T	S	T	S	T	S	T	26"
27"	-	S	T	S	T	S	T	S	T	S	T	S	T	S	T	S	T	S	T	S	T	27"
28"	-	S	T	S	T	S	T	S	T	S	T	S	T	S	T	S	T	S	T	S	T	28"
29"	-	S	T	S	T	S	T	S	T	S	T	S	T	S	T	S	T	S	T	S	T	29"
30"	-	-	-	-	-	-	-	-	-	-	-	S	T	S	T	S	T	S	T	S	T	30"
31"	-	-	-	-	-	-	-	-	-	-	-	S	T	S	T	S	T	S	T	S	T	31"
32"	-	-	-	-	-	-	-	-	-	-	-	S	T	S	T	S	T	S	T	S	T	32"

S Sport Back Q Q-Back T TRU-Comfort Back

Sport, Q-Back, and TRU-Comfort Back Availability Chart

12.19. Back Cushion Overall Back Height**Power Recline available height range:** 19"-32"H. **Odd heights include 6" gap and even heights include 7" gap** between seat pan and bottom of back cushion. See the availability chart for stock order back size ranges.19" Height
Part: SETHGHT103222" Height
Part: SETHGHT103525" Height
Part: SETHGHT103828" Height
Part: SETHGHT104131" Height
Part: SETHGHT110220" Height
Part: SETHGHT103323" Height
Part: SETHGHT103626" Height
Part: SETHGHT103929" Height
Part: SETHGHT104232" Height
Part: SETHGHT110321" Height
Part: SETHGHT103424" Height
Part: SETHGHT103727" Height
Part: SETHGHT104030" Height
Part: SETHGHT1043**12.20. Optional Back Height Extension****Optional.** Adds 2" of seat back pan height. 25"H selection includes a back height extension.Back Height Extension \$48
Part: XRF000041**13. ARMRESTS****REQUIRED.** Select the desired type of armrest below. Please read notes to ensure seating compatibility. Be sure to also select an armrest size and armpad type.**13.1a. Flip-Back, 8"-16" H Armrests**

Choose a left and right armrest and set height in the following section. Not available with waterfall armpads. Not available with Back Cane Option.

L R Flip-Back, 8"-16"H Armrests \$145 ea.
Left part: ARM1806662. Right part: ARM1806663**Armrest Support****Optional.** Available up to 2" difference between seat and back width via endomorph or mesomorph configuration. Not available less than 16" width seat or back.L R Flip-Back Armrest Downpost \$115 ea.
Left part: XRF000265. Right part: XRF000266**13.1b. Two-Post, Height Adjustable, Removable Armrests**

Select the armrest and set the left and right heights in the following section. These armrests may interfere with offset backs/ back canes.

1a. Two-Post Armrest for Power SeatingL R Full Length, Two-Post Armrest \$145 ea.
Left part: XRF000047. Right part: XRF000048

Armrests section continued on following page.

13. ARMRESTS CONTINUED**13.1c. Pediatric, Removable, Height Adjustable Armrests**Height adjustment from 7.75" to 9.75" in .5" increments. **Not available with full length armpads.**

☐ L ☐ R	Pediatric, Removable Armrests \$145 ea.
---	--

13.2. Armrest Height**REQUIRED.** Set left and right armrest height. Note the available height range in the option description.**1a. Flip-Back and Two-Post Left Armrest Height**

8" H PTOINDV3535	10" H PTOINDV3537	12" H PTOINDV3539	14" H PTOINDV3541	16" H PTOINDV3543
9" H PTOINDV3536	11" H PTOINDV3538	13" H PTOINDV3540	15" H PTOINDV3542	

1b. Pediatric, Left Armrest Height

7.75" H PTOINDV3722	8.25" H PTOINDV3724	8.75" H PTOINDV3726	9.25" H PTOINDV3728	9.75" H PTOINDV3730
------------------------	------------------------	------------------------	------------------------	------------------------

2a. Flip-Back and Two-Post Right Armrest Height

8" H PTOINDV3450	10" H PTOINDV3452	12" H PTOINDV3454	14" H PTOINDV3456	16" H PTOINDV3458
9" H PTOINDV3451	11" H PTOINDV3453	13" H PTOINDV3455	15" H PTOINDV3457	

2b. Pediatric, Right Armrest Height

7.75" H PTOINDV3721	8.25" H PTOINDV3723	8.75" H PTOINDV3725	9.25" H PTOINDV3727	9.75" H PTOINDV3729
------------------------	------------------------	------------------------	------------------------	------------------------

13.3. Armrest Size**REQUIRED.** Select a left (L) and right (R) armrest size. This also affects armpad selection.

☐ L ☐ R	10" Length L part: XRF000057. R part: XRF000123	☐ L ☐ R	14" Length L part: XRF000058. R part: XRF000124	☐ L ☐ R	17" Length L part: XRF000059. R part: XRF000125 Only available with Flip-Back, 8"-16" H Armrests and straight arm pads.
---	--	---	--	---	--

14. ARMPADS**REQUIRED.****14.1a. Omit Armpad**Select one below if aftermarket armpads or handpads are required. See [Positioning Components Order Form](#).

No Armpad Selection..... No Charge Part: XRF000060 Not available with Flip-back, 8"-16" H armrests.	Interface for Aftermarket Pad No Charge Part: XRF000308 Only available with Flip-back, 8"-16" armrests. Not available with chest bar.
---	---

14.1b. Foam Armpads

☐ L ☐ R	Straight, Armpad Standard Left part: XRF000061. Right part: XRF000126	☐ L ☐ R	Waterfall, Armpad No Charge Left part: XRF000062. Right part: XRF000127 Not available with Flip-back, 8"-16" armrests.
---	--	---	--

14.1c. Gel Armpads

Sold only in pairs. Armpads are 3" wide. Manufacturer part number for full length is 314G and desk length is 310G.

a. Falcon Gel Armpads for Flip-Back Armrests

Falcon Gel, Full Length Left and Right \$300 Part: ARM1806807	Falcon Gel, Desk Length Left, Full Length Right \$250 Part: ARM1806810
Falcon Gel, Full Length Left, Desk Length Right \$250 Part: ARM1806813	Falcon Gel, Desk Length Left and Right \$250 Part: ARM1806806

b. Falcon Gel Armpads for Tube Style Armrests

Falcon Gel, Full Length Left and Right \$300 Part: ARM136438	Falcon Gel, Desk Length Left, Full Length Right \$250 Part: ARM136439
Falcon Gel, Full Length Left, Desk Length Right \$250 Part: ARM136442	Falcon Gel, Desk Length Left and Right \$200 Part: ARM136437

15. ARMREST BRACKETS

Standard Fixed Receiver Standard Part: XRF000063	Stealth Height Adjustable Receiver \$194 Part: XRF000064. HCPCS: K0108
---	---

16. HEADRESTS**Optional.** Select a headrest pad and bracket.**16.1a. Headrest Pads, Redesign Back Assembly**

Only available with 14" to 22" back widths of Static Type or Flip Forward back, without occupied transit.

1. Headrest PadsStealth 8" Comfort Plus, Cool Core (CP08CCFPM) \$279
Part: ST-CP08CCFPM-Q. HCPCS: E0955Stealth 10" Comfort Plus, Cool Core (CP10CCFPM) \$279
Part: ST-CP10CCFPM-Q. HCPCS: E0955Stealth 14" Comfort Plus, Cool Core (CP14CCFPM) \$279
Part: ST-CP14CCFPM-Q. HCPCS: E0955Stealth 6", Comfort Plus, Cool Core (CP06CCFPM) \$279
Part: ST-CP06CCFPM-Q. HCPCS: E0955**2. Headrest Brackets**Comfort Plus Headrest Mounting Bracket, UniLink \$485
Part: ST-UHB2D33200-Q. HCPCS: E1028Comfort Plus Headrest Mounting Bracket, HMO475 \$271
Part: ST-HMO475-17M-Q.**17. LAP BELTS****REQUIRED.** Select a lap belt below.**17.1a. Standard Lap Belts**
☐ 38" Lap Belt Part: ACCBELT1006
 ☐ 50" Lap Belt Part: ACCBELT1000
 ☐ 60" Lap Belt Part: ACCBELT1002
 ☐ 70" Lap Belt Part: ACCBELT1003
 ☐ 80" Lap Belt Part: ACCBELT1004
 ☐ 90" Lap Belt Part: ACCBELT1001
17.1b. Optional Lap BeltsStealth Positioning, 60" Single Pull, Push Button, Medium \$134
(7.5") Padded, 1.5" Belt (IFPB5213)
Part: ST-IFPB5213-QStealth Positioning, 60" Single Pull, Push Button, Large \$134
Padded (9.5"), 2" Belt (IFPB5314)
Part: ST-IFPB5314-QStealth Positioning, 60" Single Pull, Push Button, \$134
Extra Large, Padded (11.5"), 2" Belt (IFPB5315)
Part: ST-IFPB5315-Q60" Sure-Lok Lap Belt (Non-retractable) \$200
Part: ACCBELT1016
Required with ISO 7176-19 Compliant Transit Option.**18. LEGRESTS**

All efforts will be made to accommodate legrest length selection, but minor changes may need to be made to meet shipping guidelines. Some adjustments may be needed upon receipt of the unit. Reference the section's compatibility matrix for legrest interference.

18.1a. Basic Center Mount Foot Platform**1. Basic Center Mount Foot Platform Assembly**Basic Center Mount \$200
Part: XRF000323**2. Basic Center Mount Extension and Angle Adjustable Footplates**Short Extension, 5"x9" Footplates
Part: XRF000312. HCPCS: K0040Medium Extension, 5"x9" Footplates
Part: XRF000313. HCPCS: K0040Long Extension, 5"x9" Footplates
Part: XRF000314. HCPCS: K0040**3. Basic Center Mount Calf Pads**Small Calf Pads, Basic Center Mount \$210
Part: RIG2104369Large Calf Pads, Basic Center Mount \$210
Part: RIG2202997
Only available with Long Extension.**4. Basic Center Mount Additional Accessory**Basic Center Mount Footplate Connector \$55
Part: RIG2104240**18.1b. Center Mount Foot Platform****1. Foot Platform Assembly**Dual Actuator, Power Articulating Foot Platform \$4,295
Part: RIG2009246. HCPCS: E1012+K0040
Positional feedback with independent control of articulation.Quick Adjustable Center Mount Foot Platform \$265
Part: XRF000311. HCPCS: K0040
Not available with power lift.**2. Lower Extension and Angle Adjustable Footplates**Medium Extension, 5"x6" Footplates
Part: RIG150206Medium Extension, 5"x8" Footplates
Part: RIG150203Medium Extension, 5"x11" Footplates
Part: RIG150212

Legrest Compatibility Matrix Articulating Foot Platform TRU Balance 4	Seat to Floor Option		
	Low	Middle	High
Medium Length Extension	11"-14.5"	11"-15"	11"-15.5"

Legrests section continued on following page.

18. LEGRESTS CONTINUED**3a. Center Mount Calf Pads**

Calf pads are standard with a power articulating foot platform., but chargeable on Quick adjustable center mount.

Small Calf Pads, AFP Only, 6"Wx6"H (pair) \$210
Part: RIGASMB7120018Large Calf Pads, AFP Only, 8" taper to 6"Wx8"H (pair) \$210
Part: RIGASMB7120019**4. Center Mount Modular Footplate Extension Components**

Optional. Assembly required. May require seating configuration adjustments to avoid front caster interference.

Width Extension for 8"D Footplates \$55
Part: RIG1525174" Outside Edge (pair) \$55
Part: RIG152510Width Extension for 11"D Footplates \$55
Part: RIG152518

Only available with 6"D footplates and 8"D footplate extension.

2" Bell-Shaped Length Extension \$55
Part: RIG1525196" Outside Edge (pair) \$55
Part: RIG152516

Not compatible with width extensions. Can be used with 6" outside edge with 11" footplates.

Available with 8"D and 11"D footplates and 11"D footplate extension.

5. Center Mount Additional AccessoriesSmall Reversed Dartex Calf Pad Cover (pair) \$150
Part: RIG1337175"x11" Tapered Gel Footplate Covers..... \$420
Part: ST-FPG5011T-Q. HCPCS: K0108Large Reversed Dartex Calf Pad Cover (pair) \$150
Part: RIG1337196.5"x8" Tapered Gel Footplate Covers for Width Extension..... \$420
Part: ST-FPG6508T-Q. HCPCS: K0108Footplate Connector \$55
Part: RIGASMB71200235"x11" Tapered Gel Footplate Cover for Width Extension \$420
Part: ST-FPG7011T-Q. HCPCS: K01085"x8" Tapered Gel Footplate Covers..... \$420
Part: ST-FPG5008T-Q. HCPCS: K0108**18.2. Legrest Accessories**Wheel Bumper..... \$15 ea.
Part: XRF000100**18.3. Legrest Length****REQUIRED with matching left and right legrest configuration. Part number is PTOINDV plus number following # below.** All efforts will be made to accommodate legrest length selection, but minor changes may need to be made to meet shipping guidelines. Some adjustments may be needed upon receipt of the unit. Reference the section's compatibility matrix for legrest interference. Non-matching or individual legrest configurations lengths are not set by the manufacturer.

11" L (#3478)	12" L (#3480)	13" L (#3482)	14" L (#3484)	15" L (#3486)	16" L (#3488)	17" L (#3490)
11.5" L (#3479)	12.5" L (#3481)	13.5" L (#3483)	14.5" L (#3485)	15.5" L (#3487)	16.5" L (#3489)	

19. POWER POSITIONING ELECTRONICS**REQUIRED.****19.1a. Q-Logic 3 Controls through Joystick**

Single Actuator Through Q-Logic 3 Controller and Interface (AM3) No Charge

2 Actuators Through Q-Logic 3 Controller and Interface (AM3).....\$1,250
Not funded by ADP3 Actuators Through Q-Logic 3 Controller and Interface (AM3)..... \$1,500
Not funded by ADPUp to 5 Actuators Through Q-Logic 3 Controller and Interface (AAM).....\$2,175
ADP: SECX15005**19.1b. iAccess Switch Control**iAccess Flat Button.....No Charge
Part: CTL144009
See next section for switch mounting.iAccess Toggle Button.....No Charge
Part: CTL144010
See next section for switch mounting.**19.2. iAccess Mount****REQUIRED with iAccess options.****1. iAccess Mount Selection**Joystick Mount No Charge
Part: MEC1907500Attendant Mount \$35
Part: MEC1906564Under Arm Mount, Flip-Back Armrest..... No Charge
Part: MEC1906531Under Arm Mount, Flip-Back Armrest, plus iAccess Guard..... \$35
Part: MEC2001392**2. Under Arm Mount Side**Left Side Under Arm Mount
Part: PTOINDV3678Right Side Under Arm Mount
Part: PTOINDV3679

20. ATTENDANT CONTROL**Optional.** Attendant control will be at no charge when ordered with Tilt and Recline

Attendant Control \$950

Part: XRF000078. HCPCS: E2331

21. ACCESSORIES**Optional.** Multiple accessories can be mounted along the various tracks but some interference may occur.

Accessory Bar \$90

Part: XRF000128

Horizontal rail attached to rear of seat. Standard with the Trilogy vent tray and oxygen holder.

Backpack Holder (Pair) (A) \$42

Part: ACC147774

Cane/Crutch Holder \$100

Part: ACC152515. HCPCS: E2207

Cell Phone/ Media Holder, X-Grip (E) \$95

Part: ACC155835

Not available with Heavy Duty armrests.

Clothing Guard (Pair) \$275

Part: ACC125007

Compact Vent Tray \$1,395

Part: ACC133889. HCPCS: E1029

Cup Holder \$20

Part: ACC125003

For use with track-style armrests.

Cup Holder (Tube-style Armrest) \$20

Part: ACCASMB2634

For use with tube-style armrests.

Flag Kit \$125

Part: ACC1700443

Flip-Back Armrest, Side Track Accessory Bracket .. \$25

Part: MEC1809178

Glove Box (C)

\$85

Part: ACC125011

Large, Black, Weather Cover \$240

Part: ACCCOVR1021

Oxygen Holder \$110

Part: ACC125006. HCPCS: E2208

Includes accessory bar, a credit will be applied PTOIN-DV3769.

Personal Item Hook \$20

Part: ACC125008

Personal Item Pouch \$175

Part: ACC2101017

Privacy Bag \$80

Part: ACC2005723

Push Handles (pair) \$140

Part: ACC2005027

Not available with oxygen holder or vent tray.

Q-Logic EX/NE Programmer \$625

Part: ELEASMB6976

Quantum Backpack \$200

Part: ACC1905853

Rear View Mirror \$35

Part: ACC155833

Reflector Kit \$35

Part: ACC124670

Stealth Large Tablet Holder \$565

Part: ACC141818

Transfer Handle (D) \$152.50 ea.

Part: MEC2205680 x2. HCPCS: E0705

Trilogy Vent Tray \$420

Part: ACC136791. HCPCS: E1029

Includes accessory bar, a credit will be applied PTOIN-DV3769. Contact custom for models other than 100 and 200.

VOCSN Vent Tray for Accessory Bar Mount \$1,000

Part:


. HCPCS: E1029

Includes accessory bar, a credit will be applied PTOIN-DV3769.

Depth Adjustable Vent Tray \$1,000

Part: ACC2302500. HCPCS: E1029

Includes accessory bar, a credit will be applied PTOIN-DV3769.

**23. POSITIONING COMPONENTS****Optional.** More options are available on the [Positioning Components Order Form](#).**23.1. Stealth Unilink Pelvic/Thigh Guide****1. UniLink Pelvic/Thigh Hardware**

UniLink, 5" Height Adjustable, \$405 ea.

Removable Mounting Hardware

Part: ST-UTTP5-Q. HCPCS: E1028

UniLink, 7" Height Adjustable, \$423 ea.

Removable Mounting Hardware

Part: ST-UTTP7-Q. HCPCS: E1028

2. UniLink Pelvic/Thigh Pads

UniLink, 4" x 4" Pad \$164 ea.

Part: ST-UTP544056100-Q. HCPCS: E0953/E0956

UniLink, 4" x 8" Pad \$164 ea.

Part: ST-UTP548056100-Q. HCPCS: E0953/E0956

UniLink, 4" x 12" Pad \$164 ea.

Part: ST-UTP54D056100-Q. HCPCS: E0953/E0956

**4. UniLink Pelvic/Thigh Guide Installation**Install on Left
Part: PTOINDV3585Install on Right
Part: PTOINDV3586Do Not Install
Part: PTOINDV3719

Positioning components section continued on following page.

23. POSITIONING COMPONENTS CONTINUED**23.2. Stealth Thoracic Laterals**

Swing-away laterals may interfere with 2-post, flip up armrests.

1. Bracket

Swing-away, UniLink, Redesign Back, pair (ULS4042000-Q)..... \$1,032
 Part: ST-ULS4042000-Q. HCPCS: E1028
 Only available with Q-Backs.

2. Pads

Curved Pads, 3"x5", pair (TWBL3X5)..... \$164
 Part: INDPART2663. HCPCS: E0956
 Flat Pads, 3"x5", pair (TWBLP3X5) \$164
 Part: INDPART2676. HCPCS: E0956
 Curved Pads, 4"x6", pair (TWBL4X6)..... \$164
 Part: INDPART2664. HCPCS: E0956
 Flat Pads, 4"x6", pair (TWBLP4X6) \$164
 Part: INDPART2675. HCPCS: E0956

3. Installation

Install Laterals
 Part: PTOINDV3499
 Do Not Install
 Part: PTOINDV3500

23.3. Stealth Pelvic/Thigh Guides

Pelvic/thigh guides may interfere with the joystick receiver block and down posts. Please provide quantity next to check box.

1. Pelvic/Thigh Guide Bracket

Flip Down Bracket (TWBTLADD-HW) \$343 ea.
 Part: INDPART2804. HCPCS: E1028
 Flip Down with LEFT Full Surface Contact \$459
 Feature (TWBADD-L & LPHW407)
 Part: POS140428. HCPCS: E1028
 Removable Bracket (LP430)..... \$383 ea.
 Part: INDPART2805. HCPCS: E1028
 Flip Down with RIGHT Full Surface Contact..... \$459
 Feature (TWBADD-R & LPHW407)
 Part: POS140429. HCPCS: E1028
 Removable Bracket with 1" Offset \$383 ea.
 (LP430 & LHW-121)
 Part: INDPART2854. HCPCS: E1028

2. Pelvic/Thigh Guide Adapter

2" Height Adapter Plate (LHW-122)..... \$35
 Part: INDPART2821. HCPCS: K0108
 Height adapter is needed with TRU-Comfort seat cushions. Not available with the flip down bracket.
 2 1/4" Adapter Plate (LHW-129) \$82
 Part: POS140430. HCPCS: K0108
 Only available with the flip down full surface bracket.

3. Pelvic/Thigh Guide Pad

Short Pad, 4"x4" (TWBADD-4X4) \$85 ea.
 Part: INDPART2801. HCPCS: E0953/E0956
 Long Pad, 4"x12" (TWBADD-4X12) \$85 ea.
 Part: INDPART2803. HCPCS: E0953/E0956
 Medium Pad, 4"x8" (TWBADD-4X8)..... \$85 ea.
 Part: INDPART2802. HCPCS: E0953/E0956

4. Pelvic/Thigh Guide Installation

Install on Left
 Part: PTOINDV3585
 Install on Right
 Part: PTOINDV3586
 Do Not Install
 Part: PTOINDV3719

23.4. Chest Bar Assembly

Chest support is required with anterior tilt 20 and 30 degree functions. Omitting a chest support selection requires the provider to facilitate aftermarket chest supports that best fits consumer needs.

Small Chest Bar, 14" Foam Pad (ST-CBATB4ATT-S-Q)..... \$300
 Part: ST-CBA-TB4-S. HCPCS: K0108
 Available with 14" to 20" back widths.
 Medium Chest Bar, 16" Foam Pad (ST-CBATB4ATT-M-Q)..... \$300
 Part: ST-CBA-TB4-M. HCPCS: K0108
 Available with 17" to 22" back widths.

23.5. Stealth Stealth Positioning Chest Harness

Upper Classic, Medium (IFPS1103) \$221
 Part: ST-IFPS1103-Q. HCPCS: E0960
 Upper Contour, Medium (IFPS2103)..... \$221
 Part: ST-IFPS2103-Q. HCPCS: E0960
 Upper Classic, Large (IFPS1104) \$221
 Part: ST-IFPS1104-Q. HCPCS: E0960
 Upper Contour, Large (IFPS2104) \$221
 Part: ST-IFPS2104-Q. HCPCS: E0960

23.6. Stealth Stealth Positioning Chest Straps

Chest Strap, Medium w/ Belt Anchor, \$140
 (IFCS4113)
 Part: ST-IFCS4113-Q. HCPCS: E0960
 Chest Strap, Large w/ Belt Anchor, \$140
 (IFCS4114)
 Part: ST-IFCS4114-Q. HCPCS: E0960

23.7. Stealth Swing-Away Headrest Mounts

These swing-away headrest mounts are only available with the Stealth Comfort Series headrests.

1. Facial Lateral Hardware

UniLink Facial Lateral, Swing Away, Right..... \$496
 Part: ST-SAR2A-R-Q. HCPCS: E1028
 UniLink Facial Lateral, Swing Away, Left \$496
 Part: ST-SAR2A-L-Q. HCPCS: E1028
 UniLink Facial Lateral, Dual Swing Away \$765
 2 Upper & 1 Lower Right
 Part: ST-SAR4A-DBU-RL-Q. HCPCS: E1028
 UniLink Facial Lateral, Dual Swing Away \$674
 2 Upper
 Part: ST-SAR4A-DBU-Q. HCPCS: E1028



Positioning components section continued on following page.

23. POSITIONING COMPONENTS CONTINUED**2. Facial Lateral Pads and Accessories**


Egg Switch Mount (ESM) (A) \$70 ea.
 QTY Part: INDPART2815

Figure 8 Bracket with 13" Rod \$366 ea.
 QTY Part: ST-F8675-13-Q
 and Egg Switch Mount (B)

Facial Pad, 4".5x2", Cool Core (900M-Q) \$96 ea.
 QTY Part: ST-900M-Q. HCPCS: K0108

Facial Pad, 6"x2", Cool Core (920M-Q) \$96 ea.
 QTY Part: ST-920M-Q. HCPCS: K0108

Swing-Away Joystick Mount, 18" Rod (C) \$100 ea.
 QTY Part: ST-SUS18-Q
 Allows mounting of facial pads and attaches to the headrest

**23.8. Stealth Elbow Support****1. Elbow Support Bracket**

Left, Swing-Away, Elbow Support \$279
 Bracket Tube Style Arms (ARMS-ES-L)
 Part: POS124064. HCPCS: E1028

Left, UniLink, Fixed, Elbow Support Bracket \$309
 for Flip-Back, 8"-16"H Armrests (UETB4021FL-Q)
 Part: ST-UETB4021FL-Q


Right, Swing-Away, Elbow Support \$279
 Bracket Tube Style Arms (ARMS-ES-R)
 Part: POS124066. HCPCS: E1028

Right, UniLink, Fixed, Elbow Support Bracket \$309
 for Flip-Back, 8"-16"H Armrests (UETB4021FR-Q)
 Part: ST-UETB4021FR-Q

2. Elbow Support Pad

3"Wx3"H Elbow Pad (ARMS-ES-3X3) \$57 ea.
 QTY Part: POS124065. HCPCS: K0108

4"Wx4"H Elbow Pad (ARMS-ES-4X4) \$57 ea.
 QTY Part: POS124067. HCPCS: K0108

**23.9. Stealth Foot Boxes**

Stealth Small Foot Boxes, 5"Wx8"Dx7"H \$725
 (SPLE107) (pair)
 Part: INDPART2786. HCPCS: E0954

Stealth Medium Foot Boxes, 6"Wx10"Dx9"H \$725
 (SPLE108) (pair)
 Part: INDPART2787. HCPCS: E0954

Stealth Large Foot Boxes, 7"Wx12"Dx11"H \$725
 (SPLE109) (pair)
 Part: INDPART2788. HCPCS: E0954

23.10. Stealth Positioning Shoe Holder

Small Shoe Holder, 3/4" Webbing (IFSH4211) \$93 ea.
 QTY Part: ST-IFSH4211-Q. HCPCS: E0951+E0952

Medium Shoe Holder, 3/4" Webbing (IFSH4212) \$93 ea.
 QTY Part: ST-IFSH4212-Q. HCPCS: E0951+E0952

Large Shoe Holder, 3/4" Webbing (IFSH4213) \$93 ea.
 QTY Part: ST-IFSH4213-Q. HCPCS: E0951+E0952

24. LAP TRAY**24.1. Quantum Lap Tray**

Mounting bracket is required with lap tray selection. These trays are compatible with Flip-Back armrests.

1a. Black Lap Tray

Black Lap Tray, 11" W \$360
 Part: MEC1905328-11B. HCPCS: E0950
 For use with 12"-13" back width.

Black Lap Tray, 13" W \$360
 Part: MEC1905328-13B. HCPCS: E0950
 For use with 14"-15" back width.

Black Lap Tray, 15" W \$360
 Part: MEC1905328-15B. HCPCS: E0950
 For use with 16"-17" back width.

Black Lap Tray, 17" W \$360
 Part: MEC1905328-17B. HCPCS: E0950
 For use with 18"-19" back width.

Black Lap Tray, 19" W \$360
 Part: MEC1905328-19B. HCPCS: E0950
 For use with 20"-22" back width.

1b. Clear Lap Tray

Clear Lap Tray, 11" W \$360
 Part: MEC1905328-11C. HCPCS: E0950
 For use with 12"-13" back width.

Clear Lap Tray, 13" W \$360
 Part: MEC1905328-13C. HCPCS: E0950
 For use with 14"-15" back width.

Clear Lap Tray, 15" W \$360
 Part: MEC1905328-15C. HCPCS: E0950
 For use with 16"-17" back width.

Clear Lap Tray, 17" W \$360
 Part: MEC1905328-17C. HCPCS: E0950
 For use with 18"-19" back width.

Clear Lap Tray, 19" W \$360
 Part: MEC1905328-19C. HCPCS: E0950
 For use with 20"-22" back width.

2. Mounting Hardware

Mounting Bracket, Chest Bar/Lap Tray No Charge
 Part: MEC2104886

Lap tray section continued on following pages.

25. LAP TRAY**24.2. Simple Solutions Lap Tray**

Mounting bracket is required with lap tray selection. See chart below for compatibility.

1. Lap TraySimple Solutions Small Tray with Cut Out \$750
Part: ST-SSLXS010001171-Q. HCPCS: E0950Simple Solutions Large Tray with Cut Out..... \$750
Part: ST-SSLXL010001171-Q. HCPCS: E0950Simple Solutions Medium Tray with Cut Out..... \$750
Part: ST-SSLXM010001171-Q. HCPCS: E0950**2. Mounting Hardware**Swing Away Mounting Hardware.....No Charge
Part: ST-SSLXHD10000060-Q

Simple Solutions Lap Tray Armrest Compatibility	Width											
	10	11	12	13	14	16	17	18	19	20	21	22
Flip-Back	S	S	S	S	S	S	S	M	M	L	L	L
Flip-Back, Cane Mounted	S	S	S	S	S	S	M	M	M	L	L	L
Drop-In, T-Arm	S	S	S	S	S	S	M	M	M	L	L	L
Two-Post, Flip-Back	S	S	S	S	S	S	M	M	M	L	L	L
Cantilever	S	S	S	S	S	M	M	M	M	L	L	L
Q-Captain, Flip-Up						S	M	M	L	L	L	L
Single Post/ Heavy Duty/ Bariatric	S	S	S	S	S	S	M	M	M	M	L	L

S:ST-SSLXS010001171-Q M:ST-SSLXM010001171-Q L:ST-SSLXL010001171-Q

S:Compatible M:Armrest configuration may not hold tray position

*Some armrests listed may not be available on this order form.

25. ALTERNATIVE CONTROLS AND ELECTRONICSOptional. Select a type of specialty controls below. More options are available on the [Alternative Controls and Electronics order form](#).**25.1. Enhanced Display**Q-Logic 3 Enhanced Display Kit..... \$2,045
Part: CTL167726

Select a color option below. Includes enhanced display, flexible gooseneck and mounting bracket. Features: 3.5" LCD display; up to 5 profiles to control drive, seat, and auxiliary functions; Bluetooth connectivity for mobile devices and mouse emulation; integrated single switch scanner; IR controls.

Q-Logic 3 Enhanced Display Shroud

Required with Q-Logic 3 Enhanced Display selection above. Reference colors shown on page 2.

Black, Soft Touch Standard
Part: DGN137728Deep Purple
Part: DGN137729-PPR5Electric Blue
Part: DGN137729-PBU9Military Digital Camo..... \$50
Part: DGN137729-WCM7Red Corvette
Part: DGN137729-PRD1Raspberry Beret
Part: DGN137729-PPK4Blue Jean
Part: DGN137729-PBU11Faded Glory \$50
Part: DGN137729-WFL2Back in Black
Part: DGN137729-PBK1Pretty in Pink
Part: DGN137729-PPK3Lemon Crush
Part: DGN137729-PYL4Carbon Fiber \$50
Part: DGN137729-WCB1Night Train
Part: DGN137729-PBK3Evolution Orange
Part: DGN137729-POR7Green Machine
Part: DGN137729-PGN11Grey Street
Part: DGN137729-PSL3Midnight Blue
Part: DGN137729-PBU10Military Green
Part: DGN137729-PGN8**25.2. SCIM Module**

Optional.

SCIM Module for Q-Logic 3 \$850
Part: CTL2001963

Allows use of joystick display for specialty devices. Available with Q-Logic 3 EX Joystick selection above. Up to 5 profiles to control drive, seat, and auxiliary functions, excluding sip-n-puff and single switch functions. A switch is required to operate the input device attached to the SCIM module.

25.3. Q-Logic 3 SAJQ-Logic 3 EX SAJ \$850
Part: CTL143980

Must be used in conjunction with Enhanced Display or Q-Logic 3 Joystick.



Enhanced Display



SCIM Module



EX SAJ

Alternative controls and electronics section continued on following page.

25. ALTERNATIVE CONTROLS AND ELECTRONICS CONTINUED**25.4. Head Array**

i-Drive 4.0 Bluetooth Pro Series,\$7,406
 3 Switch with Egg Switch
 Part: ST-IDHBT300-U-Q. HCPCS: E2330



Adult Pro-Series

25.5. Sip and Puff Drive**1. Sip and Puff Module**

Sip and Puff Module for Q-Logic 3,\$1,695
 Part: CTL168056. HCPCS: E2325
 Includes module and mounting hardware.



Sip & Puff Module

2. Sip and Puff Straw and Mount

The straw and mount do not include electrical components to operate the Sip&Puff with a head array. This component requires the selection of a Q-Logic Sip&Puff module or additional module from Stealth Products for operation.

Stealth Sip and Puff for Comfort Series,\$496
 Left Mount (SP625-CL)
 Part: ST-SP625-CL-Q. HCPCS: E2326
 Includes straw and mounting hardware. Use with Stealth Comfort Plus headrests.

Stealth Sip and Puff for Stealth i-Drive,\$496
 Head Array (SP625-U)
 Part: ST-SP625-U-Q. HCPCS: E2326
 Includes straw and mounting hardware. Use with Stealth i-Drive Head Arrays.

Stealth Sip and Puff for Comfort Series,\$496
 Right Mount (SP622-CR)
 Part: ST-SP625-CR-Q. HCPCS: E2326
 Includes straw and mounting hardware. Use with Stealth Comfort Plus headrests.



Stealth Sip & Puff

3. Sip & Puff Accessories

Sip & Puff Straws, 12 pack (SP625-S12),\$47
 Part: ST-SP625-S12-Q

25.6. Egg Switches

Red Egg Switch (ESRED),\$116
 Part: SWTMCR01077
 QTY
 Blue Egg Switch (ESBLU),\$116
 Part: SWTMCR01079
 QTY
 Black Egg Switch (ESBLK),\$116
 Part: SWTMCR01081
 QTY

Green Egg Switch (ESGRN),\$116
 Part: SWTMCR01078
 QTY
 Yellow Egg Switch (ESYEL),\$116
 Part: SWTMCR01080
 QTY



Egg Switch

25.7. Switch Adapter

2 to 1 Switch Adapter,\$135
 Part: DWR1361H037
 Connects 2 single switches to a stereo jack to act as a dual switch.

25.8. Environmental Controls

Q-Logic 3 ECU Module Kit,\$1,780
 Part: CTL168057

25.9. USB Mobile Device Charger

USB Mobile Device Charger,\$185
 Part: ELE168069
 Can charge 1 device. A multiplier harness is required and must be added when a Static Seat or Power Tilt Only seating system is chosen. The selection of more than 3 modules requires an additional multiplier.



USB Mobile Device Charger