



SCOUT EXCEL COLOR OPTIONS:		
CATALOG NO.	PLASTIC COLOR	POWDER COATING COLOR
F-SCTXL	BLUE	WHITE
F-SCTXL-T	TAN	WHITE
F-SCTXL-G	GRAY	WHITE
F-SCTXL-W	WHITE	WHITE
F-SCTXL-C	CUSTOM	CUSTOM



TOLERANCES
DIMENSIONS ARE IN INCHES
UNLESS OTHERWISE SHOWN
FRACTIONAL: ±1/16
ANGULAR: ±1°
TWO PLACE DECIMAL: ±.01
THREE PLACE DECIMAL: .. ±.005


BREAK & DEBURR ALL SHARP EDGES



SCOUT EXCEL
CATALOG NO. F-SCTXL
PART NO. INV02394
WEIGHT CAPACITY: 375lbs. [170.09kg.]

- FEATURES**
- ◆ UL CERTIFIED ◆ 375 LB. (170 kg.) WEIGHT CAPACITY ◆ SETBACK RANGE 7" TO 44" (18 - 111 cm) ◆ CLEARANCES VARY BY WALL HEIGHT & WIDTH
 - ◆ WATER DRAFT 10" ◆ 360° ROTATION ◆ FLIP UP ARMREST ◆ ALL STAINLESS STEEL CONSTRUCTION ◆ DURABLE POWDER COAT FINISH ◆ RECHARGEABLE BATTERY & WALL MOUNT CHARGER ◆ FIVE YEAR STRUCTURAL WARRANTY ◆ FIVE YEAR PRO-RATED ELECTRONICS WARRANTY

CLICK HERE TO VIEW THIS PRODUCT:
<https://aquacreekproducts.com>



MATERIAL: 304 STAINLESS STEEL
FINISH: POWDER COAT - WHITE
APPROX. WEIGHT: 196.57 LBS.

Aqua Creek Products
LEADERS IN RECREATION, FITNESS AND ABILITY

9889 GARRYMORE LANE
MISSOULA, MT 59808
(406) 549-0769
FAX (406) 549-2602
www.AquaCreek.com

© 2021 AQUA CREEK PRODUCTS. ALL RIGHTS RESERVED.

PROPRIETARY AND CONFIDENTIAL INFORMATION

THIS DRAWING AND ALL INFORMATION CONTAINED HEREIN IS THE SOLE PROPERTY OF AQUA CREEK PRODUCTS LLC. ANY REPRODUCTION IN PART OR AS A WHOLE WITHOUT WRITTEN PERMISSION FROM AQUA CREEK IS PROHIBITED.

*SPECIFICATIONS ARE SUBJECT TO CHANGE WITHOUT NOTICE

REVISIONS			
REV.	DESCRIPTION	DATE	APPROVED
P1	INITIAL RELEASE	3/30/2021	J.MAXWELL
P2	ENCLOSED LOWER CASING	2/28/2022	J.MAXWELL
P3	ADDED SOLAR PANEL ATTACHMENT	12/19/2022	J.MAXWELL