

MD Series Harnesses

PivotFit™ • Stayflex™ • Trimline

Same great harnesses, now with MD (multi-directional) Pull Straps. No need to choose in advance between Front-Pull and Rear-Pull — the versatile MD Pull Straps work both ways! Easily set up and change pull styles on the spot to adapt to individual user needs. Simplified ordering, and fewer compromises at the time of fitting.

Attach straps to fastener on harness shoulder pad. Simply reverse direction of strap assembly to switch between Front and Rear-Pull.



Front-Pull Configuration



Rear-Pull Configuration

Patented, sewn-in D-ring looks sharp and always holds its position.

Two-piece Flat-Mounts™ make it easy to adjust strap length after installation. Sized to fit M6 screws.

Updated end-fittings, too!

Our MD Series harnesses come with four of our, high-strength plastic end-fittings. Choose from plastic Flat-Mounts™ or Cinch-Mounts™. Optional Grommet Straps available separately.

End-fittings also available separately in 2-packs. See Price List for bulk offerings.

Plastic Flat-Mount™



- B5

Plastic Cinch-Mount™



- B6

Grommet Strap (sold separately)



SH025-2

206.405.4555 or 800.547.5716 | sales@bodypoint.com | bodypoint.com

View our product video library: [YouTube.com/BodypointInc](https://www.youtube.com/BodypointInc)

© 2021 Bodypoint. All Rights Reserved. BMM319 2021.6



MD Series Harness Styles & Part Numbers

PivotFit™ Shoulder Harness

Our advanced-design shoulder harness distributes pressure across a greater surface area than a traditional H-style, with pads contoured around the shoulders and ribcage to provide optimal support and comfort. Center swivel buckle equalizes tension to accommodate asymmetry. The Standard style is non-stretch to reduce movement; the Dynamic allows controlled torso movement.



PivotFit™ Standard

Size	Part number
XXS	SH290XXS
XS	SH290XS
S	SH290S
M	SH290M
L	SH290L
XL	SH290XL



PivotFit™ Dynamic

Size	Part number
XXS	SH295XXS
XS	SH295XS
S	SH295S
M	SH295M
L	SH295L
XL	SH295XL

Stayflex™ Chest Support

Provides wheelchair users firm trunk control without sacrificing comfortable shoulder movement. Patented, dual stretch zone construction eliminates upwards shift at the neckline when leaning forward, minimizing risk of unintended airway restriction. Carefully selected materials optimally balance stretch and resistance, and quick-release swivel buckles allow easy removal.



Stayflex™ Standard, with zipper

Size	Part number
S	SH350S
M	SH350M
L	SH350L
XL	SH350XL



Stayflex™ Standard, without zipper

Size	Part number
XS	SH351XS
S	SH351S
M	SH351M
L	SH351L
XL	SH351XL



Stayflex™ Narrow, with zipper

Size	Part number
M	SH360M
L	SH360L
XL	SH360XL



Stayflex™ Narrow, without zipper

Size	Part number
M	SH361M
L	SH361L
XL	SH361XL

END-FITTINGS

Add extension -B5 or -B6 to the end of the part number.

Plastic Flat-Mount™



Plastic Cinch-Mount™



Trimline Shoulder Harness

Features two separate, adjustable padded straps for more focused shoulder retraction. Comfortably accommodates asymmetry and keeps the chest clear for maximum freedom of movement.



Size	Part number
S	SH235S
M	SH235M
L	SH235L

206.405.4555 or 800.547.5716 | sales@bodypoint.com | bodypoint.com

View our product video library: [YouTube.com/BodypointInc](https://www.youtube.com/BodypointInc)

© 2021 Bodypoint. All Rights Reserved. BMM319 2021.6

