

# FREEDOM AQUA ELECTRIC PATIENT LIFT

This innovative chair not only delivers all the trusted benefits of our Freedom Patient Chair but also brings the game-changing feature of water resistance. Perfect for shower use, the Freedom Aqua ensures safety, comfort, and reliable support, making every moment worry-free.

## Best for Patient Transfer/Bathing



**300 LBS.**  
WEIGHT CAPACITY



**47 LBS.**  
UNIT WEIGHT



**200 TRANSFERS**  
PER CHARGE

### Water-Resistant Design

This feature ensures that the device remains functional and reliable, even when exposed to water during bathing.

### Convenient Transfer Design

The powered lift with built-in controls allows the seat to be adjusted from 19" to 29" to match the desired height of your destination.

### Electric Adjustable Seat Height

Engineered for durability, the Freedom Aqua features a robust construction that withstands daily use.

### Enclosing Seat

The seat opens 90° and wraps around the patient, making it easy to safely transfer from a seated position to the chair.



## PRODUCT SPECIFICATIONS

**SKU: FCOM-02**

**LEARN MORE!**

LIFT CAPACITY

300 POUNDS

BATTERY LIFE

200 LIFTS PER CHARGE

INCLUDED ACCESSORIES

BATTERY, BATTERY CHARGER,  
COMMUNE PAIL, DETACHABLE REMOTE, & TOOL KIT



[www.MobilePatientLift.com](http://www.MobilePatientLift.com)



[Info@ClimbingSteps.com](mailto:Info@ClimbingSteps.com)



(212) 366-5483

## PRODUCT DETAILS

HEIGHT	33.5 - 43.25 INCHES
DEPTH	27 INCHES
WIDTH	22.5 INCHES
ADJ. SEAT HEIGHT	19 - 29 INCHES
FRONT WHEEL HT.	5 INCHES
BACK WHEEL HT.	3 INCHES
HANDLE LENGTH	5 INCHES
HANDLE HEIGHT	10 INCHES
FOOTREST LENGTH	14 INCHES
FOOTREST WIDTH	6 INCHES
SEAT WIDTH	20 INCHES
SEAT DEPTH	17 INCHES
ARMREST LENGTH	10 INCHES
REAR FRAME HEIGHT	4 INCHES
WEIGHT	47 POUNDS

## EXTENSION DETAILS

MIN & MAX SEAT HEIGHT	16" MIN - 24" MAX
-----------------------	-------------------

## SHIPPING DETAILS

BOX LENGTH	INCHES
BOX WIDTH	17 INCHES
BOX HEIGHT	31 INCHES
BOX WEIGHT	63 POUNDS

## TECHNICAL DETAILS

FRAME MATERIAL	STEEL	BATTERY TYPE	LITHIUM ION
FRAME COLOR	WHITE	BATTERY WEIGHT	1 POUND
FRAME SEAT MATERIAL	STEEL	BATTERY SIZE (WXHXL)	8" X 3" X 1"

## ORDERING DETAILS

PRODUCT NAME	FREEDOM AQUA TRANSFER PATIENT LIFT
PRODUCT SKU	FCOM-02
PRODUCT UPC	850031392502



## FOLDED SPECIFICATIONS

LENGTH	12 INCHES
WIDTH	17 INCHES
HEIGHT	26.5 INCHES