





- Effective micro low air loss addresses moisture management.
- Interchangeable air cell design provides easy cleaning.

emotion

Features

Alternating Pressure

Provides active pressure redistribution by using 3-in-1 alternating pressure cycle. This results in a comfortable, effective therapy for patients.

Mattress Design

Mattress has an 8" air bladder with cell in cell (4.5" of therapy, 3.5" of static) to provide therapeutic alternating pressure and allowing the mattress to stay inflated in a power outage.

Advanced Pump Options

Adjustable comfort settings, variable cycle times and a max inflate option for customizable comfort to fit the patient's needs.

Weight Limit: 440 lbs. (200 kg)

Mattress: 80" x 36" x 8"

Product Weight: 20.6 lb

Cover: Top: Nylon coated with PU + Single

quilting. Bottom/Base: Polyester laminated with PVC, with adjustable safety straps

Warranty: One Year Limited Warranty on control unit

and soft goods

Medicare Code: E0277 Group 2 Support Surface

LIT-0010 | R100121



1939 S. Vandeventer Ave. St. Louis, MO 63110

T 1-800-942-7678 E info@moxiusa.com W moxiusa.com