



1. Cell with Laser Holes- 14029DP-CELL-1
(middle 9 cells)

Cell w/o Laser Holes- 14029DP-CELL

2. Pump 2 Port- [4]

Pump 3 Port- 14029XP

Pump 3 Port- 14029XPN [2]

3. Pump Hook- 14030HK

4. Intake Filter- 14029FN

Intake Filter Cover- 14029XP-FC

5. Fuse 91 amp, 250 volts)- 14530F [1]

6. Power Receptacle- 14030PR [1] [2]

Power Receptacle- 14030PR-6G [3]

7. Power Cord- 15005PC [1] [2]

8. Complete Mattress 2 Port- [4]

Complete Mattress 3 Port- 14029DPMN

9. Mattress Cover- 14029DPC

10. Mattress Base- 14029DPB

11. EZ Lock Pump Socket 2 Port- 14028S

EZ Lock Pump Socket 3 Port- 14030XP-S

"O" Rings (3 ea)- 14007OR-3

12. One Way Valve- 14028V

13. CPR Valve 2 Port- [4]

CPR Valve 3 Port- 14029DP-CPR

14. Side Bolster Left- 14029EAL

Side Bolster Right- 14029EAR

[1] s/n 4E1506016611 & later.

[2] 3 prong plug s/n 6G201210746 & before.

[3] 2 Prong Plug s/n 6G201210747 & later.

[4] This part is no longer available.