

Prime Care Bed Model P703

HCPCS:



FIGURE A




FIGURE B



FIGURE C



FIGURE D



FIGURE E

- 550 lbs. weight capacity, Resident Weight 500 lbs
- 7" to 30" height range
- Optional integrated width adjusts to 39" - 42"
- Integrated 76", 80" and 84" sleep surface length
- Two way Red/Green color coded floor lock system
- Eight (8) Function Handheld Pendant Control
- Auto Contour
- 4 DC Motors Synchronized Hi/Lo control system keeps bed level
- Four sealed bearings secured in "G" track
- IV pole holders
- Extended Heel Lift
- Eight 3" Heavy Duty casters roll at any height
- Fixed adjustable wall bumper

PRODUCTS	Description
P703	Prime Care Model P703, 1ea
SPECIFICATIONS	Description
Warranty	Lifetime on Welds
Wheels	Accommodate up to 550 lbs Maximum Resident Weight 500 lbs
Bed Height Range	7" - 30"
Warranty Limited	15 Year on Frame
Warranty on Batteries	5 years on Electrical
Weight without battery	205 lbs.
ACCESSORIES	Description
PLTCDR	LTC Deluxe Rail,1PR, PR
PLTCAR	LTC Assist Rail,1PR, PR
PLTCAB	LTC Assist Bar, 1PR, PR
P703LE	P703 Length Extension 88",1ea, EA
P703-BBK	P703 Battery Backup,1ea, EA
P703SC04	Staff Control Panel f/P703 ,1ea, EA
P703WE	P703 Width Extension 39-42",1pr, PR
35BWBE-B-SC04	35" BRNWOOD BEW/SC04,OKOAK,EA, EA