

Cruiser X4 Wheelchair

HCPCS: K0004,K0004



FIGURE A



FIGURE B



FIGURE C



FIGURE D



FIGURE E



FIGURE F



FIGURE G



FIGURE H

- Weighs 33.5 pounds or less when properly equipped in the basic configuration (excluding front riggings)
- Easy to propel and transport
- Frame style eliminates seat rail guides and allows for custom backs and accessories
- Flip-back, removable arm is height adjustable from top of seat rail 10" to 14" and features a front mounted easy-to-use lever release mechanism
- Padded armrests provide added user comfort
- Height adjustable back
- Aluminum and carbon steel blended frame with combination powder coated and silver vein finish is attractive, chip resistant and easy-to-maintain
- Durable, flame retardant nylon upholstery resists mildew and bacteria
- Patented adjustable depth seat form 16" and 18"
- Three position quick release rear axles mounted on composite mag-style wheels with composite hand rims are lightweight and maintenance free
- Polyurethane rear tires and composite casters offer superior performance, smoother ride and are lightweight
- Dual over center cross frame assists provide added rigidity
- Includes swing-away footrests or elevating leg rests with tool-free adjustable length riggings
- Composite footplates with heel loops are durable and lightweight
- Precision sealed wheel bearings ensure durability and reliability
- 8" front casters with three height adjustments and are angle adjustable
- Triple rear axle positions provides easy transition of seat height from standard to hemi position
- Comes with push-to-lock wheel locks

PRODUCTS	Description
CX416ADDA-ELR	CruiserX4, 16" Adj DeskArm, ElvLg, 1/ea
CX416ADDA-SF	CruiserX4, 16" Adj DeskArm, SwgFt, 1/ea
CX416ADFA-ELR	CruiserX4, 16" Adj FullArm, ElvLg, 1/ea
CX416ADFA-SF	CruiserX4, 16" Adj FullArm, SwgFt, 1/ea
CX418ADDA-ELR	CruiserX4, 18" Adj DeskArm, ElvLg, 1/ea
CX418ADDA-SF	CruiserX4, 18" Adj DeskArm, SwgFt, 1/ea
CX418ADFA-ELR	CruiserX4, 18" Adj FullArm, ElvLg, 1/ea
CX418ADFA-SF	CruiserX4, 18" Adj FullArm, SwgFt, 1/ea
CX420ADDA-ELR	CruiserX4, 20" Adj DeskArm, ElvLg, 1/ea
CX420ADDA-SF	CruiserX4, 20" Adj DeskArm, SwgFt, 1/ea
CX420ADFA-ELR	CruiserX4, 20" Adj FullArm, ElvLg, 1/ea
CX420ADFA-SF	CruiserX4, 20" Adj FullArm, SwgFt, 1/ea

SPECIFICATIONS	Description
Warranty	Limited Lifetime
Seat Dimensions	16" (W) x 16" - 18" (D)
Weight Capacity	300 lb
Seat To Floor Height	17" 18" 19"
Carton Shipping Weight	46 lb 48 lb 50 lb
Folded Dimensions	43" (L) w/ riggings x 12" (W) x 36"
Back Height	16" - 17.5"
ACCESSORIES	Description
STDSSHK-K3	Super Hemi Kit, Cruiser 3 & X4, 1/ea, EA
STDS851	Seat Belt, Velcro Close, 1/ea, EA
STDS855	Seat Belt, 60" Auto-type, 1/ea, EA
STDS820	Telescoping IV Pole, W/C, 1/bx, BX
STDS856	Seat Belt, 60" Velcro-type, 1/ea, EA
STDS806	Anti-Fold/Theft Device for w/c, 1/ea, BX
STDS850	Seat Belt, Clasp Type, 1/ea, EA
STDS834	Overhead Anti-Theft Bar, Single, 1/ea, EA
STDS1034	Wheelchair Cane/Crutch Holder, 1/ea, EA
WASL	Limb Support, Left, 1/cs, CS
WASR	Limb Support, Right, 1/cs, CS
STDS601	Brake Extension, 6", 1/ea, EA
STDS1002	Caster, 8in Sentra, Viper, Astaire, 1/ea, EA
STDS860	Forwrd Stablzrs/Univ. Fit, 1 pr/bx, BX