

PreserveTech Universal Raised Toilet Seat

HCPCS: E0244,E0244



FIGURE A



FIGURE B



FIGURE C



FIGURE D



FIGURE E



FIGURE F

Introducing the PreserveTech Universal Raised Toilet Seat, our first ever truly universal raised toilet seat accommodates three different seat heights (2", 4" and 6") by maneuvering just two components. This seat can fit on most standard and elongated toilets and features two hygienic cutaways for easy cleaning. The pivoting armrests can be raised or removed, depending on user needs and bathroom space and also includes a lid with a recessed grip for effortless opening.

- Accommodates three different seat heights (2", 4" and 6") by maneuvering just two components.
- Pivoting armrests can be raised out of the way or removed depending on your needs or bathroom environment.
- Fits comfortably on most standard and elongated toilets.
- Includes a lid with recessed grip for better hygiene and aesthetic and effortless opening.
- Contoured to provide a safe and comfortable elevation.
- Drive PreserveTech products are manufactured with an antimicrobial* agent that helps to minimize the growth of stain and odor-causing bacteria and mold on treated surfaces.
- *Antimicrobial properties are built in to protect the products. These products do not protect users or others against bacteria, viruses, germs or other disease organisms.

PRODUCTS

Description

RTL12C002-WH	PreserveTech Universal RTS, 1EA, 1/CV
--------------	---------------------------------------

SPECIFICATIONS

Description

Warranty	Limited Lifetime
Weight Capacity	350 lb
Width Between Arms	18.5"