

Two Wheeled Walker with Seat

HCPCS: E0143/E0156-1,E0143/E0156-1



- Features comfortable and convenient padded seat and backrest
- Quickly folds when not in use
- Rear leg tips act as brakes when pressed down
- Height adjustable handles
- Large 5" swivel wheels for added maneuverability
- Contoured handles for a comfortable grip

PRODUCTS	Description
1239RD	2 Wheeled Walker w/Seat, Red

SPECIFICATIONS	Description
Warranty	Limited Lifetime
Seat Dimensions	15.875" (W) x 8" (D)
Handle Height	30.5" - 37.8"
Weight Capacity	300 lb
Carton Shipping Weight	16.5 lb
Overall Width	24"
Weight	14.3 lb