

GRAVIS Floor Lift, Manual Base 600 lbs

HCPCS: E0635



FIGURE A



FIGURE B



FIGURE C



FIGURE D



FIGURE E



FIGURE F



FIGURE G



FIGURE H

- Weight capacity of 600 lbs
- Base option: Manual
- Elegant, contemporary design with soft lines but still heavy duty
- Oversized sealed rear casters roll easily and lock securely
- Easily lifts resident from floor
- Redundant controls allow for use even if hand control is inoperable
- Smart & Safety Technology LCD display
- Visual battery fuel indicator on Hand Pendant and Lift Control Module
- Auditable low battery alert fuel indicator
- Lift safety stop if resident or caregiver accidentally gets in the way of the boom arm
- Audible Alert if the lift detects weight overload
- Maintenance monitoring: by months or lift cycles
- Audible Alert that maintenance is due
- NEW Six-point cradle with spring clips offer a more secure and wider foot print reducing sling squeeze, shallow height
- offers greater height range
- Can be charged on, or off the lift
- Includes 1 battery and external charger stand

PRODUCTS	Description
FLNP600	Gravis Floor lift,manual base,600 lbs
SPECIFICATIONS	Description
Battery Pack Weight	7.9 lbs (5.0A), 4.9 lbs (2.9A)
Battery Warranty	1 yr
Boom Height	26" - 76"
Ground Clearance	2.9"
Warranty on Battery	1 yr
ACCESSORIES	Description
FLSCALE	Patient Lift Scale 660lbs capacity, EA
13012	Patient Sling,Canvas,1/bx, BX
13019	Chains for Patient Lft,1pr/bx, BX
13025	Sling,Polyst.r.for Lift 1/ea, BX
13026	Sling,Polyst.r.w/Cmde Open1/ea, BX
13060	BariatricSling,Canvas,1/bx, BX
13061	BariatricSlingw/ComdeOpen,Canvas,1/bx, BX

ACCESSORIES**Description**

13232	Straps for Patient Slings,1set, BX
13220L	U-Sling 4pt padded Lg,1/bx, BX
13220M	U Sling 4pt padded Med,1/bx, BX
13220S	U-Sling 4pt padded Sm,1/bx, BX
13220XL	U-Sling 4pt padded X-lg,1/bx, BX