

3-Position Recliner, Bariatric

HCPCS:



- 3 positions: upright (Figure A), deep recline and elevated leg rest/footrest (Figure B)
- Retractable lock bar secures the chair in the desired position
- Self-adjusting suspension in the seat prevents bottoming out
- Side panels "pop off" for easy cleaning and maintenance
- Large, blow-molded tray locks in 5 positions and can be stored on the side of the chair
- 5" casters (2 with locks) make it easy to move and secure the recliner
- Wide, comfortable armrests are angled to provide support in any reclined position
- Neutral-colored, sturdy steel frame is easy to clean and maintain
- Moisture barrier on seat prevents seepage
- Comfortable built-in headrest
- Foam and vinyl meet California Technical Bulletin 117

PRODUCTS	Description
D574EW-CHAR	3 Position Extra Wide Charcoal
D574EW-CHOC	3 Position Extra Wide Chocolate
D574EW-TAN	3 Position Extra Wide Tan
SPECIFICATIONS	Description
Warranty	Limited Lifetime on Frame, 1 Year on Upholstery and Moving Parts
Seat Dimensions	20" (D) x 25" (W) x 21" (H)
Weight Capacity	500 lb
Carton Shipping Weight	127 lb
Width Between Arms	24.5"