

## Transfer Tub Bench

HCPCS: E0247,E0247



FIGURE A



FIGURE B



FIGURE C



FIGURE D



FIGURE E

- Tool-free assembly of back, legs and arm
- A-frame construction provides stability
- Durable blow-molded plastic bench and backrest
- Height adjusts in ½" increments with unique "Dual Column" extension legs
- Reversible to accommodate any bathroom
- Pinch-free lever allows for push pins to be depressed without pinching fingers
- Made using recycled materials

| PRODUCTS               | Description                       |
|------------------------|-----------------------------------|
| 12011KD-1              | Transfer Bench KD,Tool Free, 1/cs |
| 12011KD-2              | Transfer Bench KD Tool Free,2/cs  |
| SPECIFICATIONS         | Description                       |
| Warranty               | Limited Lifetime                  |
| Seat Dimensions        | 26" (W) x 18.5" - 19.5" (D)       |
| Weight Capacity        | 400 lb                            |
| Carton Shipping Weight | 13.6 lb                           |
| HCPCS                  | E0247                             |
| Height Adjustment      | 17.5" - 22.5"                     |
| Weight                 | 13.6 lb                           |
| ACCESSORIES            | Description                       |
| 12011SS                | SuctionCups,Sm.12011,1pr/bx, BX   |
| 12011SL                | SuctionCups,Lg.12011,1pr/bx, BX   |