

Deluxe Bariatric Shower Chair with Cross-Frame Brace

HCPCS: E0245



- Frame has width/depth cross braces attached with aircraft-type rivets
- Durable, blow-molded plastic bench
- Aluminum frame is lightweight, durable and corrosion proof Drainage holes reduce slipping
- Euro-style clip for leg height adjustment Height adjusts in 1" increments (Figure A)
- Weight Capacity 500 lbs
- Made using recycled materials

PRODUCTS	Description
12021KD-1	Bath BenchHD,w/Bk,1/ea
12022KD-1	Bench,HD,w/o.Bk,1/ea
SPECIFICATIONS	Description
Warranty	Limited Lifetime
Seat Dimensions	20" (W) x 12" (D)
Weight Capacity	500 lb
Carton Shipping Weight	32.5 lb 8.5 lb 27 lb 8.0 lb
Weight	7.5 lb 6.0 lb
ACCESSORIES	Description
12021T	RubberTip,12021,1/ea, EA