

Bariatric Transfer Board With Hand Holes

HCPCS:




FIGURE A

- For ease of transferring to and from a wheelchair (Figure A)
- Two cut-out handles provide additional maneuverability
- Made of sturdy Baltic birch with a clear lacquer finish
- Ends tapered, corners rounded
- Slightly beveled on each end
- Weight Capacity 600 lbs
- Size: 35" x 8"

PRODUCTS

Description

RTL7047	BariatricTransBoard w/Handslots,1EA,2/CS
---------	--

SPECIFICATIONS

Description

Warranty	Limited 6 Months
Weight Capacity	600 lbs