

Shower Chair with Folding Back

HCPCS:



FIGURE A



FIGURE B



FIGURE C



FIGURE D



FIGURE E

- Adds safety and support in the bath tub and shower
- Easy, tool-free assembly and height adjustable legs
- Folding back (Figure A) and removable legs make this shower chair great for storage and transport
- Lightweight design with stable legs and non-slip feet
- Armrests and soap tray add comfort and convenience
- Retail packaged
- Made using recycled materials

PRODUCTS	Description
RTL12606	BATH SEAT WITH BACK, 1 RTL, 2 RTL/CS
SPECIFICATIONS	Description
Warranty	Limited lifetime
Seat Dimensions	21.5" (W) x 13.125" (L)
Weight Capacity	300 lb
Carton Shipping Weight	8.6 lb
Weight	8 lb