

ProBasics Quad Cane, Small Base (Black)

HCPCS:



- ProBasics Quad Canes feature a unique four-point design which allows for greater weight-bearing and stability than a single-point cane.
- The handle rotates allowing for a proper fit with right or left hand use.
- Small and large base models available in black or silver finish with a limited lifetime warranty.
- Unique four-point design allows for greater weight bearing and stability than a single-point cane
- Non-skid rubber tips on all four feet enhance durability for long term use
- Small base design allows for easier maneuverability while the large base design provides more stability
- Handle rotates for left or right hand use
- Quad canes stand up on their own when not in use
- Easy push button height adjustment in 1-inch increments with a locking, anti-rattle ring allowing for quick, safe adjustment
- Comfortable, contoured vinyl handle with sculpted grip reduces chance of slippage
- Lightweight aluminum upper shaft with welded steel base

PRODUCTS

Description

CNQSB	QUAD CANE SM BASE, BLACK, 6/CS 6EA/CV
-------	---------------------------------------