




NOMINAL SIZE CHART

DIMENSION CHART (INCHES)									
PART NO.	SIZE	A	B	C	D	E	F	MAXIMUM SAFE WORKING LOAD	SIZE TRIM
547050-10	S	35	52	32	20	8	2	625lbs/284kg	RED
547060-10	M	40	57 1/2	34 1/2	23	10	4	625lbs/284kg	YELLOW
547070-10	L	45	65	37	28	12	4	625lbs/284kg	GREEN
547080-10	XL	50	72 1/2	39 1/2	33	14	4	625lbs/284kg	BLUE
CUSTOM SIZE									

APPROVAL: _____

COMPANY: _____

DATE: _____

Scale and / or dimensions are approximate (within +/- 1" range), and are to be used as a guide only.
This drawing and the information contained is the property of Handicare and is not to be copied or issued -E & OE

REV F 21-JUN-2023