




handicare

SLING SPECIFICATION SHEET

UNIVERSAL SPACER SLING with HEAD SUPPORT



NOMINAL SIZE CHART

DIMENSION CHART (INCHES)													
PART NO.	SIZE	A	B	C	D	W1	W2	G	L	H	S	MAX SAFE WORKING LOAD	SIZE TRIM
8A4730	XXS	25	15	2	5.5	24	24	18	16	16	10.5	625lbs	BLACK
8A4630	XS	30	18	3.5	8	29.5	29.5	18	18.5	18.5	10.5	625lbs	BLACK
8A4530	SMALL	35	22	6	8	32	33.5	18	18.5	18.5	10.5	625lbs	BLACK
8A4430	MED	41	25	6	12	39	40	18	22	22	10.5	625lbs	BLACK
8A4330	LARGE	41	28.5	6	14	42	44	18	22	22	10.5	625lbs	BLACK
8A4230	X-LARGE	41	32.5	6	15	45	47	18	22	22	10.5	625lbs	BLACK
CUSTOM SIZE													

APPROVAL

COMPANY

DATE

Scale and / or dimensions are approximate (within +/- 1" range), and are to be used as a guide only.
This drawing and the information contained is the property of Handicare and is not to be copied or issued -E & OE

REV G 29-JUN-2023

10888 Metro Court, Maryland Heights, MO 63043
81 Romina Drive, Concord, ON L4K 4Z9

Handicare NA
Direct: (314) 692-9135
Direct: (905) 850-0093
www.handicare.com

Fax: (314) 692-7858
Fax: (416) 260-0256
Toll Free: (866) 891-6502
Toll Free: 1-(877) 304-5438