

## Nitro® Duet Rollator and Transport Chair

HCPCS: E0143 + E0156,E0143 + E0156




FIGURE A



FIGURE B




FIGURE C



FIGURE D



FIGURE E



FIGURE F



FIGURE G



FIGURE H

- Combines the features of a rollator and transport chair into one unit. The user can ambulate independently or be pushed safely by a caregiver
- The tool-free two-position contoured backrest is reversible depending on whether you are using the Nitro Duet as a rollator or transport chair
- Back support height easily adjusts with tool-free buttons
- Caster fork design and large 10" front casters enhances turning radius and rolling comfort
- Cross-brace design allows for side-to-side folding and added stability
- Brake cables routed inside the frame add safety and convenience

### PRODUCTS

### Description

|            |                        |
|------------|------------------------|
| RTL10266DT | Nitro Duet, Red, 1 Ea. |
|------------|------------------------|

### SPECIFICATIONS

### Description

|                      |                     |
|----------------------|---------------------|
| Warranty             | Limited Lifetime    |
| Seat Dimensions      | 17.5" (W) x 13" (D) |
| Handle Height        | 33" - 37"           |
| Weight Capacity      | 300 lb              |
| Seat To Floor Height | 20.5"               |
| Casters              | Back 6"   Front 10" |
| Weight               | 22 lb               |