

PLATFORM WALKER	PLAT1000GB	PLAT1000G	PLAT1000MB	PLAT1000M
Colour	White			
Max. load	150 kg	150 kg	110 kg	110 kg
Length	80 cm			
Width	69 - 80 cm			
Width between forearm supports	31 - 48 cm			
Height of the arm supports (adjustable)	90 - 130 cm		90 - 140 cm	
Total height (at the top of push handles)	100 - 140 cm		100 - 150 cm	
Hand brakes	2 simultaneous brake levers with parking function	n/a	2 simultaneous brake levers with parking function	n/a
Front wheels	4" turning castors with foot operated parking brake and directional locking			
Rear wheels	4" turning castors with hand operated parking brake and foot directional locking			
Gas spring	150 N	150 N	n/a	n/a
Product weight	18 kg	17,2 kg	17,4 kg	16 kg
Box size 1 x Platform disassembled	85 x 77 x 23 cm (L x W x H)			
Box weight with 1 x Platform disassembled	21,5 kg	20,7 kg	20,9 kg	19,5 kg
Pallet size 10 x Platform disassembled	80 x 100 x 200 cm (L x B x H)			
Pallet size 2 x Platform assembled	80 x 100 x 105 cm (L x W x H)			
Pallet size 4 x Platform assembled	80 x 100 x 225 cm (L x B x H)			
Pallet weight 10 x Platform disassembled	235 kg			
Pallet weight 2 x Platform assembled	60 kg			
Pallet weight 4 x Platform assembled	105 kg			