



**PCP<sup>®</sup>**  
*Personal Care Products*

*Helping you maintain your  
freedom and independence!*



## COMMITTED TO BEING YOUR PARTNER FOR HOME MEDICAL EQUIPMENT AND ACCESSORIES

PCP is one of the most well-established and major categories under the Airway Surgical Appliances umbrella, offering a sensible mix of products for home convalescence, and to help restore individuals to an active, independent lifestyle.

And PCP product designs blend with home decor – very appealing for today's care-giver shopper.

### PCP PRODUCT SEGMENTS INCLUDE:

**PCP BATH SAFETY** – provides the assurance of a safe environment with an excellent offering of in-home use benches, grab bars, rails, single and 3-in-1 frames, elevated seats, and accessories. Many bath items are portable for added safety while traveling.

**PCP HOME CARE** – delivers support and convenience during home convalescence with a value-oriented, basic assortment of bed rails, accessories and commodes.

**PCP PERSONAL CARE** – offers therapeutic support surfaces and daily living supplies such as cushions, cast shoes, pillows, reachers, home traction, simple rehab equipment and hot/cold therapy.

**PCP MOBILITY** – offers the power of independence with a best-mix assortment of canes, crutches, walkers, rollators and accessories



VALUE

You don't have to pay more to get the very best! High quality and affordability are synonymous with PCP products – a real plus in today's competitive market.



CONSISTENCY

In design and quality – with few exceptions you buy direct from the original manufacturer, which ensures consistent product quality from order to order.



NO MINIMUMS

No minimum order restrictions!



WARRANTY

PCP Homecare products are warranted to be free of material or workmanship defects. Problematic merchandise will be repaired or replaced.



FULL LINE SUPPLIER  
1 STOP SHOP

Virtually every imaginable client need can be satisfied from the Airway Surgical program! Literally the best mix of motion control, support, safety convenience and mobility products, all from one source of supply!



DEDICATED SALES/  
SERVICE ASSOCIATES

Fully trained professionals will personally visit your facility on a regular basis to assist in managing various aspects of your business.

## PRACTITIONER AND PATIENT FOCUSED



The Airway Surgical Appliances program is built on the best mix of products that are designed and engineered for therapeutic performance and patient comfort – and specifically targeted to meet the needs of the everyday consumer and the occupational injury user.

Dynamic contemporary packaging educates and sells, creating a bold retail presence to help generate faster turns. Innovative packaging design features help create planograms that are easy to shop and easy to merchandise.



## PCP CONSIDERS YOUR MERCHANDISING NEEDS!

PCP packaging is designed to be displayed in both landscape and portrait orientations to better serve your merchandising needs.



### BATH SAFETY

|                                     |    |
|-------------------------------------|----|
| Bath & shower benches .....         | 4  |
| Bath & shower safety .....          | 7  |
| Safety grab bars .....              | 8  |
| Hand showers .....                  | 9  |
| Raised toilet seats .....           | 9  |
| Safety frames and accessories ..... | 12 |

### COMMODOES

|                           |    |
|---------------------------|----|
| Bedside commodes .....    | 12 |
| Commode accessories ..... | 13 |

### IN-HOME CARE

|                             |    |
|-----------------------------|----|
| Aids for daily living ..... | 14 |
| Personal care .....         | 16 |
| Bedroom convalescence ..... | 17 |

### CUSHIONS

|                                    |    |
|------------------------------------|----|
| Pillows, cushions and wedges ..... | 18 |
| Wheelchair cushions .....          | 23 |

### EXERCISERS

|                                       |    |
|---------------------------------------|----|
| Pedal and upper body exercisers ..... | 23 |
|---------------------------------------|----|

### CANES

|                                |    |
|--------------------------------|----|
| Adjustable canes .....         | 24 |
| Folding adjustable canes ..... | 26 |
| Wood and aluminum canes .....  | 26 |
| Specialty canes .....          | 27 |
| Cane accessories .....         | 27 |

### MOBILITY

|                                     |    |
|-------------------------------------|----|
| Crutches .....                      | 30 |
| Walkers .....                       | 30 |
| Walker and crutch accessories ..... | 31 |
| Rollators .....                     | 33 |
| Mobility Chairs .....               | 34 |

*Helping you maintain  
your freedom and  
independence!*

## BATH &amp; SHOWER BENCHES



### ADJUSTABLE BATH SAFETY STEP MODEL 7100

- Height can be adjusted for individual user
- Molded step with non-slip padded finish
- Non-marring, non-slip rubber suction tips
- Overall step height adjustment: 9" – 11" (23 – 28 cm)
- Step dimensions 14" wide, 11" deep (35.5 x 28 cm)
- Base dimensions 18" x 14" (46 x 35.5 cm)
- Weight capacity: 300 lbs. (136 kg.)



### PADDED ROUND SAFETY STOOL MODEL 7101

- Adjustable seat height for true custom fit
- Molded seat with non-slip padded finish
- Non-marring, non-slip rubber suction tips
- Overall seat height adjustment: 17.5" – 21.5" (44.5 – 54.6 cm)
- Seat dimensions: 13" (33 cm) round
- Base dimensions: 14" x 14" (36 x 36 cm)
- Weight capacity: 300 lbs. (136 kg.)



### PADDED BATH SAFETY SEAT MODEL 7102

- Adjustable seat height for true custom fit
- Notched to hold handshower
- Molded seat with non-slip padded finish
- Non-marring, non-slip rubber suction tips
- Tool free installation
- Seat height adjustment: 14" – 18" (35.5 – 45.7 cm)
- Seat dimensions: 16" wide, 11.5" deep (40.6 x 27.9 cm)
- Base dimensions: 16" x 16" (41 x 41 cm)
- Weight capacity: 300 lbs. (136 kg.)



### PADDED BATH SAFETY SEAT WITH BACKREST MODEL 7103

- Adjustable seat height for true custom fit
- Notched to hold handshower
- Molded seat with non-slip padded finish
- Non-marring, non-slip rubber suction tips
- Tool free installation
- Seat height adjustment: 14" – 18" (35.5 – 45.7 cm)
- Seat dimensions: 16" wide, 11.5" deep (40.5 x 28 cm)
- Base dimensions: 16" x 16" (40.5 x 40.5 cm)
- Back rest height from seat: 13" (33 cm)
- Weight capacity: 300 lbs. (136 kg.)

## PADDED BATH SAFETY SEAT WITH HYGIENIC CUTOUT

### MODEL 7104

- Adjustable seat height for true custom fit
- Cutout allows for easy cleaning of perineal area
- Molded seat with non-slip padded finish
- Non-marring, non-slip rubber suction tips
- Tool free installation
- Overall seat height adjustment: 17" – 22" (43 – 56 cm)
- Seat dimensions: 17" wide, 14" deep (43 x 35.5 cm)
- Base dimensions: 20" x 20" (51 x 51 cm)
- Weight capacity: 300 lbs. (136 kg.)



## PADDED BATH SAFETY SEAT WITH HYGIENIC CUTOUT AND BACK REST

### MODEL 7105

- Adjustable seat height for true custom fit
- Cutout allows for easy cleaning of perineal area
- Molded seat with non-slip padded finish
- Non-marring, non-slip rubber suction tips
- Tool free installation
- Overall seat height adjustment: 17" – 22" (43 – 56 cm)
- Seat dimensions: 17" wide, 14" deep (43 x 35.5 cm)
- Base dimensions: 20" x 20" (51 x 51 cm)
- Back rest height from seat: 13" (33 cm)
- Weight capacity: 300 lbs. (136 kg.)



## FOAM PADDED BATH SAFETY SEAT WITH SWING AWAY ARMS

### MODEL 7106

- For extra stability, easy lateral transfers
- Fits in most standard tubs or showers
- Adjustable seat height for a true custom fit
- Non-slip, foam padded seat with drain holes for extra safety
- Non-marring, non-slip rubber tips
- Seat height adjustment: 16" – 20" (41 – 51 cm)
- Height from seat to top of armrests: 5" (12.7 cm)
- Seat dimension: 16" x 13" (40.6 x 33 cm)
- Base dimensions: 19" x 16.25" (48 x 41.3 cm)
- Total Weight: 2.27 kg (5 lb.)
- Weight capacity: 136 kg (300 lb)



## FOAM PADDED BATH SAFETY SEAT WITH BACK AND SWING AWAY ARMS

### MODEL 7107

- For extra stability, easy lateral transfers
- Fits in most standard tubs or showers
- Adjustable seat height for a true custom fit
- Non-slip, foam padded seat and back rest
- Drain holes on seat for extra safety
- Non-marring, non-slip rubber tips
- Seat height adjustment: 16" – 20" (41 – 51 cm)
- Height from seat to top of armrests: 5" (12.7 cm)
- Height from top of backrest to seat: 12" (30.5 cm)
- Seat dimension: 16" x 13" (40.6 x 33 cm)
- Backrest dimensions: 16.5" x 7" (41.9 x 17.8 cm)
- Base dimensions: 19" x 16.25" (48 x 41.3 cm)
- Total weight: 3.6 kg (8 lb)
- Maximum weight capacity: 136 kg (300 lb)



# BATH & SHOWER BENCHES



## BATH SAFETY SEAT WITH BACKREST

MODEL 7003

- Adjustable seat height for a true custom fit
- Durable molded seat with non-slip finish
- Drain holes on seat for extra safety
- Rigid frame with sturdy, flared legs
- Non-marring, non-slip rubber tips
- Overall seat height adjustment (in 1" increments): 14" - 22" (35.5 - 56 cm)
- Seat dimensions: 19.25" wide, 11.5" deep (48.9 x 29.2 cm)
- Base dimensions: 16" x 17" (40.6 x 43.1 cm)
- Back height from seat: 15 inches (38 cm)
- Weight capacity: 300 lbs. (136 kg)



## FOLDING TRANSFER BENCH

MODEL 7392

- Allows for an easy transfer over the tub wall for bathing
- Folds compactly for easy storage
- Durable molded non-slip seating area, easy to clean
- Drain holes on seat for extra safety
- Side hand grip for extra support
- Non-marring, non-skid rubber tips
- Handle on backrest for easy lifting
- Backrest is easily reversible to accommodate right or left facing tub
- Overall seat height adjustment (in .75" increments): 16" - 21" (40.5 - 53 cm)
- Seating area: 26" wide, 16" deep (66 x 40.5 cm)
- Seat dimensions: 29.5" x 20" (75 x 51 cm)
- Back height from seat: 13" (33 cm)
- Weight capacity: 300 lbs. (136 kg)



## ROUND SHOWER STOOL

MODEL 7001

- Adjustable seat height for a true custom fit
- Neutral colors blend perfectly with any decor
- Easy to assemble using only common household tools
- Overall seat height adjustment (in 1" increments): 14" - 22" (35.5 - 56 cm)
- Seat diameter: 12.5" (31.75 cm)
- Width between adjustable legs: 14" (35.6 cm)
- Weight capacity: 300 lbs. (136 kg)



## BATH SAFETY SEAT

MODEL 7002

- Adjustable seat height for a true custom fit
- Durable molded seat with non-slip finish
- Drain holes on seat for extra safety
- Rigid frame with sturdy, flared legs
- Non-marring, non-slip rubber tips
- Overall seat height adjustment (in 1" increments): 14" - 22" (35.5 - 56 cm)
- Seat dimensions: 19.25" wide, 11.5" deep (48.9 x 29.2 cm)
- Base dimensions: 16" x 17" (40.6 x 43.1 cm)
- Weight capacity: 300 lbs. (136 kg)

## SUCTION TIPS

MODEL: 6118 (PICTURED TO THE RIGHT)

- Suction tips for the 7100 series of bath safety seats.
- Sold in sets of four
- Suitable fit for numbers 7100, 7101, 7103, 7103, 7104 and 7105

MODEL: 7002TP

- Suction tips for the 7000 series of bath safety seats.
- Grey in colour
- Sold individually
- Suitable fit for numbers 7001, 7002, and 7003



## SURE-GRIP BATHTUB SAFETY RAIL

MODEL 7252

- Easy tool-free installation
- Fits a wide variety of tubs
- One piece steel platform with stainless steel locking mechanism
- Strong welded construction
- Attractive, easy to clean vinyl coating
- Neutral white blends perfectly with any bathroom decor
- Non-marring rubber pads protect tub wall
- Overall height: 19" (48 cm)
- Width: 8" (20 cm)
- Locking bracket adjustment: 3" - 6" (7.6 - 15.2 cm)
- Weight capacity: 300 lbs. (136 kg.)



## ECONOMY SAFETY RAIL

MODEL 7255


- Attractive, easy to clean vinyl coating.
- Neutral white blends perfectly with any bathroom decor.
- Non-marring rubber pads protect tub wall.
- Overall height: 20 inches (50.8 cm)
- Height above tub wall: 15 inches (38 cm)
- Width: 6 inches (15.2 cm)
- Locking bracket adjustment: 3 - 6 inches (7.6 - 15.2 cm)
- Product Weight: 6 lb, 10 oz
- Weight Capacity: 300 lbs. (136 kg.)



## SHOWER SAFETY MAT

MODEL 7039

- Reliable slip resistance with cushioned, comfortable feel underfoot
- Made of high-quality, allergen-free vinyl
- Holes in center over shower drain
- Mildew resistant
- Stays soft and pliable
- 20.5" x 21.5" (52 x 54.5 cm)



## BATH SAFETY MAT

MODEL 7036

- Full size for safety
- Over 200 extra large size suction cups prevent slippage
- Mildew resistant
- Reliable slip resistance with a cushioned, comfortable feel underfoot
- Made of high-quality, allergen-free vinyl
- Stays soft and pliable.
- Dimensions: 15.5" x 39.5" (39.4 x 100.3 cm)



## SAFETY GRAB BARS



## SUCTION BALANCE HAND GRIP WITH RED/GREEN SAFETY INDICATORS

MODEL 9211

- Suction cups are easy to attach, remove or relocate in seconds – no tools required
- Can be installed at any angle
- Ribbed grip area for extra help when grasping
- Works with any smooth, flat and non-porous surface
- Avoid grout lines for proper placement



## 9200 SERIES SUCTION GRAB BARS WITH RED/GREEN SAFETY INDICATORS

MODELS:

9212 – 12" (30.5 cm) | 9216 – 16" (40.6 cm) | 9226 – 180° SWIVEL

- Allows user to enter a wet shower or tub with safety and confidence
- Helps prevent dangerous slips or falls; great for travel
- Suction cups are easy to attach, remove or relocate within seconds
- Can be installed at any angle
- Ribbed grip area for extra help when grasping
- Works with any smooth, flat and non-porous surface
- Red/green safety indicators allows user to know when safely installed
- This is not intended for full body weight leverage

## 9300 SERIES BALANCE ASSIST BAR

MODELS:

9312 – 12" (30.5 cm) | 9316 – 16" (40.6 cm) | 9324 – 24" (45.7 cm)

- Allows user to enter a wet shower or tub with safety and confidence
  - Helps prevent dangerous slips or falls; great for travel
  - Suction cups are easy to attach, remove or relocate within seconds
  - Can be installed at any angle
  - Works with any smooth, flat and non-porous surface
- \* This product is not intended for full body weight leverage.

## 8400 SERIES GRAB BARS

MODELS:

8412 – 12" (30.5 cm) | 8416 – 16" (40.6 cm) | 8418 – 18" (45.7 cm)

- Flange is fixed to the bar
- Fully integrated steel brackets
- Rust resistant, thick wall chrome steel tubing
- Diamond cut finish for safe and secure gripping
- Perfectly polished finish to accent bathroom decor
- Extra long plated mounting screws enclosed
- Illustrated installation instructions enclosed
- Easy to clean – just wipe with a damp cloth
- Mounts anywhere – bathtub, shower, commode area


## 9100 SERIES GRAB BARS

MODELS:

9112K – 12" (30.5 cm) | 9116K – 16" (40.6 cm) | 9118K – 18" (45.7 cm)

9124K – 24" (61 cm) | 9132K – 32" (81.3 cm)

- Rotating flanges
- Rust resistant, thick wall chrome steel tubing
- Diamond cut finish for safe and secure gripping
- Extra long, plated mounting screws enclosed
- Illustrated installation instructions enclosed
- Easy to clean – just wipe with a damp cloth
- Mounts anywhere – bathtub, shower, commode area



## 5 SETTING HAND HELD SHOWER

MODEL 7041

- Provides a relaxing spa like setting in the palm of the hand
- Five different settings range from a soft relaxing spray to a hard massage
- Easy-turn spray dial
- 6 ft. (72" / 183 cm) tangle free flexible hose
- Attractive chrome finish
- Easy to install; requires no tools
- Hand shower holder and diverter valve both sold separately



## SUCTION HAND SHOWER HOLDER

MODEL 7047

- Extra strong suction cup sticks to walls, can be removed without damage to surfaces
- Shower head can be angled and rotated to desired position
- Easy to install and remove; requires no tools



## DIVERTER VALVE

MODEL 7043

- Attaches quickly and easily to shower arm; a turn of the valve allows use of either regular shower head or hand held shower
- Compatible with PCP model 7041 hand-held shower



# RAISED TOILET SEATS

## CONTOURED MOLDED RAISED TOILET SEAT

MODEL 7013 **FITS: ROUND BOWL**

- Lightweight and easily portable – recessed handles on sides
- Rotational molded polyethelene for comfort and cleanliness
- No drill holes where dirt and bacteria can become trapped
- Includes discreet carrying bag for use when travelling
- Seat height: 5" (12.7 cm) in front, 6" (15 cm) in back
- Seat width: 15" (38 cm)
- Seat depth: 15.5" (39.5 cm)
- Weight capacity: 300 lbs. (136 kg.)



## MOLDED RAISED TOILET SEAT WITH TIGHTENING LOCK

MODEL 7015 **FITS: ROUND BOWL**

- Adds a full 5" (12.7 cm) of seat height
- Contour molded design for comfort and cleanliness
- Twist and lock dial safely secures seat to most rims
- Front indentation allows for maximum hygiene
- Includes channels to accommodate model 7021 detachable arms
- Seat width: 15.5" (39.5 cm)
- Seat depth: 16.5" (42 cm)
- Weight capacity: 300 lbs. (136 kg.)



# RAISED TOILET SEATS



## MOLDED RAISED TOILET SEAT WITH REMOVABLE ARMS

MODEL 7016 **FITS: ROUND BOWL**

- Adds a full 5" (12.7 cm) of seat height
- Contour molded design for comfort and cleanliness
- Padded armrests for extra support when sitting, extra safety when standing
- Twist and lock dial safely secures seat to most rims
- Front indentation allows for maximum hygiene
- Seat width: 15.5" (39.5 cm)
- Seat depth: 16.5" (42 cm)
- Weight capacity: 300 lbs. (136 kg.)



## MOLDED RAISED TOILET SEAT WITH ARM SETS

MODEL 7017 **FITS: ROUND BOWL**

- Adds a full 5" (12.7 cm) of seat height
- Contour molded design for comfort and cleanliness
- Padded armrests for extra support when sitting, extra safety when standing
- Twist and lock dial safely secures seat to most rims
- Front indentation allows for maximum hygiene
- Seat width: 15.5" (39.5 cm)
- Seat depth: 16.5" (42 cm)
- Weight capacity: 300 lbs. (136 kg.)



## 4" TOILET SEAT RISER WITH REMOVABLE ARMS

MODEL 7027 **FITS: ROUND BOWL**

- Adds 4" (10 cm) to standard seat height
- Allows for use of existing toilet seat and lid
- Heavy duty molded plastic
- Durable lightweight metal arms and receiver channels
- Arms are easily attached and removed if desired
- Soft foam padding on arm rests
- Attaches easily to toilet bowl with provided hardware; no tools required
- Seat depth: 17" (43.18 cm)
- Seat width: 13.5" (34.29 cm)
- Distance between armrests: 16" (40.64 cm).
- Total weight: 1.59 kg (3.5 lb)
- Weight capacity: 158.7 kg. (350 lb.)



Recommended for use with PCP model 7021 Detachable Arms



## 3" UNIVERSAL RAISED TOILET SEAT

MODEL 7019 **FITS: ROUND & ELONGATED BOWL**

- Fits both standard and elongated toilet bowls
- Low profile: adds 3" (7.5 cm) of seat height
- Adjustable twist and lock clamps safely secure seat to most rims
- Twist and lock dial safely secures seat to most rims
- Front indentation allows for maximum hygiene
- Includes channels to accommodate model 7021 detachable arms
- Seat width: 15" (38 cm).
- Seat depth: 16.5" (42 cm).
- Weight capacity: 300 lb. (136 kg.)

## 5" UNIVERSAL RAISED TOILET SEAT

MODEL 7020 **FITS: ROUND & ELONGATED BOWL**

- Fits both standard and elongated toilet bowls
- Adds a full 5" (12.7 cm) of seat height
- Adjustable twist and lock clamps safely secure seat to most rims
- Twist and lock dial safely secures seat to most rims
- Front indentation allows for maximum hygiene
- Includes channels to accommodate model 7021 detachable arms
- Seat width: 15" (38 cm)
- Seat depth: 16.5" (42 cm)
- Weight capacity: 300 lbs. (136 kg.)



Recommended for use with PCP model 7021 Detachable Arms

## 2" MOLDED RAISED TOILET SEAT WITH LID

MODEL 7022 **FITS: ELONGATED BOWL**

- Low profile – adds 2" (5 cm) of seat height
- Convenient front and rear hygiene cutouts
- Designed for use with existing toilet seat
- Twist and lock dial safely secures seat to most rims
- Front indentation allows for maximum hygiene
- Seat width: 14" (35.5 cm)
- Seat depth: 16" (40.5 cm)
- Weight capacity: 300 lbs. (136 kg.)



## 4" MOLDED RAISED TOILET SEAT WITH LID

MODEL 7024 **FITS: ELONGATED BOWL**

- Adds a full 4" (10 cm) of seat height
- Convenient front and rear hygiene cutouts
- Designed for use with existing toilet seat
- Twist and lock dial safely secures seat to most rims
- Front indentation allows for maximum hygiene
- Seat width: 14" (35.5 cm)
- Seat depth: 16" (40.5 cm)
- Weight capacity: 300 lbs. (136 kg.)



## RAISED TOILET SEAT WITH SAFETY FRAME

MODEL 7007 **FITS: ROUND BOWL**

- Sturdy 1" (2.5 cm) frame.
- Molded plastic seat and encircling splash guard fixed to frame; easy to keep clean.
- Vinyl hand grips / arm rests and silencer rings.
- Seat height adjusts 18"-23" (46-58 cm).
- Overall width 21" (53 cm).
- Width between arms 18" (46 cm).
- Weight capacity: 300 lb. (136 kg.)



# SAFETY FRAMES & ACCESSORIES



## TOILET SEAT CUSHION

MODEL 7018

- Provides extra comfort when sitting
- Helps minimize pressure sores
- Comfortable foam padding, 2" (5 cm) thick, covered with naugahyde
- Hook-and-loop fastening straps allow for easy attachment and removal
- Fits all standard toilet seats
- Easy to clean

## TOILET SAFETY FRAME

MODEL 7009

- Provides support and balance when using toilet
- Full length waterfall armrests for comfort and support
- Sturdy, lightweight aluminum frame
- Quick and easy to assemble, requiring only common household tools
- Overall rail height adjustment (in 1" increments): 27" – 31" (68.5 – 79 cm)
- Width between armrests (in 1" increments): 17" – 22" (43 – 56 cm)
- Weight capacity: 300 lbs. (136 kg.)

## DETACHABLE ARMS FOR MOLDED RAISED TOILET SEAT

MODEL 7021

- Can easily be attached to either 7015, 7019 or 7020 raised toilet seats
- Provide extra support when sitting, extra safety when arising from seat
- Easy installation – Phillips head screwdriver is all that is required for six screws

## TOILET SEAT SAFETY RAIL

MODEL 7026

- Specially designed for use in tight bathrooms
- Durable coated steel construction with plastic armrests
- Anti-slip grip design
- Mounting bracket adjusts from 5.5" – 10.2" (14 – 26 cm) wide
- Easy to clean
- Weight capacity: 250 lb (113 kg)

# BEDSIDE COMMODES

## FOLDING STEEL COMMODE

MODEL 5523

- Opens easily and folds flat for storage or transport
- Strong steel welded construction with gray powder coat finish
- Seat assembly easily removable for cleaning
- 6 quart pail with lid and splash guard
- Width between arms: 18" (46 cm)
- Overall width: 23" (58.5 cm); overall depth 20" (51 cm)
- Seat height adjustment: 16" – 22" (40.5 – 56 cm)
- Weight capacity: 300 lbs. (136 kg.)

## HEAVY DUTY EXTRA WIDE COMMODE

MODEL 5029

- Strong welded steel frame, 7" (18 cm) wider than standard
- Molded plastic armrests for comfort
- Seat assembly easily removable for cleaning
- Includes removable 6 qt. pail with lid
- Seat height adjustment: 19" – 22" (48 – 56 cm)
- Overall width 28" (71 cm); overall depth 21" (53 cm)
- Width between arms 24" (61 cm)
- Weight capacity: 450 lbs. (204 kg.)

## COATED STEEL COMMUNE

MODEL 5526-RBE

- Performs 3 functions: free-standing commode and with backrest removed as raised toilet seat and safety frame (when placed over bathroom commode)
- Wider footprint for more stability
- Bolted gray coated steel frame
- Molded plastic armrests for added comfort
- Easily assembled using commode household tools
- Seat assembly easily removable for cleaning
- Includes removable 6 qt. pail with lid
- Replacement pails available
- Seat height adjustment: 16" - 22" (40.5 - 56 cm)
- Overall width 26" (66 cm); overall depth 20" (51 cm)
- Width between arms 18" (46 cm)
- Weight capacity: 300 lbs. (136 kg.)



## 5-IN-1 MOBILITY AND BATHROOM AID

MODEL 5024

- 5 modes: commode, folding walker, toilet seat frame, raised toilet seat, shower seat
- Lightweight aluminum frame
- Folds for easy transport and storage
- Non-slip rubber tipped legs
- Tool free assembly
- Neutral colors blend perfectly with any decor
- Seat height adjustment: 16" - 23" (41-58 cm.)
- Overall width 20.5" (48.3 cm); overall depth 19" (52.1 cm)
- Seat hole diameter: 10" (25 cm.)
- Width between adjustable legs: 19" (48 cm.)
- Product weight: 10.5lbs (4.8kg.)
- Weight capacity: 300lbs (136kg.)



## COMMUNE ACCESSORIES

### REPLACEMENT HALF PAIL WITH LID

MODELS:

5114 - WHITE

5114-GR - GREY

- 6 quart capacity
- Fits all PCP bedside commodes perfectly.



### REPLACEMENT FULL PAIL WITH LID AND HANDLE

MODELS:

5115 - WHITE

- 10 quart capacity
- Fits all PCP bedside commodes perfectly



### PADDED COMMUNE SEAT

MODEL 5024-ST

- Use this product to add extra comfort to the 5-in-1 Mobility and Bathroom Aid!





## EASY RISER SUPPORT RAIL

MODEL 6086

- Makes it easy to get in or out of a bed or seat
- Multiple grips allow use of both hands for leverage and optimum support
- Strong coated steel frame
- Quick and easy assembly with provided tool
- Wide non-skid, non-marring base of support
- Weighs less than 4 lb (1.8 kg)
- Base: 15" long x 20" wide (38.1 x 50.8 cm)
- Hand Grip: 32" high (81.28 cm)
- Weight Capacity: 250 lb (113.4 kg)



## FOOT STOOL

MODEL 5008

- Sturdy and tip proof
- Rugged non-skid 10" x 14" (25.5 x 35.5 cm) top
- Chrome steel legs, metal reinforced rubber tips
- Height 9" (23 cm).
- Weight capacity: 300 lbs.



## FOOT STOOL WITH HANDLE

MODEL 5008-H

- Sturdy and tip proof.
- Chrome steel legs, metal reinforced rubber tips.
- Rugged non-skid 11" x 14" top.
- Height 9" (23 cm)
- Handle height: 34.25" (87 cm) – does not adjust
- Weight capacity: 300 lbs.
- Tool free assembly



## INFLATABLE HAIR WASH BASIN

MODEL 6096

- Designed for comfortable and convenient in-bed shampooing
- Heavy-duty vinyl with drain stopper at the base of the flexible drain hose
- Lightweight, portable and compact for convenient storage
- Dimensions: 20" (51 cm) wide x 20" (51 cm) long x 2.5" (6 cm) deep
- Drain hose efficiently removes soap suds and water from the hair wash basin



## MOTION MEDICINE



MODELS:

6502 - 120G TUBE (SOLD IN DISPLAY PACK OF 12)

6500 - 500G JAR

- Reduces pain, bruising & swelling
- Fast acting and heals tissue
- Natural product made in Canada



## CONVALESCENT HOSPITAL GOWN

MODELS:

68033B - PRINT

68030B - BLUE

- Latex free
- Made of comfortable, durable polyester/cotton blend
- One size fits most

## METAL SHOE HORN

MODEL 513000

- Made from durable metal material
- Extra long 24" (61 cm) shaft for less strain on the lower back

## WOOD SHOE HORN

MODEL 514127

- Made from durable wooden material
- Long 21" (53 cm) shaft for easy reach and less strain on the lower back

## REACH EXTENDER

MODEL 6031

- Permits users with limited strength to get at hard-to-reach objects  
Extends average reach by 30 inches (76 cm)
- Non-slip prongs
- Easy to grasp handle

## REACH EXTENDER WITH SUCTION CUPS

MODEL 6033

- Rounded clamps make it easy to grasp objects like cans, balls, bottles and cups without losing hold
- Soft rubber grip tips help pick up small items
- Lightweight rigid aluminum shaft extends average reach by 30 inches (76 cm)
- Pistol grip handle with trigger control
- Optional suction cups attach to tips for enhanced control

## TRANSFER BOARD

MODEL 6032

- Eases the transfer from a wheelchair to bed or from the bed to a bedside commode
- Molded ribs add to the strength, minimize friction to facilitate easy transfer
- Smooth plastic surface is easy to maintain and clean
- Dimension: 8" x 24" (20 x 61 cm)
- Weight capacity: 300 lbs. (136 kg.)

## COMPACT SNOW/ICE GRABBERS

MODEL 109550

- Material is flexible to -30°
- Special non-slip studs offer excellent traction on ice and snow
- Grabbers will fit all types of footwear
- 1 size fits most
- Sold in pairs


## NON-SLIP SNOW/ICE GRABBERS

MODEL 10957

- Material is flexible to -40°
- Special non-slip studs offer excellent traction on ice and snow
- Sold in pairs

### MEASUREMENT BASED ON SHOE SIZE

|          |    | LADIES'    | MEN'S   |
|----------|----|------------|---------|
| 10957-M  | M  | 6.5 - 9.5  | 5 - 8   |
| 10957-L  | L  | 9.5 - 12.5 | 8 - 11  |
| 10957-XL | XL | 13 +       | 11 - 14 |





## FEMALE URINAL

MODEL 6088

- 28 oz. (0.79kg.) capacity
- Shaped to prevent accidental spillage
- Sure-grip handle hangs on bed rail for easy access
- Angled for use in several positions
- Easy to read calibrations for quick measurement



## MALE URINAL

MODELS:

6091 - REGULAR

6092 - AUTOCLAVABLE (NOT AVAILABLE IN RETAIL PACKAGING)

- 1 quart capacity
- Complete with hinged cover and sure-grip handle that hangs on bed rail
- Cap opens and closes easily; confines odors
- Shaped to prevent accidental spillage
- Easy to read calibrations for quick measurement



## BED PAN

MODEL 6094

- Hospital quality, generous capacity
- Contoured for comfort
- Lightweight, easy to clean
- Dimensions: 15" x 12" (38 x 30.5 cm)



## FRACTURE BED PAN

MODEL 6095

- For someone recovering from a hip fracture or replacement surgery, or for those who cannot raise their hips high enough for a regular bed pan
- Easy grip handle at the back
- Tapered front can be positioned without having to move patient
- Lightweight molded plastic, easy to clean



## SITZ BATH

MODEL 5040

- Used for treatment of hemorrhoids and other causes of rectal discomfort
- Fits and works with any standard commode
- Easy to follow instructions on solution bag
- Comfortably contoured bowl cleans easily and stays looking bright
- Designed to provide treatment around wide area surrounding the entire perineal region

## PRESSURE RELIEF MEDICAL OVERLAY

MODEL 350110

- Used to prevent decubitus ulcers (bed sores) resulting from prolonged periods of bed rest
- Made of genuine sheepskin
- Redistributes body weight to prevent decubitus ulcers
- Wicks moisture away from skin, keeping it cool and dry
- Maintains user's normal skin temperature
- Stain resistant, non-flammable, non-allergenic, non-toxic
- Remains cool in summer, warm in winter
- Machine washable and dryable
- Dimensions: 7-8 sq. ft.



## PRESSURE RELIEF PROTECTORS

MODELS:

350010 HEEL PROTECTORS

350011 ELBOW PROTECTORS

- Used to prevent decubitus ulcers (bed sores) resulting from prolonged periods of bed rest
- Made of genuine sheepskin
- Redistributes body weight to prevent decubitus ulcers
- Wicks moisture away from skin, keeping it cool and dry
- Maintains user's normal skin temperature
- Stain resistant, non-flammable, non-allergenic, non-toxic
- Remains cool in summer, warm in winter
- Machine washable and dryable
- One size fits most
- Sold in pairs



## PRESSURE RELIEF PROTECTORS

MODELS:

6206 HEEL PROTECTORS

6207 ELBOW PROTECTORS

- Redistributes body weight to prevent decubitus ulcers
- Wicks moisture away from skin, keeping it cool and dry
- Stain resistant, non-flammable, non-allergenic, non-toxic
- Made from soft polyester KODEL
- One size fits most
- Sold in pairs



## XSPAN® HEEL/ELBOW PROTECTOR PAD

MODELS:

832000 - SMALL/MEDIUM

833000 - LARGE

- Made from a comfortable knitted material
- Easy one handed application
- Features a protective foam insert
- Fits both heel and elbow
- Sold in pairs



# BEDROOM CONVALESCENCE



## SYNTHETIC SHEEPSKIN



### MODELS:

- 6261 – BED PAD, 24" X 30" (61 x 76 cm)
- 6262 – BED PAD, 30" X 40" (76 x 101.5 cm)
- 6263 – BED PAD, 30" X 60" (76 x 152.5 cm)
- 6264 – BED PAD, 60" X 80" (152.5 x 203 cm)
- 6266 – WHEELCHAIR PAD, 18" X 18" (45.7 x 45.7 cm)
- 646400 – BED PAD, 30" X 80" (76 X 203 CM)

- Provides many of the material benefits of genuine sheepskin but at a fraction of the cost
- Non flammable, non-allergenic, non-toxic – 94% polyester, 6% acrylic
- Can be cut to shape without unraveling
- Will not pack down
- Stays soft and thick through repeated use
- Cool in summer, warm in winter
- Machine wash and dry
- Stain resistant



## KODEL BED PAD

### MODELS:

- 6203 – BED PAD, 24" X 30" (61 x 76 cm)
- 6204 – BED PAD, 30" X 40" (76 x 101.5 cm)
- 6205 – BED PAD, 30" X 60" (76 x 152.5 cm)

- Stays soft and thick through repeated use
- Cool in summer, warm in winter
- Made from soft polyester KODEL
- Machine wash and dry
- Stain resistant

# PILLOWS, CUSHIONS AND WEDGES



## INFLATABLE RING CUSHION

### MODEL 6230

- Relieves pressure and fatigue troubling anyone who is required by illness, injury or occupation to lie or sit in the same position over prolonged periods
- Can be used in leg elevation and as a backrest (against a straight back chair)
- In bed, use cushion under buttocks or knees to relieve pressure and tenderness resulting from friction with bed clothes



## ROUND CERVICAL PILLOW

### MODEL 6145

- Provides correct spinal alignment and a therapeutic sleep posture
- Can also be used under the lower back, knees or ankle
- Can be used at home or when traveling
- Hypo-allergenic
- Dimensions: 18" (46 cm) long, 6" (15 cm) diameter

## MEMORY FOAM CERVICAL PILLOW

MODELS: 6152-MF FULL SIZE | 6154-MF TRAVEL SIZE

- Ergonomically designed to support the head and neck correctly and to relieve pressure on neck muscles
- Helps prevent stiffness, alleviate soreness
- Made of viscoelastic memory foam, super soft and comfortable
- Comes with a removable, washable zippered ecru terry cloth cover



## FOAM RING CUSHION

MODELS:

|                              |                               |
|------------------------------|-------------------------------|
| 6232-NV - 14" (35.6 cm) NAVY | 6232-PL - 14" (35.6 cm) PLAID |
| 6233-NV - 16" (40.6 cm) NAVY | 6233-PL - 16" (40.6 cm) PLAID |
| 6234-NV - 18" (45.7 cm) NAVY | 6234-PL - 18" (45.7 cm) PLAID |

- One piece puncture-proof molded foam conforms to the natural contours of the body and provides maximum comfort.
- Comes with a washable cover for easy maintenance.



## CONVOLUTED FOAM RING CUSHION



MODELS:

6235-CU - WITH REMOVABLE WHITE COVER

6235-NV - WITH REMOVABLE NAVY COVER

- Convoluted foam disperses surface pressure for greater sitting comfort, promotes good posture and fights pressure sores
- Open center helps ease pressure by distributing body weight evenly, reducing pressure point discomfort



## REVERSIBLE TRAVEL CUSHION

MODEL 6236-BL

- Converts from lumbar roll to neck cushion with a simple pull of a zipper
- Soft microbead fill



## TRAVEL RING CUSHION

MODEL 6237-BL

- Folds compactly for travel or storage
- Reduces pressure on lower back and tailbone; comfortably conforms to body contours
- Ideal for use after childbirth, for hemorrhoid sufferers or for those with low back pain
- Removable and machine washable cover



## MOLDED RING CUSHION

MODEL 6238-BL

- Maximum relief and comfort – reduces strain on pressure points
- Non-slip bottom keeps cushion in correct position – will not shift or slide
- Made of 3" (7.6 cm) thick premium comfort foam – maintains its shape through repeated use
- Removable and machine washable cover



## ORTHOPEDIC SEAT CUSHION WITH REMOVABLE PAD

MODEL 6239-BL

- Superior comfort for any chair, bench, car seat or wheelchair
- Coccyx pad detaches easily from unit for extra relief over pressure points
- 3" (7.6 cm) thick sculpted memory foam – will not flatten out with prolonged use
- Removable and machine washable cover





### MEMORY FOAM KNEE PILLOW

MODEL 6150-MF

- Promotes proper alignment the of spine when lying on the side
- Helps relieve discomfort and stiffness, aids in circulation and helps ease pain in the back, legs, hips and knees
- Made of viscoelastic memory foam, super soft and comfortable
- Adjustable hook and loop strap offers a custom fit and keeps pillow in place
- Removable soft cover zips off easily for machine washing



### CARVED FOAM CERVICAL PILLOW

MODEL: 6155 FULL SIZE



- Relieves neck pain for a better sleep
- Helps prevent stiffness, alleviate soreness
- Made of firm density foam for optimum support
- Comes with a removable, washable fabric cover



### FACECRADLE TRAVEL PILLOW

MODELS:

6229BL - BLACK 6229GR - GREY

- 5 modes of comfort
- Designed for the total environment that is economy class travel
- Utilize the complete 180 degrees of space in and in front of your seat



### INFLATABLE NECK CUSHION

MODEL 6231

- Comfortably supports the head and neck while reading, napping or relaxing
- Helps prevent pain and stiffness in the neck and shoulders
- Perfect for travel: folds easily and comes in a handy pouch
- Quick and easy to inflate and deflate



### MEMORY FOAM NECK CUSHION

MODEL 6151-MF | 6151-MFBL

- Helps prevent pain and stiffness in the neck and shoulders
- Helps maintain proper head and neck alignment
- Great for travel, including airplanes and long bus and car trips
- Comes with removable cover

## FOLDABLE BED WEDGES

MODELS: 6167 – 8" (20.3 cm) | 6168 – 10" (25.4 cm) | 6169 – 12" (30.5 cm)

- Provides gradually sloping support for head, foot or leg elevation
- Helps alleviate respiratory symptoms and improves sleep
- Firm density foam
- Soft terry cloth cover – removable and washable
- Completely non-allergenic
- Folding design allows for easy storage
- Latex free



## HIP ABDUCTION BLOCK WITH COVER

MODEL 614900

- Provides relief from discomfort and allows for a restful sleep
- Helps relieve discomfort and stiffness, aids in circulation and helps ease pain in the back, legs, hips and knees
- Specially designed foam block provides for proper alignment of the spine, hips and knees
- Adjustable encircling strap with hook-and-loop closure keeps block in proper position
- Removable and washable poly/cotton cover



## SACRO CUSHION WITH REMOVABLE COVER

MODELS:

6243 – BLACK    6246 – NAVY  
6244 – GREY    6247 – BURGUNDY

- Economical, comfortable and durable
- Molded polyurethane foam in a removable and washable fabric cover
- Holds lower spine in an upright, natural position for comfort during prolonged sitting or while driving
- Back strap holds cushion securely in position
- Dimensions: 12.5" x 14" (32 x 35.5 cm)



## LUMBAIRPLUS BACK SUPPORT

MODELS:

620000 – NAVY    626001 – BLACK

- Air adjustable portable Lumbar spine support
- Prevents back pain and stress
- Easy to use at home, car or travel
- Dimensions: 10.5" x 7" (26.7 cm x 17.8 cm)



## LUMBAR BACK PAD

MODELS:

625700 – NAVY    625701 – BLACK

- Extra support for low back
- Quality fabric and foam
- Handy and portable
- Dimensions: 11" x 6" (28 cm x 15.25 cm)





## PROPAIR SLEEPER CUSHION



MODEL 640400

- Air adjustable support
- Reduce pain in back and neck
- Easy to use, comfortable and adjustable
- Contains quality components like real satin
- One size fits most
- Improves sleep quality and prevents cumulative back stress
- Dimensions: 12" x 6.5" (30.5 - 16.5 cm)



## EMBRACEAIR PLUS BACK CUSHION



MODEL 6400

COLOURS: BLACK AND NAVY

- Air adjustable support
- Back cushion aligns with natural spine curvature
- Easy to adjust softness/firmness with push button valve
- Easy to fit different body heights
- Proven to prevent back pain and improve spine posture
- Dimensions: 19.5" x 20.5" (49.6 - 52.2 cm)



## EMBRACE IN-LINE SEATING



MODEL 6401

COLOURS: BLACK AND NAVY

- Stabilize body position with crafted, molded plywood
- Offers 1" (2.5 cm) Ultra cell super resilient foam and treated fabric
- Supports the posterior and evenly distributes body weight
- Complements all Embrace Products
- Dimensions: 19" x 16.5" (48.3 x 42 cm)



## EMBRACEAIR KING BACK SUPPORT



MODEL 6402

COLOURS: BLACK AND NAVY

- Air adjustable support
- Adjust softness/firmness with push button valve
- For comfort and prevention of back pain
- Vertical adjustment straps
- Dimensions: 21.5" x 24.5" (54.7 x 62.25 cm)



## EMBRACEPLUS BACK CUSHION



MODEL 6403

COLOURS: BLACK AND NAVY


- Built in additional lumbar support is movable up and down
- Relieves strain on the upper and lower back
- High quality fabric and foam built to last
- Recycled molded wood frame
- Dimensions: 19.5" x 20.5" (49.6 x 52.2 cm)

## CONVOLUTED FOAM WHEELCHAIR CUSHION

MODEL 6221



- Convoluted foam disperses surface pressure for greater sitting comfort
- Promotes good posture, fights pressure sores
- Dimensions: 16" x 18" (40.5 x 45.7 cm)



## KODEL WHEELCHAIR CUSHION

MODEL 6208




- Constructed from soft Polyester
- Dimensions: 18" x 18" (45.7 x 45.7 cm)



## GEL FOAM WHEELCHAIR CUSHION

MODEL 6226

- Combines the benefits and sitting comfort of both gel and foam cushions
- Wipes clean easily – fabric cover removable for washing
- Dimensions: 16" x 16" (40.5 x 40.5 cm)
- Anti-microbial



## 2" FOAM WHEELCHAIR CUSHION

MODEL 6227

- Designed to minimize the risk of pressure sores
- High density foam, 2" (5 cm) thick – compresses to 1" (2.5 cm)
- Dimensions: 16" x 18" (40.5 x 46 cm) – fits all standard adult wheelchairs
- Removable and washable cloth cover



# PEDAL AND UPPER BODY EXERCISER

## PEDAL EXERCISER

MODEL 6298

- For upper or lower body strengthening, improving coordination and circulation
- Adjustable tension control allows a broad range of pedal resistance
- Chrome steel frame, non-skid tips
- Straps keep feet in place while pedaling
- Compact, lightweight and very portable
- Four anti-slip rubber pads prevent sliding and protect surfaces
- Stores out of the way under a chair or bed
- Quick and easy assembly



## ARM AND UPPER BODY EXERCISER

MODEL 6299

- Used as part of a rehab program to help restore shoulder, arm and upper torso mobility
- Complete with: over-door bracket, pulley assembly, traction cord, 2 hand grips
- Easily installs over any door.



# ADJUSTABLE CANES

## CLASSIC OFFSET HANDLE CANES



212100  
BLACK



212400  
SILVER FROST

## STANDARD OFFSET HANDLE CANES

- Choice of durable, easy to clean vinyl or soft foam grip
- Fabricated using strong extruded aluminum tubing: 7/8" (2.2 cm) diameter shaft, 3/4" (1.9 cm) diameter leg extension with slip resistant rubber tip
- Wrist strap on selected models allows for use of both hands without cane falling away.
- One button height adjustment with locking silencer ring
- Handle height adjusts from 30" to 39" (76 - 99 cm)
- Replacement tip: 6105-B, 3/4" (1.9 cm) diameter black slip-resistant, metal-reinforced rubber cane tip.
- Weight capacity: 300 lbs. (136 kg.)



513672  
BRONZE  
SOFT GRIP



514165  
BLACK  
BARIATRIC  
SOFT GRIP



5185-5  
BLACK  
SOFT GRIP



5585B-5  
LIGHT BRONZE  
SOFT GRIP



5585S-5  
SILVER FROST  
SOFT GRIP

500 LBS. WEIGHT LIMIT

## PATTERN OFFSET HANDLE CANES



513540  
CANTERBURY



513605  
PAISLEY



513606  
COMPASS



513607  
BUTTERFLY



513541  
SUMMER  
GARDEN



514049  
CASINO



514156  
CATS



514157  
BLUE PEACOCK



514158  
CAMOUFLAGE



514160  
FLOWER  
GARDEN



514162  
RED MARBLE



514163  
SUN BURST



514164  
CHEETAH

## SOFT SILICONE HANDLE OFFSET ADJUSTABLE CANE

- Adjustable offset series offers deluxe styling for a cane that is both comfortable and functional for everyday use.
- Soft silicone handle for all day comfort.
- Shaft is fabricated using a strong 7/8" aluminum.
- Wrist strap allows for use of both hands without cane falling away.
- One button height adjustment.
- Handle height adjusts from 30" to 39" (76 - 99 cm).
- Slip-resistant, metal-reinforced rubber cane tip.
- Weight capacity: 300 lbs. (136 kg.)



214171  
METALLIC BLUE



214172  
METALLIC GREEN



214173  
METALLIC PINK



214174  
METALLIC RED



214175  
METALLIC BLACK



## ORTHOPEDIC HANDLE CANES

- Molded ergonomic handle evenly distributes body weight, reduces stress on user's hands and joints
- Fabricated using strong extruded aluminum tubing: 7/8" (2.2 cm) diameter shaft, 3/4" (1.9 cm) diameter leg extension with slip resistant rubber tip
- One button height adjustment with silencer ring
- Handle height adjusts from 31" to 39" (79 - 99 cm).
- Replacement tip: 6105-B, 3/4" (1.9 cm) black slip-resistant, metal-reinforced rubber cane tip.
- Weight capacity 300 lbs. (136 kg.)



5144BZ  
BRONZE  
LARGE



5145BK  
BLACK  
WITH ICE PICK



5145BZ  
BRONZE  
WITH ICE PICK



210100  
BLACK  
MEN'S



212800  
BRONZE  
LADIES'



## TOURIST HANDLE CANES

- Large round foam handle, soft to the touch and shock absorbing, easy and comfortable to grasp.
- Fabricated using strong extruded aluminum tubing: 7/8" (2.2 cm) diameter shaft, 3/4" (1.9 cm) diameter leg extension with slip resistant rubber tip
- One button height adjustment with silencer ring
- Handle height adjusts from 33" to 41" (84 - 104 cm)
- Replacement tip: 6105-B, 3/4" (1.9 cm) black slip-resistant, metal-reinforced rubber cane tip
- Weight capacity: 300 lbs. (136 kg.)



5079-B  
BLACK  
WITH ICE PICK



5079-S  
SILVER  
WITH ICE PICK



5086  
SILVER FROST  
MEN'S



5086-L  
SILVER FROST  
LADIES'



514171  
BLACK  
MEN'S



210600  
SILVER FROST  
MEN'S



211400  
SILVER FROST  
LADIES'



211800  
BRONZE  
LADIES'

\* REPLACEMENT 'ICE PICK' CANE TIPS AVAILABLE FOR MODELS 5079-B / S (5079TP)

# ADJUSTABLE CANES

## DEVON AND PALM HANDLE CANES

- Molded and textured ergonomic handle evenly distributes body weight, reduces stress on hands and joints.
- Fabricated using strong extruded aluminum tubing: 7/8" (2.2 cm) diameter shaft, 3/4" (1.9 cm) diameter leg extension with slip resistant rubber tip
- One button height adjustment with locking silencer ring
- Devon handle height adjusts from 31" to 40" (79 - 101 cm)
- Palm handle height adjusts from 28" to 36" (71 - 91 cm); specify for right or left hand use.
- Replacement tip: 6105-B, 3/4" (1.9 cm) black slip-resistant, metal-reinforced rubber cane tip.
- Weight capacity: 300 lbs. (136 kg.)



5181  
BLACK



5581B  
BRONZE



5581S  
SILVER



213601 - LEFT  
213602 - RIGHT  
BLACK

## FRITZ HANDLE

- Evenly distributes body weight, reduces stress on user's hands and joints.
- Attractive wood finish.
- 3/4" (1.9 cm) diameter wood shaft.
- Length from handle to tip: fritz handle 35" (89 cm), wooden styles can be cut to desired user height.
- Replacement tip: 6105-B, 3/4" (1.9 cm) black slip-resistant, metal-reinforced rubber cane tip.
- Weight capacity: 300 lbs. (136 kg.)



5143BK  
ALUMINUM  
BLACK



5143BZ  
ALUMINUM  
BRONZE

# NON-ADJUSTABLE CANES

## DERBY HANDLE

- Derby handle reduces stress on user's hands and joints.
- Length from handle to tip: 36" (91 cm), can be cut to desired user height
- Weight capacity: 300 lbs. (136 kg.)



513299  
MAHOGANY  
MEN'S



513308  
ROSEWOOD  
LADIES'

## FRITZ HANDLE

- Evenly distributes body weight, reduces stress on user's hands and joints.
- Attractive wood finish. 3/4" (1.9 cm) diameter wood shaft.
- Length from handle to tip: fritz handle 35" (89 cm), wooden styles can be cut to desired user height.
- Replacement tip: 6105-B, 3/4" (1.9 cm) black slip-resistant, metal-reinforced rubber cane tip.
- Weight capacity: 300 lbs. (136 kg.)



513357  
WOOD  
BLACK  
MEN'S



513412  
WOOD  
BLACK  
MEN'S

## DERBY HANDLE CANES

- Derby handle reduces stress on user's hands and joints.
- Easily collapses into four sections to fit into small spaces  
Folded dimensions: 12" x 5 1/2" x 1" (30.5 cm x 14 cm x 2.5 cm).
- Fabricated using strong extruded aluminum tubing:  
7/8" (2.2 cm) diameter shaft, 3/4" (1.9 cm) diameter leg extension with slip resistant rubber tip.
- Quick and easy one button height adjustment with silencer ring. Silencer ring prevents rattling.
- Handle height adjusts from 33" to 37" (84 - 94 cm).
- Replacement tip: 6105-B, 3/4" (1.9 cm) black slip-resistant, metal-reinforced rubber cane tip.
- Weight capacity: 300 lbs. (136 kg.)



513532  
MATTE BLACK



513533  
MATTE BRONZE



513390  
PAISLEY



## FRITZ HANDLE CANES

- Large Fritz style handle of carved wood evenly distributes body weight, reduces stress on user's hands and joints.
- Easily collapses into four sections to fit into small spaces.  
Folded dimensions: 12" x 5 1/2" x 1" (30.5 cm x 14 cm x 2.5 cm).
- Fabricated using strong extruded aluminum tubing: 7/8" (2.2 cm) diameter shaft, 3/4" (1.9 cm) diameter leg extension with slip resistant rubber tip
- Quick and easy one button height adjustment with silencer ring. Silencer ring prevents rattling.
- Handle height adjusts from 33" to 37" (84 - 94 cm).
- Replacement tip: 6105-B, 3/4" (1.9 cm) black slip-resistant, metal-reinforced rubber cane tip.
- Weight capacity: 300 lbs. (136 kg.)



513388  
BRONZE



513386  
BLACK



513528  
BUTTERFLY



514138  
CANTERBURY



514140  
CATS



514141  
BLUE  
PEACOCK



# SPECIALTY CANES



## QUAD ADJUSTABLE CANE

MODEL 529500

- Free standing – cane can be set aside so hands can be used for other activities
- Four points of floor contact provide stability when walking
- Sturdy aluminum shaft, 7/8" (2.2 cm) diameter
- Handle height adjusts from 28" to 36" (71 - 91 cm)
- Replacement tips: 6109-B, 1/2" (1.3 cm) black cane tip (pair)
- Weight capacity: 300 lbs. (136 kg.)

## TRI-SEAT ADJUSTABLE CANE

MODEL 513858

- Made from durable light weight aluminum
- Folds into round handle cane
- Weight capacity: 300 lbs. (136 kg.)

## EASY RISER ADJUSTABLE CANE

MODEL 5188

- Two specifically located handgrips allows you to apply all of your weight in the middle of the cane to relieve pressure on joints and legs when rising from a sitting position.
- Two soft ergonomic handgrips are made out of a comfortable foam material that allows the cane to fit naturally into the curvatures of the hand.
- Black aluminum shaft.
- Adjusts 34" - 42" (86.4 - 106.7 cm).
- Distance between hand grips: 8" (20.3 cm).
- Non-slip reinforced tip.
- Weight capacity: 300 lb (136 kg).

## EASY RISER QUAD ADJUSTABLE CANE

MODEL 5190

- Two specifically located handgrips allows you to apply all of your weight in the middle of the cane to relieve pressure on joints and legs when rising from a sitting position.
- Two soft ergonomic handgrips are made out of a comfortable foam material that allows the cane to fit naturally into the curvatures of the hand.
- Black aluminum shaft.
- Adjusts 30" - 39" (76.2 - 90.1 cm).
- Distance between hand grips: 8" (20.3 cm).
- Black coated steel base, 6" x 8" (15 x 20 cm).
- Non-slip reinforced tips.
- Weight capacity: 300 lb (136 kg).

# CANE ACCESSORIES



## CANE WRIST STRAP

MODEL 5131

- Easily attaches to all PCP Canes
- Made from comfortable nylon



## CANE HOLDER

MODEL 5132

- Works with all PCP Canes
- Adjustable to fit different diameters

## CANE HOLDER FOR WHEELCHAIR / ROLLATOR

MODEL 5136

- Easy to attach
- Allows canes to be easily mounted to wheelchairs and rollators



## ICE CANE ATTACHMENT

MODEL 6107-5

- For safer walking on ice and snow
- Strong five prong steel tip
- Adjustable to fit different diameters



## REPLACEMENT CANE TIPS

MODELS:

6103 - 1" (2.54 cm) | 6104 - 7/8" (2.2 cm) | 6105 - 3/4" (1.9 cm)  
| 6108 - 5/8" (1.6 cm) | 6109 - 1/2" (1.27 cm)

- Replacement cane tip, fits all PCP canes
- Set of 2
- 6103 / 6104 / 6108 / 6109 available in both Grey and Black
- 6105 available in both Grey, Black, and Brown.



## REPLACEMENT CANE TIPS FOR SOFT SILICONE HANDLE OFFSET ADJUSTABLE CANES

MODELS:

6118-171 (BLUE) | 6118-172 (GREEN) | 6118-173 (PINK)  
| 6118-174 (RED) | 6118-175 (BLACK)

- Anti-slip, first quality rubber
- Sold in pairs
- Wide range of color options (Black, Blue, Green, Pink, and Red)



## REPLACEMENT QUAD CANE TIPS

MODEL 6106 - 5/8" (1.6 cm)

- Replacement quad cane tips
- Set of 4
- Available in both Grey and Black



## TRIPOD TIP

MODEL 6114 - 3/4" (1.9 cm)

- Provides balance and stability
- Makes cane self standing
- Designed to fit 3/4 (1.9 cm) cane



## FREE STANDING CANE TIP

MODEL 6117

- Designed to increase balance, traction and stability.
- Tip with 4 rubber contact pads allows cane to be free-standing.
- Non-slip rubber base provides added stability when walking.
- Fits all PCP standard canes 3/4" (1.9 cm).
- Tool-free installation.



## CRUTCHES



### PUSH-BUTTON CRUTCHES

#### MODELS:

5091S - TALL - FOR HEIGHTS 5'10" - 6'6"

5092S - ADULT - FOR HEIGHTS 5'2" - 5'10"

5093S - JUNIOR - FOR HEIGHTS 4'6" - 5'2"

5094S - CHILD - FOR HEIGHTS 4'2" - 4'8"

- Completely outfitted including tips, handgrips, underarm pads
- Quick and easy single button height adjustment in 1" (2.54 cm) increments
- Printed height scale eliminates guesswork when sizing
- Easy push button adjustment with user height scale
- Weight capacity: 250 lbs. (114 kg.)

### PUSH-BUTTON FOREARM CRUTCHES

#### MODELS:

5090 - STANDARD - FOR HEIGHTS 5'2" - 6'2"

5090-J - JUNIOR - FOR HEIGHTS 4'2" - 5'2"

- Strong, lightweight 1" (2.5 cm) diameter aluminum shaft
- Coated cuffs with 180° pivot
- Molded vinyl hand grips on welded and coated steel handles
- Adjustment at legs and forearm cuffs
- Weight capacity: 250 lbs. (114 kg.)

## WALKERS



### SINGLE BUTTON FOLDING WALKER

#### MODELS:

5051 - STANDARD

5051W - STANDARD WITH 2 WHEELS, 2 SKIS

5051J - JUNIOR

5051JW - JUNIOR WITH 2 WHEELS, 2 SKIS

- Quick and easy folding with a single touch of a button
- Top release palm activated mechanism
- Strong riveted 1" (2.5 cm) diameter aluminum frame
- Coated steel front cross brace with internal bushing for smooth operation
- Reinforced rubber tips
- Adjustable height: 32"-39" (81-99 cm), Junior 25.5"-32.5" (65-82.5 cm).
- Overall width with wheels: 24.5" (62.2 cm), Junior 23" (58.4 cm).
- Overall width without wheels: 22.5" (57.2 cm), Junior 19.5" (49.5 cm).
- Width between handles: 18.5" (47 cm), Junior 18" (45.7 cm).
- Weight capacity: 250lb (114Kg)

### DOUBLE BUTTON FOLDING WALKER

#### MODELS:

5052 - STANDARD

5052W - STANDARD WITH 2 WHEELS, 2 SKIS

5052J - JUNIOR

5052JW - JUNIOR WITH 2 WHEELS, 2 SKIS

- Each side folds independently allowing walker to be pulled closer to assist while rising from a sitting position
- Palm activated push button folding mechanism
- Sturdy 1" (2.5 cm) diameter frame
- Coated steel front cross brace with internal bushing for smooth operation
- Double side flush cross braces for optimum stability
- Reinforced rubber tips
- Adjustable height: 32"-39" (81-99 cm), Junior 25.5"-32.5" (65-82.5 cm).
- Overall width with wheels: 24.5" (62.2 cm), Junior 23" (58.4 cm).
- Overall width without wheels: 22.5" (57.2 cm), Junior 19.5" (49.5 cm).
- Width between handles: 18.5" (47 cm), Junior 18" (45.7 cm).
- Weight capacity: 250lb (114Kg)



## HEMI FOLDING WALKER

MODEL 5063

- For users with limited or no dexterity in one arm or hand
- Strong lightweight aluminum frame, 7/8" (2.2 cm) diameter
- Folding cross brace with firm locking action
- Angled legs for full tip contact with floor
- Black foam hand grips
- Adjustable height: 33"-37" (84-94 cm).
- Weight capacity: 300 lbs.



## WALKERS AND CRUTCH ACCESSORIES

### CRUTCH TIPS

MODELS:

6099G - SMALL | 6100G - STANDARD | 6101G - LARGE

- Set of 2
- Available in Grey



### CRUTCH UNDERARM PADS

MODELS 5070G

- Extra thick cushioning to provide comfort and reduce fatigue
- Easily washable, perspiration resistant
- Does not contain natural latex rubber
- Sold in pairs
- Available in Grey



### CRUTCH HAND GRIPS

MODEL:

5071 - OPEN STYLE

5072 - CLOSED STYLE

- Extra thick cushioning to provide comfort and reduce fatigue
- Easily washable, perspiration resistant
- Does not contain natural latex rubber
- Sold in pairs



### SYNTHETIC SHEEPSKIN CRUTCH COVERS

MODELS:

6267NA - SYNTHETIC

6268NA - NATURAL

- Includes 2 piece underarm pads and 2 piece hand grips
- Adds cushioning and support to mobility aids
- Protects sensitive areas from pressure and soreness
- Helps to prevent skin ulcers and pressure sores
- Sheepskin reduces pressure at the point of body contact
- Eliminates friction and tearing of the skin
- Prevents the buildup of surface moisture, keeping the patient dry





## REPLACEMENT TIPS FOR WALKERS

MODEL:

6102-G - GREY

6102-B - BLACK

- Replacement walker tips, designed to fit all PCP walkers
- Fits 1 1/8" (2.8 cm) walker legs
- 4 per pack



## SKI GLIDE ATTACHMENTS FOR WALKER LEGS

MODELS:

5106/78 - 7/8" (2.2 cm)

5106/34 - 3/4" (1.9 cm)

5106/10 - 1" (2.54 cm)

- Fit on the rear legs of a standard walker and, with wheel attachments fitted onto the front legs.
- Makes standard walker easier to push and maneuver.
- Attach to the walker legs with a screw.



## NYLON WALKER SKIS, 1"(2.5 CM) DIAMETER

MODEL 5108

- Fit on the rear legs of a standard walker
- Turn / push directly into the tube of the legs in order to create a tight and secure fit
- Makes standard walkers easier to push and maneuver



## WHEEL ATTACHMENT FOR WALKERS

MODEL 5109

- Attachments fit on the front legs of a standard walker, and with nylon walker skis fitted into the rear legs, converts a standard walker to a wheeled walker
- Rubber wheels move smoothly over irregular surfaces
- Quick and easy to install
- Eight height adjustments
- 5" (12.7 cm) diameter



## WALKER POUCH

MODEL 5164

- Large 10" x 15" (25.5 x 38 cm) pocket with smaller dual inner pockets.
- Blue denim fabric with hook and loop fastening
- Machine washable



## WIRE BASKET FOR DUAL RELEASE WALKER

MODEL 5166

- Large 16" x 6" x 7" (40.5 x 15 x 18 cm) basket
- White vinyl coated steel wire
- Molded plastic tray at bottom

## LIGHTWEIGHT ROLLATOR

### MODEL 5305

- Adjustable height handles with molded handgrips; adjust in height from 30" - 36" (76 - 91.5 cm)
- Loop-type locking cable brakes operate on rear wheels; specially angled for those with dexterity issues
- Smooth rolling 8" (20 cm) wheels suitable for all terrains
- Rugged, high-impact plastic seat, 16" x 10" (40 x 25 cm), folds for easy transport and storage
- Curved backrest is easily adjustable
- Removable nylon and mesh carrying bag immediately beneath seat; can stay in position when unit is folded
- Distance between handle bars: 18" (45.72 cm)
- Overall width: 23" (58.4 cm.)
- Product weight: 22 lbs. (10 kg.)
- Weight capacity: 300 lbs. (136 kg.)



## ROLLATOR WITH CURVED BACKREST

### MODEL 5311

- Adjustable height handles with comfortable molded handgrips.
- Comfortably padded seat, folds for easy transport and storage.
- Curved backrest is removable, swings away to allow complete access to nylon storage pouch.
- Nylon storage pouch immediately beneath seat; can stay in position when unit is folded.
- Smooth rolling 8" (20.3 cm) wheels suitable for all terrains.
- Fits users 5'0" to 6'1 inches; Padded seat and back bar.
- Storage compartment; over all width 22.5".
- Adjustable height handles with comfortable molded handgrips; adjust in height from 30" to 36" (76 - 91 cm)
- Easy transport and folding; Collapsed dimension under 17 inches.
- Overall width: 22.5" (57.2 cm.)
- Product weight: 17 lbs. (7.7 kg.)
- Weight capacity: 300 lbs. (136 kg.)



## ROLLATOR WITH CURVED BACKREST

### MODEL 5313TI

- Easy to operate ergonomic hand brakes lock rear wheels for security
- Truly effortless folding mechanism
- Curved backrest for more comfort when sitting
- Comfortably padded seat
- Basket can be positioned either under or in front of seat
- Adjustable handle height for custom fit
- Smooth rolling 8" (20 cm) wheels suitable for all terrains
- Distance between handle bars: 18" (45.72 cm)
- Overall width: 22.5" (57.2 cm.)
- Product weight: 19 lbs. (8.6 kg.)
- Weight capacity: 300 lbs. (136 kg.)



# ROLLATORS AND MOBILITY CHAIRS



## ROLLATOR WITH CURVED BACKREST

MODELS:

5314RD - RED

5314BK - BLACK

- Easy to operate ergonomic hand brakes lock rear wheels for security
- Truly effortless folding mechanism
- 6" (16.24 cm) wheels for smoother rolling operation
- Curved backrest for more comfort when sitting
- Comfortably padded seat
- Basket can be positioned either under or in front of seat
- Adjustable handle height for custom fit
- Distance between handle bars: 18" (45.72 cm)
- Overall width: 21" (53.3 cm.)
- Product weight: 16 lbs. (7.3 kg.)
- Weight capacity: 250 lbs. (114 kg.)



## LIGHTWEIGHT TRANSPORT CHAIR

MODEL 800003

- Back folds down for easy transport and storage
- Total unit weighs less than twenty pounds (20 lbs.)
- Durable fabric upholstery, well padded
- Black aluminum frame
- 8" (20 cm) molded mag type wheels
- Comfortable padded armrests
- Swing-away and detachable footrests
- Seat belt for safety
- Distance between handle bars: 19" (48.26 cm)
- Overall width: 22.5" (57.2 cm.)
- Product weight: 22 lbs. (10 kg.)
- Weight capacity: 300 lbs. (136 kg.)



## WHEELCHAIR

MODEL 6009

- Durable fabric upholstery, well padded.
- Black aluminum frame.
- Wheelchair contains an anti-tipping support for added stability.
- Swing-away and detachable footrests.
- 8" (20 cm) hard rubber front wheels.
- Seat width: 18" (45.72 cm).
- Overall width: 22.5" (57.2 cm.)
- Product weight: 38.36 lbs. (17.4 kg.)
- Weight capacity: 350 lbs. (159 kg.)



PCP's planogram layouts make it easier for you to merchandise, for your staff to sell from, and for your more self-service clients to shop. Modular planogram components enable you to grow your business gradually, purchasing at your own pace, and getting the most return for your investment. Choose between starter and intermediate assortments depending on your business' particular needs and focus.

## PCP PRODUCT INDEX

|                |                |                 |                 |               |               |
|----------------|----------------|-----------------|-----------------|---------------|---------------|
| 10955.....15   | 5090.....30    | 514158.....24   | 6100G.....31    | 6233.....19   | 7016.....10   |
| 10957.....15   | 5090-J.....30  | 514160.....24   | 6101G.....31    | 6234.....19   | 7017.....10   |
| 210100.....25  | 5091S.....30   | 514162.....24   | 6102.....32     | 6235.....19   | 7018.....12   |
| 210600.....25  | 5092S.....30   | 514163.....24   | 6103.....29     | 6236.....19   | 7019.....10   |
| 211400.....25  | 5093S.....30   | 514164.....24   | 6104.....29     | 6237.....19   | 7020.....11   |
| 211800.....25  | 5094S.....30   | 514165.....24   | 6105.....29     | 6238.....19   | 7021.....12   |
| 212100.....24  | 5106-78.....32 | 514171.....25   | 6106.....29     | 6239.....19   | 7022.....11   |
| 212400.....24  | 5106/10.....32 | 5143BK.....26   | 6107-5.....29   | 6243.....21   | 7024.....11   |
| 212800.....25  | 5106/34.....32 | 5143BZ.....26   | 6108.....29     | 6244.....21   | 7026.....12   |
| 213601.....26  | 5108.....32    | 5144BZ.....25   | 6109.....29     | 6246.....21   | 7027.....10   |
| 213602.....26  | 5109.....32    | 5145BK.....25   | 6114.....29     | 6247.....21   | 7036.....7    |
| 214171.....25  | 5114.....13    | 5145BZ.....25   | 6117.....29     | 625700.....21 | 7039.....7    |
| 214172.....25  | 5114-GR.....13 | 5164.....32     | 6118.....7      | 625701.....21 | 7041.....9    |
| 214173.....25  | 5115.....13    | 5166.....32     | 6118-171.....29 | 62600.....21  | 7043.....9    |
| 214174.....25  | 513000.....15  | 5181.....26     | 6118-172.....29 | 62601.....21  | 7047.....9    |
| 214175.....25  | 5131.....28    | 5185-5.....24   | 6118-173.....29 | 6261.....18   | 7100.....4    |
| 350010.....17  | 5132.....28    | 5188.....28     | 6118-174.....29 | 6262.....18   | 7101.....4    |
| 350011.....17  | 513299.....26  | 5190.....28     | 6118-175.....29 | 6263.....18   | 7102.....4    |
| 350110.....17  | 513308.....26  | 529500.....28   | 6145.....18     | 6264.....18   | 7103.....4    |
| 5008.....14    | 513357.....26  | 5305.....33     | 614900.....21   | 6266.....18   | 7104.....5    |
| 5008-H.....14  | 513386.....27  | 5311.....33     | 6150-MF.....20  | 6267NA.....31 | 7105.....5    |
| 5024.....13    | 513388.....27  | 5313TI.....33   | 6151.....20     | 6268NA.....31 | 7106.....5    |
| 5024-ST.....13 | 513390.....27  | 5314BK.....34   | 6152-MF.....19  | 6298.....23   | 7107.....5    |
| 5029.....12    | 513412.....26  | 5314RD.....34   | 6154-MF.....19  | 6299.....23   | 7252.....7    |
| 5040.....16    | 513528.....27  | 5523.....12     | 6155.....20     | 6400.....22   | 7255.....7    |
| 5051.....30    | 513532.....27  | 5526-RBE.....13 | 6167.....21     | 6401.....22   | 7392.....6    |
| 5051J.....30   | 513533.....27  | 5581-B.....26   | 6168.....21     | 6402.....22   | 800003.....34 |
| 5051JW.....30  | 513540.....24  | 5581-S.....26   | 6169.....21     | 6403.....22   | 832000.....17 |
| 5051W.....30   | 513541.....24  | 5585B-5.....24  | 6203.....18     | 640400.....22 | 833000.....17 |
| 5052.....30    | 5136.....29    | 5585S-5.....24  | 6204.....18     | 646400.....18 | 8412.....8    |
| 5052J.....30   | 513605.....24  | 6009.....34     | 6205.....18     | 6500.....14   | 8416.....8    |
| 5052JW.....30  | 513606.....24  | 6031.....15     | 6206.....17     | 6502.....14   | 8418.....8    |
| 5052W.....30   | 513607.....24  | 6032.....15     | 6207.....17     | 68030B.....14 | 9112K.....8   |
| 5063.....31    | 513672.....24  | 6033.....15     | 6208.....23     | 68033B.....14 | 9116K.....8   |
| 5070G.....31   | 513858.....28  | 6086.....14     | 6221.....23     | 7001.....6    | 9118K.....8   |
| 5071.....31    | 514049.....24  | 6088.....16     | 6226.....23     | 7002.....6    | 9124K.....8   |
| 5072.....31    | 514127.....15  | 6091.....16     | 6227.....23     | 7002TP.....7  | 9132K.....8   |
| 5079-B.....25  | 514138.....27  | 6092.....16     | 6229BL.....20   | 7003.....6    | 9211.....8    |
| 5079-S.....25  | 514140.....27  | 6094.....16     | 6229GR.....20   | 7007.....11   | 9212.....8    |
| 5079TP.....25  | 514141.....27  | 6095.....16     | 6230.....18     | 7009.....12   | 9216.....8    |
| 5086.....25    | 514156.....24  | 6096.....14     | 6231.....20     | 7013.....9    | 9226.....8    |
| 5086-L.....25  | 514157.....24  | 6099G.....31    | 6232.....19     | 7015.....9    | 9312.....8    |
|                |                |                 |                 |               | 9316.....8    |
|                |                |                 |                 |               | 9324.....8    |

# PCP<sup>®</sup>

## Personal Care Products

®



**PCP**  
Personal Care Products

Sure-grip bathtub safety rail  
Baranda de seguridad para bañera de agarre firme  
Main courante à prise pour baignoire

**PCP**  
Personal Care Products

Padded bath safety seat with back rest  
Asiento de seguridad almohadillado para baño con respaldo

**PCP**  
Personal Care Products

Padded round safety stool  
Banco de seguridad redondo almohadillado  
Tabouret à siège rond rembourré

**PCP**  
Personal Care Products

Adjustable bath safety step  
Escalón ajustable de seguridad para baño  
Marche réglable pour baignoire

Do not use near open flame  
 Do not use as a weight  
 Do not use for climbing  
 100 lb weight capacity  
 Capacité de charge de 100 lb (45 kg)



**AIRWAY**  
SURGICAL APPLIANCES

PCP is a division of Airway Surgical Appliances  
Tel: 800-267-3476 • Fax: 800-267-7805  
E-mail: [custserv@airwaysurgical.ca](mailto:custserv@airwaysurgical.ca)  
[Airwaysurgical.ca](http://Airwaysurgical.ca)