

Airgo Aluminum Etched Cane

HCPCS:



FIGURE A



FIGURE B

- Our aluminum etched canes are not only lightweight they really stand out with their attractive etched tubing design
- They combine quality with the latest manufacturing technologies to create a cane that's as functional as it is fabulous
- These canes also feature a comfortable foam handle with a wrist strap for extra security and are adjustable in length
- Easy push-button adjustments with anti-rattle lock-nut
- Recommended for patients up to 6' 6"

PRODUCTS	Description
730-420	ALUM. ETCHED OFFSET CANE, BLACK,1EA,6CV
SPECIFICATIONS	Description
Weight Capacity	250 lb (113 Kg)