

Airgo Comfort-Plus Adjustable Quad Canes, Small Base

HCPCS:



- For people who need extra help with balance, our quad cane offers added support in a lightweight design
- Built for longevity from durable anodized aluminum
- The K-shaped rotating base swivels and locks into position to accommodate right-handed or left-handed users and is available in small or large size for greater stability
- Ergonomic cushion-top handle conforms to your hand and helps reduce impact shock
- Designed to fit more comfortably in your hand and is ideal for people with arthritis
- Cane strap with reflective patch for added night time safety
- 11 individual height settings offer an unparalleled range of adjustments that fits most people 5' - 6'5"
- Free Standing Conveniently stands on its own

PRODUCTS	Description
770-852	AIRGO QUAD CANE SM, BLUE, 1 EA,2CV,RET
770-854	AIRGO QUAD CANE SM, BLACK, 1 EA,2CV,RET
770-856	AIRGO QUAD CANE SM, CHAR, 1 EA,2CV,RET
SPECIFICATIONS	Description
Warranty	Limited Lifetime
Weight Capacity	300 lb (136 Kg)
Length	36" (91.4 cm)
Weight	7.0 lb (3.2 kg)