

Bellavita Dive Bath Lift

HCPCS:



FIGURE A



FIGURE B



FIGURE C




FIGURE D




FIGURE E




FIGURE F

The lightest and lowest bath lift on the market now helps users dive deeper and get cleaner!

- New removable backrest with 50° recline
- Larger seat with hygiene cutout
- Improved dual column design provides optimal stability for up to 300 pounds
- Seat height adjusts from 2.6" to 18.9"
- Easy-to-use floating remote control
- Installs easy in almost any bathtubno plumbing or electrical work required

PRODUCTS	Description
477400252	Bellavita Dive Bath Lift w/ Wh Cvr, 1/EA
SPECIFICATIONS	Description
Warranty	5yr Limited
Product Weight	21.8 lb
Material	Plastic
Weight Capacity	300 lb
Color	White
Retail Box Dimensions	15.6" (L) x 17.3" (W) X 14.4" (H)
Retail Box Weight	5.8 lb
CA Prop 65	yes DEHP chemical