

Product Portfolio

BEDS, SURFACES and POSITIONERS



About

Savaria Patient Care provides solutions that support safe patient handling, effective pressure injury prevention, and improved clinical outcomes for long-term care, hospital settings and homecare environments. The company’s research and clinically-informed product line-up includes therapeutic surfaces for pressure management and positioning, medical beds and a complete portfolio of innovative ceiling lifts and slings.

Built from a heritage of brands including Span, Handicare and Silvalea, Savaria Patient Care is a division of Savaria Corporation (TSX:SIS), a global leader in accessibility.

savariapatientcare.com

Table of Contents

BEDS	4	SPECIALTY SURFACES	44
Introduction	4	Bariatric Models	45
Key Features	6	Extra-wide and Consumer Sizes	46
Encore Bed	8	Recessed Deck Models	47
ReadyWide Bed	10	In-Home Styles	48
Advantage Bed	11		
Weigh Scale Bed	12	OVERLAYS & SEATING	49
Profiles & Finishes	13	POSITIONERS	58
Dolly	13	SEATING & SURFACE SELECTION GUIDES	64
Risk Manager	15		
PRESSUREGUARD Support Surfaces	16		
Powered Air	18		
Bariatric APM	19		
APM- <i>i</i>	20		
APM-2	22		
APM2 Safety Supreme	23		
APM	24		
Non-powered and Convertibles	25		
Custom Care Convertible LAL	26		
Custom Care Convertible	27		
Custom Care	28		
Span-Care	29		
Microclimate Management	30		
Easy Air XL	31		
Easy Air LR	32		
Easy Air	33		
Protocol	34		
Topper	35		
GEO-MATTRESS Therapeutic Foam Mattresses	36		
UltraMax	38		
UltraMax for Stretchers	39		
Max	40		
Pro	41		
Plus	40		
Wings	41		
Atlas	42		
HC	42		
Gel-80	43		

Beds


At Savaria, our bed solutions are designed with the needs of healthcare environments in mind—balancing durability, comfort, and practical functionality. Built to support both patients and caregivers, our beds are engineered for reliable performance in demanding care settings.

From dependable everyday models to more specialized configurations, Savaria beds are thoughtfully designed to enhance safety, ease of use, and patient comfort, while supporting efficient care delivery across a range of healthcare environments.



Beds


Individually powder-coated bed components



KEY FEATURES:

- Durable, anti-microbial polyester finish
- Resists chipping and flaking
- Featured on: **All models**


Readywide® integrated “two click” expansion



KEY FEATURES:

- Click once to pull expander out to 39” wide. Click again for 42” width.
- Navigates narrow doorways: expanders provide fast, hassle-free solution for in-facility bed mobility
- Caregiver advantages: Wider surface for safer, easier transfers and ADLs.
- Featured on:
 - **Encore Bed**
 - **ReadyWide Bed**
 - **Weigh Scale Bed**

Intelligent Smart Stop®



KEY FEATURES:

- Low position auto-pause improves low-bed safety and ergonomics
- Prevents falls by prohibiting the user from trying to egress from a dangerously too-low level.
- Lowering down into low height zone is only possible using the FEC.
- Featured on:
 - **Encore Bed**
 - **ReadyWide Bed**
 - **Weigh Scale Bed**

GlideAlign



KEY FEATURES:

- Keeps the user vertically aligned, prevents sliding during HOB elevation, and reduces repositioning needs.
- Featured on: **Encore Bed**



Encore® Bed

The Encore® is our premium bed, delivering enhanced comfort and care for users while improving caregiver ergonomics and reducing injury risk.



KEY FEATURES:

- GlideAlign Retractable Deck – keeps the user vertically aligned, prevents sliding during HOB elevation, and reduces repositioning needs
- Intelligent Smart Stop® low position auto-pause improves low-bed safety and ergonomics
- Advanced Clinical Positioning includes Comfort Chair, T/RT, and Cardiac Chair
- Bed travels virtually straight up and down, while also providing ample underbed clearance
- Dual-position, dual intensity under-bed lighting
- Trapezoidal slat profile offers superior strength, ventilation, and cleaning ease
- Integrated ReadyWide® expanders to 42"
- Foot-end controller with staff lockout
- Customize with a range of rails, headboard/footboard options, bed extenders, and other accessories detailed below:

ITEM CODE	DESCRIPTION
MC9RW	Encore with ReadyWide®
MC9RZ	Encore with built-in backup battery
MC9RWZ	Encore ReadyWide® with built-in backup battery
R9308H	1/2 side molded rails head, pair
R9519	R9519 Split side rail pair
R9958-L/R	Single molded Pivot Assist rail, quick width adjustable, frame mounted
R9959	Pair molded Pivot Assist rail, quick width adjustable, frame mounted
Q6731	Trapeze mounting frame
Q6730	Trapeze Pole, offset
BS79(B)	I.V. Pole (adjustable, board mounted)
Q6785	4" Bed extender kit - factory installed
Q6689	Expander kit (39"/42" quick adjustment)
Q6701	Battery Back-up kit - factory installed
Q6699	Transport Dolly Set
Q6729	Foot end equipment bar
Q6531	Roller bumper
CS82	Board mounted base for I.V. Pole
Q6760	Bed Extender Storage Kit

WARRANTY	
Electrical	5 years
Frame	15 years
Welds	15 years

SPECIFICATIONS	
Height Floor to Mattress Support	7-7/8" - 31" / 200 mm - 787 mm
Deck Angles – Back Elevation	0 – 70°
Deck Angles – Knee Elevation	0 – 30°
Bed Length – Overall Length	82", 86", 90" / 2083 mm, 2184 mm, 2286 mm
Sleep Surface – Retraction Travel	12" / 305 mm
Bed Length – Sleep Surface	76", 80", 84"* / 1930 mm, 2032 mm, 2134 mm*
Bed Width – Sleep Surface	36", 39", 42" / 889 mm, 991 mm, 1067 mm
Bed Width – With Side Rails	40-1/2", 43-1/2", 46-1/2" / 1029 mm, 1105 mm, 1181 mm
Weight of Bed (without accessories)	253 lbs / 115 kg
Safe Working Load	600 lbs / 272 kg
Mobility (at any height)	7-7/8" - 31" / 200 mm - 787 mm
Bed End Panels	Permanent or removable mounting
Ingress Protection Rating	IPX4 (Electrical components IP66)
Electric Motors	Four (4) DC self-lubricating linear actuators

* Requires the 4" Bed Extender Kit

CERTIFICATIONS

MEDICAL ELECTRICAL EQUIPMENT

ANSI/AAMI
ES60601-1:2005

★★★

MEDICAL ELECTRICAL EQUIPMENT

CAN/CSA-C22.2
NO. 60601-1-08

★★★

MEDICAL BEDS

CAN/CSA-C22.2
NO. 60601-2-52:11

★★★

EMC

CAN/CSA-C22.2
NO. 60601-1-2-08

★★★

ReadyWide Bed

ReadyWide® Integrated Expansion and smart technology that enables advanced clinical and resident positioning while boosting caregiver safety and ergonomics.



KEY FEATURES:

- Push-button ReadyWide expandability system up to 42"
- Available in 500lb or 650lb (PLUS) safe working loads
- Smart Stop low-position auto pause for under the bed safety
- Foot-end controller with staff lockout
- Trapezoidal slat profile provides unmatched strength, ventilation and ease of cleaning
- Durable, antimicrobial coating
- Vertical movement with virtually no horizontal movement
- Customize with a range of rails, headboard/footboard options, bed extenders, and other accessories detailed below:

ITEM CODE	DESCRIPTION
MC7W	ReadyWide with foot end controls
MC7WM	ReadyWide with mobility & foot end controls
MC7WL	ReadyWide with floor lock & foot end controls
MC7WML	ReadyWide with mobility, floor lock & foot end controls
MC7WPML	ReadyWide Plus (650 lbs SWL) with mobility, floor lock & foot end controls
R9959	Set of 3 position molded pivot assist rails
R9308H	Set of molded head end half side rails
R9308H-L	Left side molded head end half side rail
R9308H-R	Right side molded head end half side rail
R9519	Split side rail
Q6749	Under bed lighting - factory installed
Q6701	Battery back up (Standard model) - factory installed
Q6770	Battery back up (Plus model) - factory installed
Q6740	Floor placement arm
Q6741	Retractable power cord cable reel & placement arm
Q6785	Extender Kit (Plus model) - factory installed
Q6730	Patient helper, trapeze bar, chain and hand grip
Q6731	Trapeze bed mounting bracket
Q6703	Bed and deck angle indicator
Q6699	Bed transportation system
Q6531	Honeycomb roller bumper kit
Q6787	Underbed lighting - factory intalled
BS79	I.V. Pole

SPECIFICATIONS	
Height Floor to Mattress Support	7-7/8" – 31" / 200 mm – 787mm
Deck Angles – Back Elevation	0 – 70°
Deck Angles – Knee Elevation	0 – 30°
Bed Length – Overall Length	82", 86", 90" / 2083 mm, 2184 mm, 2286 mm
Bed Length – Sleep Surface	76", 80", 84" / 1930 mm, 2032 mm, 2134 mm
Bed Width – Sleep Surface	35", 39", 42" / 889 mm, 991 mm, 1067 mm
Weight of Bed (without accessories)	259 lbs / 118 kg (Standard model) 270 lbs / 123 kg (Plus model)
Safe Working Load	500 lb / 227 kg (Standard model) 650 lb/ 295 kg (Plus model)

WARRANTY	
Mechanical	5 years
Frame	15 years
Welds	15 years

CERTIFICATIONS			
			

Advantage Bed

Smart technology for advanced clinical and resident positioning



KEY FEATURES:


- Precision welded orthopedic grid deck for no-sag support
- Durable, antimicrobial polyester powder-coating
- Vertical movement with virtually no horizontal movement
- Arched legs allow wide stance for clearance and stability
- Baseboard bumpers protect walls
- Standard hand pendant controller
- Smooth deck surface, free from any protruding welds, for easy disinfection
- Anchored corner retainers prevent mattress shifting and dangerous gaps
- Customize with a range of rails, headboard/footboard options, bed extenders, and other accessories detailed below:

ITEM CODE	DESCRIPTION
QD1000	Advantage bed 35" x 80"
QD1000MLF	Advantage bed 35" x 80"; 360° mobility; FEC; floor lock
QD1000Y	Advantage bed 35" x 76"
QD1000YMLF	Advantage bed 35" x 76"; 360° mobility; FEC; floor lock
Q6599	Head end roller bumpers
RM930	Three position pivot assist rail - universal
RAD959	Molded Pivot Assist pair
RM950	Fixed assist rail - universal
RAD306H	Standard half head end side rail - set
RAD308H	Set of molded half head end side rails
RAD308H-L	Left side molded half head end side rail
RAD308H-R	Right side molded half head end side rail
RAD501	3/4 rail, 35" fixed width, folding, pair
Q6730	Trapeze Pole, Offset
Q6733	Mounting Frame, Trapeze
Q6587	Transport Dolly Set
CS82	IV pole adjustable
BS79B	I.V. Pole
RAD959	Pair Molded Pivot Assist Rails, quick width adjustable and frame mounted
Q6736	4" Extender Kit
Q6633	39" Fixed Expander Kit
Q6709	42" Fixed Expander Kit
Q6681	Battery Back-Up - factory Installed
Q6674	Foot end controller
Q6615	Floor Lock Kit

SPECIFICATIONS	
Height Floor to Mattress Support	7-7/8" – 30" / 200 mm – 762mm
Deck Angles – Back Elevation	0 – 70°
Deck Angles – Knee Elevation	0 – 30°
Bed Length – Overall Length	84-3/4", 88-3/4" / 2134 mm, 2235 mm
Bed Length – Sleep Surface	76", 80" / 1930 mm, 2032 mm
Bed Width – Sleep Surface	35" / 889 mm
Weight of Bed (without accessories)	160 lbs / 73 kg
Safe Working Load	500 lb / 227 kg

WARRANTY	
Electrical	5 years
Frame	15 years
Welds	15 years

CERTIFICATIONS



Weigh Scale Bed

Weigh scale low-bed design with advanced clinical positioning



KEY FEATURES:

- Built-in weigh scale with dual controllers
- Push-button ReadyWide expandability system up to 42"
- Available in 500lb or 650lb (PLUS) safe working loads
- Smart Stop low-position auto pause for under the bed safety
- Standard hand pendant and pendant bed controller
- Foot-end controller with staff lockout
- Easy to operate, ratchet-style vascular footrest
- Trapezoidal slat profile offers superior strength, ventilation, and cleaning ease
- Durable, antimicrobial coating
- Vertical movement with virtually no horizontal movement
- Under-bed lighting for safety
- Rapid head deck lowering for emergency treatment
- Customize with a range of rails, headboard/footboard options, bed extenders, and other accessories detailed below:

ITEM CODE	DESCRIPTION
MC7WSP	Weigh Scale Plus Model with ReadyWide, mobility, floor lock, and foot-end controls
MC7WSTML	Weigh Scale with ReadyWide, mobility, floor lock, and premium foot-end controls
MC7WSML	Weigh Scale with ReadyWide, mobility, floor locks, and standard foot end controls
R9958-R	3 position molded pivot assist rail with right bracket
R9958-L	3 position molded pivot assist rail with left bracket
R9959	Set of 3 position molded pivot assist rails
R9308H	Set of molded head end half side rails
R9308H-L	Left side molded head end half side rail
R9308H-R	Right side molded head end half side rail
R9519	Split side rail pair, 3/4 bed coverage
Q6770	Battery back-up – factory installed
Q6740	Floor placement arm
BS79	I.V. Pole
Q6784	Bed extender to 84"
Q6730	Patient helper, trapeze bar, chain and hand grip
Q6731	Trapeze bed mounting bracket
Q6703	Bed and deck angle indicator
Q6699	Bed transportation system
Q6531	Honeycomb roller bumper kit
Q6785	4" bed extender kit - factory installed

SPECIFICATIONS	
Height Floor to Mattress Support	8-5/8" – 32-1/4" / 219 mm – 819mm
Deck Angles – Back Elevation	0 – 70°
Deck Angles – Knee Elevation	0 – 30°
Bed Length – Overall Length	82", 86", 90" / 2083 mm, 2184 mm, 2286 mm
Bed Length – Sleep Surface	76", 80", 84"* / 1930 mm, 2032 mm, 2134 mm*
Bed Width – Sleep Surface	36", 39", 42" / 889 mm, 991 mm, 1067 mm
Weight of Bed (without accessories)	297 lb / 135 kg (Standard model) 301 lb / 137 kg (Plus model)
Safe Working Load	500 lb / 227 kg (Standard model) 650 lb/ 295 kg (Plus model)

WARRANTY	
Electrical	5 years
Frame	15 years
Welds	15 years

CERTIFICATIONS	
----------------	--

ANSI/AAMI
ES60601-1:2005
★★★

CAN/CSA-C22.2
NO. 60601-1-08
★★★

CAN/CSA-C22.2
NO. 60601-2-52:11
★★★

EMC
★★★

Head and Foot Board Profiles & Finishes

Classic

Prestige

Aria

Charleston Oak

Folkstone Grey

Ganache

Icicle

London Grey

Cocoa Pecan

Norway

Shaker Cherry

Solar Oak

Summer Flame

Talas Cherry

Wild Cherry

Williamsburg Oak

Windsor Mahogany

Ergonomic Transport Dolly

Transport dollies for easy bed relocation and maintenance are available for all Span beds.



KEY FEATURES:



- Safe, simple, and compact
 - Attaches easily to both ends of the bed
 - Bed can be safely tipped on its side for relocation ore servicing
- Easily navigates through doorways and into elevators
 - Allows compact storage of overflow bed inventory
 - Provides safe, comfortable access to under-bed components

Item #	Description
Q6587	Set of bed transport dollies for Advantage, Rexx and Rexx Fast models
Q6699	Set of bed transport dollies for Advantage ReadyWide and Encore Models

Patient Room Furniture

Create coordinated, durable patient and resident rooms with furniture solutions from Savaria Patient Care. Designed to complement our beds and support everyday use in care environments, our furniture offering includes healthcare case goods and overbed tables.


Our case goods include bedside nightstands, chests, armoires, and wardrobes, with finish and style options that can coordinate with SPAN bed headboards and footboards for a more cohesive room design. Our overbed tables are stable, height-adjustable, easy to clean, and built to support meals, charting, and bedside care.



Scan for more details online



Risk Manager® Bedside Safety Mat





KEY FEATURES:

- Durable mat made of impact-absorbing lightweight elastomeric gel surface that reduces impact from falls
- Low profile and beveled edge allow easy on and off for wheelchairs and lifts
- Provides comfort for caregivers
- Nano-silver protects mat from stains, odor-causing bacteria, mold, and mildew
- Wipe clean using most standard cleaning agents
- Slide resistant bottom surface grips floor and keeps mat in place
- Textured top provides added safety
- Rolls up for easy transport, storage and shipping
- **Warranty:** 1 year
- **Weight:** 6 lbs (24"), 9.5 lbs (36")

Item #	Description	Qty/Case
37052	24"W x 70"L x 1"H, standard – single	1
94612	24"W x 70"L x 1"H, standard – bulk pack	5
8705236	36"W x 70"L x 1"H, XL – single	1
9611836	36"W x 70"L x 1"H, XL – bulk pack	5

All components non-toxic, latex-free, and lead-free.



PressureGuard®

In the crowded world of specialty support surfaces, the PressureGuard® series truly stands apart. Packed with innovations, every model combines the effectiveness of an air flotation system with the unmatched stability and safety of a multi-component, engineered foam shell. Air and foam work together to maximize healing while promoting better comfort, rest and recovery.

PressureGuard® models are grouped into three sub-categories:


- Powered Air (pages 18-24)
- Non-powered/Convertible (pages 25-29)
- MicroClimate Management (pages 30-35)

All PressureGuard® models include these outstanding design features--

Geo-Matt®

KEY FEATURES:


Time-tested, clinically proven segmented design creates a unique, anti-shearing top surface comprised of more than 800 individually articulating cells in distinct head, torso and foot zones. Ring-of-Air® ventilation channels help eliminate heat/moisture buildup that can lead to pressure injuries and inhibit healing.



Safety Edge™ Advantage In Supine

KEY FEATURES:


On "roll-up" air only surfaces, the air system tends to collapse when the user nears the side. PressureGuard®'s two-part bolster system directs user toward center of mattress, away from fall risk and siderail hazard.



Safety Edge™ Advantage During Transfers


KEY FEATURES:

On roll-ups, edge collapses beneath user's concentrated weight. But PressureGuard® design supports seated user on mattress edge, adding safety to egress and entry.



Powered Air


PressureGuard® powered air surfaces provide pressure injury prevention and treatment thru Stage IV. Reliable, powered control units provide treatment modalities including alternating pressure (all models), plus gentle lateral rotation (APM2 and APM2 Safety Supreme) or bariatric care (Bariatric APM).



Shear Transfer Zones

Provides added shearing protection via silicone-coated, shear-minimizing fabric bands on the underside of its bi-directional stretch cover. Helps guard against the damaging effects of micro shear, macro shear, and rotational (pivot-induced) shear. Shear Transfer Zones help to “glide” the user back to the original position following HOB elevation.


Featured on: APM2 (page 22) and APM-i (page 20-21).



Smartimmerse

At the touch of a button, proprietary electronics assist in configuring the inflation profile. The embedded algorithm helps establish ideal air settings to maximize immersion, support, and redistribution. Feature is designed to complement caregiver assessment and clinical decision-making.

Featured on: APM-i (page 20-21).



Simplicity

Our powered air mattresses are designed with simplicity in mind. Intuitive control units make it easy for busy caregivers by minimizing decisions and simplifying training. The goal? Maximum positive outcomes with minimum caregiver intervention.

Featured on: all Powered Air models.

PressureGuard® Bariatric APM



KEY FEATURES:

- Simple operation
- Gentle and effective pressure changes
- Customizable firmness and time settings
- PressureGuard® design: stability, safety, anti-shearing

THERAPIES:

- Powered Pressure Redistribution
- Alternating Pressure for users up to 750 lbs

Item #	Description	Weight Limit
AP8042-29	Bariatric APM, 80"L x 42"W x 7"H	750 lbs
AP8048-29	Bariatric APM, 80"L x 48"W x 7"H	750 lbs
AP8054-29	Bariatric APM, 80"L x 54"W x 7"H	750 lbs
*AP8442-29	Bariatric APM, 84"L x 42"W x 7"H	750 lbs
AP8448-29	Bariatric APM, 84"L x 48"W x 7"H	750 lbs
*AP8454-29	Bariatric APM, 84"L x 54"W x 7"H	750 lbs
7500	Bariatric APM Control Unit Only	-

*Allow up to 10 working days lead time for this size.

QUICK FACTS:

Use For:

- Aggressive treatment (through Stage 4)
- Prevention of pressure injuries in immobile or high risk patients
- Can't or won't reposition frequently
- Can be positioned off the pressure injury

Medicare Code: E0277 (Group II)

WARRANTY:

Mattress: 3 years | **Cover:** 2 years



Control Unit

KEY FEATURES:

- Up to four times the air output of conventional units.
- Eight electronically-controlled comfort settings.
- Programmable cycle times of 20, 25, 30 or 35 minutes.
- Low pressure indicator and audible alarm.

Warranty: 18 months

PressureGuard® APM-i



KEY FEATURES:

- Smartimmerse™ Technology delivers an optimized balance of immersion, support, and redistribution- lessening the burden on caregivers.
- Proven PressureGuard® design with enhanced immersion and envelopment
- Optimized user comfort

QUICK FACTS:

Use For:


- Treatment of Stage 1- 4 pressure injuries
- Patients that can't/won't reposition frequently
- Can be positioned off pressure injury in at least two positions
- Prevention for high risk patients

Medicare Code: E0277 (Group II)

WARRANTY:

Mattress: 3 years | Cover: 2 years

Item #	Description	Weight Limit
38036	PG APM-i 80x36 with 6200 control unit	350 lbs
38436	PG APM-i 84x36 with 6200 control unit	350 lbs
38042	PG APM-i 80x42 with 6200 control unit	350 lbs
38442	PG APM-i 84x42 with 6200 control unit	350 lbs
6200	Control Unit Only	-



Control Unit

KEY FEATURES:

- Proprietary electronics within the powered control unit automatically customizes the inflation profile for each user at the press of a button. Using an embedded algorithm, the system calculates the precise air settings required to deliver an optimized balance of immersion, support and redistribution
- Multiple powered modalities, including immersion, alternating pressure, and lateral rotation

Warranty: 18 months

PressureGuard® APM2



KEY FEATURES:

- Dual-mode air therapy: Alternating Pressure and Lateral Rotation via a simple recessed toggle switch
- Geo-Matt® segmented high-density foam topper redistributes pressure and reduces heat and moisture buildup
- Safety Edge™ bolster enhances edge-of-bed stability and safer patient transfers
- Timed Auto Firm (~20 min) for transfers, ADLs, and CPR
- “Float” powered flotation mode for patient comfort and pressure redistribution
- “Low Pressure” and power failure indicators for caregiver assurance

THERAPIES:

- Powered Pressure Redistribution
- Alternating Pressure
- Lateral Rotation

Item #	Description	Weight Limit
5875LR-29	APM2 with digital control unit, 75"L x 35"W x 7"H	350 lbs
5880LR-29	APM2 with digital control unit, 80"L x 35"W x 7"H	350 lbs
5884LR-29	APM2 with digital control unit, 84"L x 35"W x 7"H	350 lbs
*AP28039-29	APM2 with digital control unit, 80"L x 39"W x 7"H	350 lbs
*AP28439-29	APM2 with digital control unit, 84"L x 39"W x 7"H	350 lbs
*3715629	APM2 with digital control unit, 80"L x 42"W x 7"H	350 lbs
5900	Digital control unit only	-

*Allow up to 10 working days lead time for this size.

QUICK FACTS:

Use For:

- Treatment of Stage 1- 4 pressure injuries
- Patients that can’t/won’t reposition frequently
- Can be positioned off pressure injury in at least two positions
- Prevention for high risk patients

Medicare Code: E0277 (Group II)

WARRANTY:

Mattress: 3 years | Cover: 2 years

PressureGuard® APM2 Safety Supreme



KEY FEATURES:

- All the features of the APM2 plus the added safety benefits of raised perimeter edges
- “Edge-to-Edge” design to help prompt user away from edge of bed without compromising skin integrity
- Dual-mode air therapy: Alternating Pressure and Lateral Rotation via a simple recessed toggle switch
- Geo-Matt® segmented high-density foam topper redistributes pressure and reduces heat and moisture buildup
- Safety Edge™ bolster enhances edge-of-bed stability and safer patient transfers
- Timed Auto Firm (~20 min) for transfers, ADLs, and CPR
- “Float” powered flotation mode for patient comfort and pressure redistribution
- “Low Pressure” and power failure indicators for caregiver assurance

THERAPIES:

- Powered Pressure Redistribution
- Alternating Pressure
- Lateral Rotation
- Raised Safety Perimeter

Item #	Description	Weight Limit
SAF5875-29	Safety Supreme, 75"L x 35"W x 7"H center, 9"H side	350 lbs
SAF5880-29	Safety Supreme, 80"L x 35"W x 7"H center, 9"H side	350 lbs
SAF5884-29	Safety Supreme, 84"L x 35"W x 7"H center, 9"H side	350 lbs
*37890	Safety Supreme, 80"L x 42"W x 7"H center, 7"H side	350 lbs
*33530	Safety Supreme, 84"L x 42"W x 7"H center, 7"H side	350 lbs
5900	Digital control unit	-

*Allow up to 10 working days lead time for this size.

QUICK FACTS:

Use For:

- Treatment of Stage 1- 4 pressure injuries
- Patients that can’t/won’t reposition frequently
- Can be positioned off pressure injury in at least two positions
- Prevention for high risk patients and patients at risk for falls

Medicare Code: E0277 (Group II)

WARRANTY:

Mattress: 3 years | Cover: 2 years



Control Unit

KEY FEATURES:

- Lightweight, quiet, dependable and easy to operate.
- Low Pressure indicator light and selectable audible alarm
- Power Failure indicator light and selectable audible alarm.
- Simple, indexed comfort adjustment.
- Timed (20-minute) “Auto Firm” mode for transfers, ADLs, & CPR.
- “Float” mode provides powered flotation therapy.
- “Alternate” mode changes loading in 10-minute cycle; controlled at mattress for choice of lateral rotation or alternating pressure therapy

Warranty: 18 months



Control Unit

KEY FEATURES:

- Lightweight, quiet, dependable and easy to operate.
- Low Pressure indicator light and selectable audible alarm
- Power Failure indicator light and selectable audible alarm.
- Simple, indexed comfort adjustment.
- Timed (20-minute) “Auto Firm” mode for transfers, ADLs, & CPR.
- “Float” mode provides powered flotation therapy.
- “Alternate” mode changes loading in 10-minute cycle; controlled at mattress for choice of lateral rotation or alternating pressure therapy

Warranty: 18 months

PressureGuard® APM



QUICK FACTS:
Use For:

- Treatment of Stage 1- 4 pressure injuries
- Patients that can't/won't reposition frequently
- Can be positioned off pressure injury in at least two positions
- Prevention for high risk patients

Medicare Code: E0277 (Group II)

WARRANTY:
Mattress: 18 months | **Cover:** 18 months

- KEY FEATURES:**
- Simple operation
 - Powered flotation mode
 - Gentle and efficient pressure changes
 - PressureGuard® design: stability, safety, anti-shearing
- THERAPIES:**
- Powered Pressure Redistribution
 - Alternating Pressure

Item #	Description	Weight Limit
5175-29	PressureGuard® APM, w/standard control unit, 75"L x 35"W x 7"H	350 lbs
5180-29	PressureGuard® APM, w/standard control unit, 80"L x 35"W x 7"H	350 lbs
5184-29	PressureGuard® APM, w/standard control unit, 84"L x 35"W x 7"H	350 lbs
5500	Standard control unit only	-



Control Unit

KEY FEATURES:

- Comfort dial provides firmness control without compromising pressure redistribution
- Choice of Modes: Alternate (cylinders inflate and deflate in a fixed 8-minute cycle) or Float
- Visual low pressure indicator

Warranty: 18 months

Non-powered and Convertibles


This popular category of specialty surfaces allows facilities to provide proven, non-powered comfort, pressure injury prevention and treatment on every bed. Convertible models facilitate convenient transformation to a powered solution when patient-specific advanced needs arise.



Star Chamber

In a Star Chamber air system, the key to ideal pressure redistribution—Air—is the star. Unlike typical foam block-filled “air” systems, this proprietary design optimizes air volume to reduce concentrated tissue load over bony prominences. Even in non-powered mode, it delivers responsive support and immersion comparable to traditional powered surfaces, enhancing comfort and clinical performance.

Featured on: Custom Care Convertible (page 27) and Custom Care Convertible LAL (page 26).



Treatment Flexibility

Easily convert to a powered surface with the simple addition of a control unit, reducing the need for costly rentals. Add-on powered control units can instantly provide active modalities including alternating pressure, lateral rotation, and microclimate management, depending on model.

Featured on: Custom Care Convertible (page 27), Custom Care Convertible LAL (page 26), Span-Care Convertible (page 29).



LifeSpan™ Fabric

Polycarbonate-fortified fabric for unsurpassed resistance to damaging effects of diluted bleach and other aggressive cleaners and disinfectants. Engineered to deliver an ideal combination of durability, immersion and patient mobility.

Featured on: Custom Care non-powered (page 28), Custom Care Convertible (page 27), Span-Care Convertible (page 29).

PressureGuard® Custom Care® Convertible LAL



QUICK FACTS:

Use For:

- Excessive skin moisture due to sweating, decreased mobility, inability or refusal to be repositioned
- Macerated skin due to any of the above
- Increased skin or body temperature due to infection, sepsis, or other conditions
- A combined score ≤ 5 on Braden Moisture, Activity and Mobility scales
- Treatment of Stage 1-4 Pressure Injuries
- Prevention of Pressure Injuries for high-risk patients
- Some patients with stabilized flaps/grafts

Medicare Code:

E0373 (Group II) (without control unit)
E0277 (Group II) (with control unit)

WARRANTY:

Mattress: 5 years | Cover: 2 years

KEY FEATURES:

- Effective pressure transfer through Geo-Matt® surface geometry and dual heel protection
- Reduced maceration through continuous airflow via Air Diffusion Matrix™ layer
- More effective powered therapies through the Star Chamber™ design
- Easy-to-use control unit with quick-connect air lines and push-button therapy selection

THERAPIES:

- Non-powered (Air) Reactive Pressure Redistribution

With control unit added:

- Microclimate Management (LAL)
- Alternating Pressure
- Lateral Rotation
- Powered Pressure Redistribution

Item #	Description	Weight Limit
CL753629	75"L x 36"W x 7"H	500 lbs (non-powered & alternating modes)
CL803629	80"L x 36"W x 7"H	
CL843629	84"L x 36"W x 7"H	
74750	80"L x 42"W x 7"H	350 lbs (lateral rotation mode)
49907	84"L x 42"W x 7"H	
*CL85V29	Custom Care Convertible LAL for recessed deck frames	350 lbs (lateral rotation mode)
*CL80CA29	Custom Care Convertible LAL 80" for CareAssist®	
*CL84CA29	Custom Care Convertible LAL 84" for CareAssist®	
8400	Add-on control unit for Custom Care Convertible LAL	-

PressureGuard® Custom Care® Convertible



QUICK FACTS:

Use For:

- Treatment of Stage 1- 4 pressure injuries
- Patient can't/won't reposition frequently
- Can be positioned off the pressure injury in at least 2 positions
- Prevention for high risk patients

Medicare Code:

E0373 (Group II) (without control unit)
E0277 (Group II) (with control unit)

WARRANTY:

Mattress: 5 years | Cover: 2 years

KEY FEATURES:

- Minimized shearing through Geo-Matt® surface geometry and Shear Transfer Zones® on bi-directional stretch cover
- Star Chamber system engineered with air support cylinders that are constantly filled with air- whether or not the surface is powered.
- Minimized risk of falls and bed entrapment through Safety Edge™ bolster design
- Flexible treatment modalities, including alternating pressure and lateral rotation
- Easy-to-use ruggedized control unit with quick connect-to-mattress port

THERAPIES:

- Non-powered (Air) Reactive Pressure Redistribution

With control unit added:

- Alternating Pressure
- Lateral Rotation
- Powered Pressure Redistribution

Item #	Description	Weight Limit
CJ753629	75"L x 36"W x 7"H	500 lbs (non-powered mode)
CJ803629	80"L x 36"W x 7"H	
CJ843629	84"L x 36"W x 7"H	
*CJ85V29	Custom Care Convertible for recessed deck frames	350 lbs (powered modes)
*CJ80CA29	Custom Care Convertible 80" for CareAssist®	
*CJ84CA29	Custom Care Convertible 84" for CareAssist®	
26130	80"L x 42"W x 7"H	350 lbs (powered modes)
95724	84"L x 42"W x 7"H	
6500	Add-on control unit for Custom Care Convertible	-



Control Unit

KEY FEATURES:

- Quiet and easy to use
- Quick-connect air lines are color-coded and gender-matched
- Anti-microbial treatment embedded in control panel
- Push-button therapy selection
- Low air loss
- Alternating pressure and lateral rotation modalities
- Simple comfort adjustment
- "Disconnect" button resets surface to ideal non-powered setting

Warranty: 2 years



Control Unit

KEY FEATURES:

- Quiet and rugged
- Quick-connect air lines are color-coded and gender-matched
- Anti-microbial treatment embedded in control panel
- Push-button therapy selection
- Alternating pressure and lateral rotation modalities
- Simple comfort adjustment
- "Disconnect" button resets surface to ideal non-powered setting

Warranty: 2 years

PressureGuard® Custom Care®



- QUICK FACTS:**
Use For:
- Prevention of single Stage 1-4 and multiple Stage 2 pressure injuries
 - Treatment of uncomplicated (non-complex) stage 3 and 4 pressure injuries with multiple turning surfaces
 - Can be positioned off the injury in a least 2 positions
 - Prevention for high risk and/or severely compromised
 - History of previously healed pressure injuries
 - Can be repositioned frequently
 - Poorly padded boney prominences

Medicare Code: E0373 (Group II)

WARRANTY:
Mattress: 5 years | **Cover:** 2 years

- KEY FEATURES:**
- Minimal shearing through Geo-Matt® surface geometry and Shear Transfer Zones® on bi-directional stretch cover
 - Maintenance-free support chambers require no adjustment or refilling
 - No initial set-up necessary
- THERAPIES:**
- Non-powered (Air) Reactive Pressure Redistribution

Item #	Description	Weight Limit
NP753629	75"L x 36"W x 7"H	500 lbs
NP803629	80"L x 36"W x 7"H	500 lbs
NP843629	84"L x 36"W x 7"H	500 lbs
*NP85V29	Custom Care for recessed deck frames	500 lbs
*NP80CA29	Custom Care 80" for CareAssist®	500 lbs
*NP84CA29	Custom Care 84" for CareAssist®	500 lbs
60787	80"L x 42"W x 7"H	500 lbs
66914	84"L x 42"W x 7"H	500 lbs

PressureGuard® Span-Care® Convertible



- QUICK FACTS:**
Use For:
- Treatment of Stage 1- 4 pressure injuries
 - Patient can't/won't reposition frequently
 - Can be positioned off the pressure injury in at least 2 positions
 - Prevention for high risk patients

WARRANTY:
Mattress: 5 years | **Cover:** 2 years

- KEY FEATURES:**
- Convertible air therapy surface – use in non-powered mode for pressure injury prevention in at-risk patients, early intervention of existing pressure injuries, or when a step-down surface is needed.
 - Add the simple, digital control unit when powered treatment modalities are necessary.
 - Minimized shearing through Geo-Matt® surface geometry and Shear Transfer Zones® on bi-directional stretch cover
 - Minimized risk of falls and bed entrapment through Safety Edge™ bolster design
 - Flexible treatment modalities, including alternating pressure and powered flotation

THERAPIES:

- Non-powered (Air) Reactive Pressure Redistribution

With control unit added:

- Alternating Pressure
- Powered Pressure Redistribution

Item #	Description	Weight Limit
SC753629	75"L x 36"W x 7"H	500 lbs (400 in powered mode)
SC803629	80"L x 36"W x 7"H	500 lbs (400 in powered mode)
SC843629	84"L x 36"W x 7"H	500 lbs (400 in powered mode)
SC803929	80"L x 39"W x 7"H	500 lbs (400 in powered mode)
SC804229	80"L x 42"W x 7"H	700 lbs (600 in powered mode)
SC844229	84"L x 42"W x 7"H	700 lbs (600 in powered mode)
*SC804829	80"L x 48"W x 7"H	700 lbs (600 in powered mode)
*SC805429	80"L x 54"W x 7"H	700 lbs (600 in powered mode)
SC85V29	Span-Care Convertible for recessed deck frames	700 lbs (600 in powered mode)
5900	Add-on control unit	-


*Allow up to 10 working days lead time for this size.



- Control Unit**
- KEY FEATURES:**
- Easy-to-use, easy-to-connect control unit
 - Push-button therapy selection with adjustable comfort settings
 - Alternating pressure and timed auto-firm
- Warranty:** 18 months

Microclimate Management


Clinicians choose Microclimate Management systems to address maceration, the weakening of the skin due to excess moisture. Savaria’s products are independently documented as industry-best in this key performance area. Remarkably, they do so with none of the compromises providers have long settled for in traditional Microclimate Management surfaces. They won’t lose inflation in a power outage. Won’t “swallow up” the user or collapse at the edge like air-only mattresses can, making transfers safer. Won’t make the user too hot or too cold, either. Ultra-simple to set up, clean and use.



Air Diffusion Matrix

Three-dimensional, non-compressible layer helps ensure uninterrupted air flow in both cover and coverlet. Coupled with a proprietary, dual-delivery powered air control unit, this design ensures an always-open pathway for a constant flow of air beneath the patient, sweeping away moisture vapor before it can reform as a liquid to protect vulnerable skin and tissue.


Featured on: All Microclimate Management models.



Add-on Microclimate Management

The outstanding performance of a complete PressureGuard® Microclimate Management surface, with the convenience and treatment flexibility of an add-on modality. Store it away until needed, then deploy it atop any prevention or treatment healthcare support surface. Convenient setup, maintenance, cleaning, and storage allow for multi-patient use.

Featured on: The Topper (page 35).



Bariatric Treatment

Patients of size benefit significantly from effective microclimate moisture management. The unique prevention and treatment challenges associated with users up to 1,000 lbs make it critical that these products perform as promised. As pioneers in the design and standardized testing of microclimate management technology, we ensure our solutions deliver reliable, proven results.

Featured on: Easy Air XL (page 31).

PressureGuard® Easy Air® XL



QUICK FACTS:

Use For:

- Patients weighing 350–1,000 lbs with excessive skin moisture (do not use for patients weighing less than 350 lbs).
- Patients with macerated skin or limited mobility
- Patients with elevated skin/body temperature due to conditions such as infection or sepsis
- Patients with low Braden Moisture, Activity and Mobility scores (≤ 5)
- Treatment of Stage 1–4 pressure injuries and prevention for high-risk patients

WARRANTY:

Mattress: 3 years | **Cover:** 2 years

KEY FEATURES:

- Bariatric Microclimate Management with Air Diffusion Matrix™ for superior moisture control
- Alternating Pressure
- Powered Pressure Redistribution (“Float” Mode)
- PressureGuard® air/foam design for comfort and support
- Safety Edge™ for added edge-of-bed stability and fall prevention
- Geo-Matt® segmented surface for improved pressure redistribution

THERAPIES:

- Microclimate Management [“Low Air Loss” (LAL)]
- Alternating Pressure
- Powered Pressure Redistribution

Item #	Description
L8042XL-29	Easy Air XL system, 80”L x 42”W x 7”H
L8048XL-29	Easy Air XL system, 80”L x 48”W x 7”H
*L8053XL-29	Easy Air XL System, 80”L x 53”W x 7”H
*L8442XL-29	Easy Air XL System, 84”L x 42”W x 7”H
*L8448XL-29	Easy Air XL System, 84”L x 48”W x 7”H
*L8453XL-29	Easy Air XL System, 84”L x 53”W x 7”H
8310	Easy Air XL Bariatric control unit only

*Allow up to 10 working days lead time for this size.



Control Unit

KEY FEATURES:

- Lightweight and quiet
- Low Pressure indicator light and selectable audible alarm
- Power Failure indicator light and selectable audible alarm
- Simple, indexed comfort adjustment
- Timed (20-minute) “Auto Firm” mode for transfers, ADLs, & CPR
- “Float” mode provides powered flotation therapy
- “Alternate” mode changes loading in 10-minute cycle

Warranty: 2 years

PressureGuard® Easy Air® LR



- QUICK FACTS:**
Use For:
- Lateral Rotation for patients who find frequent manual repositioning painful or need more frequent repositioning
 - Excessive skin moisture due to sweating, decreased mobility, inability or refusal to be repositioned
 - Macerated skin due to any of the above
 - Increased skin or body temperature due to infection, sepsis, or other conditions
 - A combined score <5 on Braden Moisture, Activity and Mobility scales
 - Treatment of Stage 1–4 pressure injuries
 - Prevention of Pressure injuries for high risk patients

Medicare Code: E0277 (Group II)

WARRANTY:
Mattress: 3 years | **Cover:** 2 years

- KEY FEATURES:**
- Air Diffusion Matrix™ technology for continuous microclimate management (True Low Air Loss)
 - Alternating pressure therapy with the additional treatment modality, “Lateral Rotation” (LR).
 - Powered pressure redistribution (“Float”)
 - Safety Edge™ for enhanced stability
 - Geo-Matt® segmented top surface for optimized pressure redistribution

- THERAPIES:**
- Microclimate Management [“Low Air Loss” (LAL)]
 - Alternating Pressure
 - Lateral Rotation
 - Powered Pressure Redistribution

Item #	Description	Weight Limit
L7535LR-29	Easy Air LR System 75”L x 35”W x 7”H	500 lbs (float & alternating modes)
L8035LR-29	Easy Air LR System 80”L x 35”W x 7”H	
L8435LR-29	Easy Air LR System 84”L x 35”W x 7”H	
*L8039LR-29	Easy Air LR System 80”L x 39”W x 7”H	
*L8439LR-29	Easy Air LR System 84”L x 39”W x 7”H	
*92726	Easy Air LR System 80”L x 42”W x 7”H	350 lbs (rotation mode)
*L8053LR-29	Easy Air LR System 80”L x 53”W x 7”H	
*L8453LR-29	Easy Air LR System 84”L x 53”W x 7”H	
L85VLR-29	Easy Air LR for recessed deck frame	-
8210	Easy Air LR control unit only	

*Allow up to 10 working days lead time for this size.



- Control Unit**
- KEY FEATURES:**
- Lightweight and quiet
 - Low Pressure indicator light and selectable audible alarm
 - Power Failure indicator light and selectable audible alarm
 - Simple, indexed comfort adjustment
 - Timed (20-minute) “Auto Firm” mode for transfers, ADLs, & CPR
 - “Float” mode provides powered flotation therapy
 - “Alternate” mode changes loading in 10-minute cycle; controlled at mattress for choice of lateral rotation or alternating pressure therapy
- Warranty:** 2 years

PressureGuard® Easy Air®



- QUICK FACTS:**
Use For:
- Excessive skin moisture due to sweating, decreased mobility, inability or refusal to be repositioned
 - Macerated skin due to any of the above
 - Increased skin or body temperature due to infection, sepsis, or other conditions
 - A combined score <5 on Braden Moisture, Activity and Mobility scales
 - Treatment of Stage 1–4 pressure injuries
 - Prevention of Pressure injuries for high risk patients

Medicare Code: E0277 (Group II)

WARRANTY:
Mattress: 3 years | **Cover:** 2 years

- KEY FEATURES:**
- Air Diffusion Matrix™ technology for continuous microclimate management (True Low Air Loss)
 - Alternating pressure therapy
 - Documented superiority in moisture removal
 - Powered pressure redistribution (“Float”)
 - Safety Edge™ for enhanced stability
 - Geo-Matt® segmented top surface for optimized pressure redistribution

- THERAPIES:**
- Microclimate Management [“Low Air Loss” (LAL)]
 - Alternating Pressure
 - Powered Pressure Redistribution

Item #	Description	Weight Limit
L7535-29	Easy Air System 75”L x 35”W x 7”H	500 lbs
*L7835-29	Easy Air System 78”L x 35”W x 7”H	500 lbs
L8035-29	Easy Air System 80”L x 35”W x 7”H	500 lbs
L8435-29	Easy Air System 84”L x 35”W x 7”H	500 lbs
*L8039-29	Easy Air System 80”L x 39”W x 7”H	500 lbs
*6470029	Easy Air System 80”L x 42”W x 7”H	500 lbs
*81300	Easy Air System 84”L x 42”W x 7”H	500 lbs
8010	Easy Air Control Unit	-

*Allow up to 10 working days lead time for this size.
Complete System includes Mattress, Control Unit, and Coverlet.
To order Coverlet only, add “CLT-” to part number.



- Control Unit**
- KEY FEATURES:**
- Lightweight and quiet
 - Low Pressure indicator light and selectable audible alarm
 - Power Failure indicator light and selectable audible alarm
 - Simple, indexed comfort adjustment
 - Timed (20-minute) “Auto Firm” mode for transfers, ADLs, & CPR
 - “Float” mode provides powered flotation therapy
 - “Alternate” mode changes loading in 10-minute cycle
- Warranty:** 2 years

PressureGuard® Protocol®



- QUICK FACTS:**
- Use For:**
- Patients requiring microclimate management to control moisture and skin heat
 - Prevention and treatment of pressure injuries (Stage 1–4) in high-risk patients
 - Patients with macerated skin due to immobility or moisture exposure

Medicare Code: E0277 (Group II)

WARRANTY:
Mattress: 3 years | **Cover:** 2 years

KEY FEATURES:

- Proprietary microclimate management with Air Diffusion Matrix™ technology to remove excess moisture vapor and support vulnerable skin
- Alternating Pressure therapy
- Powered Pressure Redistribution (“Float” mode)
- Timed Auto Firm mode for added stability during transfers and ADLs
- Geo-Matt® segmented top surface for pressure redistribution and reduced shear
- Safety Edge™ design for enhanced edge-of-bed stability and safer transfers

THERAPIES:

- Microclimate Management [“Low Air Loss” (LAL)]
- Alternating Pressure
- Powered Pressure Redistribution

Protocol Models		
Item #	Description	Weight Limit
P7535-29	Protocol, 75”L x 35”W x 7”H	400 lbs
P8035-29	Protocol, 80”L x 35”W x 7”H	400 lbs
P8435-29	Protocol, 84”L x 35”W x 7”H	400 lbs
23173	Protocol, 80”L x 42”W x 7”H	500 lbs
17108	Protocol, 84”L x 42”W x 7”H	500 lbs

Protocol models with raised perimeter (7”H, 9”H side)		
Item #	Description	Weight Limit
69212	Protocol, 80”L x 35”W with raised perimeter	400 lbs
33905	Protocol, 84”L x 35”W with raised perimeter	400 lbs
72235	Protocol, 80”L x 42”W with raised perimeter	500 lbs
81990	Protocol, 84”L x 42”W with raised perimeter	500 lbs
3400	Protocol control unit only	-

The Topper MicroEnvironment Manager™



- QUICK FACTS:**
- Use For:**
- Excessive perspiration or maceration
 - The Topper is indicated for use in combination with a pressure redistributing surface to aid in the prevention and treatment of skin breakdown and pressure injuries through Stage 4 by reducing excess heat and moisture within the microclimate, especially in patients who are obese, febrile and/or diaphoretic due to infection, antibiotics and/or analgesics

WARRANTY:
Mattress: 2 years | **Cover:** 2 years

KEY FEATURES:

- Proprietary microclimate management with Air Diffusion Matrix™ technology to remove excess moisture vapor and support vulnerable skin
- Helps keep vulnerable skin cool, dry and comfortable
- Works to reduce heat, moisture, shear, pressure, odor
- Convenient setup, maintenance and storage
- Multi-patient use – easily wiped clean and disinfected in place
- Quiet, powerful air flow
- Low profile design attaches like a fitted sheet to any standard pressure-redistributing healthcare support surface

THERAPIES:

- Microclimate Management

Item #	Description
MEM36	The Topper, system (80”L x 36”W)
MEM8436	The Topper, system (84”L x 36”W)
MEM42	The Topper, system (80”L x 42”W)
MEM8442	The Topper, system (84”L x 42”W)
MEM48	The Topper, system (80”L x 48”W)
MEM8448	The Topper, system (84”L x 48”W)
MEM7554	The Topper, system (75”L x 54”W)
MEM8054	The Topper, system (80”L x 54”W)
MEM8060	The Topper, system (80”L x 60”W)
MEM8460	The Topper, system (84”L x 60”W)
2100	Control unit only

Weight limit is equal to the support surface used with The Topper.



Control Unit

KEY FEATURES:

- Maintains therapeutic pressure during power interruption and transport
- Uses less power, generates less noise and heat than typical systems
- Heavy-duty mounting hooks
- Flush-mounted, detachable power cord

Warranty: 18 months



Control Unit

KEY FEATURES:

- Compact, quiet control unit
- Provides gentle, uninterrupted airflow for maximum control of humidity at the patient/support surface interface
- Hangs easily at footboard

Warranty: 2 years

Geo-Mattress®

SPAN Geo-Mattress® series of therapeutic foam mattresses include multiple options for standard of care, raised perimeter, premium, bariatric and home care needs.

All with proven design advantages for pressure redistribution, immersion and envelopment, shear relief, enhanced safety, lasting comfort and reliable durability.



Unique Geo-Matt® Design

KEY FEATURES:

Creates anti-shearing surface with 800 individually responsive cells in three zones.

Ring-of-Air® channels help eliminate heat and moisture build-up.



Heel Slope

KEY FEATURES:

The Heel Slope feature is designed to help reduce pressure on one of the body's most vulnerable areas.

By gently elevating and redistributing weight away from the heels, this design promotes improved pressure management while maintaining proper leg positioning and overall comfort.

Integrated seamlessly into the surface, the Heel Slope works to support circulation and reduce the risk of pressure-related concerns, providing both patients and caregivers with greater confidence in everyday care.

Geo-Mattress® UltraMax™



KEY FEATURES:

- Geo-Matt® surface geometry minimizes shearing
- LifeSpan™ stretch cover with Shear Transfer Zones® to provide protection against micro, macro, and rotational shear
- Helps reduce caregiver injury by “gliding” user back to position following head of bed elevation
- UHP 500 foam provides enhanced seat-section durability, resiliency, and comfort, helping protect against bottoming out under repeated head-of-bed elevation and concentrated patient weight.
- Viscoelastic foam and Heel Slope® provide added protection at vulnerable heels

THERAPIES:

- Non-powered Pressure Redistribution

QUICK FACTS:

Use For:

- Prevention of Stage 1-4 pressure injuries
- Early intervention or multiple Stage 1-2 and treatment of single Stage 3 pressure injuries
- Can be positioned off injury in at least 2 positions
- Reposition frequently by self or caregiver

WARRANTY:

Mattress: 8 years | **Cover:** 2 years

Item #	Description	Weight Limit
UMX7536-29	75”L x 36”W x 6”H	500 lbs
UMX8036-29	80”L x 36”W x 6”H	500 lbs
UMX8436-29	84”L x 36”W x 6”H	500 lbs
UMX85V-29	For recessed deck frames	500 lbs
UMX80CA-29	80” L for CareAssist®+	600 lbs
UMX84CA-29	84” L for CareAssist®+	600 lbs
*UMX8039-29	80”L x 39”W x 6”H	600 lbs
*UMX8042-29	80”L x 42”W x 6”H	600 lbs
*UMX8439-29	84”L x 39”W x 6”H	600 lbs
*UMX8442-29	84”L x 42”W x 6”H	600 lbs

* Allow up to 10 working days lead time for these sizes.

Geo-Mattress® UltraMax™ for Stretchers



KEY FEATURES:

- Geo-Matt® surface geometry minimizes shearing
- Engineered, multi-layer support core
- Extra-soft heel section facilitates immersion to provide added heel protection
- UHP 500 foam provides enhanced seat-section durability, resiliency, and comfort, helping protect against bottoming out under repeated head-of-bed elevation and concentrated patient weight.
- Universal hook and loop pattern to help anchor it to a wide variety of stretcher models

THERAPIES:

- Non-powered Pressure Redistribution

WARRANTY:

Mattress: 3 years | **Cover:** 2 years

All models listed are 5” high to facilitate maximum immersion

Item #	Description	(H)ead-end taper length	(F)oot-end taper length	Weight Limit
47318	75”L x 26”W	7”	0”	700 lbs
49297	75”L x 26”W	13”	4”	700 lbs
97107	76”L x 26”W	10”	4”	700 lbs
10635	76”L x 26”W	12”	5”	700 lbs
2879529	76”L x 27”W	11”	5”	700 lbs
46234	76”L x 29”W	10”	8”	700 lbs
79012	76”L x 30”W	13”	4”	700 lbs
98969	78”L x 27”W	12”	12”	700 lbs
81921	2-pc “Gynie”	n/a	n/a	700 lbs

Other dimensions also available upon request.

Geo-Mattress® Max®



QUICK FACTS:

Use For:

- Prevention of pressure injuries
- General patient population
- Moderate to high risk for skin breakdown
- Decreased mobility
- Repositions frequently by self or caregiver

WARRANTY:

Mattress: 7 years | **Cover:** 2 years

KEY FEATURES:

- Geo-Matt® anti-shearing surface
- Luxurious, visco-elastic high resiliency “memory” foam in foot section enhances the pressure reducing properties of the Heel Slope® feature.
- Firm perimeter bolsters help to provide safe transfers, more stable edge-of-bed sitting, and reduced risk of gaps at siderails.
- Extra-soft, high performance top provides maximum immersion and a premium feel against user’s skin.

THERAPIES:

- Therapeutic Foam Mattress
- Non-powered Pressure Redistribution

Institutional Sizes		
Item #	Description	Weight Limit
MX7535-29	75”L x 35”W x 6”H	400 lbs
*MX7835-29	78”L x 35”W x 6”H	400 lbs
*MX8032-29	80”L x 32”W x 6”H	400 lbs
*MX8035-29	80”L x 35”W x 6”H	400 lbs
*MX8435-29	84”L x 35”W x 6”H	400 lbs

* Allow up to 10 working days lead time for these sizes.

Consumer Sizes		
Item #	Description	Weight Limit
*MX7539-29	75”L x 39”W x 6”H Twin	500 lbs
*MX7554-29	75”L x 54”W x 6”H Full	500 lbs
*MX8060-29	80”L x 60”W x 6”H Queen	500 lbs
*MX8076-29	80”L x 76”W x 6”H King	500 lbs

Geo-Mattress® Plus®



QUICK FACTS:

Use For:

- Prevention of pressure injuries
- General patient population
- Moderate to high risk for skin breakdown
- Decreased mobility
- Repositions frequently by self or caregiver

WARRANTY:

Mattress: 6 years | **Cover:** 2 years

KEY FEATURES:

- Geo-Matt® anti-shearing surface
- Dual layer design for increased patient comfort
- Heel Slope®

THERAPIES:

- Therapeutic Foam Mattress
- Non-powered Pressure Redistribution

Institutional Sizes		
Item #	Description	Weight Limit
PL7535-29	75”L x 35”W x 6”H	400 lbs
*PL7835-29	78”L x 35”W x 6”H	400 lbs
*PL8032-29	80”L x 32”W x 6”H	400 lbs
PL8035-29	80”L x 35”W x 6”H	400 lbs
PL8435-29	84”L x 35”W x 6”H	400 lbs

* Allow up to 10 working days lead time for these sizes.

Consumer Sizes		
Item #	Description	Weight Limit
*PL7539-29	75”L x 39”W x 6”H Twin	300 lbs
*PL7554-29	75”L x 54”W x 6”H Full	300 lbs
*PL8060-29	80”L x 60”W x 6”H Queen	300 lbs
*PL8076-29	80”L x 76”W x 6”H King	300 lbs

Geo-Mattress® Pro



QUICK FACTS:

Use For:

- Prevention of pressure injuries
- General patient population
- Moderate to high risk for skin breakdown
- Decreased mobility
- Repositions frequently by self or caregiver

WARRANTY:

Mattress: 5 years | **Cover:** 2 years

KEY FEATURES:

- Geo-Matt® anti-shearing surface
- Single layer, zoned construction delivers excellent value and effectiveness
- Also available in the Geo-Mattress® Pro RP with added safety of a Raised Perimeter

THERAPIES:

- Therapeutic Foam Mattress
- Non-powered Pressure Redistribution

Institutional Sizes		
Item #	Description	Weight Limit
PR7535-29	75”L x 35”W x 6”H	300 lbs
*PR7835-29	78”L x 35”W x 6”H	300 lbs
*PR8032-29	80”L x 32”W x 6”H	300 lbs
PR8035-29	80”L x 35”W x 6”H	300 lbs
PR8435-29	84”L x 35”W x 6”H	300 lbs

* Allow up to 10 working days lead time for these sizes.

Consumer Sizes		
Item #	Description	Weight Limit
*PR7539-29	75”L x 39”W x 6”H Twin	300 lbs
*PR7554-29	75”L x 54”W x 6”H Full	300 lbs
*PR8060-29	80”L x 60”W x 6”H Queen	300 lbs
*PR8076-29	80”L x 76”W x 6”H King	300 lbs

Geo-Mattress® with Wings®



KEY FEATURES:

- Tissue-friendly Geo-Matt® surface from edge to edge
- Raised perimeter gently urges the user away from the edge of the bed and helps reduce risk of accidental falls
- Lowered, padded openings make transfers easier and safer for user and caregiver
- Heel Slope®

THERAPIES:

- Therapeutic Raised Perimeter Mattress
- Non-powered Pressure Redistribution

Item #	Description	Weight Limit
W7535-29	Geo-Mattress® with Wings, 75”L x 35”W x 6”H center, 8”H side	400 lbs
W7835-29	Geo-Mattress® with Wings, 78”L x 35”W x 6”H center, 8”H side	400 lbs
W8035-29	Geo-Mattress® with Wings, 80”L x 35”W x 6”H center, 8”H side	400 lbs
W8435-29	Geo-Mattress® with Wings, 84”L x 35”W x 6”H center, 8”H side	400 lbs
87536	80”L x 39”W x 6”H center, 8”H side	500 lbs
10312	84”L x 39”W x 6”H center, 8”H side	500 lbs
66301	80”L x 42”W x 6”H center, 8”H side	500 lbs
79243	84”L x 42”W x 6”H center, 8”H side	500 lbs

QUICK FACTS:

Use For:

- Patients at fall risk due to independent exits from bed
- General patient population
- Moderate to high risk for skin breakdown
- Decreased mobility
- Repositions frequently by self or caregiver

WARRANTY:

Mattress: 5 years | **Cover:** 2 years

Geo-Mattress® Atlas®



QUICK FACTS:

Use For:

- Prevention for high risk patients
- Patients 350-600 lbs

WARRANTY:

Mattress: 3 years | **Cover:** 2 years

KEY FEATURES:

- Bariatric Therapeutic Foam Mattress
- Progressive resistance accommodates wide range of patient weights
- Geo-Matt® anti-shearing surface
- Heel Slope®

THERAPIES:

- Non-powered Pressure Redistribution

Item #	Description	Weight Limit
A7535-29	Geo-Mattress® Atlas, 75"L x 35"W x 7"H	600 lbs
A8035-29	Geo-Mattress® Atlas, 80"L x 35"W x 7"H	600 lbs
A8435-29	Geo-Mattress® Atlas, 84"L x 35"W x 7"H	600 lbs
*A8039-29	Geo-Mattress® Atlas, 80"L x 39"W x 7"H	675 lbs
A8042-29	Geo-Mattress® Atlas, 80"L x 42"W x 7"H	750 lbs
A8048-29	Geo-Mattress® Atlas, 80"L x 48"W x 7"H	750 lbs
A8054-29	Geo-Mattress® Atlas, 80"L x 54"W x 7"H	750 lbs
*A8442-29	Geo-Mattress® Atlas, 84"L x 42"W x 7"H	750 lbs
*A8448-29	Geo-Mattress® Atlas, 84"L x 48"W x 7"H	750 lbs
*A8454-29	Geo-Mattress® Atlas, 84"L x 54"W x 7"H	750 lbs

Geo-Mattress® HC



Medicare Code: E0184 (Group 1)

WARRANTY:

Mattress: 2 years

KEY FEATURES:

- Therapeutic Home Care Foam Mattress
- Non-powered Pressure Redistribution
- Specifically designed for pressure injury prevention in the home care setting
- Folded and boxed for easy shipping and storage
- Note: This product does not include Heel Slope® feature

Item #	Description	Weight Limit
SP812-29	Geo-Mattress® HC, 75"L x 35"W x 5"H	250 lbs
SP810-29	Geo-Mattress® HC, 80"L x 35"W x 5"H	250 lbs
SP811-29	Geo-Mattress® HC, 84"L x 35"W x 5"H	250 lbs

Geo-Mattress® 350



Medicare Code: E0184 (Group 1)

WARRANTY:

Mattress: 5 years (2 years prorated)

Cover: 2 years

KEY FEATURES:

- 6" Therapeutic Foam Mattress
- Non-powered Pressure Redistribution
- Segmented top surface, 6" height, and fluid proof cover
- Economical yet effective option for long term care residents

Item #	Description	Weight Limit
67535-29	Geo-Mattress® 350, 75"L x 35"W x 6"H	350 lbs
68035-29	Geo-Mattress® 350, 80"L x 35"W x 6"H	350 lbs
68435-29	Geo-Mattress® 350, 84"L x 35"W x 6"H	350 lbs

Gel-80™



Medicare Code: E0185

WARRANTY:

Mattress: 1 year

KEY FEATURES:

- Gel/Foam Mattress Overlay
- Non-powered Pressure Redistribution
- Effective pressure redistribution from multiple gel bladders
- High-density sculpted support base
- Wipe-clean, fluid-impervious fabric top surface, elastic corner straps, and one-piece construction
- Folded and boxed for easy shipping and storage

Item #	Description	Weight Limit
GEL-80	Gel/foam overlay for hospital-sized twin (80x35) mattress	300 lbs

Specialty Surfaces


Savaria goes beyond standard therapeutic foam and air support surface configurations to equip caregivers with solutions for challenges others simply can't address. Through thoughtful design and engineering, we develop specialty surfaces that expand what's possible in patient care.

The following page highlight four categories of specialty surfaces—each extending our PressureGuard® and Geo-Mattress® designs to meet needs that might otherwise go unmet:

- Bariatric Models (page 45)
- Extra-wide and Consumer Sizes (page 46)
- Recessed Deck models (page 47)
- In-Home Styles (page 48)

Specialty Surfaces: Bariatric Models

Addressing unique needs with dignity, Savaria offers multiple bariatric support surface options to provide specialized care to patients of size. See pages referenced for complete descriptions and available sizes of each.




PressureGuard® Easy Air XL

MicroClimate Management /Alternating Pressure Mattress System.
Weight range: 350 lbs – 1000 lbs
See Page 31 for Details



PressureGuard® Bariatric APM

Alternating Pressure Mattress System.
Weight range: 150 lbs – 750 lbs
See Page 19 for Details



Geo-Mattress® Atlas

Therapeutic Foam Mattress
Weight limit: 600 lbs – 750 lbs
See Page 42 for Details

Specialty Surfaces: Extra-wide and Consumer Sizes

Several models are available in sizes wider than typical healthcare beds. These non-bariatric options provide comfortable, effective use by those who require a wider surface but do not require a bariatric weight capacity.

Options include:

PRESSUREGUARD® MICROCLIMATE MANAGEMENT MODELS:



Easy Air®
Supports up to 500 lbs
Widths: 39", 42" | See page 33



Easy Air® LR
Supports up to 500 lbs
Widths: 39", 42", 53" | See page 32



Protocol®
Supports up to 500 lbs
Widths: 42" | See page 34



Topper®
Widths: 42", 48", 54", 60"
See page 35

PRESSUREGUARD® NON-POWERED AND CONVERTIBLE MODELS:



Custom Care® Convertible
Supports up to 500 lbs
Widths: 42" | See page 27



Span-Care® Convertible
Supports up to 700 lbs
Widths: 42" | See page 29



Custom Care®
Supports up to 500 lbs
Widths: 42" | See page 28

PRESSUREGUARD® POWERED AIR MODELS:



APM2
Supports up to 350 lbs
Widths: 42" | See page 22




APM² Safety Supreme
Supports up to 350 lbs
Widths: 42" | See page 23

GEO-MATTRESS® THERAPEUTIC FOAM MODELS






UltraMax™
Supports up to 600 lbs
Widths: 39", 42" | See page 38




Max®, Plus® and Pro™ models
Supports up to 750 lbs
Widths: 39", 42", 48", 54" | See page 40-41

Specialty Surfaces: Recessed Deck models

Specially designed models with V-shaped bases engineered to accommodate the unique challenges of recessed deck beds. Featuring a modified version of proprietary "Safety Edge" two-piece safety bolster system, adapted to overcome the design restrictions inherent in recessed decks.



RECESSED VERSIONS ARE AVAILABLE FOR SIX LEADING PRESSUREGUARD® AND GEO-MATTRESS® MODELS.



Item #	Description	See Page
UMX85V-29	Geo-Mattress® UltraMax for recessed deck	38
NP85V29	PressureGuard® Custom Care for recessed deck	28
SC85V29	PressureGuard® Span-Care Convertible for recessed deck	29
CJ85V29	PressureGuard® Custom Care Convertible for recessed deck	27
CL85V29	PressureGuard® Custom Care Convertible LAL for recessed deck	26
L85VLR-29	PressureGuard® Easy Air LR for recessed deck	32

Specialty Surfaces: In-Home Styles

Extra-wide, consumer-sized versions of PressureGuard® Span-Care® convertible, with options to accommodate one user or two. Entire surface is non-powered therapy when add-on unit is not connected. Order “Center” or “Right/Left” versions depending on user’s desired location of the alternating pressure therapy when the powered control unit is connected.



“Right/Left” models:

Ports on either side allow connection of powered add-on control unit for alternating pressure on either side as desired.

DESIGNED FOR: Two users, where therapy can be added to the desired side. The remaining side maintains non-powered pressure redistribution, or a single user who cannot/prefers not to be positioned on the center of the mattress.



“Center” models:

Four alternating pressure cylinders are located at the center of the mattress, while the outside cylinders maintain non-powered pressure redistribution.

DESIGNED FOR: One user who prefers to sleep in center of mattress, with or without add-on powered control unit.


Item #	Description	Bed Size
SC7554-RL	Span Care 75" x 54" Homestyle Right/Left	Full
SC7554-CTR	Span Care 75" x 54" Homestyle Center	Full
SC8054-RL	Span Care 80" x 54" Homestyle Right/Left	Full XL
SC8054-CTR	Span Care 80" x 54" Homestyle Center	Full XL
SC8060-RL	Span Care 80" x 60" Homestyle Right/Left	Queen
SC8060-CTR	Span Care 80" x 60" Homestyle Center	Queen
SC8038-HSL	Span Care 80" x 38" Homestyle Left	Twin XL / Dual King
SC8038-HSR	Span Care 80" x 38" Homestyle Right	Twin XL / Dual King

Overlays & Seating

Savaria offers a range of overlays and seating designed to enhance comfort, support, and safety across care environments. From pressure-relief surfaces to ergonomically designed seating, each solution is built to meet diverse needs while prioritizing ease of use and durability. Our offerings provide flexibility and reliability, giving caregivers and users the confidence that comfort and support are always within reach.



Geo-Matt® Therapeutic Overlays



QUICK FACTS:

Use For:

- Pressure injury prevention in mod-high risk patients



KEY FEATURES:

- Therapeutic foam overlay
- Segmented surface redistributes pressure
- Reduces shearing and the build-up of heat and moisture
- Multiple sizes and cover options available
- Proven solution for clinical care settings

Standard & Specialty Sizes			
Item #	Description	Qty/Cs	Weight Limit
SP584-000	Geo-Matt Overlay without Soft Skin®	4	250 lbs
50960-580	Geo-Matt Overlay with Soft Skin® sleeve (fitted)	1	250 lbs
50960-581	Geo-Matt Overlay with Soft Skin® sleeve (torso)	4	250 lbs
50960-582	Geo-Matt Overlay with FR Soft Skin® sleeve (tubular)	1	250 lbs
50960-583	Geo-Matt Overlay with FR Soft Skin® sleeve (tubular)	4	250 lbs
50960-584	Geo-Matt Overlay with Soft Skin® sleeve (fitted)	4	250 lbs
50960-585	Geo-Matt Overlay, 80" with Soft Skin® sleeve	4	250 lbs
50960-455	Geo-Matt Overlay, 72" with Nylon Coverlet	4	250 lbs
50960-458	Geo-Matt Overlay, 80" with Nylon Coverlet	4	250 lbs
50960-586	Geo-Matt O.R. Size (20") Overlay with Soft Skin®	4	250 lbs
SP653-000	Geo-Matt Overlay with FR Soft Skin® sleeve	1	250 lbs
SP654-000	Geo-Matt Overlay with FR Soft Skin® sleeve	4	250 lbs

Standard size: 33" x 72" x 3.5"

Consumer Sizes (no Soft Skin®)			
Item #	Description	Qty/Cs	Weight Limit
SP760-000	Full (Double) Size 52" x 72"	1	250 lbs
SP761-000	Queen Size 58" x 75"	1	250 lbs
SP762-000	King Size (2 pads) 74" x 75"	1	250 lbs


RISK LEVEL

LOW

MODERATE

HIGH

Geo-Matt Therapeutic Overlay with Cover



KEY FEATURES:

- Vinyl/Nylon cover with built-in antibacterial agent
- Side zipper with flap provides protection from fluids
- Durable side elastic bands help secure overlay to mattress
- Easy to clean and maintain

Geo-Matt Therapeutic Overlay with Cover			
Item #	Description	Qty/Cs	Weight Limit
50960-456	Geo-Matt with Vinyl/Nylon Cover, 72"	4	250 lbs
50960-457	Geo-Matt with Vinyl/Nylon Cover, 72"	1	250 lbs
50960-459	Geo-Matt with Vinyl/Nylon Cover, 80"	4	250 lbs
50960-500	Geo-Matt with Vinyl/Nylon Cover, 80"	1	250 lbs

Geo-Matt Overlay Coverlets, Zippered Covers and Soft Skin Sleeves			
Item #	Description	Qty/Cs	Weight Limit
5000	Geo-Matt Nylon Fitted Coverlet, 72"	4	250 lbs
5001	Geo-Matt Nylon Fitted Coverlet, 80"	4	250 lbs
5006	Geo-Matt Vinyl/Nylon Zippered Cover, 72"	4	250 lbs
5009	Geo-Matt Vinyl/Nylon Zippered Cover, 72"	1	250 lbs
5007	Geo-Matt Vinyl/Nylon Zippered Cover, 80"	4	250 lbs
5010	Geo-Matt Vinyl/Nylon Zippered Cover, 80"	1	250 lbs
50960-561	Soft Skin sleeve, 72" (fitted)	12	250 lbs
50960-563	Fire Retardant Soft Skin sleeve, 72" (fitted)	12	250 lbs
SP759-000	Fire Retardant Soft Skin sleeve, 72" (tubular)	12	250 lbs
50960-590	Soft Skin sleeve, 80"	12	250 lbs


RISK LEVEL

LOW

MODERATE

HIGH

Geo-Matt® Pediatric Overlays



QUICK FACTS:

Use For:

- Skin management in infants and children



KEY FEATURES:

- Therapeutic foam overlay ideal for pediatric settings
- Segmented surface redistributes pressure

Pediatric & Neonatal Sizes			
Item #	Description	Qty/Cs	Weight Limit
50960-587	Geo-Matt Pediatric Overlay (28" x 56.5" x 2")	4	70 lbs
SP706-000	Geo-Matt Isolette Overlay (21" x 26" x 2")	4	70 lbs
SP729-000	Geo-Matt Infant Overlay (18" x 25" x 2")	4	70 lbs
SP730-000	Geo-Matt Isolette Overlay (21" x 26" x 2") w/ nylon cover	4	70 lbs
SP732-000	Geo-Matt ICU Bassinet Overlay (13" x 26" x 2")	4	70 lbs
SP769-000	Geo-Matt ICU Overlay (12" x 22.5" x 2")	4	70 lbs

All pediatric overlays (except SP730-000) include Fire Retardant Soft Skin®.

Delta Pediatric Overlays



QUICK FACTS:

Use For:

- Skin management in infants and children

KEY FEATURES:

- Economical foam overlay ideal for pediatric settings
- Proven effectiveness
- High conformance to body contours
- Protection in home or hospital

Delta Foam			
Item #	Description	Qty/Cs	Weight Limit
DELTA-C11	Overlay 32"x75"x3.5"	1	70 lbs
DELTA-120	Overlay 34"x75"x3.5"	1	70 lbs
DELTA-200	Isolette Matt Overlay 13.5"x25"x1.75"	1	70 lbs
DELTA-202	Warmer Overlay 29"x23.75"x2.25"	1	70 lbs
DELTA-204	Crib Overlay 26.5"x53"x2.75"	1	70 lbs

All Delta foam overlays include Soft Skin.

Pressure Kair™ Preventative Overlay



QUICK FACTS:


Use For:

- Pressure injury prevention in low-mod risk patients


KEY FEATURES:

- Preventative foam overlay with individually cut support cells
- Divided into distinctive head, torso, and foot zones
- Economical preventative solution

Pressure Kair		
Item #	Description	Qty/Cs
SP789-000	Pressure Kair™ 34" × 72" × 3.5" with Soft Skin®	4
SP889-000	Pressure Kair™ 34" × 80" × 3.5" with Soft Skin®	4



Convoluted Overlays



QUICK FACTS:

Use For:

- Comfort

KEY FEATURES:


- High quality polyurethane foam
- Available in premium and standard density
- Multiple thickness options for user needs
- Lightweight, economical comfort solution
- Compression rolled for efficient storage

Convoluted Overlays		
Item #	Description	Qty/Cs
SP22R-000	2" Convoluted Pads, Premium Density	12
SP22S-000	2" Convoluted Pads, Standard Density	12
SP32R-000	3" Convoluted Pads, Premium Density	6
SP35S-000	3" Convoluted Pads, Standard Density	6
SP45R-000	4" Convoluted Pads, Premium Density	4
SP45S-000	4" Convoluted Pads, Standard Density	4

All pads are 72"L x 34"W. Convoluted, or "eggcrate", overlays are recognized for comfort, not pressure management.

Caution: Convoluted, or "eggcrate", overlays are recognized for comfort, not pressure management. Call your Span-America representative when selecting the pressure management support surface most effective for your patients' needs.

The Equalizer® Positioning Cushion



QUICK FACTS:

Use For:

- Pain Relief
- Prevention of breakdown
- Assistance with positioning

Warranty: 2 years

Medicare Code: Skin Protection & Positioning – E2607

For sizes 22" or wider, use E2608.



KEY FEATURES:

- Multi-component design disperses pressure
- Viscoelastic memory foam top enhances comfort
- Contoured base supports leg positioning
- Relief zones protect bony prominences

Standard Sizes			
Item #	Description	Qty/Cs	Weight Limit
EQ1616	Equalizer cushion, 16"W × 16"L x 3.5"H	1	250 lbs
EQ1616-04	Equalizer cushion, 16"W × 16"L x 3.5"H	4	250 lbs
EQ1816	Equalizer cushion, 18"W × 16"L x 3.5"H	1	250 lbs
EQ1816-04	Equalizer cushion, 18"W × 16"L x 3.5"H	4	250 lbs
EQ1818	Equalizer cushion, 18"W × 18"L x 3.5"H	1	250 lbs
EQ1818-04	Equalizer cushion, 18"W × 18"L x 3.5"H	4	250 lbs
EQ2018	Equalizer cushion, 20"W × 18"L x 3.5"H	1	250 lbs
EQ2018-04	Equalizer cushion, 20"W × 18"Lx 3.5"H	4	250 lbs

Bariatric Sizes			
Item #	Description	Qty/Cs	Weight Limit
BEQ2218	Equalizer cushion, 22"W × 18"L x 3.5"H	1	500 lbs
BEQ2218-04	Equalizer cushion, 22"W × 18"L x 3.5"H	4	500 lbs
BEQ2418	Equalizer cushion, 24"W × 18"L x 3.5"H	1	500 lbs
BEQ2418-04	Equalizer cushion, 24"W × 18"L x 3.5"H	4	500 lbs

RISK LEVEL	LOW	MODERATE	HIGH
------------	-----	----------	------

EZ-Dish® Positioning Cushion



QUICK FACTS:

Use For:

- Pain Relief
- Prevention of breakdown
- Assistance with positioning

Warranty: 18 months

Medicare Code: Skin Protection- E2603


KEY FEATURES:

- Contoured positioning cushion designed to reduce pressure on the ischials and coccyx while supporting upright posture.
- Viscoelastic memory foam top molds to anatomy
- Pocket design reduces ischial pressure
- Anti-slip cover enhances stability

EZ-Dish® Positioning Cushion			
Item #	Description	Qty/Cs	Weight Limit
EZ1616-96	EZ-Dish w/ anti-slip cover, 16"W × 16"L x 3"H	6	250 lbs
EZ1616-96-1	EZ-Dish w/ anti-slip cover, 16"W × 16"L x 3"H	1	250 lbs
EZ1816-96	EZ-Dish w/ anti-slip cover, 18"W × 16"L x 3"H	6	250 lbs
EZ1816-96-1	EZ-Dish w/ anti-slip cover, 18"W × 16"L x 3"H	1	250 lbs
EZ1818-96	EZ-Dish w/ anti-slip cover, 18"W × 18"L x 3"H	6	250 lbs
EZ1818-96-1	EZ-Dish w/ anti-slip cover, 18"W × 18"L x 3"H	1	250 lbs
EZ2018-96	EZ-Dish w/ anti-slip cover, 20"W × 18"L x 3"H	6	250 lbs
EZ2018-96-1	EZ-Dish w/ anti-slip cover, 20"W × 18"L x 3"H	1	250 lbs

RISK LEVEL	LOW	MODERATE	HIGH
------------	-----	----------	------

Geo-Matt® Contour



QUICK FACTS:

Use For:

- Comfort
- Prevention of breakdown
- Assistance with positioning

Warranty: 18 months

Medicare Code: Skin Protection- E2603


KEY FEATURES:

- Anti-shearing Geo-Matt® top layer
- Contoured base enhances positioning
- Leg channels support alignment
- Pelvic bucketing helps reduce slouching

Geo-Matt®Contour			
Item #	Description	Qty/Cs	Weight Limit
262R	Geo-Matt Contour 16"W × 16"L x 3"H, durafilm coating, deluxe cover	1	250 lbs
262W	Geo-Matt Contour 18"W × 16"L x 3"H, durafilm coating, deluxe cover	1	250 lbs
GMC1616-96	Geo-Matt Contour, 16"W × 16"L x 3"H with anti-slip cover	6	250 lbs
GMC1616-96-1	Geo-Matt Contour, 16"W × 16"L x 3"H with anti-slip cover	1	250 lbs
GMC1816-96	Geo-Matt Contour, 18"W × 16"L x 3"H with anti-slip cover	6	250 lbs
GMC1816-96-1	Geo-Matt Contour, 18"W × 16"L x 3"H with anti-slip cover	1	250 lbs
GMC1818-96	Geo-Matt Contour, 18"W × 18"L x 3"H with anti-slip cover	6	250 lbs
GMC1818-96-1	Geo-Matt Contour, 18"W × 18"L x 3"H with anti-slip cover	1	250 lbs
GMC2018-96	Geo-Matt Contour, 20"W × 18"L x 3"H with anti-slip cover	6	250 lbs
GMC2018-96-1	Geo-Matt Contour, 20"W × 18"x 3"H L with anti-slip cover	1	250 lbs

RISK LEVEL	LOW	MODERATE	HIGH
------------	-----	----------	------

Geo-Matt® Wedge Cushion



QUICK FACTS:

Use For:

- Comfort
- Prevention of breakdown
- Assistance with positioning

Warranty: 18 months

Medicare Code: General Use-E2601

KEY FEATURES:

- High-density foam wedge design
- Tapered profile (4" front to 2" back) helps reduce sliding
- Geo-Matt® surface supports skin integrity
- Anti-slip cover and non skid bottom for secure placement

Geo-Matt® Wedge Cushion			
Item #	Description	Qty/Cs	Weight Limit
GMW1616-96	Geo-Matt Wedge w/ anti-slip cover, 16"W × 16"L	6	250 lbs
GMW1616-96-1	Geo-Matt Wedge w/ anti-slip cover, 16"W × 16"L	1	250 lbs
GMW1816-96	Geo-Matt Wedge w/ anti-slip cover, 18"W × 16"L	6	250 lbs
GMW1816-96-1	Geo-Matt Wedge w/ anti-slip cover, 18"W × 16"L	1	250 lbs
GMW1818-96	Geo-Matt Wedge w/ anti-slip cover, 18"W × 18"L	6	250 lbs
GMW1818-96-1	Geo-Matt Wedge w/ anti-slip cover, 18"W × 18"L	1	250 lbs
GMW2018-96	Geo-Matt Wedge w/ anti-slip cover, 20"W × 18"L	6	250 lbs
GMW2018-96-1	Geo-Matt Wedge w/ anti-slip cover, 20"W × 18"L	1	250 lbs

RISK LEVEL	LOW	MODERATE	HIGH
------------	-----	----------	------

Geo-Wave® Recliner Cushion



QUICK FACTS:

Use For:

- Comfort
- Minimizing shear and pressure

Warranty: 18 months

KEY FEATURES:

- Engineered, segmented fins that reduce reclining-induced shearing
- 3 zoned orientations to reduce shear in multiple directions
- Designed to also reduce pressure when user is in upright or reclined position
- Elastic hook and loop straps help anchor cushion to the chair

Geo-Wave®			
Item #	Description	Qty/Cs	Weight Limit
WAVE-04	Geo-Wave seat cushion w/cover (18"W x 68"L)	4	300 lbs
WAVE-01	Geo-Wave seat cushion w/cover (18"W x 68"L)	1	300 lbs
WAVE22-04	Geo-Wave seat cushion w/cover (22"W x 68"L)	4	300 lbs
WAVE22-01	Geo-Wave seat cushion w/cover (22"W x 68"L)	1	300 lbs
WAVE28-04	Geo-Wave seat cushion w/cover (28"W x 68"L)	4	500 lbs
WAVE28-01	Geo-Wave seat cushion w/cover (28"W x 68"L)	1	500 lbs

Gel-T®



QUICK FACTS:

Use For:

- Comfort
- Prevention of breakdown

Warranty: 18 months

Medicare Code: General Use-E2601.
For sizes 22" or wider, use E2602.

KEY FEATURES:

- Segmented gel bladder for pressure relief
- Dual foam layers cradle gel insert
- Wipe-clean cover reduces shear and friction
- Non-skid bottom for secure placement

Gel-T®			
Item #	Description	Qty/Cs	Weight Limit
8015-05	Gel-T, 15"W × 15"L x 2.5"H	1	250 lbs
8016-05	Gel-T, 16"W × 16"L x 2.5"H	1	250 lbs
801618-05	Gel-T, 16"W × 18"L x 2.5"H	1	250 lbs
8017-05	Gel-T, 17"W × 17"L x 2.5"H	1	250 lbs
801816-05	Gel-T, 18"W × 16"L x 2.5"H	1	250 lbs
8018-05	Gel-T, 18"W × 18"L x 2.5"H	1	350 lbs
8019-05	Gel-T, 19"W × 19"L x 2.5"H (fits Geri-chairs)	1	350 lbs
802018-05	Gel-T, 20"W × 18"L x 2.5"H (fits Geri-chairs)	1	350 lbs
8021-05	Gel-T, 21"W × 21"L x 2.5"H (fits Geri-chairs)	1	350 lbs

Bariatric Sizes			
Item #	Description	Qty/Cs	Weight Limit
B802218-05	Gel-T, 22"W × 18"L x 2.5"H	1	500 lbs
B802418-05	Gel-T, 24"W × 18"L x 2.5"H	1	500 lbs

RISK LEVEL	LOW	MODERATE	HIGH
------------	-----	----------	------

Geo-Matt® PRT (3” High Density)



QUICK FACTS:

Use For:

- Comfort
- Prevention of breakdown

Warranty: 18 months

Medicare Code: Skin Protection- E2603 (only versions with fabric cover). For sizes 22” or wider, use E2604.

KEY FEATURES:

- Lower, firmer base for more active wheelchair users
- Ring-of-Air® design helps manage heat and moisture
- Geo-Matt design redistributes pressure

Geo-Matt® PRT (3” high density)			
Item #	Description	Qty/Cs	Weight Limit
50746-070	Geo-Matt PRT, 18”W x 16”L, w/ Soft Skin® sleeve	6	250 lbs
50746-083	Geo-Matt PRT, 18”W x 16”L, w/ vinyl/nylon cover	6	250 lbs
PRT1616-96	Geo-Matt PRT, 16”W x 16”L, w/ anti-slip cover	6	250 lbs
PRT1616-96-1	Geo-Matt PRT, 16”W x 16”L, w/ anti-slip cover	1	250 lbs
PRT1816-96	Geo-Matt PRT, 18”W x 16”L, w/ anti-slip cover	6	250 lbs
PRT1816-96-1	Geo-Matt PRT, 18”W x 16”L, w/ anti-slip cover	1	250 lbs
PRT1818-96	Geo-Matt PRT, 18”W x 18”L, w/ anti-slip cover	6	250 lbs
PRT1818-96-1	Geo-Matt PRT, 18”W x 18”L, w/ anti-slip cover	1	250 lbs
PRT2018-96	Geo-Matt PRT, 20”W x 18”L, w/ anti-slip cover	6	250 lbs
PRT2018-96-1	Geo-Matt PRT, 20”W x 18”L, w/ anti-slip cover	1	250 lbs

Bariatric Sizes (4” high)			
Item #	Description	Qty/Cs	Weight Limit
BPRT2218-96	Geo-Matt PRT, 22”W x 18”L, w/anti-slip cover	6	500 lbs
BPRT2218-96-1	Geo-Matt PRT, 22”W x 18”L w/anti-slip cover	1	500 lbs
BPRT2418-96	Geo-Matt PRT, 24”W x 18”L w/anti-slip cover	6	500 lbs
BPRT2418-96-1	Geo-Matt PRT, 24”W x 18”L w/anti-slip cover	1	500 lbs

Geo-Matt® (4” Medium Density)



QUICK FACTS:

Use For:

- Comfort
- Prevention of breakdown

Warranty: 18 months

Medicare Code: Skin Protection- E2603 (only versions with fabric cover).


KEY FEATURES:

- Responsive cells reduce shear and cradle bony prominences
- Ring-of-Air® design helps manage heat and moisture
- Proven geometry enhances comfort and pressure relief
- 4” medium-density construction for everyday support

Geo-Matt® (4” medium density)			
Item #	Description	Qty/Cs	Weight Limit
50746-040	Geo-Matt, cushion only	6	200 lbs
50746-050	Geo-Matt with Soft Skin® sleeve	6	200 lbs
50746-051	Geo-Matt with fire-retardant Soft Skin® sleeve	6	200 lbs
50746-084	Geo-Matt with anti-slip cover	6	200 lbs

All models: 18”W x 16”L x 4”H

Cross Cut Cushion



QUICK FACTS:

Use For:

- Comfort

Warranty: 1 year


Medicare Code: General Use- E2601

KEY FEATURES:

- Economical prevention and comfort
- Segmented top surface provides pressure dispersion and comfort
- Cover provides anti-slip top and non-skid bottom

Cross Cut Cushion			
Item #	Description	Qty/Cs	Weight Limit
CC1616-96	Cross-cut Cushion w/ anti-slip cover, 16”W x 16”L x 3”H	6	200 lbs
CC1816-96	Cross-cut Cushion w/ anti-slip cover, 18”W x 16”L x 3”H	6	200 lbs
CC1818-96	Cross-cut Cushion w/ anti-slip cover, 18”W x 18”L x 3”H	6	200 lbs
CC2018-96	Cross-cut Cushion w/ anti-slip cover, 20”W x 18”L x 3”H	6	200 lbs

Convolved Comfort Cushions



QUICK FACTS:

Use For:

- Comfort

KEY FEATURES:

- High quality polyurethane foam
- Available in premium and standard density
- Multiple thickness options for user needs
- Lightweight, economical comfort solution

Convolved Comfort Cushions		
Item #	Description	Qty/Cs
50746-010	18” x 16” x 4”P (.75” base) Premium Density	12
50077-109	16” x 16” x 3”P (.5” base) Standard Density	12
50077-112	16” x 16” x 3”P (1.5” base) Standard Density	12
50746-025	18” x 16” x 4”P (.75” base) Standard Density	12
SP663-000	18” x 16” x 4”P (1.5” base) Standard Density	12

Note: Convolved cushions provide comfort only and are not intended for patients at risk of pressure ulcers.

Positioners

Positioners are an essential component of comprehensive patient care, offering targeted support and stability that enhance comfort, safety, and clinical outcomes. Designed in collaboration with healthcare professionals, our range of positioning aids—including wedges, rolls, pads, and specialty shapes—helps caregivers maintain proper anatomical alignment, facilitate movement, and reduce unnecessary pressure on vulnerable areas.

Whether used for pressure redistribution, post-operative support, or everyday positioning needs, these tools are engineered to work intuitively with therapeutic surfaces and support systems throughout the care setting. Easy to employ and adaptable to a wide range of body sites and applications, Savaria positioners help create a more comfortable and supportive care experience for patients and a more efficient workflow for caregivers



Heel Manager™ Offloading Heel Elevator



KEY FEATURES:

- Fully offloads heels with skin-friendly design as heels are the second most common pressure injury site
- Extreme comfort and support for both legs
- Shear Transfer Zones™ cover technology
- Geo-Matt® layer with Ring-of-Air® design helps to draw heat and moisture away from the skin
- Wipe-clean, anti-microbial cover with non-slip bottom
- **Warranty:** 1 year

CLINICAL NOTE - Heel Elevation and Foot Drop:

For patients at risk for, or who currently have foot drop, a support product must be used to hold the foot in the correct position. See our line of optional positioners.

Item #	Description	Qty/Case	Weight Limit
44900-1	Heel Manager	1	300 lbs
44900	Heel Manager	4	300 lbs

All components latex-free.



Cradle Boot – Original



KEY FEATURES:

Floats heel, supporting foot in neutral position, prevents external leg rotation.

Item #	Qty/Case
55555-02	*6

*Indicates items vacuum packed to save storage. All items conform to California Technical bulletin 117.

Limb Elevator




KEY FEATURES:

Elevates casted or uncasted leg. Can be used in both supine and side lying positions.

Item #	Qty/Case
50632-508	*6
50632-308	2

*Indicates items vacuum packed to save storage. All items conform to California Technical bulletin 117.

Flared Leg Extension



KEY FEATURES:

Elevates leg and knee in extended position.

Item #	Qty/Case
SP291-000	2

Economy Abduction Pillow



KEY FEATURES:

Adjustable Velfoam straps.

Item #	Description	Qty/Case
50642-320	SML (17x12x4x7)	2
50642-340	MED (22x14x4x7)	2
50642-360	LRG (25x16x7x7)	2

Abduction Pillow with Contour and Comfort Band



KEY FEATURES:

Adjustable Velfoam straps, contoured leg channels and a cotton comfort band.

Item #	Description	Qty/Case
SP122-000	SML (17x12x4x7)	6
SP119-000	MED (22x14x4x7)	6
SP197-000	WIDE (20x19x10x7)	2
SP131-000	LRG (25x16x7x7)	6

*Indicates items vacuum packed to save storage.

Abduction Pillow with Comfort Band




KEY FEATURES:

Adjustable Velfoam straps with a cotton comfort band.

Item #	Description	Qty/Case
50642-820	*SML (17x12x4x7)	6
50642-840	*MED (22x15x5x7)	6
50642-860	*LRG (25x16x7x7)	6

Abduction Pillow with Contour



KEY FEATURES:


Adjustable Velfoam straps and contoured leg channels.

Item #	Description	Qty/Case
SP357-000	*SML (17x12x4x7)	6
SP358-000	*MED (22x14x4x7)	6
SP359-000	*LRG (25x16x7x7)	6
SP204-000	SML (17x12x4x7)	2
SP205-000	MED (22x14x4x7)	2
SP206-000	LRG (25x16x7x7)	2

*Indicates items vacuum packed to save storage

Size Guide — If patient is: less than 5’ tall use SML | 5’1” to 5’11” use MED | 6’ to 6’5” use LRG

Laminectomy Arm Cradle Set




KEY FEATURES:
For use in prone position. Restricts rotation of arm, cushions axilla and provides firm support for entire limb.

Item #	Qty/Case
50670-527	*12 pairs

*Indicates items vacuum packed to save storage. All items conform to California Technical bulletin 117.


O.R. Table Pillow



KEY FEATURES:
Used during surgery in supine or lateral position to provide support for sensitive areas of the head. Segmented squares offer maximum pressure dispersion, especially during prolonged surgical procedures.

Item #	Qty/Case
50654-519	*10


Axilla Cuff



KEY FEATURES:
Supports and protects axilla and upper arm during surgery.

Item #	Qty/Case
SP010-000	12

Thoracotomy Arm Sling – Standard




KEY FEATURES:
Used with an ether screen bar, it positions hand and wrist for access to thoracic cavity without occluding IV. Adjustable hand strap allows multiple wrist positions.

Item #	Qty/Case
50668-526	*12 pairs

*Indicates items vacuum packed to save storage. All items conform to California Technical bulletin 117.


O.R. Armboard Pads



KEY FEATURES:
Comfortably supports entire arm during surgery when used on O.R. table armboards.

Item #	Qty/Case
50676-378	12 pairs

Geo-Matt® 30° Body Wedge




KEY FEATURES:
Segmented surface provides pressure dispersion while supporting trunk at a 30° angle. Also used to support arm or leg in side lying, or at foot of bed to prevent sliding in supine or gatched position.
(7 3/8" x 11-5/8" x 20")

Item #	Qty/Case
50960-300*	6
50960-301	2 (includes fire retardant Soft Skin®)
50960-300-96	6 (includes vinyl cover)

*Indicates items vacuum packed to save storage.
All items conform to California Technical bulletin 117.

45° Body Wedge




KEY FEATURES:
Supports the entire length of the torso.

Item #	Qty/Case
50630-507*	Standard 6/cs (9" x 10-1/2" x 24")
50630-307	Economy 2/cs (8" x 9-1/2" x 21")


*Indicates items vacuum packed to save storage.
All items conform to California Technical bulletin 117.

Support Surface Algorithm and Selection Guide



Seating Tissue Load Management Algorithm

Optimal Seating Device Selection Guide



SPAN PressureGuard® Air Therapy Surfaces Matrix

		Non-powered Air Therapy				Powered Air					Microclimate Management ["Low Air Loss" (LAL)]				Bariatric (Powered Air)	
		Custom Care®	Span-Care® Convertible	Custom Care® Convertible	Custom Care® Convertible LAL	APM	APM2™	APM-i	Span-Care® Convertible w/control unit	Custom Care® Convertible w/control unit	Protocol®	Easy Air®	Easy Air LR®	Custom Care® Convertible LAL w/control unit	Bariatric APM	Easy Air® XL
Therapy Modalities	Non-powered pressure redistribution	✓	✓	✓	✓											
	Convertible to powered via add-on control unit		✓	✓	✓											
	Alternating Pressure					✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
	Lateral Rotation						✓	✓		✓			✓	✓		
	Microclimate Management (LAL)									✓	✓	✓	✓	✓	✓	
	Immersion (Powered Flotation)					✓	✓	✓	✓		✓	✓	✓	✓	✓	✓
Frame Types	Versions to accommodate recessed deck frames	✓	✓	✓	✓				✓	✓			✓	✓		
Control Unit Features	SmartImmerse™ Air System Setting							✓								
	Alarm: (V)isual, (A)udible					V	V, A	V, A	V, A	V, A	V, A	V	V	V	V, A	✓
	Timed Auto Firm						✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
Mattress Construction	Raised Perimeter										Available					
	Geo-Matt® top surface	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
	Safety Edge™ bolster design	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
	Can maintain air in power outage; (A) = Automatically	✓ (A)	✓ (A)	✓ (A)	✓ (A)	✓	✓	✓	✓	✓ (A)	✓	✓ (A)	✓ (A)	✓ (A)	✓	✓ (A)
	Heel Slope™; (F) = Foam-only heel section	✓ (F)	✓ (F)	✓ (F)	✓ (F)	✓	✓	✓	✓ (F)	✓ (F)	✓	✓	✓	✓ (F)	✓	✓
	100% Latex free, 100% DHEP-free	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
Cover	Non-stretch nylon with fluid-repellent coating					✓									✓	
	Bi-directional stretch with Shear Transfer Zone™ design	✓	✓	✓			✓	✓	✓	✓						
	Zippered air flow coverlet with Air Diffusion Matrix™ layer				✓						✓	✓	✓	✓		✓
General	Medicare Code	E0373	E0373	E0373	E0373	E0277	E0277	E0277	E0277	E0277	E0277	E0277	E0277	E0277	E0277	E0277
	Therapeutic Weight Limit (lbs.)	500	500	500	500	350	350	350	400	350	400	500	500 (350 LR)	500 (350 LR)	750	1000
	Warranty: Mattress, non-pro-rated	5 yrs.	5 yrs.	5 yrs.	5 yrs.	18 mos.	3 yrs.	3 yrs.	5 yrs.	5 yrs.	3 yrs.	3 yrs.	3 yrs.	5 yrs.	3 yrs.	3 yrs.
	Warranty: Control unit, non-pro-rated					18 mos.	18 mos.	18 mos.	18 mos.	2 yrs.	18 mos.	2 yrs.	2 yrs.	2 yrs.	18 mos.	2 yrs.
	Warranty: Cover/Coverlet, non-prorated	2 yrs.	2 yrs.	2 yrs.	2 yrs.	18 mos.	2 yrs.	2 yrs.	2 yrs.	2 yrs.	2 yrs.	2 yrs.	2 yrs.	2 yrs.	2 yrs.	2 yrs.

SPAN Seating Products Matrix		Treatment w/ Positioning	Skin Protection				General Use			
		Equalizer®	Ez Dish®	Geo-Matt® Contour	Geo-Matt® PRT	Geo-Matt®	Geo-Matt® Wedge	Gel-T®	Cross Cut	Geo-Wave®
Level of Care	General Use			✓			✓	✓	✓	✓
	Skin Protection	✓	✓	✓	✓	✓				
	Positioning Feature	✓	✓	✓			✓			✓
Core	Geo-Matt®			✓	✓	✓	✓			✓
	Minimizes Shearing			✓	✓	✓	✓			
	Gel Components	✓						✓		
Cover	Non-Slip Bottom		✓	✓	✓	✓	✓	✓	✓	✓
	Two-way Stretch Fabric Top	✓								
	Wipe Clean		✓	✓	✓	✓	✓	✓	✓	✓
General	Medicare Code *Bariatric Code	E2607	E2603	E2603	E2603 E2604*	E2603	E2601	E2601 E2602*	E2601	
	Therapeutic Weight Limit (lbs.)	250	250	250	250	250	250	250	250	300
	Bariatric Weight Limit (lbs.)	500			500			500		
	Warranty	2 yrs.	18 mos.	18 mos.	18 mos.	18 mos.	18 mos.	18 mos.	1 yr.	18 mos.
	Seat and Back Cushion									✓
	Height of Cushion (Inches)	3.5"	3"	3"	3" standard 4" bariatric	4"	4" front 2" back	2.5"	3"	3"
	Viscoelastic Top	✓	✓							



Learn more about these products here:
www.savariapatientcare.com

