

AVERO Comfort Tub



**Seamlessly integrates
into wellness and bathing
areas with a homelike feel**

Height adjustable
for staff and
resident comfort

Freely accessible
by standard
bath lifts

Customizable
with a wide
range of wellness
features

science, safety, simplicity

 **savaria®**

AVERO Comfort Tub

Height adjustable care tub



Easy operation

Wired hand control with intuitive symbols for all bath functions.

Optimal Resident Support

- Generous, ergonomic inner tub body for optimal support.
- Large recessed armrests for added comfort.
- Ergonomic chrome grab handles ensure secure and comfortable bathing.




Hygienic and Safe Design

- The non-porous surface facilitates easy cleaning and maintenance, ensuring optimal hygiene at all times.
- Scalding protection for hand shower and tub spout.



Exclusive wellness options



Colored-Light System

Create a relaxing and visually stimulating atmosphere with customizable colors and an automatic sequence program.



Aroma System

Enhance sensory experiences with various scents and warm air from the Air Spa System. Easy to operate and maintenance-free.



Music-Sound System

Provide soothing and stimulating sound experience during baths. Connect via Bluetooth or USB for personalized music.



Air Spa System

Deliver soothing and calming effects for ultimate relaxation. The jets are specially protected to prevent bacterial contamination.

Create a welcoming environment

The tub can be rotated and positioned anywhere in the room, without compromising the aesthetics of your facility's bathing area. This offers an unobstructed view for the resident while bathing, greatly enhancing their wellbeing.

TECHNICAL SPECIFICATION

Technical Data Dimensions (L x W)	2.100 mm / 800 mm
Inside Dimensions (L x W x D)	1.770 mm x 575 mm x 400 mm
Lifting Height	680 mm - 1.080 mm SWL 400 kg
Power connection	230 V, 50 Hz
Filling Volume	ca. 270 l

PRODUCT INFORMATION:

- Electrical Height Adjustment
- Easy operation via hand control
- Ergonomical inner design
- Armrests and grab bars for maximum comfort

AVAILABLE OPTIONS:

- Aroma system
- Colored-light system
- Air-Relax system
- Music-sound system
- Bluetooth-music-sound system
- RAL colors available
- Automatic tub filling
- Emergency function in case of power failure
- Bathtub shortener
- Change of viewing direction
- Legionella flushing



Get more information

800.888.6752

savariapatientcare.com



**Span-America
Handicare USA**

**800.888.6752
866.891.6502**

**PO Box 5231 Greenville SC 29606
10888 Metro Court St. Louis MO 63043**

Span-America is a division of Savaria Corporation, a global leader in accessibility. Learn more about Savaria at savaria.com.

PL Beka Avero Comfort Tub CC-4631 r0