



GO BEYOND!

Who we are: passionate innovators pushing the boundaries of what is possible. **What we do:** combine leading-edge engineering with innovative materials to create movement. **What makes us proud:** having our ultralight wheelchairs echo the love of life their users demonstrate each day.

Through our unique approach and distinct attitude, our aim is to enhance lives, nourish a sense of freedom and revolutionize the way people see wheelchairs.

LEARN MORE To learn and find more information regarding our products, options and accessories, visit our website at motioncomposites.com



T 1 (866) 650-6555

F 1 (888) 966-6555

motioncomposites.com

0070838A-00-V15-OPTIONS ACCESSORIES-ENGLISH-24-08-2026

MANUAL
WHEELCHAIR

OPTIONS AND ACCESSORIES GUIDE





NO MATTER WHAT YOU NEED, WE'VE GOT YOU COVERED

Our manual wheelchair options and accessories guide lets you explore the numerous choices available to tailor your wheelchair to your specific needs, while also offering elegant enhancements that truly reflects who you are.

GO BEYOND



INDEX

06	07	10	12	13	14	15	16	18	19	21	22	23	24	26	28
CASTERS	FOOTRESTS/ FOOTPLATES	BACK CANES	BACK UPHOLSTERY	ONE ARM DRIVE	REAR WHEELS	TIRES	HANDRIMS	AXLES	WHEEL LOCKS	ARMRESTS	SIDE GUARDS	POSITIONING BELTS	ACCESSORIES	VISION PACKAGE	FRAME COLORS



CASTERS

CA01 →
Composite wheel w/PU tire (1" width)



CA10 →
Pneumatic casters 6" x 1¼"



CA11 →
Pneumatic casters 8" x 1¼"

CA14 to CA16 →
Newton UltraCaster Soft-Roll composite
4", 5", 6" x 1½"



CA17 to CA19 →
Newton UltraCaster Soft-Roll aluminum
4", 5", 6" x 1½"



OPTION

CA13 →
Carbon fiber Frog Legs
suspension forks – Phase 2



FOOTRESTS / FOOTPLATES

RIGID WHEELCHAIR FOOTRESTS



FP12 ↑
Tubular open loop



FP13 ↑
Tubular with ABS footplate



FP31 ↑
Adjustable angle –
Aluminum (one piece)



FP32 ↑
Adjustable angle – Carbon fiber with
rubber inserts (one piece)



Flip back position



FP15 ↑
Flip back, adjustable angle – Carbon with rubber inserts (one piece)
FP14 Flip back, adjustable angle – Aluminum (one piece)

FR26 ↑
High mount with adjustable angle footplate – Carbon fiber with rubber inserts (one piece)
FR25 High mount with adjustable angle footplate – Aluminum (one piece)

FOLDING WHEELCHAIR FOOTRESTS



FR03 to FR08, FR15 ↑
Swing-in/swing-out footrest
FR01, FR02 Swing away footrest



FR09 ↑
High mount (footplate attached to the hanger with clamps)



FR10 ↑
Elevating legrest with calf pad, 14"-22" length adjustable



FR27 to FR30 ↑
Residual limb support (flat cushion)



FR32, FR33 ↑
Residual limb support (contoured)

FOOTPLATES



FP01 ↑
Standard – Composite



FP02 ↑
Newton adjustable angle – Aluminum




FP03 ↑
Newton Adjustable Angle - Composite



FP04 ↑
Newton adjustable angle – Carbon fiber with rubber inserts




FP23 ↑
Newton adjustable angle extra deep – Aluminum



FP06, FP07 ↑
One piece adjustable angle flip-up

Folded position



FP28, FP29 ↑
One piece auto-folding adjustable angle

Folded Position

FOOTREST/FOOTPLATE ACCESSORIES



FP16 ↑
2" riser



FP08 ↑
Length adjustable heel loop



FP09 ↑
Calf strap



FP30 ↑
Bodypoint Aeromesh padded calf strap

BACK CANES

RIGID WHEELCHAIRS

FOLD DOWN BACKRESTS



BC20 ↑
Locking backrest



BC21 ↑
Ultralight non locking backrest
(not locking in upright or
down position)



BC56 ↑
Depth adjustable back
(provides 3" of depth adjustment)

FOLDING WHEELCHAIRS



BC09 ↑
Depth adjustable back



BC15 ↑
Angle adjustable backrest option

BACK CANE OPTIONS

BACKREST OPTIONS



← BC57
Standard 3" deep rigidizer bar –
Carbon




← BC57
Standard 3" deep rigidizer bar – Aluminum
BC55 4" deep rigidizer bar – Aluminum

PUSH HANDLE OPTIONS



Fold down position

PH05 ↑
Fold down push handles



Clamp detail

PH06 ↑
Clamp on height adjustable push handle

STROLLER HANDLE AND STABILIZER OPTIONS



← PH04
Stroller handle
Not available with the HD option




Folded position

← PH03
Newton folding stabilizer bar

For the complete back canes overview, please refer to the compatibility table on page 28.




BACK UPHOLSTERY



Front

Back


BH05 ↑
Slip-on Nylon back upholstery



Front

Back


BH06 ↑
Adjustable tension Nylon back upholstery



Front

Back


BH07 ↑
Slip-on Active-Mesh back upholstery



Front

Back


BH03 ↑
Slip-on Syntech-Suede back upholstery



Front

Back


BH04 ↑
Adjustable tension Syntech-Suede back upholstery



Front

Back

BH08 ↑
Adjustable tension Active-Mesh back upholstery



ONE ARM DRIVE



HR22 ↑
One arm drive 24"



REAR WHEELS



RW01 ↑
Mag



 RW02 ↑
Newton One – Spoke wheel



 RW03 ↑
Newton Gravity – Ultralight wheel



RW04 ↑
Spinergy Spox*



RW05 ↑
Spinergy LX*



RW07 ↑
Spinergy CLX Carbon Blade*

* Different colors are available for the Spinergy rear wheels.



TIRES




TI01 ↑
Soft Urethane 1 3/8"
Low Tread, black



TI02 ↑
Soft Urethane 1 3/8"
Medium Tread,
black



TI13 ↑
Soft Urethane 1"
No tread, black -
Sentinel



TI08 Pneumatic with
Airless insert 1 3/8"
Medium Tread
TI05 ↑
Pneumatic 1 3/8"
Medium Tread



TI06 ↑
TI06 - Pneumatic
HP 1" - Low Tread,
Puncture resistant,
black - SpeedLite



TI07 ↑
Pneumatic HP 1"
- Medium Tread,
Puncture resistant,
black - TrailBlazer




TI03 ↑
Pneumatic 2,1"
Large Tread,
black



HANDRIMS

🔗 HR07 to HR18 →
Newton Air Grip:
Black, Sunkissed Orange,
Burgundy, Sapphire Blue,
Monster Green, Charcoal,
Fuchsia, Steel Blue,
Ferrari Red, Acid Green,
Electric Purple, Caribbean Blue,
White, Pink Pearl, Plum Purple,
Sky Blue, Black Cherry,
Khaki Green, Desert Sand



HR01 ↑
Aluminum anodized



HR61 ↑
Black aluminum
hard anodized



HR02 ↑
Plastic coated



HR03 ↑
High friction coated



HR05 ↑
Surge LT
HR28 Surge



HR04 ↑
Natural Fit (shown with optional
thumb standard grip (HR23))



HR62 ↑
SIMI H - 1 3/8" oval handrim with
high friction strip



HR63 ↑
NOVA H - 1 1/4" oval handrim
with high friction strip



HR64 ↑
OPTIMUM H - 1 1/8" heart-
shaped handrim



AXLES



AX01 ↑
Quick release axle



AX02 ↑
Permanent axle



AX03 ↑
Quad release axle



AX04 ↑
Amputee adjustable axle plate



WHEEL LOCKS



WL08
 Newton wheel lock – push to lock



WL02
 Pull to lock



WL04
 Newton grade aid push to lock (suggested with pneumatic or high pressure tire)



WL13
 6" removable extension handle for grade-aid and pull to lock



WL14
 6" removable extension handle for Newton push to lock



WL05
 Scissor lock




WL15
 Unilateral wheel lock – pull to lock



WL06
 Attendant lock (additional set of wheel lock mounted to the back)



AR02
 Flip back armrest shown with 14" full length arm pad



AR04
 "T" armrest shown with 14" full length arm pad



AR14
 "T" armrest shown with 14" full length arm pad



AR14-610
 "T" armrest with 14" full length arm pad (6" - 10")



AR05
 Swing-away armrest with tubular arm pad



AR06
 Swing-away armrest with 10" desk length arm pad




AR07
 Locking flip-up armrest

* Most of these armrests are available with the 10" desk length arm pad and with 14" full length arm pad. Please check the compatibility table on page 33-34.

SIDE GUARDS




SG01 ↑
Plastic side guards
(removable)



SG03 ↑
Carbon fiber side guards
(removable)




SG06 ↑
Plastic fender side guards
(removable)




SG04 ↑
Carbon fiber fender side guards
(removable)



SG02 ↑
Carbon fiber side guards
(fixed)



SG ↑
Plastic fender side guards
straight bracket (removable)



SG ↑
Plastic fender side guards
offset bracket (removable)

POSITIONING BELTS



SB01 ↑
Velcro adjustable



SB02 ↑
Auto buckle



SB05 ↑
BodyPoint Hip Belt - 2 point padded



SB06 ↑
BodyPoint Hip Belt - 4 point padded



SB08 ↑
BodyPoint EVOFLEX® Positioning Belt Rehab latch and Band Clamps

ACCESSORIES



AC03 ↑
Anti tippers



AC27 ↑
Tip assist



AC06 ↑
Cane carrier




AC07 ↑
Oxygen tank holder – D & E type (left)



AC08 ↑
Spoke guards (pair)




AC09 ↑
Footrest neoprene impact guard (pair) – shown in red



AC29 ↑
Swing away anti-tipper



AC10 ↑
Seat pouch




AC30 ↑
Newton Multi purpose pouch



AC26 ↑
BodyPoint mobility bag



AC04 ↑
Transit tie down (unoccupied)



AC23 ↑
Newton backpack




AC13 ↑
Tool kit



AC16 ↑
Motion Composites polo shirt



AC41 ↑
Tray



VISION PACKAGE



AR02 ↑
Height adjustable flip back armrest



FR01-FR02 ↑
Swing-away footrest



WL14 ↑
6" removable extension handle for Newton push to lock



AR14 ↑
Height adjustable "T" armrest, 14" full length arm pad



FR03 ↑
Swing-in/swing-out footrest



FP03 ↑
Newton adjustable angle – Composite



FR03
Swing-in/swing-out footrest - HELIO




FR01-FR02
Swing away footrest - MOVE



FRAME COLORS



Black or High Gloss Carbon Fiber (if carbon fiber frame)



*The above frame colors are available for all wheelchair models.

*Color representation may vary from one monitor/printer to another. Therefore, the actual paint color on your wheelchair might differ slightly from what you see on your screen or your brochure.

Custom color
Motion Composites can match any desired color.

ACCENT COLORS

* Available on Apex and Veloce only.

Black



Sunkissed Orange



Fuchsia



Ferrari Red



Acid Green



White



Sky Blue



COMPATIBILITY TABLE

■ AVAILABLE ■ NOT AVAILABLE

		APEX	VELOCE	HELIO C2	HELIO A7	HELIO A6	MOVE ¹	HELIO KIDS
CASTERS								
CA01	Composite wheel w/PU tire 3", 4", 5", 6", 7", 8" x 1"	■	■	■	■	■	■	■
CA10	Pneumatic casters 6" x 1¼"	■	■	■	■	■	■	■
CA11	Pneumatic casters 8" x 1¼"	■	■	■	■	■	■	■
Ⓜ CA14 to CA16	Newton UltraCaster Soft-Roll composite 4", 5", 6" x 1½"	■	■	■	■	■	■	■
Ⓜ CA17 to CA19	Newton UltraCaster Soft-Roll aluminum 4", 5", 6" x 1½"	■	■	■	■	■	■	■
CASTER OPTIONS								
CA12	Caster pin lock	■	■	■	■	■	■	■
CA13	Carbon fiber Frog Legs suspension forks Phase 2	■	■	■	■	■	■	■
FOOTRESTS/FOOTPLATES								
Rigid wheelchair – Footrests								
FP12	Tubular open loop	■	■	■	■	■	■	■
FP13	Tubular with ABS footplate	■	■	■	■	■	■	■
FP31	Adjustable angle – Aluminum (one piece)	■	■	■	■	■	■	■
FP32	Adjustable angle – Carbon fiber with rubber inserts (one piece)	■	■	■	■	■	■	■
FP14	Flip back, adjustable angle – Aluminum (one piece)	■	■	■	■	■	■	■
FP15	Flip back, adjustable angle – Carbon fiber with rubber inserts (one piece)	■	■	■	■	■	■	■
FR25	High mount with adjustable angle footplate – Aluminum (one piece)	■	■	■	■	■	■	■
FR26	High mount with adjustable angle footplate – Carbon fiber with rubber inserts (one piece)	■	■	■	■	■	■	■

¹ Not available in USA.
Note: The availability of options and accessories may vary depending on your region. Refer to the order form.

		APEX	VELOCE	HELIO C2	HELIO A7	HELIO A6	MOVE ¹	HELIO KIDS
Folding wheelchair – Footrests								
FR01, FR02	Swing-away footrest	■	■	■	■	■	■	■
FR03 to FR08, FR15 Swing-in/swing-out footrest		■	■	■	■	■	■	■
FR09	High mount (footplate attached to the hanger with clamps)	■	■	■	■	■	■	■
FR10	Elevating legrest with calf pad, 14"-22" length adjustable	■	■	■	■	■	■	■
FR27 to FR30	Residual limb support (flat cushion)	■	■	■	■	■	■	■
FR32, FR33	Residual limb support (contour cushion)	■	■	■	■	■	■	■

Footplates

FP01	Standard – Composite	■	■	■	■	■	■	■
🔊 FP02	Newton adjustable angle – Composite	■	■	■	■	■	■	■
🔊 FP03	Newton Adjustable Angle - Composite	■	■	■	■	■	■	■
🔊 FP04	Newton adjustable angle – Carbon fiber with rubber inserts	■	■	■	■	■	■	■
🔊 FP23	Newton adjustable angle extra deep – Aluminum	■	■	■	■	■	■	■
FP06, FP07	One piece adjustable angle flip-up – Aluminum, carbon fiber with rubber inserts	■	■	■	■	■	■	■
FP28, FP29	One piece auto-folding adjustable angle – Aluminum, carbon fiber with rubber inserts	■	■	■	■	■	■	■

FOOTPLATE ACCESSORIES

FP16	2" riser	■	■	■	■	■	■	■
FP08	Length adjustable heel loop	■	■	■	■	■	■	■
FP09	Calf strap	■	■	■	■	■	■	■
FP30	Bodypoint Aeromesh padded calf strap	■	■	■	■	■	■	■

BACK CANES

Fixed angle back – with integrated fixed height push handles

BC58	Straight aluminum back cane	■	■	■	■	■	■	■
BC01	8° bend aluminum back cane	■	■	■	■	■	■	■

Adjustable angle back – with integrated fixed height push handles

BC02	8° bend aluminum back cane	■	■	■	■	■	■	■
------	----------------------------	---	---	---	---	---	---	---

Fixed angle back – fixed height without option

BC03	Straight carbon fiber back cane	■	■	■	■	■	■	■
BC04	8° bend carbon fiber back cane	■	■	■	■	■	■	■

		APEX	VELOCE	HELIO C2	HELIO A7	HELIO A6	MOVE ¹	HELIO KIDS
Fixed angle back – height adjustable with options								
BC06	Straight carbon fiber back cane height adjustable	■	■	■	■	■	■	■
BC08	8° bend carbon fiber back cane height adjustable	■	■	■	■	■	■	■
BC11	Straight aluminum back cane height adjustable	■	■	■	■	■	■	■
BC12	8° bend aluminum back cane height adjustable	■	■	■	■	■	■	■

Adjustable angle back – height adjustable with options

BC13	Adjustable angle, straight aluminum back cane	■	■	■	■	■	■	■
BC14	Adjustable angle, 8° bend aluminum back cane	■	■	■	■	■	■	■
BC09	Angle adjustable, straight carbon fiber back cane	■	■	■	■	■	■	■
BC10	Angle adjustable, 8° bend carbon fiber back cane	■	■	■	■	■	■	■

Adjustable back – depth and angle adjustable with options

BC15	Straight aluminum back cane	■	■	■	■	■	■	■
BC16	8° bend aluminum back cane	■	■	■	■	■	■	■

Fold down backrest

BC20	Locking backrest	■	■	■	■	■	■	■
BC21	Ultralight non locking backrest (not locking in upright and down position)	■	■	■	■	■	■	■

Adjustable angle fold down backrest, fixed adjustable without options

BC22	Straight back canes	■	■	■	■	■	■	■
BC30	Fixed height without push handle	■	■	■	■	■	■	■

BACK CANE OPTIONS

Backrest options

BC57	Standard 3" deep rigidizer bar – Carbon or aluminum	■	■	■	■	■	■	■
BC55	4" deep aluminum rigidizer bar	■	■	■	■	■	■	■
BC56	Depth adjustable back (provides 3" of depth adjustment)	■	■	■	■	■	■	■

Push handle options

PH01	Push handles	■	■	■	■	■	■	■
PH05	Fold down push handles	■	■	■	■	■	■	■
PH06	Clamp on height adjustable push handles	■	■	■	■	■	■	■

¹ Not available in USA.
Note: The availability of options and accessories may vary depending on your region. Refer to the order form.

		APEX	VELOCE	HELIO C2	HELIO A7	HELIO A6	MOVE ¹	HELIO KIDS
Stroller handle and stabilizer bar								
	PH04	Stroller handle (not available with the HD option)						
🔓	PH03	Newton folding stabilizer bar						

BACK UPHOLSTERY

	BH05	Slip-on Nylon back upholstery						
	BH06	Adjustable tension Nylon back upholstery						
	BH03	Slip-on Syntech-Suede back upholstery						
	BH04	Adjustable tension Syntech-Suede back upholstery						
	BH07	Slip-on Active-Mesh back upholstery						
	BH08	Adjustable tension Active-Mesh back upholstery						

ONE ARM DRIVE

	HR22	One arm drive 24"						
--	------	-------------------	--	--	--	--	--	--

REAR WHEELS

	RW01	Mag						
🔓	RW02	Newton One – spoke wheel						
🔓	RW03	Newton Gravity – ultralight wheel						
	RW04	Spinergy Spox						
	RW05	Spinergy LX						
	RW07	Spinergy CLX Carbon Blade						

TIRES

	TI01	Soft Urethane 1 3/8" Low Tread, black						
	TI02	Soft Urethane 1 3/8" Medium Tread, black						
	TI13	Soft Urethane 1" No tread, black - Sentinel						
	TI05	Pneumatic 1 3/8" Medium Tread						
	TI08	Pneumatic with Airless insert 13/8" Medium Tread						
	TI06	Pneumatic HP 1" - Low Tread, Puncture resistant, black - SpeedLite						
	TI07	Pneumatic HP 1" - Medium Tread, Puncture resistant, black - TrailBlazer						

		APEX	VELOCE	HELIO C2	HELIO A7	HELIO A6	MOVE ¹	HELIO KIDS
HANDRIMS								
	HR01	Aluminum anodized						
🔓	HR07 to HR18	Newton Air Grip <i>Black, Sunkissed Orange, Burgundy, Sapphire Blue, Monster Green, Charcoal, Fuchsia, Steel Blue, Ferrari Red, Acid Green, Electric Purple, Caribbean Blue, White, Pink Pearl, Plum Purple, Sky Blue, Black Cherry, Khaki Green, Desert Sand</i>						

	HR02	Plastic coated						
	HR03	High friction coated						
	HR28	Surge						
	HR05	Surge LT						
	HR04	Natural Fit (shown with optional thumb standard grip HR23)						
	HR61	Black aluminum hard anodized						
	HR62	SIMI H - 1 3/8" oval handrim with high friction strip						
	HR63	NOVA H - 1 1/4" oval handrim with high friction strip						
	HR64	OPTIMUM H - 1 1/8" heart-shaped handrim						

AXLES

	AX01	Quick release axle						
	AX02	Permanent axle						
	AX03	Quad release axle						
	AX04	Amputee adjustable axle plate						

WHEEL LOCKS

🔓	WL08	Newton wheel lock – push to lock						
	WL02	Pull to lock						
🔓	WL04	Newton grade aid push to lock						
		(suggested with pneumatic or high pressure tire)						
	WL13	6" removable extension handle for grade-aid and pull to lock						
	WL14	6" removable extension handle for Newton push to lock						
	WL05	Scissor Lock						
	WL06	Attendant Lock (additional set of wheel lock mounted to the back)						

¹ Not available in USA.
Note: The availability of options and accessories may vary depending on your region. Refer to the order form.

		APEX	VELOCE	HELIO C2	HELIO A7	HELIO A6	MOVE ¹	HELIO KIDS
ARMRESTS								
AR01	Height adjustable flip back armrest (convertible to "T" armrest)	■	■	■	■	■ ²	■	■
	10" desk length arm pad							
AR02	Height adjustable flip back armrest (convertible to "T" armrest)	■	■	■	■	■ ²	■	■
	14" full length arm pad							
AR03/AR13	Height adjustable "T" armrest, 10" desk length arm pad	■	■	■	■	■	■	■
AR04/AR14	Height adjustable "T" armrest, 14" full length arm pad	■	■	■	■	■	■	■
AR05	Height adjustable swing-away armrest, tubular arm pad	■	■	■	■	■	■	■
AR06	Height adjustable swing-away armrest, 10" desk length arm pad	■	■	■	■	■	■	■
AR07	Locking flip-up armrest	■	■	■	■	■	■	■
AR12	14" full length armrest gel pad	■	■	■	■	■	■	■
SIDE GUARDS								
SG01	Plastic side guards (removable)	■	■	■	■	■	■	■
SG02	Carbon fiber side guards (fixed)	■	■	■	■	■	■	■
SG03	Carbon fiber side guards (removable)	■	■	■	■	■	■	■
SG06	Plastic fender side guards (removable)	■	■	■	■	■	■	■
SG04	Carbon fiber fender side guards (removable)	■	■	■	■	■	■	■
SG07	Side guards clip	■	■	■	■	■	■	■
POSITIONING BELTS								
SB01	Velcro adjustable	■	■	■	■	■ ²	■	■
SB02	Auto buckle	■	■	■	■	■	■	■
SB04	Bodypoint Hip Belt – 2 point not padded	■	■	■	■	■ ²	■	■
SB05	BodyPoint Hip Belt – 2 point padded	■	■	■	■	■	■	■
SB06	BodyPoint Hip Belt – 4 point padded	■	■	■	■	■	■	■
SB08	BodyPoint EVOFLEX® Positioning Belt Rehab latch and Band Clamps	■	■	■	■	■	■	■
HEADRESTS								
AC01	Profiled headrest	■	■	■	■	■	■	■
AC28	Headrest mounting hardware with stabilizer bar	■	■	■	■	■	■	■

		APEX	VELOCE	HELIO C2	HELIO A7	HELIO A6	MOVE ¹	HELIO KIDS
ACCESSORIES								
AC03	Anti tippers	■	■	■	■	■	■	■
AC29	Swing away anti-tipper	■	■	■	■	■	■	■
AC27	Tip assist	■	■	■	■	■	■	■
AC04	Transit tie down (unoccupied)	■	■	■	■	■	■	■
AC05	WC-19 transit tie down w/Q'straint belt	■	■	■	■	■	■	■
AC06	Cane carrier	■	■	■	■	■	■	■
AC07	Oxygen tank holder – D & E type (left)	■	■	■	■	■	■	■
AC08	Spoke guards (pair)	■	■	■	■	■	■	■
AC09	Footrest neoprene impact guard (pair) – Available with various APEX A colour frames	■	■	■	■	■	■	■
AC10	Seat pouch	■	■	■	■	■	■	■
AC30	Newton multi-purpose pouch	■	■	■	■	■	■	■
AC26	BodyPoint mobility bag	■	■	■	■	■	■	■
AC12	Velcro closing strap	■	■	■	■	■	■	■
AC41	Tray	■	■	■	■	■	■	■

VISION PACKAGE								
VP01	Vision package	■	■	■	■	■	■	■

FRAME COLORS								
FC01	Black or carbon fiber black (if carbon fiber frame)	■	■	■	■	■	■	■
FC03	Sunkissed Orange, Burgundy, Sapphire Blue, Monster Green, Charcoal, Fuchsia, Steel Blue, Ferrari Red, Acid Green, Electric Purple, Caribbean Blue, White, Pink Pearl, Plum Purple, Sky Blue, Black Cherry, Khaki Green, Desert Sand (not available on Apex A)	■	■	■	■	■	■	■
FC10	Highlight color Ferrari Red	■	■	■	■	■	■	■
FC11	Highlight color White	■	■	■	■	■	■	■
FC12	Highlight color Black	■	■	■	■	■	■	■
FC13	Highlight color Sunkissed Orange	■	■	■	■	■	■	■
FC14	Highlight color Sky Blue	■	■	■	■	■	■	■
FC15	Highlight color Fuchsia	■	■	■	■	■	■	■
FC16	Highlight color Acid Green	■	■	■	■	■	■	■
FC99	Custom color	■	■	■	■	■	■	■

¹ Not available in USA.
Note: The availability of options and accessories may vary depending on your region. Refer to the order form.

		APEX	VELOCE	HELIO C2	HELIO A7	HELIO A6	MOVE ¹	HELIO KIDS
ARMRESTS								
AR01	Height adjustable flip back armrest (convertible to "T" armrest)	■	■	■	■	■ ²	■	■
	10" desk length arm pad							
AR02	Height adjustable flip back armrest (convertible to "T" armrest)	■	■	■	■	■ ²	■	■
	14" full length arm pad							
AR03/AR13	Height adjustable "T" armrest, 10" desk length arm pad	■	■	■	■	■	■	■
AR04/AR14	Height adjustable "T" armrest, 14" full length arm pad	■	■	■	■	■	■	■
AR05	Height adjustable swing-away armrest, tubular arm pad	■	■	■	■	■	■	■
AR06	Height adjustable swing-away armrest, 10" desk length arm pad	■	■	■	■	■	■	■
AR07	Locking flip-up armrest	■	■	■	■	■	■	■
AR12	14" full length armrest gel pad	■	■	■	■	■	■	■
SIDE GUARDS								
SG01	Plastic side guards (removable)	■	■	■	■	■	■	■
SG02	Carbon fiber side guards (fixed)	■	■	■	■	■	■	■
SG03	Carbon fiber side guards (removable)	■	■	■	■	■	■	■
SG06	Plastic fender side guards (removable)	■	■	■	■	■	■	■
SG04	Carbon fiber fender side guards (removable)	■	■	■	■	■	■	■
SG07	Side guards clip	■	■	■	■	■	■	■
POSITIONING BELTS								
SB01	Velcro adjustable	■	■	■	■	■ ²	■	■
SB02	Auto buckle	■	■	■	■	■	■	■
SB04	Bodypoint Hip Belt – 2 point not padded	■	■	■	■	■ ²	■	■
SB05	BodyPoint Hip Belt – 2 point padded	■	■	■	■	■	■	■
SB06	BodyPoint Hip Belt – 4 point padded	■	■	■	■	■	■	■
SB08	BodyPoint EVOFLEX® Positioning Belt Rehab latch and Band Clamps	■	■	■	■	■	■	■
HEADRESTS								
AC01	Profiled headrest	■	■	■	■	■	■	■
AC28	Headrest mounting hardware with stabilizer bar	■	■	■	■	■	■	■

		APEX	VELOCE	HELIO C2	HELIO A7	HELIO A6	MOVE ¹	HELIO KIDS
ACCESSORIES								
AC03	Anti tippers	■	■	■	■	■	■	■
AC29	Swing away anti-tipper	■	■	■	■	■	■	■
AC27	Tip assist	■	■	■	■	■	■	■
AC04	Transit tie down (unoccupied)	■	■	■	■	■	■	■
AC05	WC-19 transit tie down w/Q'straint belt	■	■	■	■	■	■	■
AC06	Cane carrier	■	■	■	■	■	■	■
AC07	Oxygen tank holder – D & E type (left)	■	■	■	■	■	■	■
AC08	Spoke guards (pair)	■	■	■	■	■	■	■
AC09	Footrest neoprene impact guard (pair) – Available with various APEX A colour frames	■	■	■	■	■	■	■
AC10	Seat pouch	■	■	■	■	■	■	■
AC30	Newton multi-purpose pouch	■	■	■	■	■	■	■
AC26	BodyPoint mobility bag	■	■	■	■	■	■	■
AC12	Velcro closing strap	■	■	■	■	■	■	■
AC41	Tray	■	■	■	■	■	■	■

VISION PACKAGE								
VP01	Vision package	■	■	■	■	■	■	■

FRAME COLORS								
FC01	Black or carbon fiber black (if carbon fiber frame)	■	■	■	■	■	■	■
FC03	Sunkissed Orange, Burgundy, Sapphire Blue, Monster Green, Charcoal, Fuchsia, Steel Blue, Ferrari Red, Acid Green, Electric Purple, Caribbean Blue, White, Pink Pearl, Plum Purple, Sky Blue, Black Cherry, Khaki Green, Desert Sand (not available on Apex A)	■	■	■	■	■	■	■
FC10	Highlight color Ferrari Red	■	■	■	■	■	■	■
FC11	Highlight color White	■	■	■	■	■	■	■
FC12	Highlight color Black	■	■	■	■	■	■	■
FC13	Highlight color Sunkissed Orange	■	■	■	■	■	■	■
FC14	Highlight color Sky Blue	■	■	■	■	■	■	■
FC15	Highlight color Fuchsia	■	■	■	■	■	■	■
FC16	Highlight color Acid Green	■	■	■	■	■	■	■
FC99	Custom color	■	■	■	■	■	■	■

¹ Not available in USA.
Note: The availability of options and accessories may vary depending on your region. Refer to the order form.