

DASH3

FLIGHTWEIGHT
BY TRAVEL BUGGY



Less
weight
More
value



Upgraded
comfort
seating



Caregiver
Remote Pairing
Optional Add-On



OUR LIGHTEST, SMARTEST DASH YET

Meet the
DASH3

The DASH has always been a fan favourite, and now it's better than ever. Introducing the DASH 3: **lighter, smarter, and ready to go where life takes you.**

KEY FEATURES

4 LBS LIGHTER. BUILT STRONGER.

Featuring magnesium components for improved durability and performance.

ENHANCED PRECISION & CONTROL

Smoother, more responsive control with an upgraded joystick and CPU.



A FRESH NEW LOOK

Sleek dark grey.
Sharp performance.



ROLLS LIKE A SUITCASE



POWER YOU CAN FEEL, NOT HEAR



BUILT-IN STORAGE



WIDER, CONTOURED ARMRESTS



STRONGER, EXTENDABLE FOOTPLATE



YOUR DRINK, WITHIN REACH

SPECIFICATIONS

MODEL SKU	16MNA0060100
FRAME COLOUR	DARK GREY METALLIC
NET WEIGHT	42 lbs
GROSS WEIGHT	46 lbs
WEIGHT CAPACITY	265 lbs
DRIVING RANGE	15 km (9.5 mi)
INCLINE CAPABILITY	10°
TURNING RADIUS	24"
MAXIMUM SPEED	6 km/h (3.7 mph)
MOTOR POWER	2 x 170 W Brushless
UNFOLDED DIMENSIONS	38"L x 24"W x 37"H
FOLDED DIMENSIONS	30"L x 24"W x 13.5"H
FRONT / REAR WHEEL DIAMETER	8" / 12"

SEATING DIMENSIONS

SEAT SIZE	17.5"W x 17"D x 2"H
BACKREST HEIGHT	19"
BETWEEN ARMRESTS	18.5"
SEAT TO FLOOR HEIGHT	20"

BATTERY

BATTERY TYPE	Li-Ion 24 V / 10 Ah
BATTERY INCLUDED	One (1)

POPULAR ADD-ONS



CAREGIVER REMOTE

SKU: 25PET0150100



EXPLORER BACKPACK

SKU: 99ASG0900100



MULTI-POCKET ORGANIZER

SKU: 99ASG0652100



SIDEKICK POUCH

SKU: 99ASG2060100



360° SWIVEL SAFETY LIGHT

SKU: 18AET04502F0



BUMPER BLANKET

SKU: 99ASG1500100

FLIGHTWEIGHT
BY TRAVEL BUGGY

Our newest product line, built for people who want the freedom to go further without paying more. Every **FLIGHTWEIGHT** model is designed, tested, and built to the same quality and reliability standards you already know and trust.