

Wedge Cushion with Stretch Cover

HCPCS:



FIGURE A



FIGURE B

- Wedge design helps prevent patient from sliding forward
- Firm foundation reduces pressure caused by sling-seat hammocks
- The top cover is fluid resistant, low shear stretch fabric
- The bottom cover is waterproof
- High-density foam increases stability and comfort

PRODUCTS	Description
8134	WedgeWheelchairCushion18x16x3-1 1ea,4/CS
SPECIFICATIONS	Description
Warranty	1 Year Limited
Weight Capacity	250 lbs
Cover Material	Fluid-Resistant Stretch Top Waterproof Vinyl Bottom