

Deluxe Folding Travel Walker with 5" Wheels

HCPCS: E0143



FIGURE A



FIGURE B



FIGURE C



FIGURE D



FIGURE E



FIGURE F

- Easily folds tool-free to 1/4 the size of a traditional walker (Figure A - B) and can be stored in a "Ready Set Go" Pouch (Figure C)
- Legs are attached with strong elastic cords so legs can be folded (Figure D)
- Comes with 5" front wheels and rear glide skis that are ideal for most surfaces
- Two side-storage pouches with large pockets for personal items are held in place by hook and loop fastener tabs to both side frames

PRODUCTS

PRODUCTS	Description
RTL10263KDR	Walker w/ 5" Whls, Travel, Red, 1 EA

Description

SPECIFICATIONS

SPECIFICATIONS	Description
Warranty	Limited Lifetime
Weight Capacity	300 lb
Height	29.5" - 36.5"
Carton Shipping Weight	8.5 lb
Overall Width	26.5"
Weight	6 lb
Wheel Size	5"

Description