

									
	Q42	EDGE 3 STRETTO®	EDGE® 3	EDGE® 4	R-TRAK	4FRONT® 2	EDGE® HD	4FRONT® 2 HD	Q1450
ADP LISTING #	Power BaseType 2: Adult: WP0001541	Power Base Type 3: Adult: WP0001806 Paediatric: WP0001807	Power BaseType 3: Adult: WP0001575 Paediatric: WP0001576	Power BaseType 3: Adult: WP0002382 Paediatric: ADP WP0002383	Power BaseType 3: Adult: WP0001987 Paediatric: WP0001988	Power BaseType 3: Adult: WP0001441 Paediatric: WP0001442	Power BaseType 3: Adult: WP0001058 Paediatric: WP0001059	Power BaseType 3: Adult: WP0001900	Power BaseType 3: Adult: WP0001021
DRIVE WHEEL CONFIGURATION	Mid-wheel Drive	Mid-wheel Drive	Mid-wheel Drive	Mid-wheel Drive	Rear-wheel Drive	Front-wheel Drive	Mid-wheel Drive	Front-wheel Drive	Front-wheel Drive
WEIGHT CAPACITY ⁷	300 lb / 136 kg	300 lb / 136 kg	300 lb / 136 kg	300 lb / 136 kg	300 lb / 136 kg	300 lb / 136 kg	450 lb / 204 kg ³	450 lb / 204 kg	600 lb / 272 kg
MAXIMUM SPEED ⁴	Up to 6.25 mph / 10 km/h	Up to 6.25 mph / 10 km/h (3.5 mph / 5.6 km/h with iLevel®)	Up to 6.25 mph / 10 km/h (4.5 mph / 7.2 km/h with iLevel®) Optional: Up to 7 mph / 11.26 km/h ¹	Up to 7 mph / 11.26 km/h (4.5 mph / 7.2 km/h with iLevel®)	Up to 7 mph / 11.26 km/h (3.5 mph / 5.6 km/h when elevated)	Up to 6.25 mph / 10 km/h (3.5 mph / 5.6 km/h when elevated)	Up to 4.7 mph /7.5 km/h (3.5 mph / 5.6 km/h with iLevel®) ³	Up to 5 mph / 8 km/h (2.7 mph / 4.3 km/h when elevated)	Up to 4.5 mph 7.2 km/h
GROUND CLEARANCE ⁸	3.03" (battery tray)	2.9"	2.7" (battery tray)	2.87"	3.97" (battery tray)	3.14" (battery tray)	3" (base frame)	3.14" (battery tray)	1.5" (anti-tip wheels)
TURNING RADIUS ^{5,8}	20.5"	21"	20.5"	22.8"	25.7"	25.9"	19.8"	25.9"	26.5"
LENGTH ^{5,8}	35.5"	36"	35.5"	38.27"	34.8"	39.9"	35.6"	39.9"	38.3"
WIDTH ⁸	24.41"	21.5" (14" drive wheels) 20.5" (12.5" drive wheels)	24.4"	24.5"	24.2"	24"	26.3"	24"	29.3"
DRIVE WHEELS	14" Flat Free Optional: 14" Pneumatic	14" or 12.5" Pneumatic or Flat Free tires Optional: 14" Knobby	14" Flat Free Optional: 14" Pneumatic; 14" Knobby	14" Flat Free Optional: 14" Pneumatic; 14" Knobby	14" Flat Free Optional: 14" Pneumatic; 14" Knobby	14" Flat Free Optional: 14" Pneumatic; 14" Knobby	14" Flat Free Optional: 14" Pneumatic; 14" Knobby	14" Flat Free Optional: 14" Pneumatic; 14" Knobby	14" Flat Free
CASTER WHEELS	Front: 6" Solid Rear: 6" Solid	Front: 5" Solid Rear: 6" Solid	Front: 6" Solid Rear: 6" Solid	Front: 7" Solid Rear: 7" Solid	Front: 9" Flat Free	Rear: 9" Flat Free	Front: 6" Solid Rear: 6" Solid	Rear: 9" Flat Free	Rear: 9" Flat Free
ANTI-TIP WHEELS	N/A	N/A	N/A	N/A	Rear: 4" Solid	Front: 4" Solid	N/A	Front: 4" Solid	Front: 4" Solid
BASE WEIGHT	138.6 lb / 63 kg	126.4 lb / 57.3 kg	138.6 lb / 62.8 kg	154 lb / 70 kg	128 lb / 58.1 kg	130 lb / 59 kg	143.5 lb / 65 kg	130 lb / 59 kg	195 lb / 88.5 kg
BATTERY WEIGHT ¹¹	38 lb / 17.24 kg each	40 AH: 31 lb / 14.1 kg each NF-22: 40 lb / 23.6 kg each	40 lb / 23.6 kg each	47.8 lb / 21.68 kg each	51.5 lb / 23.4 kg each	Group 24: 51.5 lb / 23.4 kg each Group 34: 44.8 lb / 20.3 kg each	51.5 lb / 23.4 kg each	Group 24: 51.5 lb / 23.4 kg each	51.5 lb / 23.4 kg each
MOTOR PACKAGES	2-pole	2-pole	Standard: 2-pole Optional: 4-pole	4-pole	4-pole	4-pole	4-pole	4-pole	4-pole
BATTERIES ⁶	NF-22 (2)	40AH (2) or NF-22 (2)	NF-22 (2)	Group 24 (2)	Group 24 (2)	Group 24 (2) or Group 34 (2)	Group 24 (2)	Group 24 (2)	Group 24 (2)
BATTERY CHARGER	8A Off-board	8A Off-board	8A Off-board	8A Off-board	8A Off-board	8A Off-board	8A Off-board	8A Off-board	8A Off-board
AVAILABLE ELECTRONICS ¹⁰	Q-Logic 3e QLNES	Q-Logic 3EX Q-Logic 3e QLNES	Q-Logic 3EX Q-Logic 3e QLNES	Q-Logic 3EX Q-Logic 3e QLNES	Q-Logic 3EX Q-Logic 3e QLNES	Q-Logic 3EX Q-Logic 3e QLNES	Q-Logic 3EX Q-Logic 3e QLNES	Q-Logic 3EX Q-Logic 3e QLNES	Q-Logic 3EX Q-Logic 3e VR2 6-key
SEATING SIZES ⁸									
TRU-BALANCE® WITH ANTERIOR TILT ⁴	N/A	W: 16-22" D: 14-22"	W: 16-22" D: 14-22"	W: 16-22" D: 14-22"	N/A	W: 16-22" D: 14-22"	N/A	N/A	N/A
TRU-BALANCE® 4 POWER TILT	N/A	W: 12-22" D: 12-22"	W: 12-22" D: 12-22"	W: 12-22" D: 12-22"	N/A	W: 12-22" D: 12-22"	N/A	N/A	N/A
TRU-BALANCE® STATIC	W: 12-22" D: 12-22"	W: 12-22" D: 12-22"	W: 12-22" D: 12-22"	N/A	W: 12-22" D: 12-22"	W: 12-22" D: 12-22"	W: 20-32" D: 20-25"	W: 20-24" D: 20-24"	W: 20-32" D: 20-28"
TRU-BALANCE® 3 POWER TILT	W: 12-22" D: 12-22"	W: 12-22" D: 12-22"	W: 12-22" D: 12-22"	W: 12-22" D: 12-22"	W: 12-22" D: 12-22"	W: 12-22" D: 12-22"	W: 20-32" D: 20-25"	W: 20-24" D: 20-24"	N/A
TRU-BALANCE® 3 LIFT & TILT	N/A	W: 12-22" D: 12-22"	W: 12-22" D: 12-22"	W: 12-22" D: 12-22"	W: 12-22" D: 12-22"	W: 12-22" D: 12-22"	W: 20-24" D: 20-25" ⁹	W: 20-24" D: 20-24"	N/A
TRU-BALANCE® BARIATRIC POWER TILT	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	W: 20-32" D: 20-28"
TRU-BALANCE® FLEX SEAT	W: 15-22" D: 15-20	W: 15-22" D: 15-20"	W: 15-22" D: 15-20"	W: 10-22" D: 10-20"	W: 15-22" D: 15-20"	W: 15-22" D: 15-20"	N/A	N/A	N/A
Q CAPTAIN SEAT	W: 16-22" D: 16-22"	W:16-20" D: 16-20"	W: 16-22" D: 16-22"	W: 16-22" D: 16-22"	W: 16-22" D: 16-22"	W: 16-22" D: 16-22"	W: 20-26" D: 18-24"	N/A	N/A
SEAT-TO-FLOOR HEIGHTS ^{2, 8}									
TRU-BALANCE® 4 WITH ANNTERIOR TILT	N/A	18-19"	17.5-19.25"	18.25-21.25"	N/A	17.75-18.75"	N/A	N/A	N/A
TRU-BALANCE® 4 TILT	N/A	17.5-19.5"	17.5-19.5"	17.75-19.75"	N/A	17.75-19.75"	N/A	N/A	N/A
STRU-BALANCE® STATIC	16.5 - 18.5"	16-18"	16.5 - 18.5"	N/A	16-18"	15 - 17"	16.5 - 18.5"	15-17"	16.5-18.5"
TRU-BALANCE® 3 POWER TILT	17.5-19.5"	17.5-19.5"	17.5 - 19.5"	17.75 - 19.75"	18.5-19.5"	17.75-19.75"	18.5-20.5"	17.75-19.75"	19.5 - 20.5"
TRU-BALANCE® 3 LIFT & TILT	N/A	18-20"	17.5 - 19.5"	17.5 - 19.5"	17.25-18.25"	17.75-19.75"	18.5 - 20.5"	17.75-19.75"	N/A
TRU-BALANCE® FLEX SEAT	16.5 - 18.25	16.5-18.5"	16.5" - 18.25"	16.25" - 18.25"	16.5-17.5"	16.5 -18.5"	N/A	N/A	N/A
Q-CAPTAIN SEAT	SSP: 16.5" - 18.25" (static) Captain: 20.5"- 22.5" (static)	SSP: 16.5" - 18.5" (static) 18"-20" (lift) Captain: 20.5"- 22.5" (static) 22"- 24" (lift)	SSP: 16.5" - 18.25" (static) 17.5"-19.5" (lift) Captain: 20.5"- 22.5" (static) 21.5"- 23.25" (lift)	SSP: 16.5" - 18.25" (static) 17.5"-19.5" (lift) Captain: 20.5"- 22.5" (static) 21.5"- 23.25" (lift)	SSP: 15.5-17.5" (static) 17.25-18.25 (lift) Captain: 19.5"- 21.5" (static) 21.25"-22.25" (lift)	SSP: 17.5-22.5" (static) 17.75-19.75(lift) Captain: 21.75"- 22.75" (static) 21.75"-23.75" (lift)	SSP: 16.5" - 18.5" (static) Captain: 20.5"- 22.5" (static)	N/A	N/A

FOOTNOTES:

- 7 mph can only be achieved with optional 4-pole motors.
- Seat-to-floor heights are measured with the unit unloaded. Due to the 4Front's suspension system, seat-to-floor heights are approximately 0.25" lower when the unit is loaded. Tire type and inflation can also impact the seat-to-floor heights.
- Reduced to 400 lbs. with iLevel® power adjustable seat height
- Speed varies with user weight, terrain type, battery charge, battery condition and tire pressure.

- Without front riggings
- Battery capacity may change when adding power positioning systems.
- Weight capacity is dictated by the selection of seating.
- All measurements +/- 1/8" with 18" x 18" seat, fixed towers, 300 lbs. weight capacity and solid tires (add 3/8" with pneumatic tires). Depths may increase seat-to-floor height. Contact Customer Service for details. Listed measurements are in the lowest position; does not include elevating system unless specified.

- A minimum seating size of 20" x 20" is required with iLevel® on the Q6 Edge® HD; maximum seating size is 24" x 24"
- All Quantum units can be ordered with specialty controls
- Battery weight may vary based on manufacturer.