


TITAN AXS



1. Shroud- C08-084-00100

2. Shroud Knob- VENT-86

3. Decoration Set- S08-084-00118

4. Front Cover & Hardware- TITANMID-01

5. Motor, Gear Brake Assy. Left- C09-084-00201

6. Motor, Gear Brake Assy. Right- C09-084-00101

7. Joystick w/o Extension Cable- LINX-JOYSTICK
Joystick Tip Linx- DC GPL55470

8. Joystick Arm- LDM400001-04
Joystick Mtg. Screw- LDM400001-01

9. Joystick Mtg. Bracket- TITAN-24

10. Joystick Knob- LDM400001-07

11. Joystick Extension Cable- GLM-BUS150-A

12. Joystick Velcro Cable Tie- C34-072-00100

13. Controller- DLX-PM40-A-AXS
Controller to Shroud Cable- GLM-EXT064-A

14. Wire Harness- TITAN-12

15. Battery 35 ah (2 needed)- LRA402208

16. Charger 3 Amp- LRA402209

17. Rear Caster Axle & Nut- TITAN-04

18. Rear Caster- TITAN-03
Rear Caster Bearing- STDS1011

19. Rear Fork- LDM400013

20. Fork Bearing Set Assembly- LRB402348

21. Dust Cap- LDR4110

22. Drive Wheel- TITAN-11

23. Axle Nut- ECI6018

24. Hub Cap- LRM412202-HC

25. Front Fork- LDM400013
Front Fork- TITANMID-09 [2]

26. Caster Axle & Nut- TITAN-04
Caster Axle & Nut- TITAN-04N [2]

27. Front Caster - TITAN-03

27. Front Caster- TITAN-03N

27. Front Caster Bearing- STDS1011

28. Main Frame- TITANMID-05

29. Frame Bushing Assembly- TITANMID-06

30. Side Frame Assembly Left- TITANMID-07

31. Side Frame Assembly Right- TITANMID-08

32. Seat Post - TITAN-39

33. Seat Post Bolt Assembly- TITAN-06

34. Seat Plate- C07-072-01001

35. Captain Seat 16 x 16- TITANAXS-01

35. Captain Seat 16 x 18- TITANAXS-02

35. Captain Seat 18 x 16- TITANAXS-03

35. Captain Seat 18 x 18- TITANAXS-04

35. Captain Seat 20 x 18- TITANAXS-05

35. Captain Seat 20 x 20- TITANAXS-06

35. Captain Seat 22 x 20- TITANAXS-07

35. Pan Seat 18 x 18- ST501-181818-AXS

35. Pan Seat 20 x 18- ST502-201820-AXS

36. Headrest- ST-HEADREST

37. Armrest Assembly Left- LDR410011L

37. Armrest Assembly Right- LDR410011R

37. Arm Pad- LDR410110

37. Arm Pad Screw (pr)- LDR420018

38. Armrest Height Adj Knob- S31052

39. Armrest Width Adj Knob- S31052

40. Seatbelt 48"- STDS850-PW

40. Seatbelt 60"- STDS855-PW

41. Footrest Hanger Assembly- TITANMID-02

42. Footrest Assembly- TITANMID-03

43. Footrest Mtg Bracket Assembly- TITANMID-04
not shown ELR Adapter- AE2500
not shown ELR - LK3JELR

[1] When current chrome stock is gone the part will no longer be chrome plated. It will be black.
[2] Serial number (21)2W2006200536 and after.