

ProBasics Full Body Sling With Commode Opening Mesh, Medium

HCPCS:

- Mesh sling offers total head, neck and partial under-thigh support.
- Easy-to-use sling aids in transferring users from a bed to a wheelchair, toilet or shower chair or from the floor to a bed.
- Use with the ProBasics Hydraulic Lift (Item # PB0087X).
- Four point strap hookup
- Mesh polyester sling is washer safe (see instructions on sling label)
- Can be used in wet or dry environments
- Brushed lining helps prevent sliding and protects skin
- Size: Medium (41.5" x 53.3")
- Commode opening: 11" L x 8.5" W
- Weight Capacity: 450 lb
- Color: Blue

PRODUCTS	Description
PBSL114	MESH SLING W/COMMODE - MED 20EA/CV