



GRAND TOURER **G**

**Built-In
loading
wheels**



3
**Driving
modes**



**Built-In
safety
lock**




KEY FEATURES



**HEIGHT
ADJUSTABLE ARMRESTS**



**HEIGHT
ADJUSTABLE
FOOTREST**



**FULL CLEAR SPAN
ARMRESTS**

Enjoy complete barrier-free side entry access for easy transfers on and off the chair.



ERGONOMIC JOYSTICK



**CONTOURED
COMFORT**

Breathable padding for lasting comfort and support.



**BUILT-IN
STORAGE**



**HEIGHT ADJUSTABLE
BACKREST**

SPECIFICATIONS

SKU (MATTE BLACK)	18MNA0050109
SKU (MATTE WHITE)	18MNA0050108
NET WEIGHT	49 lbs
GROSS WEIGHT	53 lbs
WEIGHT CAPACITY	350 lbs
DRIVING RANGE	25 km / 15.5 mi
INCLINE CAPABILITY	12°
TURNING RADIUS	25"
MAXIMUM SPEED	6 kph / 3.7 mph
MOTOR POWER	2 x 250 W Brushless
UNFOLDED DIMENSIONS	39"L x 23.6"W x 40"H
FOLDED DIMENSIONS	29.5"L x 23.6"W x 13.8"H
FRONT / REAR WHEEL DIAMETER	8" / 10"

SEATING DIMENSIONS

SEAT SIZE	18"W x 18"D x 3"H
BACKREST HEIGHT	17.5" - 22"
BETWEEN ARMRESTS	19.5"
SEAT TO FLOOR HEIGHT	20"


BATTERY

BATTERY TYPE	Li-Ion 24 V / 12 Ah
BATTERY INCLUDED	One (1)

POPULAR ADD-ONS



**GT -
FLIGHT KIT**
SKU: 18ASG1850200



**GT -
SWING-AWAY
JOYSTICK MOUNT**
SKU: 18AJS0650100



**GT -
ADJUSTABLE
HEADREST**
SKU: 18ACM0250100



**GT - CANE
HOLDER SET**
SKU: 18ASG1400100



**GT -
CUP HOLDER**
SKU: 18ASG0150100



**GT -
BATTERY**
SKU: 18PBA0150100