

Folding Exercise Peddler with Electronic Display

HCPCS:



FIGURE A



FIGURE B



FIGURE C



FIGURE D

- Five function display indicates exercise time, revolution count, revolutions per minute (rpm) and calories burned (Figure A)
- Scan feature alternately displays all measurements automatically
- Four anti-slip rubber pads prevent sliding and protect surfaces
- Adjustable knob to help fine-tune tension
- Attractive silver finish
- Button cell battery included
- Completely assembled, foldable design (Figure B - C)
- Retail Packaged

PRODUCTS	Description
RTL10273	Electronic Pedal Exerciser
SPECIFICATIONS	Description
Warranty	Limited Lifetime