

Offset Handle Aluminum Cane

HCPCS:



- Made from heavy-gauge anodized aluminum with an offset hypalon sponge handle
- Featuring push-button adjust with anti-rattle lock-nut, nylon wrist strap, and cane tip with wide suction grip

PRODUCTS

Description

730-320	ALU. CANE OFFSET HAND. SILVER, 1EA, 6CV
730-325	588. CANE OFFSET HAND. BRONZE, 1EA, 6CV

SPECIFICATIONS

Description

Warranty	Limited Lifetime
Weight Capacity	250 lb (113 Kg)