

22" Bariatric Aluminum Transport Chair

HCPCS: E1039



FIGURE A



FIGURE B



FIGURE C



FIGURE D




FIGURE E




FIGURE F



FIGURE G

- Weighs only 33 lbs
- Supports individuals up to 450 lbs
- Attachable push bar with padded foam grip
- 12" rear flat-free wheels
- Back folds down for convenient storage
- Padded armrests provide additional comfort
- Heavy-duty, nylon reinforced upholstery
- Comes with hand wheel lock levers and companion activated wheel locks
- Comes with seat belt
- Comes with swing-away footrests and anti-tippers

PRODUCTS

Description

ATC22-R	Aluminum HD 22" Red Transport Chair,1/ea
---------	--

SPECIFICATIONS

Description

Warranty	Limited Lifetime
Product Dimensions	41" (L) w/ riggings x 30" (W) x 37.5" (H)
Seat Dimensions	22" (W) x 18" (D)
Weight Capacity	450 lb
Seat To Floor Height	21"
Carton Shipping Weight	39 lb
Overall Width	30" (Open) 12" (Closed)
Weight	33 lb
Back Height	16"