

SOLAX
easy life

**FULLY
ASSEMBLED**

- No Tools Required
- Fold & Unfold Easily
- Fits in Any Vehicle
- Light Weight, Easy Carry
- Ideal for Travel
- Lithium Battery



REMOTE CONTROL
AUTOMATIC FUNCTION



SOLAX Specifications	
Model:	Transformer
Overall Length:	36.9 Inches
Overall Width:	18.1 Inches
Tires (Solid Wheels):	Front 6x1.5 Inches Rear 7 x 2.36 Inches
Maximum Speed:	6 kmh/3.75 mph
Safe Gradient:	Up to 15 Degrees
Range:	Up to 20 km/13.7 miles
Turning Circle:	55 Inches
Weight (with battery):	53 lbs
Color:	Red, Blue, Yellow, Black
Brake:	Intelligent Regenerative Electromagnetic Brake
Drive System:	Rear Wheel Drive
Handlebar:	PU Foam
Maximum Weight Capacity:	300 lbs
Seat Width:	17 Inches
Seat Height (from the ground):	22 Inches
Battery:	25.9V/10Ah (Lithium Battery)
Motor:	120W
Battery Charger:	DC24V/2A Off-Board (Automatic Protection Charger)



Enhance Your Life