

Med-Aire Plus 10" Bariatric Alternating Pressure and Low Air Loss Mattress Replacement System



FIGURE A

- Active combination bariatric therapy system provides both alternating pressure and low air loss to optimize pressure redistribution, shear/friction reduction, and microclimate control
- Designed to prevent, treat and heal all stages of pressure ulcers in the home or long term care setting
- Digital control unit offers 4 cycle times and adjustable comfort weight settings for optimal immersion and patient satisfaction
- 3 diagnostic audio and visual alarms alert for issues with air pressure, the system, or power
- System includes 20 individual air cells to maximize immersion and envelopment with pillow feature for head stability & comfort
- Advanced Clinical and Care Giver Features: Auto-Firm, Seat Inflate, auto-recovery, microprocessing feedback, lock-out and more!
- Cell-on-cell designed air cells prevent "bottoming out" and provide up to 24 hours of power outage protection
- 10" deep air cells are easily removed with quick disconnects for convenient servicing
- Removable quilted Nylon cover is fluid resistant, low shear and vapor permeable helping protect the skin from friction and moisture
- CPR valve allows for rapid deflation
- Tubing quickly and easily disconnects with EZ Lock/Quick Release design

PRODUCTS	Description
14030	MedAireBariatricLowAirLoss,42",1/cs
14048	MedAireBariatricLowAirLoss,48",1/cs
14054	MedAireBariatricLowAirLoss,54",1/cs
14060	MedAireBariatricLowAirLoss,60",1/cs

SPECIFICATIONS	Description
Control Unit Airflow LPM	12 LPM
Control Unit Alarms	Visual / Audible
Control Unit Cycle Time	10, 15, 20, 25 Minutes
Control Unit Type	Compressor
Mattress Cover Material	Fluid-Resistant, 4-way Stretch
Mattress Dimensions	42", 48", 54", 60" (W) X 80" (L) X 10" (H)
Warranty	18 Month Limited
Weight Capacity	600 lbs. - 1000 lbs.

ACCESSORIES	Description
14030XP	3-Outlet Pump,14030/48/54/60,1ea, EA
14060MN	Mattress Only,14060,3Tube,1/ea, EA
14054-CELL	Replacement Bladder,14054,1/ea, EA

ACCESSORIES

Description

14060-CELL	Replacement Bladder,14060,1/ea, EA
14048-CELL	Replacement Bladder,14048,1/ea, EA
14030MN	Mattress Only, 14030, 3-Port, 1/ea, EA
14030C	Cover Only,14030,1/ea, EA
14054MN	Mattress Only/14054/3Tube,1/ea, EA
14048C	Cover,14048,1/ea, EA
14054C	Cover for 14054C, 1/ea, EA
14030-CELL	Bladder Only,14030,1/ea", EA
14060C	Cover for 14060, 1/ea, EA
14048MN	Mattress Only,14048,3Tube,1/ea, EA