

Half Length Bed Rails

HCPCS: E0305



FIGURE A



FIGURE B

- Unique design prevents individuals from being lodged between bars
- Ideal for maximum protection when bed is elevated
- Fits most hospital beds with steel frame (Figure A)
- Attractive brown-vein finish is easy-to-maintain and less institutional looking (Figure B)
- Constructed of 1" steel
- Not for use with "Pan-Style" bed surface
- Easy-to-install

| PRODUCTS | Description |
|----------------|--|
| 15201BV | Bed Rail, No-Gap, Half, BV, 1 pair= 1/cs |
| SPECIFICATIONS | Description |
| Warranty | Limited Lifetime |
| ACCESSORIES | Description |
| 14198BV | Clamp,15201BV No Gap Rail, 1/ea, BX |
| PARTS | Description |
| 15302P1005 | CompositeHeadPanels,15302,1/pr/cs, CS |
| 15302P1006 | CompositeFootPanel,15302,1pr/ea, EA |
| 15302P1051 | PanelFrame,HeadBoard,15302,1/ea, EA |
| 15302P1052 | PanelFrame,FootBoard,15302,1/ea, EA |
| 15309N | HeadFrame,15302,BariatricBed,1/cs, CS |
| 15320N | FootFrame,15302,BariatricBed,1/cs, CS |