

Bariatric Foam Mattresses

HCPCS: E0184-1,E0184



FIGURE A

- Nylon top cover is fluid-resistant and easy to clean for extended life and service
- Vinyl bottom is waterproof
- Fire Retardant: Mattress meets the requirements of 16 CFR Part 1632 Standard for the Flammability of Mattresses and Mattress Pads and 16 CFR Part 1633 (Federal Flammability (open flame) standard for mattress sets) when used without a foundation
- High-Density foam mattress for the bariatric patient
- Solid foam inner core provides comfort, support and durability in a one-piece construction
- Latex free

| PRODUCTS | Description |
|----------|---------------------------------------|
| 15301 | Bariatric Mattress, 42in,1/ea |
| 15301-84 | Bariatric Foam Matt,84"x42"x6.5",1/ea |

| SPECIFICATIONS | Description |
|-------------------------|---|
| Cover Material | Fluid-Resistant Nylon Top Waterproof Vinyl Bottom |
| Dedicated Heel Section | No |
| FireRating | CFR 16 Part 1633 Federal Fire Code |
| Limited Warranty | 1 Year |
| Mattress Layers of Foam | 1 |
| Product Dimensions | 15301: 80(L) x 42(W) x 6.5(H) 15301-84: 84(L) x 42(W) x 6.5(H) 15310-84: 84(L) x 48(W) x 6.5(H) |
| Side Rails | No |
| Warranty | 1 Year |
| Weight Capacity | 15301: 585 lbs 15301-84: 585 lbs 15310-84: 665 lbs |

| ACCESSORIES | Description |
|-------------|--|
| 15301-84CVR | Cover Bari Foam Matt,84"x42"x6.5",1/ea, EA |
| 15310-84CVR | Cover Bari Foam Matt,84"x48"x6.5",1/ea, EA |
| 15312-84CVR | Cover Bari Foam Matt,84"x54"x6.5",1/ea, EA |
| 15301C | Bariatric Mattress Cover, 42", 1/ea, EA |
| 15310C | Mattress Cover 48" 1/ea, CS |