

## Bariatric Mattress

HCPCS: E0184,E0184



FIGURE A

- Mattress Top Cover Unitek 11 Nylon

PRODUCTS	Description
15312	Bariatric Mattress,54",1/ea
SPECIFICATIONS	Description
Mattress Base	4 lbs. Density Foam
Mattress Bottom Cover	Vinyl
Mattress Fire retardant	Compliance with 16 CFR 1633
Mattress Middle Layer	2.5 lbs. Density Foam
Mattress Top Cover	Unitek 11 Nylon
Mattress Warranty	1 year
Mattress Top Layer	2.5 lbs. Density Foam
Product Dimensions	80(L) x 54(W) x 6.5(H)
Warranty	Mattress: 1 Year   Mattress Cover: 60 Days
Weight Capacity	750 lbs