



Pride Mobility Products Co.

5096 South Service Rd, Beamsville, ON, L3J 1V4

Phone: 888-570-1113 | Fax: 866-514-1303 | Email: productorders@pridemobility.com

Q6 Edge® 3 Stretto Order Form with TRU Balance® Flex Seating

Account Number: _____ Date: _____

Provider Name: _____

Contact: _____

Phone: _____ Fax: _____

Email: _____

PO Number: _____

Marked for: _____

Ship to Address: _____

City: _____ Prov: _____ PC: _____



INTRODUCTION

This form is interactive when viewed with [Adobe Acrobat Reader](#) and may not function correctly if opened with applications other than Acrobat. Complete the form by placing checks in the desired boxes and provide information in the interactive fields. Buttons shown at the bottom of the form may be utilized to print or submit the order form through a desktop email application. To email via a web-based application, please 'Save As' and attach the PDF to your email.


Thank you for choosing Quantum!

PATIENT INFORMATION

Weight: _____ Height: _____ A. Shoulder Width: _____ D. Max Sitting Height: _____ G. Thigh Depth: _____

B. Chest Width: _____ E. Shoulder Height: _____ H. Lower Leg Length: _____

C. Hip Width: _____ F. Axilla Height: _____ I. Elbow Height: _____



1. BASE MODEL

REQUIRED. Select a group 3 power base. 300 lb. weight capacity with standard 6mph programmable motors and off board charger. Dimensions of base without legrests: 20.5"Wx35.5"L

Quantum Q6 Edge 3 Stretto\$8,450

Part: Q6 EDGE 3 STRETTO

2. DRIVE MOTORS

Optional. Requires Q-Logic 3 electronics. The standard 6 mph motor can be upgraded to Accu-Trac motors that delivers enhanced traction and tracking, especially with switched drive controls.

Accu-Trac Motor Technology\$995

Part: XRF000268. HCPCS: K0108

3. BASE COLOR**REQUIRED.** Select one shroud color. Colors shown here may differ from actual product depending on monitor/printer calibration used.**3.1a. Standard Shrouds**

 Red Corvette Part: DGN1900199-PRD1	 Deep Purple Part: DGN1900199-PPR5	 Midnight Blue Part: DGN1900199-PBU10	 Green Machine Part: DGN1900199-PGN11
 Back in Black Part: DGN1900199-PBK1	 Raspberry Beret Part: DGN1900199-PPK4	 Electric Blue Part: DGN1900199-PBU9	 Military Green Part: DGN1900199-PGN8
 Night Train Part: DGN1900199-PBK3	 Pretty in Pink Part: DGN1900199-PPK3	 Blue Jean Part: DGN1900199-PBU11	
 Grey Street Part: DGN1900199-PSL3	 Evolution Orange Part: DGN1900199-POR7	 Lemon Crush Part: DGN1900199-PYL4	

3.1b. Hydrographic Shrouds

 Military Digital Camo\$325 Part: DGN1801101-WCM7	 Faded Glory\$325 Part: DGN1801101-WFL2
--	--

4. BASE MOUNTED LIGHTS**Optional.** Select one option.

Switch Fender Lights..... No Charge
Part: DGN1810041
Fender lights operate by pressing included single switch button. Fender lighting is not available with Single Actuator Control through Toggle Box.

**5. QUANTUM BACKUP CAMERA****Optional.** The backup camera may interfere with armrest mounted accessories such as media device holder.

Quantum Backup Camera \$475
Part: MEC1910365

1. Optional Camera Mounting Bracket

Backup Camera Headrest Bracket \$35 Part: MEC1910373 Works with .5" x .5" vertical supports only.	Backup Camera Track Mount Bracket \$35 Part: MEC2104227
--	---

2. Optional Display Mounting Bracket

Backup Camera Q-Logic Joystick Mount Bracket \$145
Part: MEC2106123

**6. DRIVE WHEELS****REQUIRED.** Select a type of drive wheel. Accent colors are available for 14" tires.**6.1. Drive Wheels**

3"x14" Pneumatic Tires No Charge Part: WHL1713925	2"x12.5" Pneumatic Tires No Charge Part: WHL1810666
3"x14" Flat Free Tires No Charge Part: WHL1704058	2"x12.5" Flat Free Tires No Charge Part: WHL1810665
4"x14" Knobby, Flat Free Black Tires \$120 Part: XRF000368	4"x14" Knobby, Pneumatic Black Tires \$120 Part: XRF000367

6.2. Drive Wheel Accent Color for 14" Tires

Not available with 12.5" tires.

Back in Black No Charge Part: DGN164245-PBK1 x2	Deep Purple No Charge Part: DGN164245-PPR5 x2	Midnight Blue No Charge Part: DGN164245-PBU10 x2	Green Machine No Charge Part: DGN164245-PGN11 x2
Red Corvette No Charge Part: DGN164245-PRD1 x2	Raspberry Beret No Charge Part: DGN164245-PPK4 x2	Electric Blue No Charge Part: DGN164245-PBU9 x2	Military Green No Charge Part: DGN164245-PGN8 x2
Night Train No Charge Part: DGN164245-PBK3 x2	Pretty in Pink No Charge Part: DGN164245-PPK3 x2	Blue Jean No Charge Part: DGN164245-PBU11 x2	Military Digital Camo \$50 ea. Part: DGN164245-WCM7 x2
Grey Street No Charge Part: DGN164245-PSL3 x2	Evolution Orange No Charge Part: DGN164245-POR7 x2	Lemon Crush No Charge Part: DGN164245-PYL4 x2	Faded Glory \$50 ea. Part: DGN164245-WFL2 x2

7. BATTERIES**Optional.** The powerchair requires 2 batteries and will be installed if ordered.

NF-22, Gel Cell, Interceptor 2 Battery \$300 ea.
 Part: BAT2004376. HCPCS: E2361
 Uses 8 AMP charger

Gel Cell, 40AH, Interceptor 2 Battery \$200 ea.
 Part: BAT2100316. HCPCS: K0108
 Uses 5 AMP charger

8. PATIENT'S WEIGHT**REQUIRED.** The patient's weight is required to quote and order the power chair.

Up to 150 lb User Weight
 Part: PTOINDV3977

151-250 lb User Weight
 Part: PTOINDV3985

250-300 lb User Weight
 Part: PTOINDV3984

9. ELECTRONICS**REQUIRED.** Select one type of drive control from the following sections. The No Joystick option is used to omit the joystick and utilize another control, such as a specialty control found on the [Alternative Controls and Electronics Order Form](#).**9.1a. Non-Expandable Controls**

Select a joystick below or omit selection if another controller type is desired. Non-expandable does not allow supplementary drive controllers.

QLNES 6 Key Joystick Standard
 Part: CTL2205379
 Controls up to 2 actuators.

QLNES 10 Key Joystick No Charge
 Part: CTL2205380
 Controls up to 2 actuators through joystick and seat mounted, LED, full light package.

**9.1b. Expandable Controls**

Select expandable controls or omit selection if another controller type is desired. Expandable controls allow multiple or specialized controllers.

2. Q-Logic 3 Joystick

Q-Logic 3 EX Joystick includes XLR USB charger.

Q-Logic 3 EX Joystick \$900
 Part: CTL167724

Q-Logic 3 EX Joystick for Full Lights \$900
 Part: CTL167725

Q-Logic 3e Joystick \$620
 Part: CTL1910600
 Can operate full light package. Not available with alternative or specialty controls.

No Joystick, Setup for Q-Logic 3
 Part: XRF000010
 Requires selection of Enhanced Display

**9.2. Q-Logic Magnetic Charger**

Available with QLNES. Included with Q-Logic 3e Joystick.

Magnetic to XLR Charger Adapter with Wall Holder \$100
 Part: ELE2103No Charge5

9.3. Q-Logic 3 Joystick Armor**REQUIRED with Q-Logic 3 EX Joystick selections above. Not available with Q-Logic 3e Joystick.** Reference colors shown on page 2. *Night train joystick armor will be the black, soft touch.

Black, Soft Touch Standard
 Part: DGN124307
 Also matches Night Train shrouds.

Red Corvette No Charge
 Part: DGN124336-PRD1

Back in Black No Charge
 Part: DGN124336-PBK1

Grey Street No Charge
 Part: DGN124336-PSL3

Deep Purple No Charge
 Part: DGN124336-PPR5

Raspberry Beret No Charge
 Part: DGN124336-PPK4

Pretty in Pink No Charge
 Part: DGN124336-PPK3

Evolution Orange No Charge
 Part: DGN124336-POR7

Midnight Blue No Charge
 Part: DGN124336-PBU10

Electric Blue No Charge
 Part: DGN124336-PBU9

Blue Jean No Charge
 Part: DGN124336-PBU11

Lemon Crush No Charge
 Part: DGN124336-PYL4

Green Machine No Charge
 Part: DGN124336-PGN11

Military Green No Charge
 Part: DGN124336-PGN8

Military Digital Camo No Charge
 Part: DGN124336-WCM7

Faded Glory No Charge
 Part: DGN124336-WFL2

10. JOYSTICK HANDLES**Optional.** The optional handle will be mounted to the joystick and the standard handle will be shipped with the unit.Mushroom Handle (B) \$120
Part: ACCASMB1860.Chin Cup (C) \$85
Part: ACCASMB186190T-Bar (D) \$160
Part: ACCASMB1858.Medium Ball Handle (E) \$70
Part: ACCASMB1859.

Soft Ball Handle (F) \$120

Part: ACCASMB1856.

Extended I-Handle (G) \$140
Part: ACCASMB1857.Bodypoint 3" Goal Post Handle \$180
Part: PLSKNCP1057.Bodypoint 4" Goal Post Handle \$180
Part: PLSKNCP1058.**11. JOYSTICK MOUNTING BRACKETS****Optional.** Select a mount to be used with a joystick selected in previous section.**11.1a. Fixed, In-line Joystick Bracket**Fixed, In-line, Bracket for Q-Logic 3 Standard
Part: FRMASMB10617**11.1b. Swing-Away (Retractable) Joystick Bracket**

Allows the joystick to retract to the side of the armrest while maintaining the same joystick orientation.


Swing-away, Q-Logic 3 \$265
Part: ACC1908434. HCPCS: E1032**11.1c. Basic Swing Away Joystick Mount**Basic Swing Away, Q-Logic, Left Side \$265
Part: ACC2201041. HCPCS: E1032Basic Swing Away, Q-Logic, Right Side \$265
Part: ACC2201042. HCPCS: E1032**11.1d. Offset Swing Away Joystick Mount**Offset Swing Away, Q-Logic, Left Side \$265
Part: ACC142742. HCPCS: E1028Offset Swing Away, Q-Logic, Right Side \$265
Part: ACC142743. HCPCS: E1028**11.1e. Specialty Joystick Mounts**

Select a bracket and side to mount it.


Retract4 Mount \$450
Part: XRF000003. HCPCS: E1032Multi-Axis Mount \$100
Part: XRF000004

Full 180° of movement that translates into multi-position use, straight, right offset, left offset, or fold away. Joystick orientation varies as bracket is adjusted.

Allows height adjustment of joystick.



Swing Away Joystick



Retract4



Multi-Axis Joystick Mount

11.2. Joystick Mount LocationLeft Side Mount
Part: PTOINDV1225Right Side Mount
Part: PTOINDV1226**12. TRU BALANCE FLEX SEATING SYSTEM**

Seat width and depth can be adjusted in 1" increments from 15"W-20"W and 16"D-20"D. Maximum seat size 20"W x 20"D.

12.1. TRU Balance Flex Frame Selection**REQUIRED.**TRU Balance Flex Seating System \$650
Part: PTOINDV3853**12.2. Occupied Transit Kit****Optional.** The power base is standard with Unoccupied Transit Loops. 60" Sure-Lok Lap Belt is required with Occupied Transit Option.WC-19 Occupied Transit Option \$400
Part: KIT152047To determine if your chair is WC-19 compliant (occupied), please refer to this website: <https://wc-transportation-safety.umtri.umich.edu/>**12.3. Seat Width****REQUIRED.** Maximum seat width is 20". Part number is SETWDTH plus number following # below.

10" W (#1000)	11" W (#1001)	12" W (#1002)	13" W (#1003)	14" W (#1004)	15" W (#1005)	16" W (#1006)	17" W (#1007)	18" W (#1008)	19" W (#1009)	20" W (#1010)
------------------	------------------	------------------	------------------	------------------	------------------	------------------	------------------	------------------	------------------	------------------

12.4. Seat Depth**REQUIRED.** Maximum seat depth is 20". Part number is SETDPTH plus number following # below.

10" D (#1000)	12" D (#1002)	14" D (#1004)	16" D (#1006)	18" D (#1008)	20" D (#1010)
11" D (#1001)	13" D (#1003)	15" D (#1005)	17" D (#1007)	19" D (#1009)	


TRU Balance Flex seating section continued on following pages.

12. TRU BALANCE FLEX SEATING SYSTEM CONTINUED**12.5. Seat to Floor Height****Lowest Seat to Floor**
Part: PTOINDV3558**Middle Seat to Floor**
Part: PTOINDV3559**Highest Seat to Floor**
Part: PTOINDV3560

High	18.5" at front with 4" posterior angle
Med	17.5" at front with 2" posterior angle
TRU Balance Flex Seat to Floor Heights	

12.6. Stealth Seat Cushion

Optional. Specific cushion/component sizes available, part numbers, weight capacities, MSRP, HCPCS code, and special order cushions will be found on the respective complete Stealth Cushion Order Forms. Part numbers in this section are for reference only. This section allows a quick selection of a single cushion that matches the seat pan width and depth to alleviate listing all part numbers and codes on this power base form. The complete cushion part number and MSRP is a product of the cushion model, top foam, liner, cover, width, and depth as shown below. Reference the subsequent chart for available top foam, liner, and cover options for each cushion type. See specific cushion order form for all sizes, top foam, liners, covers, parts numbers, and MSRP.

**1. Cushion Type****TRU-Comfort 2 SPP**
Reference part: XRF000116.
Part series: TC2-###XXYY.**Essence SPP**
Reference part: XRF000118.
Part series: ESS-###XXYY.**Essence Plus SPP**
Part series: ESSP-###XXYY.**Solution SPP**
Reference part: XRF000113.
Part series: SOL-###XXYY.**Spectrum Foam SPP**
Reference part: XRF000112.
Part series: SPCF-###XXYY.**Spectrum Gel SPP**
Reference part: XRF000117.
Part series: SPCG-###XXYY.**Premiere P**
Reference part: XRF000122.
Part series: PREM-###XXYY.**Glacial SP**
Reference part: XRF000114.
Part series: GLA-###XXYY.**Zen SP**
Reference part: XRF000115.
Part series: ZEN-###XXYY.**Simplicity G**
Reference part: XRF000111.
Part series: SIM-###XXYY.**2. Essence Cushion Wedge Kit**

Only available with Essence cushion selection above.

Small Wedge Positioning Kit\$45
Part: POP1909006-S
14"W and smaller**Small Wedge Positioning Kit**\$45
Part: POP1712340-S
15"W and larger**Medium Wedge Positioning Kit**\$85
Part: POP1909006-M
14"W and smaller**Medium Wedge Positioning Kit**\$85
Part: POP1712340-M
15"W and larger**Large Wedge Positioning Kit**\$90
Part: POP1909006-L
14"W and smaller**Large Wedge Positioning Kit**\$90
Part: POP1712340-L
15"W and larger**3. Cushion Top Foam****Visco Top Foam**
Reference part: XRF000380**FRG Top Foam**
Reference part: XRF000381**No Top Foam**
Reference part: XRF000382**4. Cushion Liner****Stealth Urethane Liner**
Reference part: XRF000383**Stealth Polyurethane Liner**
Reference part: XRF000384**No Liner**
Reference part: XRF000385**5. Cushion Cover****CoolCore Cover**
Reference part: XRF000386**Mesh Cover**
Reference part: XRF000387**Stealth Tek Cover**
Reference part: XRF000388

	None	Visco/V.Elastic	FRG	None	Urethane	Polyurethane	CoolCore	Mesh	Stealth Tek
	Top Foam			Liner				Cover	
TRU-Comfort 2	—	✓	—	✓	—	✓	✓	—	—
Essence	—	✓	—	✓	—	✓	✓	✓	—
Solution	—	✓	—	—	✓	✓	✓	✓	✓
Spectrum Foam	—	✓	✓	—	✓	✓	✓	✓	✓
Spectrum Gel	—	✓	✓	—	✓	✓	✓	✓	✓
Premiere	—	✓	✓	✓	✓	✓	✓	—	—
Glacial	✓	—	—	✓	✓	✓	✓	✓	✓
Zen	—	✓	✓	✓	✓	✓	✓	✓	✓
Simplicity	✓	—	—	✓	—	—	—	—	✓

TRU Balance Flex seating section continued on following pages.

12. TRU BALANCE FLEX SEATING SYSTEM CONTINUED**6. Optional Base Insert**

Option selection available with Glacial, Premiere, Simplicity, Solution, Spectrum Foam, Spectrum Gel, and Zen cushions.

ABS Flat Rigidizer

Reference part: XRF000389

7. Additional Cushion Liner

Up to 5 may be selected below. See respective cushion order form for MSRP and part numbers.

Additional Stealth Urethane Liner

Reference Part: XRF000390

Additional Stealth Polyurethane Liner

Reference Part: XRF000391

8. Additional Cushion Cover

Up to 5 may be selected below. See respective cushion order form for MSRP and part numbers.

Additional CoolCore Cover

Reference Part: XRF000392

Additional Stealth Tek Cover

Reference Part: XRF000394

Additional Mesh Cover

Reference Part: XRF000393

12.7. TRU Balance Flex Back Cane Selection**REQUIRED.**

Static Back Strap, 24" Quick A/A Back Canes No Charge

Angle Adjustable Back Cane \$195

Part: MEC140266

Flip Forward Back Cane \$195

Part: XRF000274

Not available with oxygen holder, or cane mounted and cantilever armrests.

2" Mesomorph (1" per side) Back Cane \$205

Part: MEC1804736

2" Endomorph (1" per side) Back Cane \$205

Part: XRF000270

4" Endomorph (2" per side) Back Cane \$205

Part: XRF000272

Manual Recline Back \$395

Part: FRMASMB12910. HCPCS: E1226

Back canes will be 24" H. Not available with flip back armrests, cantilever armrests, TRU-Comfort Back, Sport Back, rear accessories, or articulating foot platform. Pivot point is 2.5" above the seat pan.

12.8. TRU Balance Flex Back Cane Height Selection**REQUIRED.** 16" and 17" back canes are not available with a Synergy Back and cane mounted, flip up armrests.

14" High Canes

Part: FRMASMB15832

17" High Canes

Part: FRMASMB15835

20" High Canes

Part: FRMASMB15838

23" High Canes

Part: FRMASMB15841

15" High Canes

Part: FRMASMB15833

18" High Canes

Part: FRMASMB15836

21" High Canes

Part: FRMASMB15839

24" High Canes

Part: FRMASMB15842

16" High Canes

Part: FRMASMB15834

19" High Canes

Part: FRMASMB15837

22" High Canes

Part: FRMASMB15840

25" High Canes

Part: FRMASMB15843

12.9. TRU Balance Flex Back Cushion Options

Select one type of back. See the seat back availability chart for stock order back size ranges. Other size backs may be available with special order and extended lead time. Select a back type and overall back height for configured part number.

1b. Sport Back Selection

Sport Back Cushion, 21" Wide or Less No Charge

Part: XRF000361

Sport Back Cushion, 22" Wide or More No Charge

Part: XRF000362

1c. TRU-Comfort Back Selection

Back size is the product of seat width, back cane, and overall back height selections. See availability chart for sizes.

TRU-Comfort Back, Back Cane Option, 16" to 20" Width \$835

Part: TC4V10W21WBC. HCPCS: E2620

TRU-Comfort Back, Back Cane Option, 22" Width \$835

Part: TC4V22W32WBC. HCPCS: E2621

1h. Omit Back Selection

Omit Back No Charge

Part: PTOINDV1287

		Back Width										
		12"	13"	14"	15"	16"	17"	18"	19"	20"	21"	22"
Overall Height	13"	-	-	S	-	-	-	-	-	-	-	-
	14"	P	T	-	S	P	-	-	-	-	-	-
	15"	S	P	T	-	S	P	-	-	-	-	-
	16"	S	P	T	-	S	P	T	P	T	P	T
	17"	P	T	S	S	P	T	S	P	T	S	P
	18"	-	S	S	S	S	P	T	S	P	T	S
	19"	-	-	-	S	S	P	T	S	P	T	S
	20"	-	-	-	S	S	P	T	S	P	T	S
	21"	-	-	-	S	S	P	T	S	P	T	S
	22"	-	-	S	S	S	P	T	S	P	T	S
	23"	-	-	S	S	S	P	T	S	P	T	S
	24"	-	-	S	S	S	P	T	S	P	T	S
	25"	-	-	S	S	S	P	T	S	P	T	S

S Synergy Back P Sport Back T TRU-Comfort Back

TRU Balance Flex Seat Back Availability Chart

12. TRU BALANCE FLEX SEATING SYSTEM CONTINUED**12.10. Back Width Package****REQUIRED with Sport and TRU-Comfort back selections.**16"-20" Back Width Package
Part: XRF00003818"-22" Back Width Package
Part: XRF000039**12.11. Overall Back Height****REQUIRED with back selection.** Refer to the Synergy and TRU-Comfort back availability chart for stock order back size ranges.

14" Height Part: SETHGHT1027

Part: SETHGHT1027

15" Height Part: SETHGHT1028

Part: SETHGHT1028

16" Height Part: SETHGHT1029

Part: SETHGHT1029

17" Height Part: SETHGHT1030

Part: SETHGHT1030

18" Height Part: SETHGHT1031

Part: SETHGHT1031

19" Height Part: SETHGHT1037

Part: SETHGHT1032

20" Height

Part: SETHGHT1033

21" Height

Part: SETHGHT1034

22" Height

Part: SETHGHT1035

23" Height

Part: SETHGHT1036

24" Height

Part: SETHGHT1037

25" Height

Part: SETHGHT1038

26" Height

Part: SETHGHT1039

27" Height

Part: SETHGHT1040

28" Height

Part: SETHGHT1041

12. TRU BALANCE FLEX SEATING SYSTEM CONTINUED**1c. ADI2.0 Backrest**Specify back type, width, height, and contour in the part number configurator below. Previously selected back cane height must be equal to the greater of overall height range. **Available SOF for Provider installation only.**For more information please see [ADI2.0 Backrest order form](#)

ADI2.0 Backrest 10"W - 20"W \$1,495

For use with back canes only.

Indicate part number below.

ADP: SE0002152

ADI2.0 Backrest 22"W - 24"W \$1,695


For use with back canes only.

Indicate part number below.

ADP: SE0002153

Specify Part Number for Back Selection

OPTIONS INDEX	ADI2	Back Type: AL AFT	Width: 10"-20" (1" increments), 22", 24"	Height: N = Nino L = Low M = Medium S = Standard T = Tall X = Extra Tall	Contour: A = Active D = Deep E = Extra Deep	2 = FRG Foam	P = Quick-Release
ENTER SELECTION	ADI2					2	P

BACK CONTOUR*:**BACK SERIES TYPES*:**

(AL) Aluminum Series

(*ALL BACKS SHOWN ARE 18" WIDE)

(AFT) Aluminum Flat Top Series

**BACK HEIGHTS AND WIDTH:**

(Available for AL Back Shell Only)

Nino: 8" inch (10"W - 14"W)

Low: 10 inch (10"W - 24"W)

Medium: 13 inch (12"W - 24"W)

(Available for AL and AFT Shell Types)

Standard: 16 inch (13"W - 24"W)

Tall: 18 inch (16"W - 24"W)

Extra Tall: 20 inch (16"W - 24"W)

12. TRU BALANCE FLEX SEATING SYSTEM CONTINUED**12.12. Armrests****REQUIRED.** Left (L) and right (R) boxes may be checked to configure different size arms within one section arm type.**1a. TRU Balance Flex Seat, Drop-In Armrests**

Includes full length, flat armpads. Available with 14" deep seats and greater.

☐ ☐ Full Length, Drop-In Armrests No Charge
☐ ☐ ☐ ☐ Left part: ARM1803253. Right part: ARM156742

☐ ☐ Desk Length, Drop-In Armrests No Charge
☐ ☐ ☐ ☐ Left part: ARM1803254. Right part: ARM156743
1b. Pediatric, Height Adjustable, Removable Armrests

Available with straight armpads only.

☐ ☐ Desk Length, Pediatric Armrests \$145 ea.
☐ ☐ ☐ ☐ Left part: FRMASMB10608. Right part: FRMASMB10609. HCPCS: E0973
1c. Flip Back, Cane Mounted Armrests

May interfere with swing-away laterals.

☐ ☐ Full Length, Flip Back Armrests No Charge
☐ ☐ ☐ ☐ Left part: XRF000371. Right part: XRF000372

☐ ☐ Desk Length, Flip Back Armrests No Charge
☐ ☐ ☐ ☐ Left part: XRF000373. Right part: XRF000374
1c. Two-Post, Flip Back, Height Adjustable, Removable Armrests

Not available with a depth less than 16". Required with flip forward back and laptrays.

☐ ☐ Full Length, Two-Post Armrests \$145 ea.
☐ ☐ ☐ ☐ Left part: XRF000051. Right part: XRF000052. HCPCS: E0973

☐ ☐ Desk Length, Two-Post Armrests \$145 ea.
☐ ☐ ☐ ☐ Left part: XRF000053. Right part: XRF000054. HCPCS: E0973
1d. Single Post, Height Adjustable, Removable Armrests

Height adjustment from 11"-17" in half-inch increments.

☐ ☐ Full Length, Single Post Armrests \$145 ea.
☐ ☐ ☐ ☐ Left part: ARM2000906. Right part: ARM2000907. HCPCS: E0973

☐ ☐ Desk Length, Single Post Armrests \$145 ea.
☐ ☐ ☐ ☐ Left part: ARM2000908. Right part: ARM2000909. HCPCS: E0973
12.12. TRU Balance Flex Armpads

Waterfall armpads may cause interference with the seat back. Armpads are configured as pairs.

1a. Straight Armpads
 Straight, Full Length Left and Right No Charge
 Part: ARM134275

 Straight, Full Length Left, Desk Length Right No Charge
 Part: ARM134281

 Straight, Desk Length Left, Full Length Right No Charge
 Part: ARM134280

 Straight, Desk Length Left and Right No Charge
 Part: ARM134254
1b. Waterfall Armpads
 Waterfall, Full Length Left and Right No Charge
 Part: ARM134285

 Waterfall, Full Length Left, Desk Length Right No Charge
 Part: ARM134289

 Waterfall, Desk Length Left, Full Length Right No Charge
 Part: ARM134288

 Waterfall, Desk Length Left and Right No Charge
 Part: ARM134282
1c. Falcon Gel Armpads
 Falcon Gel, Full Length Left and Right \$300
 Part: ARM136438

 Falcon Gel, Full Length Left, Desk Length Right \$250
 Part: ARM136442

 Falcon Gel, Desk Length Left, Full Length Right \$250
 Part: ARM136439

 Falcon Gel, Desk Length Left and Right \$200
 Part: ARM136437

Account Number: _____ Date: _____ Quote Name/Marked For: _____

13. HEADRESTS

Optional. Select a headrest pad and bracket.

13.1. Headrests, Comfort Plus with Cool Core

1. Headrest Pads

Stealth 8" Comfort Plus, Cool Core (CP08CCFPM) \$279
Part: ST-CP08CCFPM-Q. HCPCS: E0955

Stealth 10" Comfort Plus, Cool Core (CP10CCFPM)..... \$279
Part: ST-CP10CCFPM-Q. HCPCS: E0955

Stealth 14" Comfort Plus, Cool Core (CP14CCFPM) \$279
Part: ST-CP14CCFPM-Q. HCPCS: E0955

Stealth 6", Comfort Plus, Cool Core (CP06CCFPM)..... \$279
Part: ST-CP06CCFPM-Q. HCPCS: E0955

2. Headrest Brackets

Comfort Plus Headrest Mounting Bracket, UniLink..... \$485
Part: ST-UHB2D33200-Q. HCPCS: E1033

Multi-Axis Removable Headrest Bracket \$271
Part: ST-HMO475-17M-Q. HCPCS: E1033

13.2. Q Cam Lock

Optional.

Q Cam Lock..... No Charge
Part: XRF000376

Q Cam Lock Availability (Overall height)	TRU-Comfort/Sport		Q-Back 2	
	Standard	\$50	Standard	\$50
Tilt/Static/Back Canes	14"-21"H	22"-29"H	19"-24"H	25"-26"H
Power Recline Back	19"-21"H	22"-32"H	22"-25"H	26"-29"H

14. LAP BELTS

REQUIRED. Select one lap belt below.

☐ 38" Lap Belt ☐ 50" Lap Belt ☐ 60" Lap Belt ☐ 70" Lap Belt ☐ 80" Lap Belt ☐ 90" Lap Belt
Part: ACCBELT1006 Part: ACCBELT1000 Part: ACCBELT1002 Part: ACCBELT1003 Part: ACCBELT1004 Part: ACCBELT1001

15. LEGRESTS

Optional. All efforts will be made to accommodate legrest length selection, but minor changes may need to be made to meet shipping guidelines. Some adjustments may be needed upon receipt of the unit. Reference the section's compatibility matrix for legrest interference.

15.1a. Basic Center Mount Foot Platform

Select legrest length in following section.

1. Center Mount Foot Platform Type

Basic Center Mount..... No Charge
Part: XRF000323

Short	Medium
8"-10"	10"-14"
Basic Center Mount Extension Range	

2. Basic Center Mount Foot Platform Extension and Angle Adjustable Footplates

Short Extension, 5"x9" FootplatesNo Charge
Part: MEC1906675

Medium Extension, 5"x9" FootplatesNo Charge
Part: MEC1906680

3. Basic Center Mount Calf Pads

Calf pads are standard with a power articulating foot platform.

Small Calf Pads, Basic Center Mount.....\$210
Part: RIG2104369

4. Basic Center Mount Additional Accessories

Basic Center Mount Footplate Connector \$55
Part: RIG2104240

Legrests seating section continued on following pages.

15. LEGRESTS CONTINUED**15.1b. Center Mount Foot Platform**

Select legrest length in following section.

1. Foot Platform AssemblyQuick Adjustable Center Mount Foot Platform..... \$265
Part: XRF000311. HCPCS: K0040Power Articulating Foot Platform \$3,145
Part: PTOINDV3754. HCPCS: E1012+K0040**2. Extension and Angle Adjustable Footplates**

Pediatric extensions are not compatible with Power Articulating Foot Platform.

Pediatric Extension, 5"x6" Footplates
Part: RIG152577Medium Extension, 5"x11" Footplates
Part: RIG2404018Pediatric Extension, 5"x8" Footplates
Part: RIG2404015Long Extension, 5"x6" Footplates
Part: RIG151505Medium Extension, 5"x6" Footplates
Part: RIG150206Long Extension, 5"x8" Footplates
Part: RIG2404017Medium Extension, 5"x8" Footplates
Part: RIG2404016Long Extension, 5"x11" Footplates
Part: RIG2404019

Medium	Long
11"-17"	14"-20"
Power AFP Extension Range	

Pediatric	Medium	Long
8.5"-12.5"	11"-16.5"	14"-20"
Quick Adj Center Mount Extension Range		

3a. Center Mount Calf Pads

Calf pads are standard and required with a power articulating foot platform.

Small Calf Pads, AFP Only, 6"Wx6"H (pair) Standard
Part: RIGASMB7120018Large Calf Pads, AFP Only, 8" taper to 6"Wx8"H (pair) Standard
Part: RIGASMB7120019Small Calf Pads, Quick Adjustable Center Mount Only, \$210
6"Wx6"H (pair)
Part: XRF000156Large Calf Pads, Quick Adjustable Center Mount Only, \$210
8" taper to 6"Wx8"H (pair)
Part: XRF000157**3b. Center Mount Calf Panel**Small Calf Panel, 18"Wx8"H..... \$210
Part: FRMASMB7110012Large Calf Panel, 24"Wx8"H \$265
Part: RIG123825
Only available on widths of 20" or greater.**4. Center Mount Modular Footplate Extension Components**

Optional. Assembly required. May require seating configuration adjustments to avoid front caster interference.

2" Bell-Shaped Length Extension \$55
Part: RIG152519

Only available with the 5x6 footplates.

Width Extension for 8" \$55
Part: RIG2404919Width Extension for 11" \$55
Part: RIG2404920Outside Edge \$55
Part: RIG2404921

Available with 5"x8" and 5"x11" footplates.

Front and Rear Cap \$16
Part: RIG2404923

Available with 5"x8" and 5"x11" footplates.

5. Center Mount Additional AccessoriesSmall Foam Calf Pad Cover (pair)..... \$150
Part: ST-CSCS0775-QSmall Foam Panel Cover (pair) \$150
Part: ST-CSCP0920-QFootplate Connector \$150
Part: RIGASMB7120023Footplate Riser \$150
Part: XRF000348

Select legrest length between 2" to 9" at 1" increments.

Footplate Riser Calf Panel No Charge
Part: FRMASMB7110011

Only available and required with Footplate Riser selection and extension lengths of 6" to 9".

6. Stealth Gel Footplate Covers5"x8" Tapered Gel Footplate Covers \$275
Part: ST-FPG5008T-Q. HCPCS: K01085"x8" Tapered Gel Footplate Covers for Width Extension \$275
Part: ST-AFPG5008E-Q. HCPCS: K01085"x11" Tapered Gel Footplate Covers..... \$275
Part: ST-AFPG5011T-Q. HCPCS: K01085"x11" Tapered Gel Footplate Cover for Width Extension \$275
Part: ST-AFPG5011E-Q. HCPCS: K0108**15.1c. Swing-Away Legrests**

Select a legrest upper, a legrest extension, and type of footplate. Pair (P), left (L), and right (R) configurations are selected by checking the appropriate boxes.

1a. 60°-90° Swing-Away Legrest Upper

Adjustable in 5° increments.

☐	☐	☐	60°-90° Swing-Away Legrest, Static Seat \$142.50 ea. Pair part: XRF000082. Left: XRF000083. Right: XRF000084 Not available with 11"-13" extensions.
P	L	R	

☐	☐	☐	60°-90° Swing-Away Legrest, Short, Static Seat \$142.50 ea. Pair part: XRF000085. Left: XRF000086. Right: XRF000087 For use only with 11"-13" lower extensions.
P	L	R	

2. Lower Extension

☐	☐	☐	11"-13" Length Pair part: XRF000088. Left: XRF000089. Right: XRF000090.
P	L	R	

☐	☐	☐	13"-16.5" Length Pair part: XRF000094. Left: XRF000095. Right: XRF000096.
P	L	R	

☐	☐	☐	15.5"-19.5" Length Pair part: XRF000097. Left: XRF000098. Right: XRF000099.
P	L	R	

Legrests section continued on following pages.

15. LEGRESTS CONTINUED**3. Angle Adjustable Footplate(s)**

¹The 6"x11" footplates may cause interference with seat widths less than 18". ²The 8"x11" footplates may cause interference with seat widths less than 20".

☐ ☐ ☐	5"Wx8"D with Heel Loops..... \$87.50 ea. Pair part: FRMASMB4735. Left: RIG137791. Right: RIG137792 HCPCS: K0040+E0951	☐ ☐ ☐	6"Wx11"D¹ with Heel Loops..... \$87.50 ea. Pair part: FRMASMB4741. Left: RIG137796. Right: RIG137797 HCPCS: K0040+E0951
☐ ☐ ☐	6"Wx8"D¹ with Heel Loops..... \$87.50 ea. Pair part: FRMASMB4738. Left: RIG137793. Right: RIG137795 HCPCS: K0040+E0951	☐ ☐ ☐	8"Wx11"D² with Heel Loops..... \$87.50 ea. Pair part: FRMASMB4744. Left: RIG137798. Right: RIG137799 HCPCS: K0040+E0951

15.1d. Elevating Legrests

Select a legrest upper, a legrest extension, and type of footplate. Pair (P), left (L), and right (R) configurations are selected by checking the appropriate boxes. Select legrest length in following section.

1. Elevating Legrests Upper

☐ ☐ ☐	Elevating Legrest Assembly \$275 ea. Pair part: FRMASMB6483. Left: RIG138012. Right: RIG138013. HCPCS: E0990 Select lower extension and footplates
--	---

2a. Lower Extension

☐ ☐ ☐	11"-14" Length Pair: FRMASMB4755. Left: RIG138026. Right: RIG138030 For reinforced only	☐ ☐ ☐	14"-17" Length Pair: FRMASMB4756. Left: RIG138033. Right: RIG138035 For reinforced only
☐ ☐ ☐	8.5"-12" HD Drop-Ins/12.5"-16.5" ELR, Ext Length Pair part: XRF000101. Left: XRF000102. Right: XRF000103 For angle adjustable only	☐ ☐ ☐	13"-18" HD Drop-Ins/15.5"-19.5" ELR, Ext Length Pair part: XRF000104. Left: XRF000105. Right: XRF000106 For angle adjustable only

3a. Reinforced Cast Footplate(s)

☐ ☐ ☐	7.25"Wx6"D Reinforced Footplate..... Standard Pair part: FRMASMB4748. Left: RIG138055. Right: RIG138056	☐ ☐ ☐	8.5"Wx6"D Reinforced Footplate Standard Pair part: FRMASMB4747. Left: RIG138060. Right: RIG138062 May cause interference on seat widths less than 18."
--	---	--	---

3b. Angle Adjustable Footplate(s)

¹The 6"x11" footplates may cause interference with seat widths less than 18". ²The 8"x11" footplates may cause interference with seat widths less than 20".

☐ ☐ ☐	5"Wx8"D with Heel Loops..... \$87.50 ea. Pair part: FRMASMB4735. Left: RIG137791. Right: RIG137792 HCPCS: K0040+E0951	☐ ☐ ☐	6"Wx11"D¹ with Heel Loops..... \$87.50 ea. Pair part: FRMASMB4741. Left: RIG137796. Right: RIG137797 HCPCS: K0040+E0951
☐ ☐ ☐	6"Wx8"D¹ with Heel Loops..... \$87.50 ea. Pair part: FRMASMB4738. Left: RIG137793. Right: RIG137795 HCPCS: K0040+E0951	☐ ☐ ☐	8"Wx11"D² with Heel Loops..... \$87.50 ea. Pair part: FRMASMB4744. Left: RIG137798. Right: RIG137799 HCPCS: K0040+E0951

15.2. Legrest Accessories

Heel Loops (pair)..... \$35 Part: FRMASMB7873. HCPCS: E0951 Standard on angle adjustable and multi axis footplates when noted.	Wheel Bumper..... \$15 ea. Part: XRF000100
---	--

15.3. Set Legrest Extension Length

REQUIRED with center mounts and matching left and right rigging configurations. Part number is PTOINDV plus number following # below.

1a. Centermount Extension Length

2" L (#3989)	4" L (#3991)	6" L (#3993)	8" L (#3680)	9" L (#3682)	10" L (#3684)	11" L (#3686)	12" L (#3688)	13" L (#3690)	14" L (#3692)	15" L (#3694)
3" L (#3990)	5" L (#3992)	7" L (#3994)	8.5" L (#3681)	9.5" L (#3683)	10.5" L (#3685)	11.5" L (#3687)	12.5" L (#3689)	13.5" L (#3691)	14.5" L (#3693)	

1b. Front Rigging Extension Length

Non-matching or individual legrest configurations lengths are not set by the manufacturer. See previous section for center mount extension lengths.

7" L (#3470)	8" L (#3472)	9" L (#3474)	10" L (#3476)	11" L (#3478)	12" L (#3480)	13" L (#3482)	14" L (#3484)	15" L (#3486)
7.5" L (#3471)	8.5" L (#3473)	9.5" L (#3475)	10.5" L (#3477)	11.5" L (#3479)	12.5" L (#3481)	13.5" L (#3483)	14.5" L (#3485)	

Account Number: _____ Date: _____ Quote Name/Marked For: _____

16. POWER OPTION ELECTRONICS

REQUIRED for articulating foot platform.

Single Actuator Through Controller Interface No Charge
Part: XRF000335. HCPCS: E2310

Single Actuator Through Controller and Interface (AAM)..... \$1,675
Part: XRF000068. HCPCS: E2310

17. SEAT MOUNTED LED LIGHTS

Optional. Requires Q-Logic 3 Joystick for Full Light Package or Q-Logic 3e Joystick.

Seat Mounted, LED, Full Light Package\$770
Part: XRF000130

Not available with Base Mounted, LED, Full Light Package.

18. ACCESSORIES

Optional. Multiple accessories can be mounted along the various tracks but some interference may occur. Accessories may need to be installed by the provider.

Flex Seat Cup Holder \$20
Part: ACCASMB2634

Not available on TRU Balance Flex seating without two-post armrests.

Oxygen Holder for Flex Seat\$110
Part: ACC125006. HCPCS: E2208

Not available with Manual Recline.

Cane/Crutch Holder\$195
Flex part: ACC152515. HCPCS: E2207

Flag Holder\$125
Flex part: ACCASMB2415

Saddle Bag \$15
Part: DWR1020D010

Stealth Cell Phone/Media Device Holder \$400
Part: ACC7120006

Stealth Large Tablet Holder \$565
Part: ACC141818

Quantum Backpack \$200
Part: ACC1905853

Grip Lap Board (20005)..... \$100
Part: ACC152992

Grip Drink Holder (20085) \$100
Part: ACC152993

Accessory Bar\$90
Part: XRF000128

Large, Black, Weather Cover \$240
Part: ACCCOVR1021

Privacy Bag\$80
Part: ACC2005723

Transfer Handle \$152.50 ea.
Part: MEC2205680 x2. HCPCS: E0705

Aftermarket Accessory Mounting Bracket

Select a quantity up to 5.

Aftermarket Accessory Mounting Bracket, Slotted \$100 ea.
Part: ACC2407046

QTY

Aftermarket Accessory Mounting Bracket, Non-slotted \$100 ea.
Part: ACC2405496

QTY