



Pride Mobility Products Co.

5096 South Service Rd, Beamsville, ON, L3J 1V4

Phone: 888-570-1113 | Fax: 866-514-1303 | Email: productorders@pridemobility.com

Q6 Edge® 3 Order Form with TRU Balance® 3 Positioning

Account Number: _____ Date: _____

Provider Name: _____

Contact: _____

Phone: _____ Fax: _____

Email: _____

PO Number: _____

Marked for: _____

Ship to Address: _____

City: _____ Prov: _____ PC: _____

EDGE 3



INTRODUCTION

This form is interactive when viewed with [Adobe Acrobat Reader](#) and may not function correctly if opened with applications other than Acrobat. Complete the form by placing checks in the desired boxes and provide information in the interactive fields. Buttons shown at the bottom of the form may be utilized to print or submit the order form through a desktop email application. To email via a web-based application, please 'Save As' and attach the PDF to your email.

Thank you for choosing Quantum!


PATIENT INFORMATION

Weight: _____ Height: _____

A. Shoulder Width: _____ D. Max Sitting Height: _____ G. Thigh Depth: _____

B. Chest Width: _____ E. Shoulder Height: _____ H. Lower Leg Length: _____

C. Hip Width: _____ F. Axilla Height: _____ I. Elbow Height: _____



1. BASE MODEL

REQUIRED. Select a group 3 power base. 300 lb. weight capacity with standard 6mph programmable motors and 8 amp off board charger. Dimensions of base without legrests: 24"Wx35.5"L

Quantum Q6 Edge 3\$8,450

Part: Q6 EDGE 3. HCPCS: K0848

2. DRIVE MOTORS

REQUIRED.

6 mph Drive Motors Standard
Part: MOT147433

6 mph, Accu-Trac Motor Technology \$995
Part: MOT147881. HCPCS: K0108
Requires Q-Logic 3 electronics. Delivers enhanced traction and tracking, especially with switched drive controls.

7 mph, 4 Pole Drive Motors..... \$700
Part: MOT1902455

Requires Q-Logic 3 electronics.

7 mph, 4 Pole, Accu-Trac Motor Technology..... \$1,695
Part: MOT2008501. HCPCS: K0108
Requires Q-Logic 3 electronics. Delivers enhanced traction and tracking, especially with switched drive controls.

3. BASE COLOR

REQUIRED. Select one shroud color. Colors shown here may differ from actual product depending on monitor/printer calibration used.

3.1a. Standard Shrouds

 Red Corvette
Part: DGN1801101-PRD1

 Deep Purple
Part: DGN1801101-PPR5

 Midnight Blue
Part: DGN1801101-PBU10

 Green Machine
Part: DGN1801101-PGN11

 Back in Black
Part: DGN1801101-PBK1

 Raspberry Beret
Part: DGN1801101-PPK4

 Electric Blue
Part: DGN1801101-PBU9

 Military Green
Part: DGN1801101-PGN8

 Night Train
Part: DGN1801101-PBK3

 Pretty in Pink
Part: DGN1801101-PPK3

 Blue Jean
Part: DGN1801101-PBU11

 Grey Street
Part: DGN1801101-PSL3

 Evolution Orange
Part: DGN1801101-POR7

 Lemon Crush
Part: DGN1801101-PYL4

3.1b. Hydrographic Shrouds

☐  Military Digital Camo..... \$325
Part: DGN1801101-WCM7

☐  Faded Glory..... \$325
Part: DGN1801101-WFL2

☐  Carbon Fiber..... \$325
Part: DGN1801101-WCB1

4. BASE MOUNTED LIGHTS

Optional. Select one option. Seat Mounted, LED, Full Light Package option can be found later on this form.

Switch Fender Lights..... No Charge
Part: DGN154117


Fender lights operate by pressing included single switch button. Fender lighting is not available with Single Actuator Control through Toggle Box.

Base Mounted, LED, Full Light Package \$770
Part: XRF000129

Requires Q-Logic 3 EX Joystick for Full Lights or Q-Logic 3e Joystick.



Base Mounted Lights



Fender Lights

5. QUANTUM BACKUP CAMERA

Optional. Requires Q-Logic electronics. The backup camera may interfere with armrest mounted accessories such as media device holder.

Quantum Backup Camera \$475
Part: MEC1910475

1. Optional Camera Mounting Bracket

Backup Camera Headrest Bracket..... \$35
Part: MEC1910373

Works with .5" x .5" vertical supports only.

Backup Camera Track Mount Bracket \$35
Part: MEC2104227

2. Optional Display Mounting Bracket

Backup Camera Q-Logic Joystick Mount Bracket \$145
Part: MEC2106123



Backup Camera

6. DRIVE WHEELS

REQUIRED. Select a type of drive wheels and accent color. Knobby tires increase overall base width by 1.5".

6.1. Drive Wheels

3"x14" Flat Free Black Tires Standard
Part: PTOINDV3525 (pair)

3"x14" Pneumatic Black Tires No Charge
Part: WHL1713925 x2

4"x14" Knobby, Flat Free Black Tires \$120
Part: XRF000368

4"x14" Knobby, Pneumatic Black Tires..... \$120
Part: XRF000367

Drive wheels section continued on following page.

6. DRIVE WHEELS CONTINUED**6.2. Drive Wheel Accent Color****Optional.** Not available with knobby tires.Back in Black Standard
Part: DGN164245-PBK1 x2Red Corvette No Charge
Part: DGN164245-PRD1 x2Night Train No Charge
Part: DGN164245-PBK3 x2Grey Street No Charge
Part: DGN164245-PSL3 x2Deep Purple No Charge
Part: DGN164245-PPR5 x2Raspberry Beret No Charge
Part: DGN164245-PPK4 x2Pretty in Pink No Charge
Part: DGN164245-PPK3 x2Evolution Orange No Charge
Part: DGN164245-POR7 x2Midnight Blue No Charge
Part: DGN164245-PBU10 x2Electric Blue No Charge
Part: DGN164245-PBU9 x2Blue Jean No Charge
Part: DGN164245-PBU11 x2Lemon Crush No Charge
Part: DGN164245-PYL4 x2Green Machine No Charge
Part: DGN164245-PGN11 x2Military Green No Charge
Part: DGN164245-PGN8 x2Military Digital Camo \$50 ea.
Part: DGN164245-WCM7 x2Faded Glory \$50 ea.
Part: DGN164245-WFL2 x2Carbon Fiber \$50 ea.
Part: DGN164245-WCB1 x2**7. BATTERIES****Optional.** The powerchair requires 2 batteries and will be installed if ordered.NF-22, Gel Cell, Interceptor 2 Battery \$300 ea.
Part: BAT2004376. HCPCS: E2361NF-22, Gel Cell, MK Battery \$400 ea.
Part: BATGEL1008. HCPCS: E2361**8. PATIENT'S WEIGHT****REQUIRED.** The patient's weight is required to quote and order the power chair.Up to 200 lb User Weight
Part: PTOINDV3837201 lb-250 lb User Weight
Part: PTOINDV3861251 lb-300 lb User Weight
Part: PTOINDV3984**9. EXTENDED SEAT TOWERS****Optional.** Only available with TRU Balance Static Seat or Power Recline Only. Not available with base mounted foot platform.Extended Seat Towers No Charge
Part: MEC2107219**10. ELECTRONICS****REQUIRED.** Select a type of drive control from the following sections. The No Joystick option is used to omit the joystick and utilize another control, such as a specialty control at the end of this form or found on the [Alternative Controls and Electronics Order Form](#).**10.1a. Non-Expandable Controls**

Select a joystick below or omit selection if another controller type is desired. Non-expandable does not allow supplementary drive controllers.

QLNES 6 Key Joystick Standard
Part: CTL2205379
Controls up to 2 actuators.QLNES 10 Key Joystick No Charge
Part: CTL2205380
Controls up to 2 actuators through joystick and seat mounted, LED, full light package.**10.1b. Expandable Controls**

Select expandable controls or omit selection if another controller type is desired.

2. Q-Logic 3 Joystick

Q-Logic 3 EX Joystick includes XLR USB charger.

Q-Logic 3 EX Joystick \$900
Part: CTL167724Q-Logic 3e Joystick \$600
Part: CTL1910600
Can operate full light package. Not available with alternative or specialty controls.Q-Logic 3 EX Joystick for Full Lights \$900
Part: CTL167725No Joystick, Setup for Q-Logic
Part: XRF000010
Requires selection of Enhanced Display.**10.2. Q-Logic Magnetic Charger**

Available with QLNES. Included with Q-Logic 3e Joystick.

Magnetic to XLR Charger Adapter with Wall Holder \$100
Part: ELE2103635

Electronics section continued on following pages.

10. ELECTRONICS CONTINUED**10.3. Q-Logic 3 EX Joystick Armor****REQUIRED with Q-Logic 3 EX Joystick selections above. Not available with Q-Logic 3e Joystick.** Reference colors shown on page 2.Black, Soft TouchStandard
Part: DGN124307

Also matches Night Train shrouds.

Red Corvette No Charge
Part: DGN124336-PRD1Back in Black No Charge
Part: DGN124336-PBK1Grey Street No Charge
Part: DGN124336-PSL3Deep Purple No Charge
Part: DGN124336-PPR5Raspberry Beret No Charge
Part: DGN124336-PPK4Pretty in Pink No Charge
Part: DGN124336-PPK3Evolution Orange No Charge
Part: DGN124336-POR7Midnight Blue No Charge
Part: DGN124336-PBU10Electric Blue No Charge
Part: DGN124336-PBU9Blue Jean No Charge
Part: DGN124336-PBU11Lemon Crush No Charge
Part: DGN124336-PYL4Green Machine No Charge
Part: DGN124336-PGN11Military Green No Charge
Part: DGN124336-PGN8Military Digital Camo No Charge
Part: DGN124336-WCM7Faded Glory No Charge
Part: DGN124336-WFL2Carbon Fiber No Charge
Part: DGN124336-WCB1**11. JOYSTICK HANDLES****Optional.** The optional handle will be mounted to the joystick and the standard handle will be shipped with the unit.Mushroom Handle (B) \$120
Part: ACCASMB1860.Chin Cup (C) \$85
Part: ACCASMB186190T-Bar (D) \$160
Part: ACCASMB1858.Medium Ball Handle (E) \$70
Part: ACCASMB1859.Soft Ball Handle (F) \$120
Part: ACCASMB1856.Extended I-Handle (G) \$140
Part: ACCASMB1857.Bodypoint 3" Goal Post Handle \$180
Part: PLSKNCP1057.Bodypoint 4" Goal Post Handle \$180
Part: PLSKNCP1058.**12. JOYSTICK MOUNTING BRACKETS****Optional.** Select a mount to be used with a joystick selected in previous section.**12.1a. Fixed, In-line Joystick Bracket**Fixed, In-line, Bracket for Q-Logic Standard
Part: FRMASMB10617**12.1b. Swing-Away (Retractable) Joystick Bracket**

Allows the joystick to retract to the side of the armrest while maintaining the same joystick orientation.

Swing-away, Q-Logic \$265
Part: ACC1908434. HCPCS: E1032**12.1c. Basic Swing Away Joystick Mount**Basic Swing Away, Q-Logic, Left Side \$265
Part: ACC2201041. HCPCS: E1032Basic Swing Away, Q-Logic, Right Side \$265
Part: ACC2201042. HCPCS: E1032**12.1d. Offset Swing Away Joystick Mount**Offset Swing Away, Q-Logic, Left Side \$265
Part: ACC142742. HCPCS: E1028Offset Swing Away, Q-Logic, Right Side \$265
Part: ACC142743. HCPCS: E1028**12.1e. Specialty Joystick Mounts**


Select a bracket and side to mount it.

Retract4 Mount \$450
Part: XRF000003. HCPCS: E1032

Full 180° of movement that translates into multi-position use, straight, right offset, left offset, or fold away. Joystick orientation varies as bracket is adjusted.

Multi-Axis Mount \$100
Part: XRF000004

Allows height adjustment of joystick.

12.2. Joystick Mount Location☐ Left Side Mount
Part: PTOINDV1225☐ Right Side Mount
Part: PTOINDV1226

13. TRU BALANCE® 3 SEATING**REQUIRED.** Maximum seat size is 22"W x 22"D.**13.1. TRU Balance 3 Configuration**

TRU Balance 3 Static Seat.....\$1,045
 Part: PTOINDV3442
 20° of limited recline adjustment (-5/ +15) via gas strut within back cover.

TRU Balance 3 Power Tilt Only\$6,245
 Part: PTOINDV3443. HCPCS: E1002
 50° Center of Gravity (CG) power tilt seating system.

TRU Balance 3 Power Recline Only\$7,875
 Part: PTOINDV3444. HCPCS: E1004
 175° of power recline.

TRU Balance 3 Power Tilt & Recline\$13,125
 Part: PTOINDV3446. HCPCS: E1007
 50° CG power tilt with 175° of power recline. Full drive lockout at 20°.

13.2. 12" Power Lift**Optional.** ^Please verify the appropriate billing code to be used for this item with all third-party payors.

TRU Balance 3, 12" Power Lift with iLevel®\$2,995
 Part: PTOINDV3773. HCPCS: E2298/E2300^
 12" of scissor lift with iLevel System.

**13.3. Occupied Transit Kit****Optional.** The power base is standard with Unoccupied Transit Loops. 60" Sure-Lok Lap Belt is required with Occupied Transit Option.

WC-19 Occupied Transit Option.....\$400
 Part: KIT152047

To determine if your chair is WC-19 compliant (occupied), please refer to this website: <https://wc-transportation-safety.umtri.umich.edu/>

13.4. Seat Width**REQUIRED with seat selection.** Maximum seat width is 22". Part number is SETWDTH plus number following # below. 12" Width is not available with 16" and 17" Depth. 13" and 15" widths are only available with non-matching back width.

12" W (#1002)	13" W (#1003)	14" W (#1004)	15" W (#1005)	16" W (#1006)	17" W (#1007)	18" W (#1008)	19" W (#1009)	20" W (#1010)	21" W (#1011)	22" W (#1012)
------------------	------------------	------------------	------------------	------------------	------------------	------------------	------------------	------------------	------------------	------------------

13.5. Seat Width Package**REQUIRED with seat selection.**

12"-16" Width Package
 Part: XRF000006

14"-16" Width Package
 Part: XRF000007

16"-20" Width Package
 Part: XRF000008

20"-24" Width Package
 Part: XRF000158
 Maximum width is 22".

13.6. Seat Depth**REQUIRED with seat selection.** Maximum seat depth is 22". Part number is SETDPTH plus number following # below.

12" D (#1002)	13" D (#1003)	14" D (#1004)	15" D (#1005)	16" D (#1006)	17" D (#1007)	18" D (#1008)	19" D (#1009)	20" D (#1010)	21" D (#1011)	22" D (#1012)
------------------	------------------	------------------	------------------	------------------	------------------	------------------	------------------	------------------	------------------	------------------

13.7. Seat Depth Package**REQUIRED with seat selection.** Maximum depth is 22".

12"-16" Depth Package
 Part: XRF000011

14"-16" Depth Package
 Part: XRF000012

16"-20" Depth Package
 Part: XRF000013

20"-24" Depth Package
 Part: XRF000159
 Maximum depth is 22".

Low	Med	High
Static Seat, Recline		
16.5"	17.5"	18.5"
Tilt, Tilt & Recline, Lift		
17.5"	18.5"	19.5"
TRU Balance 3 Seat to Floor Heights		

13.8. Seat to Floor Height**REQUIRED.**


Lowest Seat to Floor
 Part: PTOINDV3558

Middle Seat to Floor
 Part: PTOINDV3559

Highest Seat to Floor
 Part: PTOINDV3560

13.9. Stealth Seat Cushion

Optional. Specific cushion/component sizes available, part numbers, weight capacities, MSRP, HCPCS code, and special order cushions will be found on the respective complete order forms. Part numbers in this section are for reference only. This section allows a quick selection of a single cushion that matches the seat pan width and depth to alleviate listing all part numbers and codes on this power base form. The complete cushion part number and MSRP is a product of the cushion model, top foam, liner, cover, width, and depth as shown below. Reference the subsequent chart for available top foam, liner, and cover options for each cushion type. See specific cushion order form for all sizes, top foam, liners, covers, parts numbers, and MSRP.



TRU Balance 3 Seating section continued on following pages.

13. TRU BALANCE 3 SEATING CONTINUED**1. Cushion Type**

TRU-Comfort 2 SPP
Reference part: XRF000116.
Part series: TC2-###XXYY.

Essence SPP
Reference part: XRF000118.
Part series: ESS-###XXYY.

Essence Plus SPP
Part series: ESSP-###XXYY.

Solution SPP
Reference part: XRF000113.
Part series: SOL-###XXYY.

Spectrum Foam SPP
Reference part: XRF000112.
Part series: SPCF-###XXYY.

Spectrum Gel SPP
Reference part: XRF000117.
Part series: SPCG-###XXYY.

Premiere P
Reference part: XRF000122.
Part series: PREM-###XXYY.

Glacial SP
Reference part: XRF000114.
Part series: GLA-###XXYY.

Zen SP
Reference part: XRF000115.
Part series: ZEN-###XXYY.

Simplicity G
Reference part: XRF000111.
Part series: SIM-###XXYY.

2. Essence Cushion Wedge Kit

Only available with Essence cushion selection above.

Small Wedge Positioning Kit\$45
Part: POP1909006-S
14"W and smaller

Small Wedge Positioning Kit\$45
Part: POP1712340-S
15"W and larger

Medium Wedge Positioning Kit\$85
Part: POP1909006-M
14"W and smaller

Medium Wedge Positioning Kit\$85
Part: POP1712340-M
15"W and larger

Large Wedge Positioning Kit\$90
Part: POP1909006-L
14"W and smaller

Large Wedge Positioning Kit\$90
Part: POP1712340-L
15"W and larger

3. Cushion Top Foam

Visco Top Foam
Reference part: XRF000380

FRG Top Foam
Reference part: XRF000381

No Top Foam
Reference part: XRF000382

4. Cushion Liner

Stealth Polyurethane Liner
Reference part: XRF000384

No Liner
Reference part: XRF000385

5. Cushion Cover

CoolCore Cover
Reference part: XRF000386

Mesh Cover
Reference part: XRF000387

Stealth Tek Cover
Reference part: XRF000388

6. Optional Base Insert

Option selection available with Glacial, Premiere, Simplicity, Solution, Spectrum Foam, Spectrum Gel, and Zen cushions.

ABS Flat Rigidizer
Reference part: XRF000389

7. Additional Cushion Liner

Up to 5 may be selected below. See respective cushion order form for MSRP and part numbers.

Additional Stealth Urethane Liner
Reference Part: XRF000390
QTY _____

Additional Stealth Polyurethane Liner
Reference Part: XRF000391
QTY _____

8. Additional Cushion Cover

Up to 5 may be selected below. See respective cushion order form for MSRP and part numbers.

Additional CoolCore Cover
Reference Part: XRF000392
QTY _____

Additional Mesh Cover
Reference Part: XRF000393
QTY _____

Additional Stealth Tek Cover
Reference Part: XRF000394
QTY _____

13.10. Back Type Selection**REQUIRED.** Match a choice in this section with the options in the following back cushion configuration section.

Redesigned, Static Type Back \$735
Part: XRF000014
12" and 24" back width will use the TRU Balance 3 Back assembly. Select a back cushion configuration in the subsequent section.

Back Cane Option
Part: XRF000015
Available only with TRU Balance 3 Static Seat and Power Tilt (1 actuator) selections. Select cane type, back cane height, and back cushion configuration below.

Redesigned, Flip Forward Back \$735
Part: XRF000016
12" and 24" back width will use the TRU Balance 3 Back assembly. Not available with power recline systems. Only available with 2-post, flip back or pediatric armrests. Not available with 12"-14" width. Select a back cushion type in the subsequent section.

	None	Visco/V Elastic	FRG	None	Urethane	Polyurethane	CoolCore	Mesh	Stealth Tek
	Top Foam			Liner			Cover		
TRU-Comfort 2	—	✓	—	✓	—	✓	✓	—	—
Essence	—	✓	—	✓	—	✓	✓	✓	—
Solution	—	✓	—	—	✓	✓	✓	✓	✓
Spectrum Foam	—	✓	✓	—	✓	✓	✓	✓	✓
Spectrum Gel	—	✓	✓	—	✓	✓	✓	✓	✓
Premiere	—	✓	✓	✓	✓	✓	✓	—	—
Glacial	✓	—	—	✓	✓	✓	✓	✓	✓
Zen	—	✓	✓	✓	✓	✓	✓	✓	✓
Simplicity	✓	—	—	✓	—	—	—	—	✓

TRU Balance 3 Seating section continued on following pages.

13. TRU BALANCE 3 SEATING CONTINUED**13.11. Back Cane Selection (Back Cane Option Only)****REQUIRED when Back Cane Option is selected above.** Select type and height of back canes.**1. Back Cane Type**Angle Adjustable Back Cane No Charge
Part: MEC140266Flip Forward Back Cane No Charge
Part: XRF0002732" Endomorph (1" per side) Back Cane No Charge
Part: XRF0002694" Endomorph (2" per side) Back Cane No Charge
Part: XRF000271**2. Back Cane Height**16"H Cane
Part: FRMASMB1583420"H Cane
Part: FRMASMB1583824"H Cane
Part: FRMASMB1584217"H Cane
Part: FRMASMB1583521"H Cane
Part: FRMASMB1583925"H Cane
Part: FRMASMB15843**13.12. Back Cushion Configuration****REQUIRED.** Sport Back and TRU-Comfort Back selections up to 22" back width will be the Sport 4 or TRU-Comfort 4 back cushions. Sport back and TRU-Comfort back selections of 24" back width will be the Sport or TRU-Comfort 2 back cushions. Select a configuration type and applicable dimensions. See the Sport, Q-Back, and TRU-Comfort back availability chart for stock order back size ranges. Other size backs may be available with special order and extended lead time.**1a. Sport Back/Omit/After-Market Interface**Sport Back No Charge
Part: XRF000017

Complete following width package and size.

Omit Back Pan/Cushion No Charge
Part: PTOINDV1287

Select a following width package.

After-market Back Interface No Charge
Part: ST-SSAPTB4-Q. Select a set width in following section, range of 12"-21". Overall seat depth may need to be adjusted to accommodate certain backs. May interfere with gel/waterfall armpads and arm through. Does not include brackets to mount after-market backs.**1b. TRU-Comfort Back**TRU-Comfort Back (Static/Flip Forward), 12"-15" Width \$835
Part: SET123622. HCPCS: E2620

Complete following width package and size.

TRU-Comfort Back (Static/Flip Forward), 16" to 20" Width ... \$1,035
Part: TC4V10W21W. HCPCS: E2620

Complete following width package and size.

TRU-Comfort Back, Back Cane Option \$835
Part: TC4V10W21WBC. HCPCS: E2620

Complete following width package and size.

TRU-Comfort Back (Static/Flip Forward), 22" Width \$1,235
Part: TC4V22W32W. HCPCS: E2621

Complete following width package and size.

Overall Height	Tilt Only Width												Recline Width											
	12"	14"	15"	16"	17"	18"	19"	20"	21"	22"	12"	14"	15"	16"	17"	18"	19"	20"	21"	22"				
14"	S	T	S	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	14"			
15"	S	T	S	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	15"			
16"	S	S	T	S	T	S	T	S	T	S	T	S	T	-	-	-	-	-	-	-	16"			
17"	S	S	T	S	T	S	T	S	T	S	T	S	T	-	-	-	-	-	-	-	17"			
18"	-	S	T	S	T	S	T	S	T	S	T	S	T	-	-	-	-	-	-	-	18"			
19"	-	S	T	S	T	S	T	S	T	S	T	S	T	S	S	T	S	T	S	T	19"			
20"	-	S	T	S	T	S	T	S	T	S	T	S	T	S	S	T	S	T	S	T	20"			
21"	-	S	T	S	T	S	T	S	T	S	T	S	T	S	S	T	S	T	S	T	21"			
22"	-	S	T	S	T	S	T	S	T	S	T	S	T	-	S	T	S	T	S	T	22"			
23"	-	S	T	S	T	S	T	S	T	S	T	S	T	-	S	T	S	T	S	T	23"			
24"	-	S	T	S	T	S	T	S	T	S	T	S	T	-	S	T	S	T	S	T	24"			
25"	-	S	T	S	T	S	T	S	T	S	T	S	T	-	S	T	S	T	S	T	25"			
26"	-	S	T	S	T	S	T	S	T	S	T	S	T	-	S	T	S	T	S	T	26"			
27"	-	S	T	S	T	S	T	S	T	S	T	S	T	-	S	T	S	T	S	T	27"			
28"	-	S	T	S	T	S	T	S	T	S	T	S	T	-	S	T	S	T	S	T	28"			
29"	-	S	T	S	T	S	T	S	T	S	T	S	T	-	S	T	S	T	S	T	29"			
30"	-	-	-	-	-	-	-	-	-	-	-	-	-	S	T	S	T	S	T	30"				
31"	-	-	-	-	-	-	-	-	-	-	-	-	-	S	T	S	T	S	T	31"				
32"	-	-	-	-	-	-	-	-	-	-	-	-	-	S	T	S	T	S	T	32"				

 Sport Back

 Q-Back

 TRU-Comfort Back

Sport, Q-Back, and TRU-Comfort Back Availability Chart

S Sport Back
Q Q-Back
T TRU-Comfort Back

Sport, Q-Back, and TRU-Comfort Back Availability Chart

TRU Balance 3 Seating section continued on following pages.

13. TRU-BALANCE 3 SEATING CONTINUED

1c. ADI2.0 Backrest

Specify back type, width, height, and contour in the part number configurator below. Previously selected back cane height must be equal to the greater of overall height range. **Available SOF for Provider installation only.**

For more information please see [ADI2.0 Backrest order form](#)

ADI2.0 Backrest 10"W - 20"W \$1,495
For use with back canes only.
Indicate part number below.
ADP: SE0002152

ADI2.0 Backrest 22"W - 24"W \$1,695
For use with back canes only.
Indicate part number below.
ADP: SE0002153

Specify Part Number for Back Selection

OPTIONS INDEX	ADI2	Back Type: AL AFT	Width: 10"-20" (1" increments), 22", 24"	Height: N = Nino L = Low M = Medium S = Standard T = Tall X = Extra Tall	Contour: A = Active D = Deep E = Extra Deep	2 = FRG Foam	P = Quick-Release
ENTER SELECTION	ADI2					2	P

BACK CONTOUR*:



BACK SERIES TYPES*:

(AL) Aluminum Series

(*ALL BACKS SHOWN ARE 18" WIDE)

(AFT) Aluminum Flat Top Series



BACK HEIGHTS AND WIDTH:

(Available for AL Back Shell Only)

Nino: 8" inch (10"W - 14"W)

Low: 10 inch (10"W - 24"W)

Medium: 13 inch (12"W - 24"W)

(Available for AL and AFT Shell Types)

Standard: 16 inch (13"W - 24"W)

Tall: 18 inch (16"W - 24"W)

Extra Tall: 20 inch (16"W - 24"W)

13. TRU BALANCE 3 SEATING CONTINUED

1f. Q-Back 2.0, AL Series, Active Contour

Not available with Back Cane Option. Select a back type and the specific overall back height in subsequent section.

Active Contour, 14"Wx19-20"H Tilt/22"-23" H Recline \$9,85
Part: QB2AL14AS. HCPCS: E2613

Active Contour, 14"Wx21-22"H Tilt/24"-25" H Recline \$985
Part: QB2AL14AT. HCPCS: E2613

Active Contour, 16"Wx19-20"H Tilt/22"-23" H Recline \$1,145
Part: QB2AL16AS. HCPCS: E2613

Active Contour, 16"Wx21-22"H Tilt/24"-25" H Recline \$1,145
Part: QB2AL16AT. HCPCS: E2613

Active Contour, 16"Wx23-24"H Tilt/26"-27" H Recline \$1,145
Part: QB2AL16AX. HCPCS: E2613

Active Contour, 16"Wx25-26"H Tilt/28"-29" H Recline \$1,145
Part: QB2AL16AXX. HCPCS: E2613

Active Contour, 18"Wx19-20"H Tilt/22"-23" H Recline \$1,145
Part: QB2AL18AS. HCPCS: E2613

Active Contour, 18"Wx21-22"H Tilt/24"-25" H Recline \$1,145
Part: QB2AL18AT. HCPCS: E2613

Active Contour, 18"Wx23-24"H Tilt/26"-27" H Recline \$1,145
Part: QB2AL18AX. HCPCS: E2613

Active Contour, 18"Wx25-26"H Tilt/28"-29" H Recline \$1,145
Part: QB2AL18AXX. HCPCS: E2613

Active Contour, 20"Wx19-20"H Tilt/22"-23" H Recline \$1,145
Part: QB2AL20AS. HCPCS: E2613

Active Contour, 20"Wx21-22"H Tilt/24"-25" H Recline \$1,145
Part: QB2AL20AT. HCPCS: E2613

Active Contour, 20"Wx23-24"H Tilt/26"-27" H Recline \$1,145
Part: QB2AL20AX. HCPCS: E2613

Active Contour, 20"Wx25-26"H Tilt/28"-29" H Recline \$1,145
Part: QB2AL20AXX. HCPCS: E2613

TRU Balance 3 Seating section continued on following pages.

13. TRU BALANCE 3 SEATING CONTINUED**1g. Q-Back 2.0, AL Series, Deep Contour****Not available with Back Cane Option.** Select a back type and the specific overall back height in subsequent section.

Deep Contour, 14"Wx19-20"H Tilt/22"-23" H Recline \$1,145 Part: QB2AL14DS. HCPCS: E2620	Deep Contour, 18"Wx21-22"H Tilt/24"-25" H Recline \$1,325 Part: QB2AL18DT. HCPCS: E2620
Deep Contour, 14"Wx21-22"H Tilt/24"-25" H Recline \$1,145 Part: QB2AL14DT. HCPCS: E2620	Deep Contour, 18"Wx23-24"H Tilt/26"-27" H Recline \$1,325 Part: QB2AL18DX. HCPCS: E2620
Deep Contour, 16"Wx19-20"H Tilt/22"-23" H Recline \$1,325 Part: QB2AL16DS. HCPCS: E2620	Deep Contour, 18"Wx25-26"H Tilt/28"-29" H Recline \$1,325 Part: QB2AL18DXX. HCPCS: E2620
Deep Contour, 16"Wx21-22"H Tilt/24"-25" H Recline \$1,325 Part: QB2AL16DT. HCPCS: E2620	Deep Contour, 20"Wx19-20"H Tilt/22"-23" H Recline \$1,325 Part: QB2AL20DS. HCPCS: E2620
Deep Contour, 16"Wx23-24"H Tilt/26"-27" H Recline \$1,325 Part: QB2AL16DX. HCPCS: E2620	Deep Contour, 20"Wx21-22"H Tilt/24"-25" H Recline \$1,325 Part: QB2AL20DT. HCPCS: E2620
Deep Contour, 16"Wx25-26"H Tilt/28"-29" H Recline \$1,325 Part: QB2AL16DXX. HCPCS: E2620	Deep Contour, 20"Wx23-24"H Tilt/26"-27" H Recline \$1,325 Part: QB2AL20DX. HCPCS: E2620
Deep Contour, 18"Wx19-20"H Tilt/22"-23" H Recline \$1,325 Part: QB2AL18DS. HCPCS: E2620	Deep Contour, 20"Wx25-26"H Tilt/28"-29" H Recline \$1,325 Part: QB2AL20DXX. HCPCS: E2620

1h. Q-Back 2.0, AL Series, Extra Deep Contour**Not available with Back Cane Option.** Select a back type and the specific overall back height in subsequent section.

Extra Deep Contour, 14"Wx19-20"H Tilt/22"-23" H Recline \$1,225 Part: QB2AL14ES. HCPCS: E2620	Extra Deep Contour, 18"Wx21-22"H Tilt/24"-25" H Recline \$1,410 Part: QB2AL18ET. HCPCS: E2620
Extra Deep Contour, 14"Wx21-22"H Tilt/24"-25" H Recline \$1,225 Part: QB2AL14ET. HCPCS: E2620	Extra Deep Contour, 18"Wx23-24"H Tilt/26"-27" H Recline \$1,410 Part: QB2AL18EX. HCPCS: E2620
Extra Deep Contour, 16"Wx19-20"H Tilt/22"-23" H Recline..... \$1,410 Part: QB2AL16ES. HCPCS: E2620	Extra Deep Contour, 18"Wx25-26"H Tilt/28"-29" H Recline \$1,410 Part: QB2AL18EXX. HCPCS: E2620
Extra Deep Contour, 16"Wx21-22"H Tilt/24"-25" H Recline..... \$1,410 Part: QB2AL16ET. HCPCS: E2620	Extra Deep Contour, 20"Wx19-20"H Tilt/22"-23" H Recline \$1,410 Part: QB2AL20ES. HCPCS: E2620
Extra Deep Contour, 16"Wx23-24"H Tilt/26"-27" H Recline \$1,410 Part: QB2AL16EX. HCPCS: E2620	Extra Deep Contour, 20"Wx21-22"H Tilt/24"-25" H Recline \$1,410 Part: QB2AL20ET. HCPCS: E2620
Extra Deep Contour, 16"Wx25-26"H Tilt/28"-29" H Recline \$1,410 Part: QB2AL16EXX. HCPCS: E2620	Extra Deep Contour, 20"Wx23-24"H Tilt/26"-27" H Recline \$1,410 Part: QB2AL20EX. HCPCS: E2620
Extra Deep Contour, 18"Wx19-20"H Tilt/22"-23" H Recline \$1,410 Part: QB2AL18ES. HCPCS: E2620	Extra Deep Contour, 20"Wx25-26"H Tilt/28"-29" H Recline \$1,410 Part: QB2AL20EXX. HCPCS: E2620

1i. Q-Back 2.0, AFT Series, Active Contour**Not available with Back Cane Option.** Select a back type and the specific overall back height in subsequent section.

AFT, Active, 16"Wx19-20"H Tilt/22"-23" H Recline \$1,145 Part: QB2AFT16AS. HCPCS: E2613	AFT, Active, 18"Wx19-20"H Tilt/22"-23" H Recline \$1,145 Part: QB2AFT18AS. HCPCS: E2613
--	--

1j. Q-Back 2.0, AFT Series, Deep Contour**Not available with Back Cane Option.** Select a back type and the specific overall back height in subsequent section.

AFT, Deep Contour, 16"Wx19-20"H Tilt/22"-23" H Recline..... \$1,325 Part: QB2AFT16DS. HCPCS: E2620	AFT, Deep Contour, 18"Wx21-22"H Tilt/24"-25" H Recline..... \$1,325 Part: QB2AFT18DT. HCPCS: E2620
AFT, Deep Contour, 16"Wx21-22"H Tilt/24"-25" H Recline \$1,325 Part: QB2AFT16DT. HCPCS: E2620	AFT, Deep Contour, 20"Wx21-22"H Tilt/24"-25" H Recline \$1,325 Part: QB2AFT20DT. HCPCS: E2620
AFT, Deep Contour, 18"Wx19-20"H Tilt/22"-23" H Recline \$1,325 Part: QB2AFT18DS. HCPCS: E2620	

1k. Q-Back 2.0, AFT Series, Extra Deep Contour**Not available with Back Cane Option.** Select a back type and the specific overall back height in subsequent section.

AFT, Extra Deep Contour, 16"Wx25-26"H Tilt/28"-29" H Recline..... \$1,410 Part: QB2AFT16EXX. HCPCS: E2620	
--	--

13.13. TRU Balance 3 Back Shroud**REQUIRED.** The back shroud is a feature of the TRU Balance 3 seating system and is available in matching shroud colors. The back shroud is not available when Omit Back Pan/Cushion is selected. Reference colors shown on page 2.

Red Corvette Part: XRF000019	Grey Street Part: XRF000023	Pretty in Pink Part: XRF000027	Electric Blue Part: XRF000030	Green Machine Part: XRF000033
Back in Black Part: XRF000021	Deep Purple Part: XRF000025	Evolution Orange Part: XRF000028	Blue Jean Part: XRF000031	Military Green Part: XRF000034
Night Train Part: XRF000022	Raspberry Beret Part: XRF000026	Midnight Blue Part: XRF000029	Lemon Crush Part: XRF000032	

TRU Balance 3 Seating section continued on following pages.

13. TRU BALANCE 3 SEATING CONTINUED**13.14. Lateral Positioning Wedge****Optional.**

Lateral Positioning Wedge \$15 ea.

Part: SET123465 x2

Available with Sport and TRU-Comfort backs. Not available on power tilt backs 16"-17"H with limited adjustability on 18" H. Not available on power recline back 19"-20"H with limited adjustability on 21"H.

13.15. Back Width Package**REQUIRED** except when Omit Back Pan/Cushion is selected. Available with +/-4" of seat width.

Match Seat Width Selections

Part: PTOINDV3468

14"W-16"W Package

Part: XRF000037

18"W-22"W Package

Part: XRF000039

10"W-14"W Package

Part: XRF000036

16"W-20"W Package

Part: XRF000038

20"W-24"W Package

Part: XRF000040

13.16. Back Width**REQUIRED** except when Omit Back Pan/Cushion is selected. Available with +/-4" of seat width.

12" Back Width

Part: PTOINDV3562

16" Back Width

Part: PTOINDV3564

19" Back Width

Part: PTOINDV3567

22" Back Width

Part: PTOINDV3570

14" Back Width

Part: PTOINDV3563

17" Back Width

Part: PTOINDV3565

20" Back Width

Part: PTOINDV3568

15" Back Width

Part: PTOINDV3856

18" Back Width

Part: PTOINDV3566

21" Back Width

Part: PTOINDV3569

13.17. Overall Back Height**Static Seat and Power Tilt available height range:** 14"-25"H. **Odd heights include 4" gap and even heights include 3" gap** between seat pan and bottom of back cushion. **Power Recline available height range:** 19"-28"H. **Odd heights include 6" gap and even heights include 7" gap.** See the Sport and TRU-Comfort back availability chart for stock order back size ranges.

14" Height

Part: SETHGHT1027

17" Height

Part: SETHGHT1030

20" Height

Part: SETHGHT1033

23" Height

Part: SETHGHT1036

26" Height

Part: SETHGHT1039

29" Height

Part: SETHGHT1042

15" Height

Part: SETHGHT1028

18" Height

Part: SETHGHT1031

21" Height

Part: SETHGHT1034

24" Height

Part: SETHGHT1037

27" Height

Part: SETHGHT1040

16" Height

Part: SETHGHT1029

19" Height

Part: SETHGHT1032

22" Height

Part: SETHGHT1035

25" Height

Part: SETHGHT1038

28" Height

Part: SETHGHT1041

13.18. Power Tilt Privacy Flap**Optional.** A privacy flap is standard with power recline systems. This selection is to add a flap for power tilt only systems.

Privacy Flap \$40

Part: XRF000351

14. TRU BALANCE 3 ARMRESTS**REQUIRED.** Select the desired type of armrest below. Please read notes to ensure seating compatibility. Be sure to also select an armrest size and armpad type.**14.1a. Flip-Back, 8"-16" H Armrests**

Choose a left and right armrest and set height in the following section. Not available with waterfall armpads. Not available with Back Cane Option.

☐ ☐	Flip-Back, 8"-16"H Armrests\$145 ea.	☐ ☐	Flip-Back, 8"-16"H Armrests for Static Seat\$145 ea.
☐ ☐	Left part: ARM1806662. Right part: ARM1806663	☐ ☐	Left part: XRF000320. Right part: XRF000321

Not available with TRU Balance 3 Static Seat.

Armrest Support**Optional.** Available up to 2" difference between seat and back width via endomorph or mesomorph configuration. Not available less than 16" width seat or back.

☐ ☐	Flip-Back Armrest Downpost \$115 ea.
☐ ☐	Left part: XRF000265. Right part: XRF000266

14.1b. Drop-In, Height Adjustable, Removable, T-Arm Armrests

Includes flat armpads.

☐ ☐	Full Length, Drop-In Armrests..... No Charge	☐ ☐	Desk Length, Drop-In Armrests No Charge
☐ ☐	Left part: ARM1803253. Right part: ARM156742	☐ ☐	Left part: ARM1803254. Right part: ARM156743

14.1c. Two-Post, Height Adjustable, Removable Armrests

Select the armrest and set the left and right heights in the following section. These armrests may interfere with offset backs/ back canes.

1a. Two-Post Armrest for Power Seating

☐ ☐	Full Length, Two-Post Armrest.....\$145 ea.	☐ ☐	Desk Length, Two-Post Armrest\$145 ea.
☐ ☐	Left part: XRF000047. Right part: XRF000048	☐ ☐	Left part: XRF000049. Right part: XRF000050

TRU Balance 3 armrests section continued on following pages.

14. TRU BALANCE 3 ARMRESTS CONTINUED**14.1d. Height Adjustable, Removable, Heavy Duty Armrests**

Select a type of HD removable armrests and the set height in the following section. Only one height package is selectable as a normal configuration. Payer coding, coverage, and bundling guidelines may apply.

8"-12"H, Single Post, Height Adjustable, Removable,\$130 ea.
Heavy Duty (Sold in pairs)

Part: XRF000055. HCPCS w/ Static Seat: E0973

10"-14"H, Single Post, Height Adjustable, Removable,\$130 ea.
Heavy Duty (Sold in pairs)

Part: XRF000056. HCPCS w/ Static Seat: E0973

14.2. Armrest Height

REQUIRED. Set left and right armrest height. Note the available height range in the option description.

1. Left Armrest Height

8" H

PTOINDV3535

10" H

PTOINDV3537

12" H

PTOINDV3539

14" H

PTOINDV3541

16" H

PTOINDV3543

9" H

PTOINDV3536

11" H

PTOINDV3538

13" H

PTOINDV3540

15" H

PTOINDV3542

2. Right Armrest Height

8" H

PTOINDV3450

10" H

PTOINDV3452

12" H

PTOINDV3454

14" H

PTOINDV3456

16" H

PTOINDV3458

9" H

PTOINDV3451

11" H

PTOINDV3453

13" H

PTOINDV3455

15" H

PTOINDV3457

14.3. Armrest Size

REQUIRED. Select a left (L) and right (R) armrest size. This also affects armpad selection.



10" Length

L part: XRF000057. R part: XRF000123



14" Length

L part: XRF000058. R part: XRF000124



17" Length

L part: XRF000059. R part: XRF000125
Only available with Flip-Back, 8"-16" H
Armrests and straight arm pads.

15. TRU BALANCE 3 ARPADS

REQUIRED. Please note that waterfall and gel armpads may cause interference with the seat back.

15.1a. Omit Armpad

Select one below if aftermarket armpads or handpads are required.

No Armpad Selection..... No Charge

Part: XRF000060

Not available with Flip-back, 8"-16" H armrests.

Interface for Aftermarket Pad No Charge

Part: XRF000308

Only available with Flip-back, 8"-16" armrests.

15.1b. Foam Armpads

Straight, Armpad Standard

Left part: XRF000061. Right part: XRF000126



Waterfall, Armpad No Charge

Left part: XRF000062. Right part: XRF000127

Not available with Flip-back, 8"-16" armrests.

15.1c. Gel Armpads

Sold only in pairs. Armpads are 3" wide. Manufacturer part number for full length is 314G and desk length is 310G.

a. Falcon Gel Armpads for Flip-Back Armrests

Falcon Gel, Full Length Left and Right \$300

Part: ARM1806807

Falcon Gel, Desk Length Left, Full Length Right \$250

Part: ARM1806810

Falcon Gel, Full Length Left, Desk Length Right \$250

Part: ARM1806813

Falcon Gel, Desk Length Left and Right \$250

Part: ARM1806806

b. Falcon Gel Armpads for Tube Style Armrests

Falcon Gel, Full Length Left and Right \$300

Part: ARM136438

Falcon Gel, Desk Length Left, Full Length Right \$250

Part: ARM136439

Falcon Gel, Full Length Left, Desk Length Right \$250

Part: ARM136442

Falcon Gel, Desk Length Left and Right \$200

Part: ARM136437

16. TRU BALANCE 3 JOYSTICK BRACKET

REQUIRED.

Standard Fixed Receiver Standard

Part: XRF000063

Stealth Height Adjustable Receiver \$194

Part: XRF000064. HCPCS: K0108

17. TRU BALANCE 3 HEADRESTS**Optional.** Select a headrest pad and bracket.**17.1. Headrests, Comfort Plus with Cool Core****1. Headrest Pads**Stealth 8" Comfort Plus, Cool Core (CP08CCFPM) \$279
Part: ST-CP08CCFPM-Q. HCPCS: E0955Stealth 10" Comfort Plus, Cool Core (CP10CCFPM) \$279
Part: ST-CP10CCFPM-Q. HCPCS: E0955Stealth 14" Comfort Plus, Cool Core (CP14CCFPM) \$279
Part: ST-CP14CCFPM-Q. HCPCS: E0955Stealth 6", Comfort Plus, Cool Core (CP06CCFPM) \$279
Part: ST-CP06CCFPM-Q. HCPCS: E0955**2. Headrest Brackets**Comfort Plus Headrest Mounting Bracket, UniLink \$485
Part: ST-UHB2D33200-Q. HCPCS: E1033Multi-Axis Removable Headrest Bracket \$271
Part: ST-HMO475-17M-Q. HCPCS: E1033**17.2. Q Cam Lock****Optional.**Q Cam Lock No Charge
Part: XRF000376

Q Cam Lock Availability (Overall height)	TRU-Comfort/Sport		Q-Back 2	
	Standard	\$50	Standard	\$50
Tilt/Static/Back Canes	14"-21"H	22"-29"H	19"-24"H	25"-26"H
Power Recline Back	19"-21"H	22"-32"H	22"-25"H	26"-29"H

18. LAP BELTS**REQUIRED.** Select a lap belt below.**18.1a. Standard Lap Belts**38" Lap Belt
Part: ACCBELT100650" Lap Belt
Part: ACCBELT100060" Lap Belt
Part: ACCBELT100270" Lap Belt
Part: ACCBELT100380" Lap Belt
Part: ACCBELT100490" Lap Belt
Part: ACCBELT1001**18.1b. Optional Lap Belts**Stealth Retractable Lap Belt (IFRB5350) \$175
Part: ST-IFRB5350Stealth Positioning, 60" Single Pull, Push Button, Medium \$134
(7.5") Padded, 1.5" Belt (IFPB5213)
Part: ST-IFPB5213-QStealth Positioning, 60" Single Pull, Push Button, Large \$134
Padded (9.5"), 2" Belt (IFPB5314)
Part: ST-IFPB5314-QStealth Positioning, 60" Single Pull, Push Button, Extra Large. \$134
Padded (11.5"), 2" Belt (IFPB5315)
Part: ST-IFPB5315-Q60" Sure-Lok Lap Belt (Non-retractable) \$200
Required with ISO 7176-19 Compliant Transit Option.
Part: ACCBELT1016**19. TRU BALANCE 3 LEGREST CONTROLS****Optional.** Specific power legrest options are selected later in this form. Select the control type here, if power legrests are desired, to assist the selection of the power positioning electronics.Independent Legrests (adds 2 actuator control requirement)
Part: XRF000108Articulating Foot Platform (adds 1 actuator control requirement)
Part: XRF000110Combined Legrests (adds 2 actuator control requirement)
Part: XRF000109**20. LEGRESTS****Optional.** Please note the set lower leg length may vary based on your seat type and seat to floor height selections. All efforts will be made to accommodate legrest length selection, but minor changes may need to be made to meet shipping guidelines. Some adjustments may be needed upon receipt of the unit. Reference the section's compatibility matrix for legrest interference.**20.1a. Base Mounted Foot Platform**

May interfere with seat mounted legrests. Not available with center mount foot platforms.

Base Mounted Foot Platform Standard
Part: FRMASMB16077Omit base mounted foot platform No Charge
Part: PTOINDV3674Bell Shaped, Base Mounted Foot Platform \$235
Part: RIGASMB7120041**20.1b. Basic Center Mount Foot Platform****1. Basic Center Mount Foot Platform Assembly**Basic Center Mount \$200
Part: XRF000323

Short	Medium	Long
8"-10"	10"-14"	14"-19"
Basic Center Mount Extension Range		

2. Basic Center Mount Extension and Angle Adjustable FootplatesShort Extension, 5"x9" Footplates
Part: XRF000312. HCPCS: K0040Long Extension, 5"x9" Footplates
Part: XRF000314. HCPCS: K0040Medium Extension, 5"x9" Footplates
Part: XRF000313. HCPCS: K0040**3. Basic Center Mount Calf Pads**Small Calf Pads, Basic Center Mount \$210
Part: RIG2104369Large Calf Pads, Basic Center Mount \$210
Part: RIG2202997
Only available with Long Extension.**4. Basic Center Mount Additional Accessory**Basic Center Mount Footplate Connector \$55
Part: RIG2104240

Legrests section continued on following pages.

20. LEGRESTS CONTINUED**20.1c. Center Mount Foot Platform**

If power articulating foot platform is selected, be sure power legrest control has been selected in later section. Select legrest length in following section.

1. Center Mount Foot Platform TypePower Articulating Foot Platform \$3,145
Part: PTOINDV3754. HCPCS: E1012+K0040Quick Adjustable Center Mount Foot Platform \$265
Part: XRF000311. HCPCS: K0040**2. Extension and Angle Adjustable Footplates**

Pediatric extensions are not compatible with Power Articulating Foot Platform.

Pediatric Extension, 5"x6" Footplates
Part: RIG152577Medium Extension, 5"x11" Footplates
Part: RIG2404018Pediatric Extension, 5"x8" Footplates
Part: RIG2404015Long Extension, 5"x6" Footplates
Part: RIG151505Medium Extension, 5"x6" Footplates
Part: RIG150206Long Extension, 5"x8" Footplates
Part: RIG2404017Medium Extension, 5"x8" Footplates
Part: RIG2404016Long Extension, 5"x11" Footplates
Part: RIG2404019

Medium	Long
11"-17"	14"-20"
Power AFP Extension Range	

Pediatric	Medium	Long
8.5"-12.5"	11"-16.5"	14"-20"
Quick Adj Center Mount Extension Range		

3a. Center Mount Calf Pads

Calf pads are standard and required with a power articulating foot platform.

Small Calf Pads, AFP Only, 6"Wx6"H (pair) Standard
Part: RIGASMB7120018Large Calf Pads, AFP Only, 8" taper to 6"Wx8"H (pair) Standard
Part: RIGASMB7120019Small Calf Pads, Manually Adjustable Center Mount Only, \$210
6"Wx6"H (pair)
Part: XRF000156Large Calf Pads, Manually Adjustable Center Mount Only, \$210
8" taper to 6"Wx8"H (pair)
Part: XRF000157**3b. Center Mount Calf Panel**Small Calf Panel, 18"Wx8"H \$210
Part: FRMASMB7110012Large Calf Panel, 24"Wx8"H \$265
Part: RIG123825
Only available on widths of 20" or greater.**4. Center Mount Modular Footplate Extension Components**

Optional. Assembly required. May require seating configuration adjustments to avoid front caster interference.

2" Bell-Shaped Length Extension \$55
Part: RIG152519
Only available with the 5x6 footplates.Outside Edge \$55
Part: RIG2404921
Available with 5"x8" or 5"x11" footplates.Width Extension for 8" \$55
Part: RIG2404919Front and Rear Cap \$55
Part: RIG2404923Width Extension for 11" \$55
Part: RIG2404920

Available with 5"x8" or 5"x11" footplates.

5. Center Mount Additional AccessoriesSmall Foam Calf Pad Cover (pair) \$150
Part: ST-CSCS0775-QFootplate Riser \$150
Part: XRF000348Large Foam Calf Pad Cover (pair) \$150
Part: ST-CSCS0985-Q

Not available with Long Extension. Select legrest length between 2" to 9" at 1" increments.

Small Foam Panel Cover (pair) \$150
Part: ST-CSCP0920-QFootplate Riser Calf Panel No Charge
Part: FRMASMB7110011Large Foam Panel Cover (pair) \$150
Part: ST-CSCP0927-Q

Only available and required with Footplate Riser selection and extension lengths of 6" to 9".

Footplate Connector \$55
Part: RIGASMB7120023**6. Stealth Gel Footplate Covers**5"x8" Tapered Gel Footplate Covers \$275
Part: ST-FPG5008T-Q. HCPCS: K01085"x8" Tapered Gel Footplate Covers for Width Extension \$275
Part: ST-AFPG5008E-Q. HCPCS: K01085"x11" Tapered Gel Footplate Covers \$275
Part: ST-AFPG5011T-Q. HCPCS: K01085"x11" Tapered Gel Footplate Cover for Width Extension \$275
Part: ST-AFPG5011E-Q. HCPCS: K0108

Legrests section continued on following pages.

20. LEGRESTS CONTINUED**20.1d. Swing-Away Legrests**

Select a legrest upper, a legrest extension, and type of footplate. Pair (P), left (L), and right (R) configurations are selected by checking the appropriate boxes.

1a. 60°-90° Swing-Away Legrest Upper

Adjustable in 5° increments.

60°-90° Swing-Away Legrest \$285
 Pair part: RIG120062. Left: FRMASMB17558.
 Right: FRMASMB17559
 Not available with 11"-11" or 11"-16.5" extensions.

60°-90° Swing-Away Legrest, Short \$285
 Pair part: RIG136481. Left: RIG132607.
 Right: RIG132609.
 For use only with 11"-13" or 11"-16.5" lower extensions.

60°-90° Swing-Away Legrest Compatibility Matrix	Power Recline or Static Seat			Power Tilt			Power Lift		
	Low	Med.	High	Low	Med.	High	Low	Med.	High
Extruded, 11"-13"	!	✓	✓	✓	✓	✓	✓	✓	✓
Extruded, 13"-16.5"	!	✓	✓	✓	✓	✓	✓	✓	✓
Extruded, 15.5"-19.5"	!	!	!	!	!	!	!	!	✓
Angle Adjustable, 13"-16.5"	!	✓	✓	✓	✓	✓	✓	✓	✓
Angle Adjustable, 15.5"-19.5"	!	!	!	!	!	!	!	!	✓
Multi-Axis, 11"-16.5"	!	✓	✓	✓	✓	✓	✓	✓	✓
Multi-Axis, 15.5"-19.5"	!	!	!	!	!	!	!	!	✓
11"-13" extension causes no interference at all seat-to-floor heights.									

☒ This rigging configuration causes no interference issues at this seat-to-floor height.

☐ This rigging configuration may have less than 1.75" clearance when fully extended and may interfere with the ground or reduce obstacle climbing ability.

2. Lower Extension

11"-13" Length
 Pair part: XRF000088.
 Left: XRF000089.
 Right: XRF000090.

13"-16.5" Length
 Pair part: XRF000094.
 Left: XRF000095.
 Right: XRF000096.

15.5"-19.5" Length
 Pair part: XRF000097.
 Left: XRF000098.
 Right: XRF000099.

3. Angle Adjustable Footplate(s)

*The 6"x11" footplates may cause interference with seat widths less than 18". *The 8"x11" footplates may cause interference with seat widths less than 20".

5"Wx8"D Footplates with Heel Loops \$87.50 ea.
 Pair part: FRMASMB4735. Left: RIG137791. Right: RIG137792
 HCPCS: E0951

6"Wx11"D¹ Footplates with Heel Loops \$87.50 ea.
 Pair part: FRMASMB4741. Left: RIG137796. Right: RIG137797
 HCPCS: E0951

6"Wx8"D¹ Footplates with Heel Loops \$87.50 ea.
 Pair part: FRMASMB4738. Left: RIG137793. Right: RIG137795
 HCPCS: E0951

8"Wx11"D² Footplates with Heel Loops \$87.50 ea.
 Pair part: FRMASMB4744. Left: RIG137798. Right: RIG137799
 HCPCS: E0951

20.1e. Power Elevating Legrests

Be sure power legrest control has been selected in previous section. Select an upper extension type and footplate. Pair (P), left (L), and right (R) configurations are selected by checking the appropriate boxes. Select legrest length in following section.

1. Power Elevating Legrest Upper Assembly

Includes calf pads.

Power Elevating Legrests, Pair \$2,865
 Part: XRF000341. HCPCS: E1010+K0040+E0951

Legrest Length Matrix Center Mount	Seat to Floor Option			Legrest Length Matrix Center Mount	Seat to Floor Option		
	Low	Middle	High		Low	Middle	High
Power ELR: Static Seat, Recline Only				Power ELR: Tilt Only, Tilt and Recline, Any Power Lift			
8.5"-14.5" Lower Length	8.5"-14"	8.5"-14.5"	8.5"-14.5"	8.5"-14.5" Lower Length	8.5"-14.5"	8.5"-14.5"	8.5"-14.5"
14.5"-21" Lower Length	12.5"-14"	12.5"-15"	12.5"-16"	14.5"-21" Lower Length	12.5"-15"	12.5"-16"	12.5"-17"

2. Lower Extension

8.5" to 14.5" Lower Length
 Part: XRF000344

12.5" to 21" Lower Length
 Part: XRF000347

3. Angle Adjustable Footplate with Heel Loop and Wheel Bumper

*The 6"W footplates may cause interference with seat widths less than 16". *The 8"W footplates may cause interference with seat widths less than 19".

4"Wx6"D, Pair
 Part: FRMASMB7448

5"Wx8"D, Pair
 Part: FRMASMB7450

6"Wx11"D¹, Pair
 Part: FRMASMB7452

5"Wx6"D, Pair
 Part: FRMASMB7449

6"Wx8"D¹, Pair
 Part: FRMASMB7451

8"Wx11"D², Pair
 Part: FRMASMB7453

Legrests section continued on following pages.

20. LEGRESTS CONTINUED**20.1f. Elevating Legrests**

Select a legrest upper, a legrest extension, and type of footplate. Pair (P), left (L), and right (R) configurations are selected by checking the appropriate boxes. Select legrest length in following section.

1. Elevating Legrests Upper

P L R Elevating Legrest Assembly \$275 ea.
 Pair part: FRMASMB483. Left: RIG138012. Right: RIG138013. HCPCS: E0990
 Select lower extension and footplates

2a. Lower Extension

P L R 11"-14" Length
 Pair: FRMASMB4755. Left: RIG138026. Right: RIG138030
 For reinforced only

P L R 14"-17" Length
 Pair: FRMASMB4756. Left: RIG138033. Right: RIG138035
 For reinforced only

P L R 8.5"-12" HD Drop-Ins/12.5"-16.5" ELR, Ext Length
 Pair part: XRF000101. Left: XRF000102. Right: XRF000103
 For angle adjustable only

P L R 13"-18" HD Drop-Ins/15.5"-19.5" ELR, Ext Length
 Pair part: XRF000104. Left: XRF000105. Right: XRF000106
 For angle adjustable only

Elevating Legrest Compatibility Matrix	Power Recline or Static Seat			Power Tilt			Power Lift		
	Low	Med.	High	Low	Med.	High	Low	Med.	High
Reinforced, 11"-14"	✓	✓	✓	✓	✓	✓	✓	✓	✓
Reinforced, 14"-17"	!	!	✓	!	✓	✓	✓	✓	✓
Angle Adjustable, 12.5"-16.5"	✓	✓	✓	✓	✓	✓	✓	✓	✓
Angle Adjustable, 15.5"-19.5"	!	!	!	!	!	✓	!	✓	✓

☒ This rigging configuration causes no interference issues at this seat-to-floor height.

☐ This rigging configuration may have less than 1.75" clearance when fully extended and may interfere with the ground or reduce obstacle climbing ability.

3a. Reinforced Cast Footplate(s)

P L R 7.25"Wx6"D Reinforced Footplate..... Standard
 Pair part: FRMASMB4748. Left: RIG138055. Right: RIG138056

P L R 8.5"Wx6"D Reinforced Footplate Standard
 Pair part: FRMASMB4747. Left: RIG138060. Right: RIG138062
 May cause interference on seat widths less than 18."

3b. Angle Adjustable Footplate(s)

*The 6"x11" footplates may cause interference with seat widths less than 18". *The 8"x11" footplates may cause interference with seat widths less than 20". 4"x6" and 5"x6" footplates are only available in pairs.

P L R 5"Wx8"D Footplates with Heel Loops \$87.50 ea.
 Pair part: FRMASMB4735. Left: RIG137791. Right: RIG137792
 HCPCS: K0040+E0951

P L R 6"Wx11"D¹ Footplates with Heel Loops \$87.50 ea.
 Pair part: FRMASMB4741. Left: RIG137796. Right: RIG137797
 HCPCS: K0040+E0951

P L R 6"Wx8"D¹ Footplates with Heel Loops \$87.50 ea.
 Pair part: FRMASMB4738. Left: RIG137793. Right: RIG137795
 HCPCS: K0040+E0951

P L R 8"Wx11"D² Footplates with Heel Loops \$87.50 ea.
 Pair part: FRMASMB4744. Left: RIG137798. Right: RIG137799
 HCPCS: K0040+E0951

20.2. Legrest Accessories

Heel Loops (pair)..... \$35
 Part: FRMASMB7873. HCPCS: E0951
 Standard on angle adjustable and multi axis footplates when noted.

Wheel Bumper..... \$15 ea.
 Part: XRF000100

20.3. Set Legrest Extension Length

REQUIRED with center mounts and matching left and right rigging configurations. Part number is PTOINDV plus number following # below.

1a. Centermount Extension Length

2" L (#3989)	5" L (#3992)	8" L (#3680)	9.5" L (#3683)	11" L (#3686)	12.5" L (#3689)	14" L (#3692)	15.5" L (#3695)	17" L (#3698)	18.5" L (#3701)
3" L (#3990)	6" L (#3993)	8.5" L (#3681)	10" L (#3684)	11.5" L (#3687)	13" L (#3690)	14.5" L (#3693)	16" L (#3696)	17.5" L (#3699)	19" L (#3702)
4" L (#3991)	7" L (#3994)	9" L (#3682)	10.5" L (#3685)	12" L (#3688)	13.5" L (#3691)	15" L (#3694)	16.5" L (#3697)	18" L (#3700)	

1b. Front Rigging Extension Length

Non-matching or individual legrest configurations lengths are not set by the manufacturer. See previous section for center mount extension lengths.

7" L (#3470)	8.5" L (#3473)	10" L (#3476)	11.5" L (#3479)	13" L (#3482)	14.5" L (#3485)	16" L (#3488)	17.5" L (#3491)	19" L (#3494)	20.5" L (#3497)	22" L (#3669)
7.5" L (#3471)	9" L (#3474)	10.5" L (#3477)	12" L (#3480)	13.5" L (#3483)	15" L (#3486)	16.5" L (#3489)	18" L (#3492)	19.5" L (#3495)	21" L (#3498)	22.5" L (#3670)
8" L (#3472)	9.5" L (#3475)	11" L (#3478)	12.5" L (#3481)	14" L (#3484)	15.5" L (#3487)	17" L (#3490)	18.5" L (#3493)	20" L (#3496)	21.5" L (#3600)	

21. TRU BALANCE 3 POWER POSITIONING ELECTRONICS**REQUIRED** with power seating or legrest options. Select one option from section 21.1 or 21.3.**21.1a. Tilt Only through Toggle**

Single Actuator Control, Power Tilt Only, Through Toggle No Charge

Part: ELE1905907

Available with QLNES. Select toggle side. May require extended lead time when selected with Fender Lights.

1. Toggle Mount Side Selection

Left Side Toggle

Part: PTOINDV1839

Right Side Toggle

Part: PTOINDV1838

21.1b. Actuator Control through Joystick**1. Actuator Control Selection**

Single Actuator Through Q-Logic 3 Controller and Interface (AM3) No Charge

2 Actuators Through Q-Logic 3 Controller and Interface (AM3) \$1,250

Not funded by ADP

3 Actuators Through Q-Logic 3 Controller and Interface (AM3) \$1,500

Not funded by ADP

Up to 5 Actuators Through Q-Logic 3 Controller and Interface (AAM) \$1,675

ADP: SECX15005

2. Q-Logic 3 Optional Switch Control

Only available with Q-Logic 3 AAM, multiple actuator configurations.

4-Function Push Button \$225

Part: CTL123218

See next section for switch mounting.

Only available with AAM.

4-Way Toggle \$575

Part: ELEASMB5925

Only available with AAM.

4-Way Button Toggle \$575

Part: ELEASMB5926

Only available with AAM.

3. User Controls Side Selection**Required with 4-Way Toggle and 4-Way Button Toggle.**

Mount User Controls on Left

Part: PTOINDV3996

Mount User Controls on Right

Part: PTOINDV3997

21.1c. Actuator Control with iAccess Switch

iAccess Flat Button \$575

Part: CTL144009

See next section for switch mounting.

iAccess Toggle Button \$575

Part: CTL144010

See next section for switch mounting.

21.2. 4-Function Push Button/iAccess Mount**REQUIRED** with 4-Function Push Button or iAccess options. Under arm mount is not available with TRU Balance 3 HD armrests.

Joystick Mount No Charge

Part: MEC1907500

Attendant Mount \$35

Part: MEC1906564

Under Arm Mount, Flip-Back Armrest No Charge

Part: MEC1906531

Under Arm Mount, Flip-Back Armrest, plus iAccess Guard \$35

Part: MEC2001392

Under Arm Mount, Tube Style Armrest No Charge

Part: MEC1906530

Under Arm Mount, Tube Style Armrest, plus iAccess Guard \$35

Part: MEC2001391

1. Under Arm Mount Side Selection

Left Side Under Arm Mount

Part: PTOINDV3678

Right Side Under Arm Mount

Part: PTOINDV3679

**22. TRU BALANCE 3 ATTENDANT CONTROL****Optional.** *Attendant control will be at no charge when ordered with Tilt and Recline. Mounting bracket is included.

Attendant Control \$950

Part: XRF000078. HCPCS: E2331

23. SEAT MOUNTED LED LIGHTS**Optional.** Requires Q-Logic 3 Joystick for Full Light Package or Q-Logic 3e Joystick.


Seat Mounted, LED, Full Light Package \$770

Part: XRF000130

Not available with Base Mounted, LED, Full Light Package.

24. ACCESSORIES**Optional.** Multiple accessories can be mounted along the various tracks but some interference may occur.**Accessory Bar** \$90
Part: XRF000128

Horizontal rail attached to rear of seat. Standard with the Trilogy vent tray and oxygen holder.

Backpack Holder (Pair) (A)..... \$42
Part: ACC147774**Cane/Crutch Holder** \$195
Part: ACC152515. HCPCS: E2207**Cell Phone/ Media Holder, X-Grip (E)**..... \$95
Part: ACC155835
Not available with Heavy Duty armrests.**Clothing Guard (Pair)**..... \$275
Part: ACC125007**Compact Vent Tray** \$1,395
Part: ACC133889. HCPCS: E1029**Cup Holder** \$20
Part: ACC125003
For use with track-style armrests.**Cup Holder (Tube-style Armrest)**..... \$20
Part: ACCASMB2634
For use with tube-style armrests.**Flag Kit** \$125
Part: ACC1700443**Flip-Back Armrest, Side Track Accessory Bracket ..** \$25
Part: MEC1809178**Glove Box (C)**..... \$85
Part: ACC125011**Grip Drink Holder (20085)**..... \$100
Part: ACC152993**Grip Lap Board (20005)**..... \$100
Part: ACC152992**Hydration System (B)**..... \$328
Part: ACC125010**Large, Black, Weather Cover** \$240
Part: ACCCOVR1021**Oxygen Holder**..... \$110
Part: ACC125006. HCPCS: E2208
Includes accessory bar, a credit will be applied PTOIN-DV3769.**Personal Item Hook**..... \$20
Part: ACC125008**Personal Item Pouch** \$175
Part: ACC2101017**Privacy Bag**..... \$80
Part: ACC2005723**Push Handles (pair)**..... \$140
Part: ACC2005027
Not available with oxygen holder or vent tray.**Q-Logic EX/NE Programmer** \$625
Part: ELEASMB6976**Quantum Backpack** \$200
Part: ACC1905853**Rear View Mirror** \$35
Part: ACC155833**Reflector Kit** \$35
Part: ACC124670**Stealth Large Tablet Holder** \$565
Part: ACC141818**Transfer Handle (D)**..... \$152.50 ea.
Part: MEC2205680 x2. HCPCS: E0705**Depth Adjustable Vent Tray** \$1,000
Part: ACC2302500. HCPCS: E1029
Includes accessory bar, a credit will be applied PTOIN-DV3769.**Remote Stop Switch**..... \$1,495
Part: XRF000379**Aftermarket Accessory Mounting Bracket**

Select a quantity up to 5.

Aftermarket Accessory Mounting Bracket, Slotted \$100 ea.
Part: ACC2407046**Aftermarket Accessory Mounting Bracket, Non-slotted** \$100 ea.
Part: ACC2405496

26. POSITIONING COMPONENTS**Optional.** More options are available on the [Positioning Components Order Form](#). Please provide quantity next to check box where applicable.**26.1. Stealth Thoracic Laterals**

Swing-away laterals may interfere with TRU Balance 3 armrests and 2-post, flip up armrests. These laterals are not available with the back cane option.

1. BracketSwing-away, UniLink, Redesign Back, pair (ULS4042000-Q)... \$1032
Part: ST-ULS4042000-Q. HCPCS: E1034

Only available with Redesigned back assemblies.

Swing-away, UniLink, 125° w/Track, pair (ULPK404108B00) ... \$966
Part: ST-ULPK404108B00-Q. HCPCS: E1034

Only available with Sport or TRU Comfort backs.

Swing-away, UniLink, Back Cane, pair (ULPK404104B00) \$966
Part: ST-ULPK404104B00-Q. HCPCS: E1034

Only available with back cane and omit back configurations.

Swing-away, UniLink, Q-Back 2, Left (ULPK404103L00) \$483
Part: ST-ULPK404103L00-Q. HCPCS: E1034

Only available with Active Contour Q-Backs.

Swing-away, UniLink, Q-Back 2, Right (ULPK404103R00) \$483
Part: ST-ULPK404103R00-Q. HCPCS: E1034

Only available with Active Contour Q-Backs.

2. PadsCurved Pads, 3"x5", pair (TWBL3X5)..... \$164
Part: INDPART2663. HCPCS: E0956Curved Pads, 4"x6", pair (TWBL4X6)..... \$
Part: INDPART2664. HCPCS: E0956**3. Installation**Install Laterals
Part: PTOINDV3499Do Not Install
Part: PTOINDV3500**26.2. Stealth UniLink Pelvic/Thigh Guide****1. UniLink Pelvic/Thigh Hardware**

^ Please verify the appropriate billing code depending upon the placement of the support pad.

UniLink, 5" Height Adjustable, \$405 ea.
Removable Mounting Hardware
QTY Part: ST-UTTP5-Q. HCPCS: E1028/E1034^UniLink, 7" Height Adjustable, \$423 ea.
Removable Mounting Hardware
QTY Part: ST-UTTP7-Q. HCPCS: E1028/E1034^UniLink Multi-Angle Thigh Guide, Left \$585 ea.
Part: ST-UTT8A11602-L-Q. HCPCS: E1028/E1034^UniLink Multi-Angle Thigh Guide, Right..... \$585 ea.
Part: ST-UTT8A11602-R-Q. HCPCS: E1028/E1034^UniLink Flip Down Hip/Thigh Support, Left..... \$483 ea.
Part: ST-UTS25100-L-Q. HCPCS: E1028/E1034^UniLink Flip Down Hip/Thigh Support, Right \$483 ea.
Part: ST-UTS25100-R-Q. HCPCS: E1028/E1034^Flip Down with LEFT Full Surface Contact \$459 ea.
Feature (TWBADD-L & LPHW407)
Part: POS140428. HCPCS: E1028/E1034^Flip Down with RIGHT Full Surface Contact..... \$459 ea.
Feature (TWBADD-R & LPHW407)
Part: POS140429. HCPCS: E1028/E1034^**2. UniLink Pelvic/Thigh Pads**

^ Please verify the appropriate billing code to be used with all third-party payors.

UniLink, 4" x 4" Pad \$164 ea.
Part: ST-UTP544056100-Q. HCPCS: E0953/E0956^UniLink, 4" x 8" Pad..... \$164 ea.
Part: ST-UTP548056100-Q. HCPCS: E0953/E0956^UniLink, 4" x 12" Pad \$164 ea.
Part: ST-UTP54D056100-Q. HCPCS: E0953/E0956^**3. UniLink Pelvic/Thigh Additional Links**UniLink, 1.5" Link \$121 ea.
Part: ST-UIL001150A000-QUniLink, 2" Link..... \$121 ea.
Part: ST-UIL001200A000-Q**4. UniLink Pelvic/Thigh Guide Installation**Install on Left
Part: PTOINDV3585Install on Right
Part: PTOINDV3586Do Not Install
Part: PTOINDV3719**26.3. Stealth Positioning Chest Harness**Upper Classic, Medium (IFPS1103) \$221
Part: ST-IFPS1103-Q. HCPCS: E0960Upper Classic, Large (IFPS1104) \$221
Part: ST-IFPS1104-Q. HCPCS: E0960Upper Contour, Medium (IFPS2103)..... \$221
Part: ST-IFPS2103-Q. HCPCS: E0960Upper Contour, Large (IFPS2104) \$221
Part: ST-IFPS2104-Q. HCPCS: E0960**26.4. Stealth Positioning Chest Straps**Chest Strap, Medium w/ Belt Anchor, \$140
(IFCS4113)
Part: ST-IFCS4113-Q. HCPCS: E0960Chest Strap, Large w/ Belt Anchor, \$140
(IFCS4114)
Part: ST-IFCS4114-Q. HCPCS: E0960



Positioning components section continued on following page.

26. POSITIONING COMPONENTS CONTINUED**26.5. Stealth Swing-Away Headrest Mounts**

These swing-away headrest mounts are only available with the Stealth Comfort Series headrests.

1. Facial Lateral HardwareUniLink Facial Lateral, Swing Away, Right.....\$496
Part: ST-SAR2A-R-Q. HCPCS: E1028UniLink Facial Lateral, Swing Away, Left\$496
Part: ST-SAR2A-L-Q. HCPCS: E1028UniLink Facial Lateral, Dual Swing Away \$765
2 Upper & 1 Lower Right
Part: ST-SAR4A-DBU-RL-Q. HCPCS: E1028UniLink Facial Lateral, Dual Swing Away \$765
2 Upper & 1 Lower Left

Part: ST-SAR4A-DBU-LL-Q. HCPCS: E1028

UniLink Facial Lateral, Dual Swing Away \$674
2 Upper
Part: ST-SAR4A-DBU-Q. HCPCS: E1028**2. Facial Lateral Pads and Accessories**Egg Switch Mount (ESM) (A) \$70 ea.
Part: INDPART2815Figure 8 Bracket with 13" Rod.....\$366 ea.
and Egg Switch Mount (B)
Part: ST-F8675-13-QFacial Pad, 4"x2", Cool Core (900M-Q) \$96 ea.
Part: ST-900M-Q. HCPCS: K0108Facial Pad, 6"x2", Cool Core (920M-Q)..... \$96 ea.
Part: ST-920M-Q. HCPCS: K0108Swing-Away Joystick Mount, 18" Rod (C).....\$100 ea.
Part: ST-SUS18-Q
Allows mounting of facial pads and attaches to the headrest**26.6. Stealth Elbow Support****1. Elbow Support Bracket**Left, Swing-Away, Elbow Support\$279
Bracket Tube Style Arms (ARMS-ES-L)
Part: POS124064. HCPCS: E1028Left, UniLink, Fixed, Elbow Support Bracket..... \$309
for Flip-Back, 8"-16"H Armrests (UETB4021FL-Q)
Part: ST-UETB4021FL-QRight, Swing-Away, Elbow Support.....\$279
Bracket Tube Style Arms (ARMS-ES-R)
Part: POS124066. HCPCS: E1028Right, UniLink, Fixed, Elbow Support Bracket..... \$309
for Flip-Back, 8"-16"H Armrests (UETB4021FR-Q)
Part: ST-UETB4021FR-Q**2. Elbow Support Pad**3"Wx3"H Elbow Pad (ARMS-ES-3X3) \$85 ea.
Part: POS124065. HCPCS: K01084"Wx4"H Elbow Pad (ARMS-ES-4X4)..... \$85 ea.
Part: POS124067. HCPCS: K0108**26.7. Stealth Foot Boxes**Stealth Small Foot Boxes, 5"Wx8"Dx7"H..... \$725
(SPLE107) (pair)
Part: INDPART2786. HCPCS: E0954Stealth Medium Foot Boxes, 6"Wx10"Dx9"H \$725
(SPLE108) (pair)
Part: INDPART2787. HCPCS: E0954Stealth Large Foot Boxes, 7"Wx12"Dx11"H..... \$725
(SPLE109) (pair)
Part: INDPART2788. HCPCS: E0954**26.8. Stealth Positioning Shoe Holder**Small Shoe Holder, 3/4" Webbing (IFSH4211)..... \$93 ea.
Part: ST-IFSH4211-Q. HCPCS: E0951+E0952Medium Shoe Holder, 3/4" Webbing (IFSH4212) \$93 ea.
Part: ST-IFSH4212-Q. HCPCS: E0951+E0952Large Shoe Holder, 3/4" Webbing (IFSH4213) \$93 ea.
Part: ST-IFSH4213-Q. HCPCS: E0951+E0952

27. LAP TRAY**27.1. Quantum Lap Tray**

Mounting bracket is required with lap tray selection. These trays are compatible with Flip-Back armrests.

1a. Black Lap Tray

Black Lap Tray, 11" W \$360
 Part: MEC1905328-11B. HCPCS: E0950
 For use with 12"-13" back width.

Black Lap Tray, 13" W \$360
 Part: MEC1905328-13B. HCPCS: E0950
 For use with 14"-15" back width.

Black Lap Tray, 15" W \$360
 Part: MEC1905328-15B. HCPCS: E0950
 For use with 16"-17" back width.

Black Lap Tray, 17" W \$360
 Part: MEC1905328-17B. HCPCS: E0950
 For use with 18"-19" back width.

Black Lap Tray, 19" W \$360
 Part: MEC1905328-19B. HCPCS: E0950
 For use with 20"-22" back width.

1b. Clear Lap Tray

Clear Lap Tray, 11" W \$360
 Part: MEC1905328-11C. HCPCS: E0950
 For use with 12"-13" back width.

Clear Lap Tray, 13" W \$360
 Part: MEC1905328-13C. HCPCS: E0950
 For use with 14"-15" back width.

Clear Lap Tray, 15" W \$360
 Part: MEC1905328-15C. HCPCS: E0950
 For use with 16"-17" back width.

Clear Lap Tray, 17" W \$360
 Part: MEC1905328-17C. HCPCS: E0950
 For use with 18"-19" back width.

Clear Lap Tray, 19" W \$360
 Part: MEC1905328-19C. HCPCS: E0950
 For use with 20"-22" back width.

2. Mounting Hardware

Mounting Bracket, Chest Bar/Lap Tray No Charge
 Part: MEC2104886

27.2. Simple Solutions Lap Tray

Mounting bracket is required with lap tray selection. See chart below for compatibility.

1. Lap Tray

Simple Solutions Small Tray with Cut Out \$750
 Part: ST-SSLXS010001171-Q. HCPCS: E0950

Simple Solutions Medium Tray with Cut Out \$750
 Part: ST-SSLXM010001171-Q. HCPCS: E0950

Simple Solutions Large Tray with Cut Out \$750
 Part: ST-SSLXL010001171-Q. HCPCS: E0950

2. Mounting Hardware

Swing Away Mounting Hardware No Charge
 Part: ST-SSLXHD10000060-Q

Simple Solutions Lap Tray Armrest Compatibility	Width											
	10	11	12	13	14	16	17	18	19	20	21	22
Flip-Back	S	S	S	S	S	S	SM	ML	ML	ML	L	L
Flip-Back, Cane Mounted	S	S	S	S	S	SM	SM	ML	ML	L	L	
Drop-In, T-Arm	S	S	S	S	S	SM	SM	ML	ML	L	L	
Two-Post, Flip-Back	S	S	S	S	S	SM	SM	ML	ML	L	L	
Cantilever	S	S	S	S	S	M	M	ML	ML	L	L	
Q-Captain, Flip-Up						SM		ML		L		
Single Post/ Heavy Duty/ Bariatric	S	S	S	S	S	SM	SM	ML	ML	ML	L	L

S:ST-SSLXS010001171-Q M:ST-SSLXM010001171-Q L:ST-SSLXL010001171-Q

:Compatible

:Armrest configuration may not hold tray position

*Some armrests listed may not be available on this order form.

28. ALTERNATIVE CONTROLS AND ELECTRONICS**Optional.** Select a type of specialty controls below. More options are available on the [Alternative Controls and Electronics order form](#).**28.1. Enhanced Display****Q-Logic 3 Enhanced Display Kit**..... \$2,045

Part: CTL167726

Select a color option below. Includes enhanced display, flexible gooseneck and mounting bracket. Features: 3.5" LCD display; up to 5 profiles to control drive, seat, and auxiliary functions; Bluetooth connectivity for mobile devices and mouse emulation; integrated single switch scanner; IR controls.

Q-Logic 3 Enhanced Display Shroud

Required with Q-Logic 3 Enhanced Display selection above. Reference colors shown on page 2.

Black, Soft Touch Standard

Part: DGN137728

Deep Purple

Part: DGN137729-PPR5

Electric Blue

Part: DGN137729-PBU9

Military Digital Camo \$50

Part: DGN137729-WCM7

Red Corvette

Part: DGN137729-PRD1

Raspberry Beret

Part: DGN137729-PPK4

Blue Jean

Part: DGN137729-PBU11

Faded Glory \$50

Part: DGN137729-WFL2

Back in Black

Part: DGN137729-PBK1

Pretty in Pink

Part: DGN137729-PPK3

Lemon Crush

Part: DGN137729-PYL4

Carbon Fiber \$50

Part: DGN137729-WCB1

Night Train

Part: DGN137729-PBK3

Evolution Orange

Part: DGN137729-POR7

Green Machine

Part: DGN137729-PGN11

Grey Street

Part: DGN137729-PSL3

Midnight Blue

Part: DGN137729-PBU10

Military Green

Part: DGN137729-PGN8

28.2. SCIM Module**Optional.**☐ **SCIM Module for Q-Logic 3** \$850

Part: CTL2001963

Allows use of joystick display for specialty devices. Available with Q-Logic 3 EX Joystick selection above. Up to 5 profiles to control drive, seat, and auxiliary functions, excluding sip-n-puff and single switch functions. A switch is required to operate the input device attached to the SCIM module.

28.3. Q-Logic 3 SAJ**Q-Logic 3 EX SAJ** \$850

Part: CTL143980

Must be used in conjunction with Enhanced Display or Q-Logic 3 Joystick.

28.4. Head Array**i-Drive 4.0 Bluetooth Tri-array Series**..... \$7,406**3 Switch with Egg Switch**

Part: ST-IDHBT200-U-Q. HCPCS: E2330

Adult, Curved, i-Drive Tri-Array with Mounting Hardware \$7,406

Part: ST-IDHBT400-U-Q. HCPCS: E2330

i-Drive 4.0 Bluetooth Pro Series..... \$7,486**3 Switch with Egg Switch**

Part: ST-IDHBT300-U-Q. HCPCS: E2330

Junior, Curved, i-Drive Tri-Array with Mounting Hardware \$7,406

Part: ST-IDHBT450-U-Q. HCPCS: E2330

28.5. Sip and Puff Drive**1. Sip and Puff Module****Sip and Puff Module for Q-Logic 3** \$1,695

Part: ST-SP625M-U-Q. HCPCS: E2325

Includes module and mounting hardware.

2. Sip and Puff Straw and Mount**The straw and mount do not include electrical components to operate the Sip&Puff with a head array. This component requires the selection of a Q-Logic Sip&Puff module or additional module from Stealth Products for operation.****Stealth Sip and Puff for Comfort Series, Left (SP625M-CL)** \$496

Part: ST-SP625M-CL-Q. HCPCS: E2326

Includes straw and mounting hardware. Use with Stealth Comfort Plus headrests.

Stealth Sip and Puff for Comfort Series, Right (SP625M-CR) \$496

Part: ST-SP625M-CR-Q. HCPCS: E2326

Includes straw and mounting hardware. Use with Stealth Comfort Plus headrests.

Part: ST-SP625-U-Q. HCPCS: E2326

Includes straw and mounting hardware. Use with Stealth i-Drive Head Arrays.

3. Sip & Puff Accessories**Sip & Puff Straws, 12 pack (SP625-S12)**..... \$47

Part: ST-SP625-S12-Q



Enhanced Display



SCIM Module



EX SAJ



Adult Pro-Series



Sip & Puff Module



Stealth Sip & Puff

Alternative controls and electronics section continued on following page.

28. ALTERNATIVE CONTROLS AND ELECTRONICS CONTINUED**28.6. Egg Switches**

Red Egg Switch (ESRED)..... \$116
Part: SWTMCRO1077

Blue Egg Switch (ESBLU)..... \$116
Part: SWTMCRO1079

Black Egg Switch (ESBLK)..... \$116
Part: SWTMCRO1081

Green Egg Switch (ESGRN) \$116
Part: SWTMCRO1078

Yellow Egg Switch (ESYEL) \$116
Part: SWTMCRO1080

**28.7. Switch Adapter**

2 to 1 Switch Adapter..... \$135
Part: DWR1361H037
Connects 2 single switches to a stereo jack to act as a dual switch.

28.8. Environmental Controls

Q-Logic 3 ECU Module Kit \$1,780
Part: CTL168057

28.9. USB Mobile Device Charger

USB Mobile Device Charger \$185
Part: ELE168069

Can charge 1 device. A multiplier harness is required and must be added when a Static Seat or Power Tilt Only seating system is chosen. The selection of more than 3 modules requires an additional multiplier.

