



Pride Mobility Products Co.

5096 South Service Rd, Beamsville, ON, L3J 1V4

Phone: 888-570-1113 | Fax: 866-514-1303 | Email: productorders@pridemobility.com

Q6 Edge® 3 Order Form with TRU Balance® 4 Positioning

Account Number: _____ Date: _____

Provider Name: _____

Contact: _____

Phone: _____ Fax: _____

Email: _____


PO Number: _____

Marked for: _____

Ship to Address: _____

City: _____ State: _____ Zip: _____

EDGE 3



INTRODUCTION

This form is interactive when viewed with [Adobe Acrobat Reader](#) and may not function correctly if opened with applications other than Acrobat. Complete the form by placing checks in the desired boxes and provide information in the interactive fields. Buttons shown at the bottom of the form may be utilized to print or submit the order form through a desktop email application. To email via a web-based application, please 'Save As' and attach the PDF to your email.

Thank you for choosing Quantum!


PATIENT INFORMATION

Weight: _____ Height: _____

A. Shoulder Width: _____ D. Max Sitting Height: _____ G. Thigh Depth: _____

B. Chest Width: _____ E. Shoulder Height: _____ H. Lower Leg Length: _____

C. Hip Width: _____ F. Axilla Height: _____ I. Elbow Height: _____



1. BASE MODEL

REQUIRED. Select a group 3 power base. 300 lb. weight capacity with standard 6mph programmable motors and 8 amp off board charger. Dimensions of base without legrests: 24"Wx35.5"L

Quantum Q6 Edge 3\$8,450

Part: Q6 EDGE 3. HCPCS: K0848

2. DRIVE MOTORS**REQUIRED.**

6 mph Drive Motors Standard
Part: MOT147433

6 mph, Accu-Trac Motor Technology \$995
Part: MOT147881. HCPCS: K0108
Requires Q-Logic 3 electronics. Delivers enhanced traction and tracking, especially with switched drive controls.

3. BASE COLOR

REQUIRED. Select one shroud color. Colors shown here may differ from actual product depending on monitor/printer calibration used.

3.1a. Standard Shrouds

 Red Corvette Part: DGN1801101-PRD1	 Deep Purple Part: DGN1801101-PPR5	 Midnight Blue Part: DGN1801101-PBU10	 Green Machine Part: DGN1801101-PGN11
 Back in Black Part: DGN1801101-PBK1	 Raspberry Beret Part: DGN1801101-PPK4	 Electric Blue Part: DGN1801101-PBU9	 Military Green Part: DGN1801101-PGN8
 Night Train Part: DGN1801101-PBK3	 Pretty in Pink Part: DGN1801101-PPK3	 Blue Jean Part: DGN1801101-PBU11	
 Grey Street Part: DGN1801101-PSL3	 Evolution Orange Part: DGN1801101-POR7	 Lemon Crush Part: DGN1801101-PYL4	

3.1b. Hydrographic Shrouds

☐  Military Digital Camo \$325 Part: DGN1801101-WCM7	☐  Faded Glory \$325 Part: DGN1801101-WFL2	☐  Carbon Fiber \$325 Part: DGN1801101-WCB1
--	--	---

4. BASE MOUNTED LIGHTS


Optional. Select one option. Seat Mounted, LED, Full Light Package option can be found later on this form.

Switch Fender Lights..... No Charge
Part: DGN154117
Fender lights operate by pressing included single switch button. Fender lighting is not available with Single Actuator Control through Toggle Box.

Base Mounted, LED, Full Light Package \$770
Part: XRF000129
Requires Q-Logic 3 EX Joystick for Full Lights or Q-Logic 3e Joystick.



Base Mounted Lights



Fender Lights

5. QUANTUM BACKUP CAMERA

Optional. The backup camera may interfere with armrest mounted accessories such as media device holder.

Quantum Backup Camera \$475
Part: MEC1910365

1. Optional Mounting Bracket

Backup Camera Headrest Bracket..... \$35
Part: MEC1910373
Works with .5" x .5" vertical supports only.

Backup Camera Track Mount Bracket \$35
Part: MEC2104227

2. Optional Display Mounting Bracket

Backup Camera Q-Logic Joystick Mount Bracket \$145
Part: MEC2106123



Backup Camera

6. DRIVE WHEELS

REQUIRED. Select a type of drive wheels and accent color. Knobby tires increase overall base width by 1.5".

6.1. Drive Wheels

3"x14" Flat Free Black Tires Standard Part: PTOINDV3525 (pair)	3"x14" Pneumatic Black Tires No Charge Part: WHL1713925 x2
4"x14" Knobby, Flat Free Black Tires \$120 Part: XRF000368	4"x14" Knobby, Pneumatic Black Tires \$120 Part: XRF000367

6.2. Drive Wheel Accent Color

Optional. Not available with knobby tires.

Back in Black Standard Part: DGN164245-PBK1 x2	Pretty in Pink No Charge Part: DGN164245-PPK3 x2	Green Machine No Charge Part: DGN164245-PGN11 x2
Red Corvette No Charge Part: DGN164245-PRD1 x2	Evolution Orange No Charge Part: DGN164245-POR7 x2	Military Green No Charge Part: DGN164245-PGN8 x2
Night Train No Charge Part: DGN164245-PBK3 x2	Midnight Blue No Charge Part: DGN164245-PBU10 x2	Military Digital Camo \$50 ea. Part: DGN164245-WCM7 x2
Grey Street No Charge Part: DGN164245-PSL3 x2	Electric Blue No Charge Part: DGN164245-PBU9 x2	Faded Glory \$50 ea. Part: DGN164245-WFL2 x2
Deep Purple No Charge Part: DGN164245-PPR5 x2	Blue Jean No Charge Part: DGN164245-PBU11 x2	Carbon Fiber \$50 ea. Part: DGN164245-WCB1 x2
Raspberry Beret No Charge Part: DGN164245-PPK4 x2	Lemon Crush No Charge Part: DGN164245-PYL4 x2	

7. BATTERIES**Optional.** The powerchair requires 2 batteries and will be installed if ordered.NF-22, Gel Cell, Interceptor 2 Battery \$300 ea.
Part: BAT2004376. HCPCS: E2361NF-22, Gel Cell, MK Battery \$400 ea.
Part: BATGEL1008. HCPCS: E2361**8. PATIENT'S WEIGHT****REQUIRED.** The patient's weight is required to quote and order the power chair.Up to 200 lb User Weight
Part: PTOINDV3837201 lb-250 lb User Weight
Part: PTOINDV3861251 lb-300 lb User Weight
Part: PTOINDV3984**9. ELECTRONICS****REQUIRED.** Select a type of drive control from the following sections. The No Joystick option is used to omit the joystick and utilize another control, such as a specialty control at the end of this form or found on the [Alternative Controls and Electronics Order Form](#).**9.1. Expandable Controls - Quantum Intelligent Systems™**

Select expandable controls or omit selection if another controller type is desired.


2. Q-Logic 3 Joystick

Q-Logic 3 EX Joystick includes XLR USB charger.


Q-Logic 3 EX Joystick \$900
Part: CTL167724Q-Logic 3 EX Joystick for Full Lights \$900
Part: CTL167725Q-Logic 3e Joystick \$600
Part: CTL1910600
Can operate full light package. Not available with alternative or specialty controls.No Joystick, Setup for Q-Logic
Part: XRF000010
Requires selection of Enhanced Display.**9.2. Q-Logic 3 EX Joystick Armor****REQUIRED with Q-Logic 3 EX Joystick selections above. Not available with Q-Logic 3e Joystick.** Reference colors shown on page 2.Black, Soft Touch Standard
Part: DGN124307

Also matches Night Train shrouds.

Red Corvette No Charge
Part: DGN124336-PRD1Back in Black No Charge
Part: DGN124336-PBK1Grey Street No Charge
Part: DGN124336-PSL3Deep Purple No Charge
Part: DGN124336-PPR5Raspberry Beret No Charge
Part: DGN124336-PPK4Pretty in Pink No Charge
Part: DGN124336-PPK3Evolution Orange No Charge
Part: DGN124336-POR7Midnight Blue No Charge
Part: DGN124336-PBU10Electric Blue No Charge
Part: DGN124336-PBU9Blue Jean No Charge
Part: DGN124336-PBU11Lemon Crush No Charge
Part: DGN124336-PYL4Green Machine No Charge
Part: DGN124336-PGN11Military Green No Charge
Part: DGN124336-PGN8Military Digital Camo No Charge
Part: DGN124336-WCM7Faded Glory No Charge
Part: DGN124336-WFL2Carbon Fiber No Charge
Part: DGN124336-WCB1**10. JOYSTICK HANDLES****Optional.** The optional handle will be mounted to the joystick and the standard handle will be shipped with the unit.Mushroom Handle (B) \$120
Part: ACCASMB1860.Chin Cup (C) \$85
Part: ACCASMB186190T-Bar (D) \$160
Part: ACCASMB1858.Medium Ball Handle (E) \$70
Part: ACCASMB1859.Soft Ball Handle (F) \$120
Part: ACCASMB1856.Extended I-Handle (G) \$140
Part: ACCASMB1857.Bodypoint 3" Goal Post Handle \$180
Part: PLSKNCP1057.Bodypoint 4" Goal Post Handle \$180
Part: PLSKNCP1058.

11. JOYSTICK MOUNTING BRACKETS**Optional.** Select a mount to be used with a joystick selected in previous section.**11.1a. Fixed, In-line Joystick Bracket**Fixed, In-line, Bracket for Q-Logic 3 Standard
Part: FRMASMB10617**11.1b. Swing-Away (Retractable) Joystick Bracket**

Allows the joystick to retract to the side of the armrest while maintaining the same joystick orientation.

Swing-away, Q-Logic 3 \$265
Part: ACC1908434. HCPCS: E1032**11.1c. Basic Swing Away Joystick Mount**Basic Swing Away, Q-Logic, Left Side \$265
Part: ACC2201041. HCPCS: E1032Basic Swing Away, Q-Logic, Right Side \$265
Part: ACC2201042. HCPCS: E1032**11.1d. Offset Swing Away Joystick Mount**Offset Swing Away, Q-Logic, Left Side \$265
Part: ACC142742. HCPCS: E1028Offset Swing Away, Q-Logic, Right Side \$265
Part: ACC142743. HCPCS: E1028**11.1e. Specialty Joystick Mounts**Retract4 Mount \$450
Part: XRF000003. HCPCS: E1032
Full 180° of movement that translates into multi-position use, straight, right offset, left offset, or fold away. Joystick orientation varies as bracket is adjusted.Multi-Axis Mount \$100
Part: XRF000004
Allows height adjustment of joystick.**11.2. Joystick Mount Location****Required except with basic swing away joystick mount.**Left Side Mount
Part: PTOINDV1225Right Side Mount
Part: PTOINDV1226**12. TRU BALANCE® 4 SEATING****REQUIRED.** Maximum seat size is 22"W x 22"D.**12.1. TRU Balance 4 Configuration****REQUIRED.**TRU-Balance 4 Power Tilt Only \$8,495
Part: PTOINDV4014. HCPCS: E1002
50° Center of Gravity (CG) power tilt seating system.TRU-Balance 4 Power Tilt and Recline \$14,995
Part: PTOINDV3954. HCPCS: E1007
50° of posterior tilt and 165° of recline.**12.2. 12" Power Lift****Requires combined selection of TRU Balance 4 Power Tilt and Recline and Articulating Foot Platform. Requires selection in 12.3 Anterior Tilt.**^Please verify the appropriate billing code to be used for this item with all third-party payors.Power Lift with iLevel \$2,995
Part: PTOINDV3956. HCPCS: E2298/E2300^
12" scissor lift with positional feedback.**12.3. Anterior Tilt****Chest support is required with anterior tilt 20 and 30 degree functions. See chest support section. Omitting a chest support selection requires the provider to facilitate aftermarket chest supports that best fits consumer needs.****1. Anterior Tilt Angle****Only available with Power Lift.**10° Anterior Tilt Angle No Charge
Part: MEC200935630° Anterior Tilt Angle \$3,195
Part: MEC2009400. HCPCS: K0108
Requires chest support and knee block.20° Anterior Tilt Angle \$2,695
Part: MEC2009399. HCPCS: K0108
Requires chest support and knee block.0° Anterior Tilt No Charge
Part: PTOINDV3980**12.4. Knee Block Assembly****REQUIRED with 20° or 30° anterior tilt. Optional with 10° and 0° anterior tilt.**Anterior Knee Block Assembly (ST-KBATB4ATC-1-Q) \$700
Part: ST-KBATB4ATC-1-Q. HCPCS: K0108Anterior Knee Block Assembly, Curved (ST-KBATB4ATC-U-Q) .. \$700
Part: ST-KBATB4ATC-U-Q. HCPCS: K0108**12.5. Occupied Transit Kit****Optional.** The power base is standard with Unoccupied Transit Loops. 60" Sure-Lok Lap Belt is required with Occupied Transit Option.WC-19 Occupied Transit Option \$400
Part: KIT152047To determine if your chair is WC-19 compliant (occupied), please refer to this website: <https://wc-transportation-safety.umtri.umich.edu/>

TRU Balance 4 Seating section continued on following pages.

12. TRU BALANCE 4 SEATING CONTINUED**12.6. Seat Width****REQUIRED with seat selection. 12"-15" widths not available with Power Lift.** Maximum seat width is 22". Part number is SETWDTH plus number following # below.

12" W (#1002)	13" W (#1003)	14" W (#1004)	15" W (#1005)	16" W (#1006)	17" W (#1007)	18" W (#1008)	19" W (#1009)	20" W (#1010)	21" W (#1011)	22" W (#1012)
------------------	------------------	------------------	------------------	------------------	------------------	------------------	------------------	------------------	------------------	------------------

12.7. Seat Width Package**REQUIRED with seat selection.**

12"-16" Width Package

Part: XRF000006

Available with Power Tilt Only.

16"-20" Width Package

Part: XRF000008

20"-24" Width Package

Part: XRF000158

Maximum width is 22".

12.8. Seat Depth**REQUIRED with seat selection. 12"-13" depth not available with Power Lift.** Maximum seat depth is 22". Part number is SETDPTH plus number following # below.

12" D (#1002)	13" D (#1003)	14" D (#1004)	15" D (#1005)	16" D (#1006)	17" D (#1007)	18" D (#1008)	19" D (#1009)	20" D (#1010)	21" D (#1011)	22" D (#1012)
------------------	------------------	------------------	------------------	------------------	------------------	------------------	------------------	------------------	------------------	------------------

12.9. Seat Depth Package**REQUIRED with seat selection.****1a. Power Tilt Depth Package**

12"-16" Depth Package

Part: XRF000011

16"-20" Depth Package

Part: XRF000013

20"-24" Depth Package

Part: XRF000159

Maximum depth is 22".

1b. Power Lift Depth Package

16"-22" Depth Package

Part: XRF000309

12.10. Seat to Floor Height**REQUIRED.**

Lowest Seat to Floor

Part: PTOINDV3558

Low Seat to Floor

Part: PTOINDV3976

Middle Seat to Floor

Part: PTOINDV3559


Highest Seat to Floor

Part: PTOINDV3560

Lowest	Low	Med	High
Q6 Edge 3			
17.50"	18.25"	18.75"	19.25"
TRU Balance 4 Seat to Floor Heights			

12.11. Stealth Seat Cushion

Optional. Specific cushion/component sizes available, part numbers, weight capacities, MSRP, HCPCS code, and special order cushions will be found on the respective complete [Stealth Cushion Order Forms](#). Part numbers in this section are for reference only. This section allows a quick selection of a single cushion that matches the seat pan width and depth to alleviate listing all part numbers and codes on this power base form. The complete cushion part number and MSRP is a product of the cushion model, top foam, liner, cover, width, and depth as shown below. Reference the subsequent chart for available top foam, liner, and cover options for each cushion type. See specific cushion order form for all sizes, top foam, liners, covers, parts numbers, and MSRP.

**1. Cushion Type**

TRU-Comfort 2 SPP
Reference part: XRF000116.
Part series: TC2-###XXYY.

Essence SPP
Reference part: XRF000118.
Part series: ESS-###XXYY.

Essence Plus SPP
Part series: ESSP-###XXYY.

Solution SPP
Reference part: XRF000113.
Part series: SOL-###XXYY.

Spectrum Foam SPP
Reference part: XRF000112.
Part series: SPCF-###XXYY.

Spectrum Gel SPP
Reference part: XRF000117.
Part series: SPCG-###XXYY.

Premiere P
Reference part: XRF000122.
Part series: PREM-###XXYY.

Glacial SP
Reference part: XRF000114.
Part series: GLA-###XXYY.

Zen SP
Reference part: XRF000115.
Part series: ZEN-###XXYY.

Simplicity G
Reference part: XRF000111.
Part series: SIM-###XXYY.

TRU Balance 4 Seating section continued on following pages.

12. TRU BALANCE 4 SEATING CONTINUED**2. Essence Cushion Wedge Kit**

Only available with Essence cushion selection above.

Small Wedge Positioning Kit\$45
 Part: POP1909006-S
 14"W and smaller

Small Wedge Positioning Kit\$45
 Part: POP1712340-S
 15"W and larger

Medium Wedge Positioning Kit.....\$85
 Part: POP1909006-M
 14"W and smaller

Medium Wedge Positioning Kit.....\$85
 Part: POP1712340-M
 15"W and larger

Large Wedge Positioning Kit.....\$90
 Part: POP1909006-L
 14"W and smaller

Large Wedge Positioning Kit.....\$90
 Part: POP1712340-L
 15"W and larger

3. Cushion Top Foam

Visco Top Foam
 Reference part: XRF000380

FRG Top Foam
 Reference part: XRF000381

No Top Foam
 Reference part: XRF000382

4. Cushion Liner

Stealth Urethane Liner
 Reference part: XRF000383

Stealth Polyurethane Liner
 Reference part: XRF000384

No Liner
 Reference part: XRF000385

5. Cushion Cover

CoolCore Cover
 Reference part: XRF000386

Mesh Cover
 Reference part: XRF000387

Stealth Tek Cover
 Reference part: XRF000388

	None	Visco/V.Elastic	FRG	None	Urethane	Polyurethane	CoolCore	Mesh	Stealth Tek
	Top Foam			Liner			Cover		
TRU-Comfort 2	—	✓	—	✓	—	✓	✓	—	—
Essence	—	✓	—	✓	—	✓	✓	—	—
Solution	—	✓	—	—	✓	✓	✓	✓	✓
Spectrum Foam	—	✓	✓	—	✓	✓	✓	✓	✓
Spectrum Gel	—	✓	✓	—	✓	✓	✓	✓	✓
Premiere	—	✓	✓	✓	✓	✓	✓	—	—
Glacial	✓	—	—	✓	✓	✓	✓	✓	✓
Zen	—	✓	✓	✓	✓	✓	✓	✓	✓
Simplicity	✓	—	—	✓	—	—	—	—	✓

6. Optional Base Insert

Option selection available with Glacial, Premiere, Simplicity, Solution, Spectrum Foam, Spectrum Gel, and Zen cushions.

ABS Flat Rigidizer
 Reference part: XRF000389

7. Additional Cushion Liner

Up to 5 may be selected below. See respective cushion order form for MSRP and part numbers.

Additional Stealth Urethane Liner
 QTY Reference Part: XRF000390

Additional Stealth Polyurethane Liner
 QTY Reference Part: XRF000391

8. Additional Cushion Cover

Up to 5 may be selected below. See respective cushion order form for MSRP and part numbers.

Additional CoolCore Cover
 QTY Reference Part: XRF000392

Additional Stealth Tek Cover
 QTY Reference Part: XRF000394

Additional Mesh Cover
 QTY Reference Part: XRF000393

12.12. Back Type Selection**REQUIRED.** Match a choice in this section with the options in the following back cushion configuration section.

Redesigned, Static Type Back..... \$735
 Part: XRF000014
 Select a back cushion configuration in the subsequent section.

Back Cane Option
 Part: XRF000015
 Available only with TRU-Balance 3 Power Tilt (1 actuator) selections.
 Select cane type, back cane height, and back cushion configuration below.

Redesigned, Flip Forward Back..... \$735
 Part: XRF000016
 Not available with power recline systems. Only available with 2-post, flip back or pediatric armrests. Not available with 12"-14" width. Select a back cushion type in the subsequent section.

12.13. Back Cane Selection (Back Cane Option Only)**REQUIRED when Back Cane Option is selected above.** Select type and height of back canes.**1. Back Cane Type**

Angle Adjustable Back Cane..... No Charge
 Part: MEC140266

Flip Forward Back Cane No Charge
 Part: XRF000273

2" Endomorph (1" per side) Back Cane No Charge
 Part: XRF000269

4" Endomorph (2" per side) Back Cane..... No Charge
 Part: XRF000271

2. Back Cane Height

16"H Cane
 Part: FRMASMB15834

17"H Cane
 Part: FRMASMB15835

20"H Cane
 Part: FRMASMB15838

21"H Cane
 Part: FRMASMB15839

24"H Cane
 Part: FRMASMB15842

25"H Cane
 Part: FRMASMB15843

TRU Balance 4 Seating section continued on following page.

12. TRU BALANCE 4 SEATING CONTINUED**12.14. Back Cushion Configuration**

REQUIRED. Select a configuration type and applicable dimensions. See the Sport, Q-Back, and TRU-Comfort back availability chart for stock order back size ranges. Other size backs may be available with special order and extended lead time.

1a. Sport Back/Omit/After-Market Interface

Sport Back No Charge
Part: XRF000017

Complete following width package and size.

Omit Back Pan/Cushion No Charge
Part: PTOINDV1287

Select a following width package.

After-market Back Interface No Charge
Part: ST-SSAPT4-Q. Select a set width in following section, range of 12"-21". Overall seat depth may need to be adjusted to accommodate certain backs. May interfere with gel/waterfall armpads and arm through. Does not include brackets to mount after-market backs.

1b. TRU-Comfort Back

TRU-Comfort Back (Static/Flip Forward), Up to 15" Width \$835
Part: SET123622. HCPCS: E2620

Complete following width package and size.

TRU-Comfort Back (Static/Flip Forward), 16" to 20" Width ... \$1,035
Part: TC4V10W21W. HCPCS: E2620

Complete following width package and size.

TRU-Comfort Back (Static/Flip Forward), 22" Width \$1,235
Part: TC4V22W32W. HCPCS: E2621

Complete following width package and size.

TRU-Comfort Back, Back Cane Option \$835
Part: TC4V10W21WBC. HCPCS: E2620

Complete following width package and size.

1c. Q-Back 2.0, AL Series, Active Contour

Select a back type and the specific overall back height in subsequent section.

Active Contour, 14"Wx19-20"H Tilt/22"-23" H Recline \$9,85
Part: QB2AL14AS. HCPCS: E2613

Active Contour, 14"Wx21-22"H Tilt/24"-25" H Recline \$985
Part: QB2AL14AT. HCPCS: E2613

Active Contour, 16"Wx19-20"H Tilt/22"-23" H Recline \$1,145
Part: QB2AL16AS. HCPCS: E2613

Active Contour, 16"Wx21-22"H Tilt/24"-25" H Recline \$1,145
Part: QB2AL16AT. HCPCS: E2613

Active Contour, 16"Wx23-24"H Tilt/26"-27" H Recline \$1,145
Part: QB2AL16AX. HCPCS: E2613

Active Contour, 16"Wx25-26"H Tilt/28"-29" H Recline \$1,145
Part: QB2AL16AXX. HCPCS: E2613

Active Contour, 18"Wx19-20"H Tilt/22"-23" H Recline \$1,145
Part: QB2AL18AS. HCPCS: E2613

Active Contour, 18"Wx21-22"H Tilt/24"-25" H Recline \$1,145
Part: QB2AL18AT. HCPCS: E2613

Active Contour, 18"Wx23-24"H Tilt/26"-27" H Recline \$1,145
Part: QB2AL18AX. HCPCS: E2613

Active Contour, 18"Wx25-26"H Tilt/28"-29" H Recline \$1,145
Part: QB2AL18AXX. HCPCS: E2613

Active Contour, 20"Wx19-20"H Tilt/22"-23" H Recline \$1,145
Part: QB2AL20AS. HCPCS: E2613

Active Contour, 20"Wx21-22"H Tilt/24"-25" H Recline \$1,145
Part: QB2AL20AT. HCPCS: E2613

Active Contour, 20"Wx23-24"H Tilt/26"-27" H Recline \$1,145
Part: QB2AL20AX. HCPCS: E2613

Active Contour, 20"Wx25-26"H Tilt/28"-29" H Recline \$1,145
Part: QB2AL20AXX. HCPCS: E2613

1d. Q-Back 2.0, AL Series, Deep Contour

Select a back type and the specific overall back height in subsequent section.

Deep Contour, 14"Wx19-20"H Tilt/22"-23" H Recline \$1,145
Part: QB2AL14DS. HCPCS: E2620

Deep Contour, 14"Wx21-22"H Tilt/24"-25" H Recline \$1,145
Part: QB2AL14DT. HCPCS: E2620

Deep Contour, 16"Wx19-20"H Tilt/22"-23" H Recline \$1,325
Part: QB2AL16DS. HCPCS: E2620

Deep Contour, 16"Wx21-22"H Tilt/24"-25" H Recline \$1,325
Part: QB2AL16DT. HCPCS: E2620

Deep Contour, 16"Wx23-24"H Tilt/26"-27" H Recline \$1,325
Part: QB2AL16DX. HCPCS: E2620

Deep Contour, 16"Wx25-26"H Tilt/28"-29" H Recline \$1,325
Part: QB2AL16DXX. HCPCS: E2620

Deep Contour, 18"Wx19-20"H Tilt/22"-23" H Recline \$1,325
Part: QB2AL18DS. HCPCS: E2620

Deep Contour, 18"Wx21-22"H Tilt/24"-25" H Recline \$1,325
Part: QB2AL18DT. HCPCS: E2620

Deep Contour, 18"Wx23-24"H Tilt/26"-27" H Recline \$1,325
Part: QB2AL18DX. HCPCS: E2620

Deep Contour, 18"Wx25-26"H Tilt/28"-29" H Recline \$1,325
Part: QB2AL18DXX. HCPCS: E2620

Deep Contour, 20"Wx19-20"H Tilt/22"-23" H Recline \$1,325
Part: QB2AL20DS. HCPCS: E2620

Deep Contour, 20"Wx21-22"H Tilt/24"-25" H Recline \$1,325
Part: QB2AL20DT. HCPCS: E2620

Deep Contour, 20"Wx23-24"H Tilt/26"-27" H Recline \$1,325
Part: QB2AL20DX. HCPCS: E2620

Deep Contour, 20"Wx25-26"H Tilt/28"-29" H Recline \$1,325
Part: QB2AL20DXX. HCPCS: E2620

1e. Q-Back 2.0, AL Series, Extra Deep Contour

Select a back type and the specific overall back height in subsequent section.

Extra Deep Contour, 14"Wx19-20"H Tilt/22"-23" H Recline \$1,225
Part: QB2AL14ES. HCPCS: E2620

Extra Deep Contour, 14"Wx21-22"H Tilt/24"-25" H Recline \$1,225
Part: QB2AL14ET. HCPCS: E2620

Extra Deep Contour, 16"Wx19-20"H Tilt/22"-23" H Recline \$1,410
Part: QB2AL16ES. HCPCS: E2620

Extra Deep Contour, 16"Wx21-22"H Tilt/24"-25" H Recline \$1,410
Part: QB2AL16ET. HCPCS: E2620

Extra Deep Contour, 16"Wx23-24"H Tilt/26"-27" H Recline \$1,410
Part: QB2AL16EX. HCPCS: E2620

Extra Deep Contour, 16"Wx25-26"H Tilt/28"-29" H Recline \$1,410
Part: QB2AL16EXX. HCPCS: E2620

Extra Deep Contour, 18"Wx19-20"H Tilt/22"-23" H Recline \$1,410

Part: QB2AL18ES. HCPCS: E2620

Extra Deep Contour, 18"Wx21-22"H Tilt/24"-25" H Recline \$1,410
Part: QB2AL18ET. HCPCS: E2620

Extra Deep Contour, 18"Wx23-24"H Tilt/26"-27" H Recline \$1,410
Part: QB2AL18EX. HCPCS: E2620

Extra Deep Contour, 18"Wx25-26"H Tilt/28"-29" H Recline \$1,410
Part: QB2AL18EXX. HCPCS: E2620

Extra Deep Contour, 20"Wx19-20"H Tilt/22"-23" H Recline \$1,410
Part: QB2AL20ES. HCPCS: E2620

Extra Deep Contour, 20"Wx21-22"H Tilt/24"-25" H Recline \$1,410
Part: QB2AL20ET. HCPCS: E2620

Extra Deep Contour, 20"Wx23-24"H Tilt/26"-27" H Recline \$1,410
Part: QB2AL20EX. HCPCS: E2620

TRU Balance 4 Seating section continued on following page.

13. TRU-BALANCE 3 SEATING CONTINUED

1c. ADI2.0 Backrest

Specify back type, width, height, and contour in the part number configurator below. Previously selected back cane height must be equal to the greater of overall height range. **Available SOF for Provider installation only.**

For more information please see [ADI2.0 Backrest order form](#)

ADI2.0 Backrest 10"W - 20"W \$1,495
 For use with back canes only.
 Indicate part number below.
 ADP: SE0002152

ADI2.0 Backrest 22"W - 24"W \$1,695
 For use with back canes only.
 Indicate part number below.
 ADP: SE0002153

Specify Part Number for Back Selection

OPTIONS INDEX	ADI2	Back Type: AL AFT	Width: 10"-20" (1" increments), 22", 24"	Height: N = Nino L = Low M = Medium S = Standard T = Tall X = Extra Tall	Contour: A = Active D = Deep E = Extra Deep	2 = FRG Foam	P = Quick-Release
ENTER SELECTION	ADI2					2	P

BACK CONTOUR*:



BACK SERIES TYPES*:

(AL) Aluminum Series

(*ALL BACKS SHOWN ARE 18" WIDE)

(AFT) Aluminum Flat Top Series



BACK HEIGHTS AND WIDTH:

(Available for AL Back Shell Only)

Nino: 8" inch (10"W - 14"W)

Low: 10 inch (10"W - 24"W)

Medium: 13 inch (12"W - 24"W)

(Available for AL and AFT Shell Types)

Standard: 16 inch (13"W - 24"W)

Tall: 18 inch (16"W - 24"W)

Extra Tall: 20 inch (16"W - 24"W)

12. TRU BALANCE 4 SEATING CONTINUED**12.15. TRU Balance 4 Back Shroud**

REQUIRED. The back shroud is a feature of the TRU Balance 4 seating system and is available in matching shroud colors. The back shroud is not available when Omit Back Pan/Cushion is selected. Reference colors shown on page 2.

Red Corvette
Part: XRF000019

Deep Purple
Part: XRF000025

Midnight Blue
Part: XRF000029

Green Machine
Part: XRF000033

Back in Black
Part: XRF000021

Raspberry Beret
Part: XRF000026

Electric Blue
Part: XRF000030

Military Green
Part: XRF000034

Night Train
Part: XRF000022

Pretty in Pink
Part: XRF000027

Blue Jean
Part: XRF000031

Faded Glory
Part: DGN1714085-WFL2

Grey Street
Part: XRF000023

Evolution Orange
Part: XRF000028

Lemon Crush
Part: XRF000032

Digital Camo
Part: DGN1714085-WCM7

12.16. Lateral Positioning Wedge

Optional. Available with Sport and TRU-Comfort backs. Not available on power recline back 19"-20"H with limited adjustability on 21"H.

Lateral Positioning Wedge \$16
Part: SET123465 x2

12.17. Back Width Package

REQUIRED except when Omit Back Pan/Cushion is selected. Available with +/-4" of seat width.

Match Seat Width Selections
Part: PTOINDV3468

14"W-16"W Package
Part: XRF000037

16"W-20"W Package
Part: XRF000038

18"W-22"W Package
Part: XRF000039

12.18. Back Width

REQUIRED except when Omit Back Pan/Cushion is selected. Available with +/-4" of seat width.

12" Back Width
Part: PTOINDV3562

16" Back Width
Part: PTOINDV3564

18" Back Width
Part: PTOINDV3566

20" Back Width
Part: PTOINDV3568

22" Back Width
Part: PTOINDV3570

14" Back Width
Part: PTOINDV3563

17" Back Width
Part: PTOINDV3565

19" Back Width
Part: PTOINDV3567

21" Back Width
Part: PTOINDV3569

	Tilt Only Width											Recline Width										
	12"	14"	15"	16"	17"	18"	19"	20"	21"	22"		12"	14"	15"	16"	17"	18"	19"	20"	21"	22"	
14"	S	T	S	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	14"
15"	S	T	S	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	15"
16"	S	S	T	S	T	S	T	S	T	S	T	S	T	S	T	S	T	S	T	S	T	16"
17"	S	S	T	S	T	S	T	S	T	S	T	S	T	S	T	S	T	S	T	S	T	17"
18"	-	S	T	S	T	S	T	S	T	S	T	S	T	S	T	S	T	S	T	S	T	18"
19"	-	S	Q	T	S	Q	T	S	Q	T	S	S	S	T	S	T	S	T	S	T	S	19"
20"	-	S	Q	T	S	Q	T	S	Q	T	S	S	S	T	S	T	S	T	S	T	S	20"
21"	-	S	Q	T	S	Q	T	S	Q	T	S	-	S	T	S	T	S	T	S	T	S	21"
22"	-	S	Q	T	S	Q	T	S	Q	T	S	-	S	Q	T	S	Q	T	S	Q	T	22"
23"	-	S	T	S	T	S	Q	T	S	Q	T	-	S	Q	T	S	Q	T	S	Q	T	23"
24"	-	S	T	S	T	S	Q	T	S	Q	T	-	S	Q	T	S	Q	T	S	Q	T	24"
25"	-	S	T	S	Q	T	S	Q	T	S	S	-	S	Q	T	S	Q	T	S	Q	T	25"
26"	-	S	T	S	Q	T	S	Q	T	S	S	-	S	T	S	Q	T	S	Q	T	S	26"
27"	-	S	T	S	T	S	T	S	T	S	S	-	S	T	S	Q	T	S	Q	T	S	27"
28"	-	S	T	S	T	S	T	S	T	S	S	-	S	T	S	Q	T	S	Q	T	S	28"
29"	-	S	T	S	T	S	T	S	T	S	S	-	S	T	S	Q	T	S	Q	T	S	29"
30"	-	-	-	-	-	-	-	-	-	-	-	-	S	T	S	T	S	T	S	T	S	30"
31"	-	-	-	-	-	-	-	-	-	-	-	-	S	T	S	T	S	T	S	T	S	31"
32"	-	-	-	-	-	-	-	-	-	-	-	-	S	T	S	T	S	T	S	T	S	32"

S Sport Back **Q** Q-Back **T** TRU-Comfort Back

Sport, Q-Back, and TRU-Comfort Back Availability Chart

12.19. Overall Back Height

Power Recline available height range: 19"-32"H. Odd heights include 6" gap and even heights include 7" gap between seat pan and bottom of back cushion. See the availability chart for stock order back size ranges.

19" Height
Part: SETHGHT1032

22" Height
Part: SETHGHT1035

25" Height
Part: SETHGHT1038

28" Height
Part: SETHGHT1041

31" Height
Part: SETHGHT1102

20" Height
Part: SETHGHT1033

23" Height
Part: SETHGHT1036

26" Height
Part: SETHGHT1039

29" Height
Part: SETHGHT1042

32" Height
Part: SETHGHT1103

21" Height
Part: SETHGHT1034

24" Height
Part: SETHGHT1037

27" Height
Part: SETHGHT1040

30" Height
Part: SETHGHT1043

13. ARMRESTS**REQUIRED.** Select the desired type of armrest below. Please read notes to ensure seating compatibility. Be sure to also select an armrest size and armpad type.**13.1a. Flip-Back, 8"-16" H Armrests**

Choose a left and right armrest and set height in the following section. Not available with waterfall armpads. Not available with Back Cane Option.

☐ L	☐ R	Flip-Back, 8"-16"H Armrests \$145 ea.
		Left part: ARM1806662. Right part: ARM1806663

Armrest Support**Optional.** Available up to 2" difference between seat and back width via endomorph or mesomorph configuration. Not available less than 16" width seat or back.

☐ L	☐ R	Flip-Back Armrest Downpost \$115 ea.
		Left part: XRF000265. Right part: XRF000266

13.1b. Two-Post, Height Adjustable, Removable Armrests

Select the armrest and set the left and right heights in the following section. These armrests may interfere with offset backs/ back canes.

1a. Two-Post Armrest for Power Seating

☐ L	☐ R	Full Length, Two-Post Armrest for Power Seating \$145 ea.	☐ L	☐ R	Desk Length, Two-Post Armrest for Power Seating \$145 ea.
		Left part: XRF000047. Right part: XRF000048			Left part: XRF000049. Right part: XRF000050

13.1c. Pediatric, Removable, Height Adjustable ArmrestsHeight adjustment from 7.75" to 9.75" in .5" increments. **Not available with full length armpads.**

☐ L	☐ R	Pediatric, Removable Armrests \$145 ea.
		Left Part: ARM142876. Right part: ARM142877

13.2. Armrest Height**REQUIRED.** Set left and right armrest height. Note the available height range in the option description.**1a. Flip-Back and Two-Post Left Armrest Height**

8" H PTOINDV3535	10" H PTOINDV3537	12" H PTOINDV3539	14" H PTOINDV3541	16" H PTOINDV3543
9" H PTOINDV3536	11" H PTOINDV3538	13" H PTOINDV3540	15" H PTOINDV3542	

1b. Pediatric, Left Armrest Height

7.75" H PTOINDV3722	8.25" H PTOINDV3724	8.75" H PTOINDV3726	9.25" H PTOINDV3728	9.75" H PTOINDV3730
------------------------	------------------------	------------------------	------------------------	------------------------

2a. Flip-Back and Two-Post Right Armrest Height

8" H PTOINDV3450	10" H PTOINDV3452	12" H PTOINDV3454	14" H PTOINDV3456	16" H PTOINDV3458
9" H PTOINDV3451	11" H PTOINDV3453	13" H PTOINDV3455	15" H PTOINDV3457	

2b. Pediatric, Right Armrest Height

7.75" H PTOINDV3721	8.25" H PTOINDV3723	8.75" H PTOINDV3725	9.25" H PTOINDV3727	9.75" H PTOINDV3729
------------------------	------------------------	------------------------	------------------------	------------------------

13.3. Armrest Size**REQUIRED.** Select a left (L) and right (R) armrest size. This also affects armpad selection.

☐ L	☐ R	10" Length L part: XRF000057. R part: XRF000123	☐ L	☐ R	14" Length L part: XRF000058. R part: XRF000124	☐ L	☐ R	17" Length L part: XRF000059. R part: XRF000125 Only available with Flip-Back, 8"-16" H Armrests and straight arm pads.
----------------------------	----------------------------	--	----------------------------	----------------------------	--	----------------------------	----------------------------	---

14. ARMPADS**REQUIRED.****14.1a. Omit Armpad**Select one below if aftermarket armpads or handpads are required. See [Positioning Components Order Form](#).

No Armpad Selection No Charge Part: XRF000060 Not available with Flip-back, 8"-16" H armrests.	Interface for Aftermarket Pad No Charge Part: XRF000308 Only available with Flip-back, 8"-16" armrests. Not available with chest bar.
--	---

14.1b. Foam Armpads

☐ L	☐ R	Straight, Armpad Standard Left part: XRF000061. Right part: XRF000126	☐ L	☐ R	Waterfall, Armpad No Charge Left part: XRF000062. Right part: XRF000127 Not available with Flip-back, 8"-16" armrests.
----------------------------	----------------------------	--	----------------------------	----------------------------	--

14.1c. Gel Armpads

Sold only in pairs. Armpads are 3" wide. Manufacturer part number for full length is 314G and desk length is 310G.

a. Falcon Gel Armpads for Flip-Back Armrests

Falcon Gel, Full Length Left and Right \$300 Part: ARM1806807	Falcon Gel, Desk Length Left, Full Length Right \$250 Part: ARM1806810
Falcon Gel, Full Length Left, Desk Length Right \$250 Part: ARM1806813	Falcon Gel, Desk Length Left and Right \$250 Part: ARM1806806

Armpads section continued on following page.

Account Number: _____ Date: _____ Marked for: _____

14. ARMPADS CONTINUED

b. Falcon Gel Armpads for Tube Style Armrests

Falcon Gel, Full Length Left and Right \$300
 Part: ARM136438
 Falcon Gel, Full Length Left, Desk Length Right..... \$250
 Part: ARM136442

Falcon Gel, Desk Length Left, Full Length Right..... \$250
 Part: ARM136439
 Falcon Gel, Desk Length Left and Right \$200
 Part: ARM136437

15. ARMREST BRACKETS

REQUIRED.

Standard Fixed Receiver Standard
 Part: XRF000063

Stealth Height Adjustable Receiver..... \$194
 Part: XRF000064

16. HEADRESTS

Optional. Select a headrest pad and bracket.

16.1a. Headrests, Comfort Plus with Cool Core

1. Headrest Pads

Stealth 8" Comfort Plus, Cool Core (CP08CCFPM) \$279
 Part: ST-CP08CCFPM-Q. HCPCS: E0955
 Stealth 10" Comfort Plus, Cool Core (CP10CCFPM)..... \$279
 Part: ST-CP10CCFPM-Q. HCPCS: E0955

Stealth 14" Comfort Plus, Cool Core (CP14CCFPM) \$279
 Part: ST-CP14CCFPM-Q. HCPCS: E0955
 Stealth 6", Comfort Plus, Cool Core (CP06CCFPM)..... \$279
 Part: ST-CP06CCFPM-Q. HCPCS: E0955

2. Headrest Brackets

Comfort Plus Headrest Mounting Bracket, UniLink..... \$485
 Part: ST-UHB2D33200-Q. HCPCS: E1028

Comfort Plus Headrest Mounting Bracket, HMO475..... \$271
 Part: ST-HMO475-17M-Q.

Q Cam Lock..... No Charge
 Part: XRF000376

Q Cam Lock Availability (Overall height)	TRU-Comfort/Sport		Q-Back 2	
	Standard	\$50	Standard	\$50
Tilt/Static/Back Canes	14"-21"H	22"-29"H	19"-24"H	25"-26"H
Power Recline Back	19"-21"H	22"-32"H	22"-25"H	26"-29"H

17. LAP BELTS

REQUIRED. Select a lap belt below.

17.1a. Standard Lap Belts

☐ 38" Lap Belt Part: ACCBELT1006 ☐ 50" Lap Belt Part: ACCBELT1000 ☐ 60" Lap Belt Part: ACCBELT1002 ☐ 70" Lap Belt Part: ACCBELT1003 ☐ 80" Lap Belt Part: ACCBELT1004 ☐ 90" Lap Belt Part: ACCBELT1001

17.1b. Optional Lap Belts

Stealth Positioning, 60" Single Pull, Push Button, Medium..... \$134
 (7.5") Padded, 1.5" Belt (IFPB5213)
 Part: ST-IFPB5213-Q

Stealth Positioning, 60" Single Pull, Push Button, \$134
 Extra Large, Padded (11.5"), 2" Belt (IFPB5315)
 Part: ST-IFPB5315-Q

Stealth Positioning, 60" Single Pull, Push Button, Large \$134
 Padded (9.5"), 2" Belt (IFPB5314)
 Part: ST-IFPB5314-Q

60" Sure-Lok Lap Belt (Non-retractable)..... \$200
 Part: ACCBELT1016
 Required with ISO 7176-19 Compliant Transit Option.

18. LEGRESTS

All efforts will be made to accommodate legrest length selection, but minor changes may need to be made to meet shipping guidelines. Some adjustments may be needed upon receipt of the unit.

18.1a. Basic Center Mount Foot Platform

1. Basic Center Mount Foot Platform Assembly

Basic Center Mount..... \$200
 Part: XRF000323

2. Basic Center Mount Extension and Angle Adjustable Footplates

Short Extension, 5"x9" Footplates
 Part: XRF000312. HCPCS: K0040

Long Extension, 5"x9" Footplates
 Part: XRF000314. HCPCS: K0040

Medium Extension, 5"x9" Footplates
 Part: XRF000313. HCPCS: K0040

3. Basic Center Mount Calf Pads

Small Calf Pads, Basic Center Mount..... \$210
 Part: RIG2104369

Large Calf Pads, Basic Center Mount \$210
 Part: RIG2202997
 Only available with Long Extension.

4. Basic Center Mount Additional Accessory

Basic Center Mount Footplate Connector \$55
 Part: RIG2104240

Short	Medium	Long
8"-10"	10"-14"	14"-19"
Basic Center Mount Extension Range		

Legrests section continued on following page.

18. LEGRESTS CONTINUED**18.1b. Center Mount Foot Platform****1. Foot Platform Assembly**

Dual Actuator, Power Articulating Foot Platform.....\$4,295
 Part: RIG2009246. HCPCS: E1012+K0040
 Positional feedback with independent control of articulation.

Quick Adjustable Center Mount Foot Platform..... \$265
 Part: XRF000311. HCPCS: K0040
 Not available with power lift.

2. Lower Extension and Angle Adjustable Footplates

Medium length extension range is 11"-17" and long extension range is 14"-20". See compatibility matrix for possible interference.

Pediatric Extension, 5"x8" Footplates
 Part: RIG2404015

Medium Extension, 5"x11" Footplates
 Part: RIG2404018

Long Extension, 5"x11" Footplates
 Part: RIG2404019

Medium Extension, 5"x6" Footplates
 Part: RIG150206

Long Extension, 5"x6" Footplates
 Part: RIG151505

Medium Extension, 5"x8" Footplates
 Part: RIG2404016

Long Extension, 5"x8" Footplates
 Part: RIG2404017

Legrest Compatibility Articulating Foot Platform TRU Balance 4	Seat to Floor Option		
	Low	Middle	High
Medium Length Extension	11.5"-14.5"	11.5"-15.5"	11.5"-16.5"
Long Length Extension	N/A	14.5"-15.5"	14.5"-16.5"

Legrest Compatibility Quick Adjustable CM TRU Balance 4	Seat to Floor Option		
	Low	Middle	High
Pediatric Length Extension	8.5"-12.5"	8.5"-12.5"	8.5"-12.5"
Medium Length Extension	11"-15"	11"-16"	11"-17"
Long Length Extension	14"-15"	14"-16"	14"-17"

3a. Center Mount Calf Pads

Calf pads are standard with a power articulating foot platform.

Small Calf Pads, AFP Only, 6"Wx6"H (pair) \$210
 Part: RIGASMB7120018

Large Calf Pads, AFP Only, 8" taper to 6"Wx8"H (pair) \$210
 Part: RIGASMB7120019

4. Center Mount Modular Footplate Extension Components

Optional. Assembly required. May require seating configuration adjustments to avoid front caster interference.

2" Bell-Shaped Length Extension \$55
 Part: RIG152519
 Only available with the 5x6 footplates.

Outside Edge \$55
 Part: RIG2404921
 Available with 5"x8" footplates.

Width Extension for 8" \$55
 Part: RIG2404919

Front and Rear Cap \$55
 Part: RIG2404923
 Available with 5"x8" footplates.

Width Extension for 11" \$55
 Part: RIG2404920

5. Center Mount Additional Accessories

Small Reversed Dartex Calf Pad Cover (pair) \$150
 Part: RIG133717

Footplate Connector \$55
 Part: RIGASMB7120023

6. Stealth Gel Footplate Covers

5"x8" Tapered Gel Footplate Covers \$275
 Part: ST-FPG5008T-Q. HCPCS: K0108

5"x8" Tapered Gel Footplate Covers for Width Extension \$275
 Part: ST-AFPG5008E-Q. HCPCS: K0108

5"x11" Tapered Gel Footplate Covers \$275
 Part: ST-AFPG5011T-Q. HCPCS: K0108

5"x11" Tapered Gel Footplate Cover for Width Extension \$275
 Part: ST-AFPG5011E-Q. HCPCS: K0108

18.2. Legrest Accessories

Wheel Bumper \$15 ea.
 Part: XRF000100

18.3. Legrest Length

REQUIRED with matching left and right legrest configuration. Part number is PTOINDV plus number following # below. All efforts will be made to accommodate legrest length selection, but minor changes may need to be made to meet shipping guidelines. Some adjustments may be needed upon receipt of the unit. Reference the section's compatibility matrix for legrest interference. Non-matching or individual legrest configurations lengths are not set by the manufacturer.

11" L (#3478)	11.5" L (#3479)	12" L (#3480)	12.5" L (#3481)	13" L (#3482)	13.5" L (#3483)	14" L (#3484)	14.5" L (#3485)	15" L (#3486)	15.5" L (#3487)	16" L (#3488)	16.5" L (#3489)	17" L (#3490)
------------------	--------------------	------------------	--------------------	------------------	--------------------	------------------	--------------------	------------------	--------------------	------------------	--------------------	------------------

19. POWER POSITIONING ELECTRONICS**REQUIRED.****19.1a. Q-Logic 3 Controls through Joystick**

Single Actuator Through Q-Logic 3 Controller and Interface (AM3) No Charge

2 Actuators Through Q-Logic 3 Controller and Interface (AM3).....\$1,250

Not funded by ADP

3 Actuators Through Q-Logic 3 Controller and Interface (AM3)..... \$1,500

Not funded by ADP

Up to 5 Actuators Through Q-Logic 3 Controller and Interface (AAM).....\$2,175

*ADP: SECX15005***19.1b. iAccess Switch Control**

iAccess Flat Button.....No Charge

*Part: CTL144009**See next section for switch mounting.*

iAccess Toggle Button.....No Charge

*Part: CTL144010**See next section for switch mounting.***19.2. iAccess Mount****REQUIRED with iAccess options.****1. iAccess Mount Selection**

Joystick Mount No Charge

Part: MEC1907500

Attendant Mount \$35

Part: MEC1906564

Under Arm Mount, Flip-Back Armrest..... No Charge

Part: MEC1906531

Under Arm Mount, Flip-Back Armrest, plus iAccess Guard \$35

*Part: MEC2001392***2. Under Arm Mount Side**

Left Side Under Arm Mount

Part: PTOINDV3678

Right Side Under Arm Mount

*Part: PTOINDV3679***20. ATTENDANT CONTROL****Optional.** Attendant control will be at no charge when ordered with Tilt and Recline

Attendant Control \$950

*Part: XRF000078. HCPCS: E2331***21. SEAT MOUNTED LED LIGHTS****Optional.** Requires Q-Logic 3 Joystick for Full Light Package or Q-Logic 3e Joystick.

Seat Mounted, LED, Full Light Package \$1,045

*Part: XRF000130**Not available with Base Mounted, LED, Full Light Package.*

22. ACCESSORIES**Optional.** Multiple accessories can be mounted along the various tracks but some interference may occur.

Accessory Bar \$90
 Part: XRF000128

Horizontal rail attached to rear of seat. Standard with the Trilogy vent tray and oxygen holder.

Backpack Holder (Pair) (A)..... \$42
 Part: ACC147774

Cane/Crutch Holder \$195
 Part: ACC152515. HCPCS: E2207

Cell Phone/ Media Holder, X-Grip (E)..... \$95
 Part: ACC155835

Not available with Heavy Duty armrests.

Clothing Guard (Pair)..... \$275
 Part: ACC125007

Compact Vent Tray \$1,395
 Part: ACC133889. HCPCS: E1029

Cup Holder \$20
 Part: ACC125003

For use with track-style armrests.

Cup Holder (Tube-style Armrest)..... \$20
 Part: ACCASMB2634

For use with tube-style armrests.

Flag Kit \$125
 Part: ACC1700443

Flip-Back Armrest, Side Track Accessory Bracket .. \$25
 Part: MEC1809178

Glove Box (C)..... \$85
 Part: ACC125011

Grip Drink Holder (20085)..... \$100
 Part: ACC152993

Grip Lap Board (20005)..... \$100
 Part: ACC152992

Hydration System (B)..... \$328
 Part: ACC125010

Large, Black, Weather Cover \$240
 Part: ACCCOVR1021

Oxygen Holder..... \$110
 Part: ACC125006. HCPCS: E2208
 Includes accessory bar, a credit will be applied PTOIN-DV3769.

Personal Item Hook..... \$20
 Part: ACC125008

Personal Item Pouch \$175
 Part: ACC2101017

Privacy Bag..... \$80
 Part: ACC2005723

Push Handles (pair)..... \$140
 Part: ACC2005027

Not available with oxygen holder or vent tray.

Q-Logic EX/NE Programmer \$625
 Part: ELEASMB6976

Quantum Backpack \$200
 Part: ACC1905853

Rear View Mirror \$35
 Part: ACC155833

Reflector Kit \$35
 Part: ACC124670


Stealth Large Tablet Holder \$565
 Part: ACC141818

Transfer Handle (D)..... \$152.50 ea.
 Part: MEC2205680 x2. HCPCS: E0705

Depth Adjustable Vent Tray \$1,000
 Part: ACC2302500. HCPCS: E1029

Includes accessory bar, a credit will be applied PTOIN-DV3769.

Remote Stop Switch \$1,495
 Part: XRF000379

**Aftermarket Accessory Mounting Bracket**

Select a quantity up to 5.

Aftermarket Accessory Mounting Bracket, Slotted \$100 ea.
 Part: ACC2407046

QTY

Aftermarket Accessory Mounting Bracket, Non-slotted \$100 ea.
 Part: ACC2405496

QTY

24. POSITIONING COMPONENTS**Optional.** More options are available on the [Positioning Components Order Form](#).**24.1. Stealth Unilink Pelvic/Thigh Guide****1. UniLink Pelvic/Thigh Hardware**

^ Please verify the appropriate billing code depending upon the placement of the support pad.

_____ Qty UniLink, 5" Height Adjustable, \$405 ea.
Removable Mounting Hardware

Part: ST-UTTP5-Q. HCPCS: E1028/E1034^

_____ Qty UniLink, 7" Height Adjustable, \$423 ea.

Removable Mounting Hardware

Part: ST-UTTP7-Q. HCPCS: E1028/E1034^

_____ Qty UniLink Multi-Angle Thigh Guide, Left \$585 ea.

Part: ST-UTT8A11602-L-Q. HCPCS: E1028/E1034^

_____ Qty UniLink Multi-Angle Thigh Guide, Right \$585 ea.

Part: ST-UTT8A11602-R-Q. HCPCS: E1028/E1034^

_____ Qty UniLink Flip Down Hip/Thigh Support, Left \$483 ea.
Part: ST-UTS25100-L-Q. HCPCS: E1028/E1034^

_____ Qty UniLink Flip Down Hip/Thigh Support, Right \$483 ea.

Part: ST-UTS25100-R-Q. HCPCS: E1028/E1034^

_____ Qty Flip Down with LEFT Full Surface Contact \$459 ea.

Feature (TWBADD-L & LPHW407)

Part: POS140428. HCPCS: E1028/E1034^

_____ Qty Flip Down with RIGHT Full Surface Contact \$459 ea.

Feature (TWBADD-R & LPHW407)

Part: POS140429. HCPCS: E1028/E1034^

2. UniLink Pelvic/Thigh Pads

^ Please verify the appropriate billing code to be used with all third-party payors.

_____ Qty UniLink, 4" x 4" Pad \$164 ea.

Part: ST-UTP544056100-Q. HCPCS: E0953/E0956^

_____ Qty UniLink, 4" x 8" Pad \$164 ea.

Part: ST-UTP548056100-Q. HCPCS: E0953/E0956^

_____ Qty UniLink, 4" x 12" Pad \$164 ea.

Part: ST-UTP54D056100-Q. HCPCS: E0953/E0956^

3. UniLink Pelvic/Thigh Guide Installation

Install on Left

Part: PTOINDV3585

Install on Right

Part: PTOINDV3586

Do Not Install

Part: PTOINDV3719

24.2. Stealth Thoracic Laterals

Swing-away laterals may interfere with flip back armrests. These laterals are not available with the back cane option.

1. Bracket

Swing-away, UniLink, Redesign Back, pair (ULS4042000-Q)... \$1032

Part: ST-ULS4042000-Q. HCPCS: E1034

Only available with Redesigned back assemblies.

Swing-away, UniLink, 125° w/Track, pair (ULPK404108B00) ... \$966

Part: ST-ULPK404108B00-Q. HCPCS: E1034

Only available with Sport or TRU Comfort backs.

Swing-away, UniLink, Back Cane, pair (ULPK404104B00) \$966

Part: ST-ULPK404104B00-Q. HCPCS: E1034

Only available with back cane and omit back configurations.

Swing-away, UniLink, Q-Back 2, Left (ULPK404103L00) \$483

Part: ST-ULPK404103L00-Q. HCPCS: E1034

Only available with Active Contour Q-Backs.

Swing-away, UniLink, Q-Back 2, Right (ULPK404103R00) \$483

Part: ST-ULPK404103R00-Q. HCPCS: E1034

Only available with Active Contour Q-Backs.

Positioning components section continued on following page.

24. POSITIONING COMPONENTS CONTINUED**2. Pads**

Curved Pads, 3"x5", pair (TWBL3X5)..... \$164
 Part: INDPART2663. HCPCS: E0956

Curved Pads, 4"x6", pair (TWBL4X6)..... \$164
 Part: INDPART2664. HCPCS: E0956

3. Installation

Install Laterals
 Part: PTOINDV3499

Do Not Install
 Part: PTOINDV3500

24.3. Chest Bar Assembly

Chest support is required with anterior tilt 20 and 30 degree functions. Omitting a chest support selection requires the provider to facilitate aftermarket chest supports that best fits consumer needs.

Small Chest Bar, 14" Foam Pad (ST-CBATB4ATT-S-Q)..... \$300
 Part: ST-CBA-TB4-S. HCPCS: K0108
 Available with 14" to 20" back widths.

Medium Chest Bar, 16" Foam Pad (ST-CBATB4ATT-M-Q)..... \$300
 Part: ST-CBA-TB4-M. HCPCS: K0108
 Available with 17" to 22" back widths.

24.4. Stealth Positioning Chest Harness

Upper Classic, Medium (IFPS1103) \$221
 Part: ST-IFPS1103-Q. HCPCS: E0960

Upper Contour, Medium (IFPS2103)..... \$221
 Part: ST-IFPS2103-Q. HCPCS: E0960

Upper Classic, Large (IFPS1104) \$221
 Part: ST-IFPS1104-Q. HCPCS: E0960

Upper Contour, Large (IFPS2104) \$221
 Part: ST-IFPS2104-Q. HCPCS: E0960

24.5. Stealth Positioning Chest Straps

Chest Strap, Medium w/ Belt Anchor, \$140
 (IFCS4113)
 Part: ST-IFCS4113-Q. HCPCS: E0960

Chest Strap, Large w/ Belt Anchor, \$140
 (IFCS4114)
 Part: ST-IFCS4114-Q. HCPCS: E0960

24.6. Stealth Swing-Away Headrest Mounts

These swing-away headrest mounts are only available with the Stealth Comfort Series headrests.

1. Facial Lateral Hardware

UniLink Facial Lateral, Swing Away, Right.....\$496
 Part: ST-SAR2A-R-Q. HCPCS: E1028

UniLink Facial Lateral, Swing Away, Left\$496
 Part: ST-SAR2A-L-Q. HCPCS: E1028

UniLink Facial Lateral, Dual Swing Away \$765
 2 Upper & 1 Lower Right
 Part: ST-SAR4A-DBU-RL-Q. HCPCS: E1028

UniLink Facial Lateral, Dual Swing Away \$765
 2 Upper & 1 Lower Left
 Part: ST-SAR4A-DBU-LL-Q. HCPCS: E1028

UniLink Facial Lateral, Dual Swing Away \$764
 2 Upper
 Part: ST-SAR4A-DBU-Q. HCPCS: E1028

2. Facial Lateral Pads and Accessories

Egg Switch Mount (ESM) (A) \$70 ea.
 Part: INDPART2815

Figure 8 Bracket with 13" Rod..... \$366 ea.
 and Egg Switch Mount (B)
 Part: ST-F8675-13-Q

Facial Pad, 4"x5x2", Cool Core (900M-Q) \$96 ea.
 Part: ST-900M-Q. HCPCS: K0108

Facial Pad, 6"x2", Cool Core (920M-Q)..... \$96 ea.
 Part: ST-920M-Q. HCPCS: K0108

Swing-Away Joystick Mount, 18" Rod (C)..... \$100 ea.
 Part: ST-SUS18-Q
 Allows mounting of facial pads and attaches to the headrest




Contoured and Classic Chest Harness



Chest Strap



Left and Right Facial Lateral with Pads



Headrest Accessories

24.7. Stealth Elbow Support**1. Elbow Support Bracket**

Left, Swing-Away, Elbow Support \$279
 Bracket Tube Style Arms (ARMS-ES-L)
 Part: POS124064. HCPCS: E1028

Left, UniLink, Fixed, Elbow Support Bracket..... \$309
 for Flip-Back, 8"-16"H Armrests (UETB4021FL-Q)
 Part: ST-UETB4021FL-Q


Right, Swing-Away, Elbow Support..... \$279
 Bracket Tube Style Arms (ARMS-ES-R)
 Part: POS124066. HCPCS: E1028

Right, UniLink, Fixed, Elbow Support Bracket..... \$309
 for Flip-Back, 8"-16"H Armrests (UETB4021FR-Q)
 Part: ST-UETB4021FR-Q

2. Elbow Support Pad

3"Wx3"H Elbow Pad (ARMS-ES-3X3)\$57 ea.
 Part: POS124065. HCPCS: K0108

4"Wx4"H Elbow Pad (ARMS-ES-4X4).....\$57 ea.
 Part: POS124067. HCPCS: K0108



Elbow Supports

24.8. Stealth Foot Boxes

Stealth Small Foot Boxes, 5"Wx8"Dx7"H..... \$725
 (SPLE107) (pair)
 Part: INDPART2786. HCPCS: E0954

Stealth Medium Foot Boxes, 6"Wx10"Dx9"H \$725
 (SPLE108) (pair)
 Part: INDPART2787. HCPCS: E0954

Stealth Large Foot Boxes, 7"Wx12"Dx11"H..... \$725
 (SPLE109) (pair)
 Part: INDPART2788. HCPCS: E0954

24.9. Stealth Positioning Shoe Holder

Small Shoe Holder, 3/4" Webbing (IFSH4211)..... \$93 ea.
 Part: ST-IFSH4211-Q. HCPCS: E0951+E0952

Medium Shoe Holder, 3/4" Webbing (IFSH4212) \$93 ea.
 Part: ST-IFSH4212-Q. HCPCS: E0951+E0952

Large Shoe Holder, 3/4" Webbing (IFSH4213) \$93 ea.
 Part: ST-IFSH4213-Q. HCPCS: E0951+E0952

Positioning components section continued on following page.

25. LAP TRAY**25.1. Quantum Lap Tray**

Mounting bracket is required with lap tray selection. These trays are compatible with Flip-Back armrests.

1a. Black Lap Tray

Black Lap Tray, 11" W \$360

Part: MEC1905328-11B. HCPCS: E0950

For use with 12"-13" back width.

Black Lap Tray, 13" W \$360

Part: MEC1905328-13B. HCPCS: E0950

For use with 14"-15" back width.

Black Lap Tray, 15" W \$360

Part: MEC1905328-15B. HCPCS: E0950

For use with 16"-17" back width.

Black Lap Tray, 17" W \$360

Part: MEC1905328-17B. HCPCS: E0950

For use with 18"-19" back width.

Black Lap Tray, 19" W \$360

Part: MEC1905328-19B. HCPCS: E0950

For use with 20"-22" back width.

1b. Clear Lap Tray

Clear Lap Tray, 11" W \$360

Part: MEC1905328-11C. HCPCS: E0950

For use with 12"-13" back width.

Clear Lap Tray, 13" W \$360

Part: MEC1905328-13C. HCPCS: E0950

For use with 14"-15" back width.

Clear Lap Tray, 15" W \$360

Part: MEC1905328-15C. HCPCS: E0950

For use with 16"-17" back width.

Clear Lap Tray, 17" W \$360

Part: MEC1905328-17C. HCPCS: E0950

For use with 18"-19" back width.

Clear Lap Tray, 19" W \$360

Part: MEC1905328-19C. HCPCS: E0950

For use with 20"-22" back width.

2. Mounting Hardware

Mounting Bracket, Chest Bar/Lap Tray No Charge

Part: MEC2104886

25.2. Simple Solutions Lap Tray

Mounting bracket is required with lap tray selection. See chart below for compatibility.

1. Lap Tray

Simple Solutions Small Tray with Cut Out \$750

Part: ST-SSLXS010001171-Q. HCPCS: E0950

Simple Solutions Large Tray with Cut Out \$750

Part: ST-SSLXL010001171-Q. HCPCS: E0950

Simple Solutions Medium Tray with Cut Out \$750

Part: ST-SSLXM010001171-Q. HCPCS: E0950

2. Mounting Hardware

Swing Away Mounting Hardware No Charge

Part: ST-SSLXHD10000060-Q

Simple Solutions Lap Tray Armrest Compatibility	Width											
	10	11	12	13	14	16	17	18	19	20	21	22
Flip-Back	S	S	S	S	S	S	M	ML	ML	ML	L	L
Flip-Back, Cane Mounted	S	S	S	S	S	S	M	ML	ML	L	L	
Drop-In, T-Arm	S	S	S	S	S	S	M	ML	ML	L	L	
Two-Post, Flip-Back	S	S	S	S	S	S	M	ML	ML	L	L	
Cantilever	S	S	S	S	S	M	M	ML	ML	L	L	
Q-Captain, Flip-Up						S		ML		L		
Single Post/ Heavy Duty/ Bariatric	S	S	S	S	S	S	S	ML	ML	ML	L	L

S:ST-SSLXS010001171-Q M:ST-SSLXM010001171-Q L:ST-SSLXL010001171-Q

S:Compatible

M:Armrest configuration may not hold tray position

*Some armrests listed may not be available on this order form.

26. ALTERNATIVE CONTROLS AND ELECTRONICS**Optional.** Select a type of specialty controls below. More options are available on the [Alternative Controls and Electronics order form](#).**26.1. Enhanced Display**

Q-Logic 3 Enhanced Display Kit..... \$2,045

Part: CTL167726

Select a color option below. Includes enhanced display, flexible gooseneck and mounting bracket. Features: 3.5" LCD display; up to 5 profiles to control drive, seat, and auxiliary functions; Bluetooth connectivity for mobile devices and mouse emulation; integrated single switch scanner; IR controls.

Q-Logic 3 Enhanced Display Shroud

Required with Q-Logic 3 Enhanced Display selection above. Reference colors shown on page 2.

Black, Soft Touch Standard

Part: DGN137728

Deep Purple

Part: DGN137729-PPR5

Electric Blue

Part: DGN137729-PBU9

Military Digital Camo..... \$50

Part: DGN137729-WCM7

Red Corvette

Part: DGN137729-PRD1

Raspberry Beret

Part: DGN137729-PPK4

Blue Jean

Part: DGN137729-PBU11

Faded Glory \$50

Part: DGN137729-WFL2

Back in Black

Part: DGN137729-PBK1

Pretty in Pink

Part: DGN137729-PPK3

Lemon Crush

Part: DGN137729-PYL4

Carbon Fiber \$50

Part: DGN137729-WCB1

Night Train

Part: DGN137729-PBK3

Evolution Orange

Part: DGN137729-POR7

Green Machine

Part: DGN137729-PGN11

Grey Street

Part: DGN137729-PSL3

Midnight Blue

Part: DGN137729-PBU10

Military Green

Part: DGN137729-PGN8



Enhanced Display

26.2. SCIM Module**Optional.**

SCIM Module for Q-Logic 3 \$850

Part: CTL2001963

Allows use of joystick display for specialty devices. Available with Q-Logic 3 EX Joystick selection above. Up to 5 profiles to control drive, seat, and auxiliary functions, excluding sip-n-puff and single switch functions. A switch is required to operate the input device attached to the SCIM module.



SCIM Module

26.3. Q-Logic 3 SAJ

Q-Logic 3 EX SAJ \$850

Part: CTL143980

Must be used in conjunction with Enhanced Display or Q-Logic 3 Joystick.



EX SAJ

26.4. Head Array

i-Drive 4.0 Bluetooth Pro Series, \$7,406

3 Switch with Egg Switch

Part: ST-IDHBT300-U-Q. HCPCS: E2330



Adult Pro-Series

26.5. Sip and Puff Drive**1. Sip and Puff Module**

Sip and Puff Module for Q-Logic 3 \$1,695

Part: CTL168056. HCPCS: E2325

Includes module and mounting hardware.

2. Sip and Puff Straw and Mount

The straw and mount do not include electrical components to operate the Sip&Puff with a head array. This component requires the selection of a Q-Logic Sip&Puff module or additional module from Stealth Products for operation.

Stealth Sip and Puff for Comfort Series, \$496

Left Mount (SP625-CL)

Part: ST-SP625-CL-Q. HCPCS: E2326

Includes straw and mounting hardware. Use with Stealth Comfort Plus headrests.

Stealth Sip and Puff for Stealth i-Drive, \$496

Head Array (SP625-U)

Part: ST-SP625-U-Q. HCPCS: E2326

Includes straw and mounting hardware. Use with Stealth i-Drive Head Arrays.

Stealth Sip and Puff for Comfort Series, \$496

Right Mount (SP622-CR)

Part: ST-SP625-CR-Q. HCPCS: E2326

Includes straw and mounting hardware. Use with Stealth Comfort Plus headrests.



Sip & Puff Module



Stealth Sip & Puff

3. Sip & Puff Accessories

Sip & Puff Straws, 12 pack (SP625-S12)..... \$35

Part: ST-SP625-S12-Q

26.6. Egg Switches

Red Egg Switch (ESRED)..... \$116

Part: SWTMCRO1077

Blue Egg Switch (ESBLU)..... \$116

Part: SWTMCRO1079

Black Egg Switch (ESBLK)..... \$116

Part: SWTMCRO1081

Green Egg Switch (ESGRN) \$116

Part: SWTMCRO1078

Yellow Egg Switch (ESYEL) \$116

Part: SWTMCRO1080



Egg Switch

Alternative controls and electronics section continued on following page.

26. ALTERNATIVE CONTROLS AND ELECTRONICS CONTINUED

26.7. Switch Adapter

2 to 1 Switch Adapter \$135
Part: DWR1361H037
Connects 2 single switches to a stereo jack to act as a dual switch.

26.8. Environmental Controls

Q-Logic 3 ECU Module Kit \$1,780
Part: CTL168057

26.9. USB Mobile Device Charger

USB Mobile Device Charger \$185
Part: ELE168069
Can charge 1 device. A multiplier harness is required and must be added when a Static Seat or Power Tilt Only seating system is chosen. The selection of more than 3 modules requires an additional multiplier.

